



ANTECH




2024

DIRECTORY OF SERVICES

[ANTECHDIAGNOSTICS.COM](https://www.antechdiagnostics.com)

UNITED STATES



BETTER DIAGNOSTICS. BETTER MEDICINE.

You deserve a choice when making diagnostic decisions for your patients.

Antech is your diagnostics partner with North America's largest reference laboratory network, best-in-class in-house diagnostics from Heska, the industry's most-trusted imaging equipment from Sound, and breakthrough telemedicine from Antech Imaging Services.

**The right diagnostics at the right time
— all from Antech.**



NORTH AMERICA'S LARGEST REFERENCE LAB NETWORK

- Core diagnostics, including hematology, chemistry, parasitology, microbiology, and clinical and anatomic pathology
- Novel innovations like KeyScreen® GI Parasite PCR — the most advanced GI parasite screening solution in veterinary medicine
- Medical consulting services provided by board-certified specialists across all key specialties



HESKA®

IN-HOUSE DIAGNOSTICS BY HESKA, AN ANTECH COMPANY

- Most comprehensive in-clinic diagnostics portfolio on the market
- Innovative, accurate, and cost-effective core diagnostics including hematology, chemistry, immunodiagnostics, and ova and parasites
- Emergency and critical care diagnostics including blood, gases, and acid-base
- Best-in-class serological allergy assessment
- Rapid infectious disease testing with expanding menu (heartworm, FIV/FelV)



HESKA'S IN-HOUSE DIAGNOSTIC SUITE, NOW OFFERED BY ANTECH



DIAGNOSTIC IMAGING AND EDUCATION FROM SOUND

- Digital radiology for superior image quality and diagnostic capabilities
- A comprehensive, accessible ultrasound portfolio
- Advanced therapeutic laser solutions for impactful treatment
- Efficient mobile and fixed CT solutions
- Innovative PACS to view your images confidently from any modality
- Gold standard for training and education



TELEMEDICINE AND AI-POWERED SOLUTIONS

- Expert interpretations of digital radiographs
- On-the-spot interpretations of ultrasounds
- Consultations by the world's largest team of board-certified veterinary radiologists
- AI-powered interpretations of routine digital radiographs in just minutes in many circumstances



IN-HOUSE DIAGNOSTICS BY HESKA, AN ANTECH COMPANY



The world's only fully automated
fecal and urine point-of-care lab

Element AIM™ works to perfectly prepare, automate, and take images of feces and urine. It finds important cells, bacteria, crystals, casts, and parasites (including roundworm, hookworm, tapeworm, whipworm, *Giardia*, and more) using a fast and sanitary load-and-go workflow. The Element AIM provides quality results and images to you and your clients — faster than ever.

Learn more

heska.com/element-aim

- Sealed test cassettes eliminate the mess and smell of handling slides and cover slips
- Superior lab-quality optics for crystal clear images in just minutes
- Digital images are easily stored in medical records and shared with clients and specialists
- Automation removes user error and wasted time at the microscope and during sample prep



DIAGNOSTIC IMAGING AND EDUCATION FROM SOUND



Sound SmartDR®
An unparalleled digital
radiography system



Sound HD Dental
Superior diagnostic
intraoral and
periodontal images



GE® Ultrasound
Helping veterinarians
incorporate ultrasound
into practice since 1996



**doctorVet Plus
Therapy Laser**
Everything you need
in a therapy laser and
nothing that you don't



CereTom® CT
A portable, wireless,
self-shielded CT
scanner to roll into
your surgical suite



OEC C-Arm
High-quality images from
a proven, dependable,
classic C-Arm



Sound Assurance®
Protecting your
investment and
your practice

ACADEMY *of* VETERINARY IMAGING

Education
Possess the training
to perform specific
ultrasound exams

Sound imaging equipment and services are not only high-tech, but they are also intuitive. We are here for you and your patients, from customer service and tech support to customer education.

Speak with a Sound Representative today

1-800-268-5354

soundvet.com/request-demo



TELEMEDICINE AND AI-POWERED SOLUTIONS

Antech Imaging Services (AIS®) is here to provide your clients and their pets with the best care possible by offering board-certified specialists in radiology, cardiology, internal medicine, avian/exotics, oncology, neurology, and surgery — all by using cutting-edge digital technology partnered with the best practices of traditional veterinary medicine.

AIS PennHIP — the most effective hip screening tool available for dogs.

All dogs can benefit from PennHIP testing. For dogs at risk of developing hip osteoarthritis (OA), early intervention can help prevent or lessen the severity of canine hip dysplasia (CHD).

It gives an estimate of the risk for painful osteoarthritis (OA) of canine hip dysplasia (CHD) later in life. With this information, preventive and palliative strategies can be recommended by the PennHIP-trained veterinarian.

AIS RapidRead™ — accurate radiology reports at the speed of light.

AIS Radiology RapidRead combines an ever-expanding, continuously trained neural network and a world-class team of radiologists and data scientists that delivers:

- 63+ findings and growing
- Contextual assessments for deeper insight so you can diagnose and treat your patient with confidence, sooner
- We continue to add new findings and will soon add ECG, Dental, CT, and MRI to our AIS RapidRead platform
- AIS online picture archiving and communication system (PACS) service is available 24 hours per day, 7 days a week, 365 days a year to securely store and digitally transmit electronic images and clinically relevant reports

Find a specialist that's right for you

1-877-727-6800

antechimaging.com/
[antechweb/contact](https://antechweb.com/contact)

Our goal is to build close working relationships between practitioners and specialists. AIS personalization allows you to pick a preferred specialist, or a group of specialists, based on your preferences.






COMPREHENSIVE GENETIC SCREENING TO IDENTIFY POTENTIAL HEALTH CONDITIONS EARLY — FOR DOGS AND CATS

Wisdom Panel™ Premium enables better risk management in pre-anesthetic and preventative care:

- Detect before clinical signs develop — and proactively counsel pet owners on risk, diet, exercise, and breeding
- Detects MDR1 Medication Sensitivity (cat and dog), von Willebrand's Disease Type 1 (dog), Intervertebral Disc Disease Risk (dog), Polycystic Kidney Disease (cat), Hypertrophic Cardiomyopathy (cat), and many more
- Canine panel includes 265+ disorder-associated genetic health variants
- Feline panel includes 45+ health conditions
- Backed by 35+ published papers
- Includes pet owner-friendly breed and trait identifications



Wisdom's veterinary consultants support you in interpreting genetic results so you can build proactive, tailored healthcare plans for your patients.

Canine
Wisdom
Panel™
Premium


\$14497
page 133



Feline
Wisdom
Panel™
Complete

\$14498
page 133





KEYSCREEN® GI PARASITE PCR

The new standard in
parasite diagnostics

T991

page 129

KeyScreen GI Parasite PCR brings the power of PCR to routine parasite screening. Use KeyScreen to find more GI parasites and treat them with greater precision, speed, and confidence — all at an affordable price.

- Screens for 20 intestinal parasites
- Detects benzimidazole treatment' resistance in hookworms
- Determines the zoonotic potential of *Giardia*

**START FINDING MORE GI
PARASITES, FASTER, TODAY**


antechdiagnostics.com/parasites



Better pet health in one test



Detects
benzimidazole-resistant canine
hookworm infections¹⁻⁵



KeyScreen GI Parasite
PCR found parasites in

1 in 4

samples⁴

Reference lab O&P
only found parasites in

1 in 10

samples⁷



**Treatment
resistant
Hookworms**
found in more dog breeds and
across the US and Canada¹⁻⁵



Giardia is NOT
typically zoonotic
less than **4%** are⁴



Parasites
are evolving
so is KeyScreen GI Parasite PCR^{2, 3, 5, 6}

1. Comparative Study of KeyScreen and traditional fecal flotation methods: <https://parasitesandvectors.biomedcentral.com/articles/10.1186/s13071-023-05904-z>
2. Leutenegger CM, et al. Emergence of Ancylostoma caninum parasites with the benzimidazole resistance F167Y polymorphism in the US dog population. Int. J. Parasitol. Drugs Drug Resist. 2023;14:131-140. <https://www.sciencedirect.com/science/article/pii/S2211320723000015?via%3Dihub>
3. Evason, MD, et al. Emergence of canine hookworm treatment resistance: Novel detection of Ancylostoma caninum anthelmintic resistance markers by fecal PCR in 11 dogs from Canada, Am J Vet Res. 2023 July: <https://doi.org/10.2460/ajvr.23.05.0116>
4. Leutenegger CM, et al. Frequency of intestinal parasites in dogs and cats identified by molecular diagnostics. ACVIM, Philadelphia, June 2023.
5. Leutenegger CM, et al. Association of the novel benzimidazole resistance marker Q134H with F167Y in dogs with Ancylostoma caninum. ACVIM, Philadelphia June 2023.
6. Venkatesan A, et al. Molecular evidence of widespread benzimidazole drug resistance in Ancylostoma caninum from domestic dogs throughout the USA and discovery of a novel -tubulin benzimidazole resistance mutation. PLoS Pathog. Mar 2023;19:e1011146. <https://journals.plos.org/plospathogens/article?id=10.1371/journal.ppat.1011146>
7. KeyScreen Whitepaper <https://www.antechediagnostics.com/keyscreen/#keyscreen-whitepaper>

ACCUPLEX® WITH C6

Screen for canine vector-borne
disease with confidence

AC100
page 20

Accuplex, now with C6, tests for Lyme disease, heartworm, *Ehrlichia*, and *Anaplasma*.

- Lyme disease detection with Accuplex includes screening for antibodies to C6, and is in alignment with the 2018 ACVIM consensus statement for Lyme borreliosis.*
- The C6 peptide is only expressed following natural infection or exposure of a host. Therefore, a positive C6 result will not be reported from cross-reactivity with vaccinal antibodies.

**EXPERIENCE IT NOW AND MOVE YOUR CANINE
VECTOR-BORNE DISEASE SCREENING FROM
IN-HOUSE TO OUR HOUSE**

antechdiagnostics.com/accuplex

*Meryl P Littman, Bernhard Gerber, Richard E Goldstein, Mary Anna Labato, Michael R Lappin, George E Moore. ACVIM consensus update on Lyme borreliosis in dogs and cats. J Vet Intern Med. 2018 May;32(3):887-903. <https://onlinelibrary.wiley.com/doi/full/10.1111/jvim.15085>



WELCOME TO

ANTECH'S DIRECTORY OF SERVICES



Quickly find the test you need, when you need it

- See the Table of Contents to explore all that is available at your fingertips.
- Head to our Test Index ([page 185](#)) to find the exact location of the test or profile you need.
- Find our frequently used phone numbers in the Support Services section ([page 14](#)).
- The Quick Reference section ([page 16](#)) contains all the information you need for scheduling pickups and ordering supplies.

Reference Laboratories Terms of Service

Unless otherwise agreed in an agreement signed by an authorized representative of Antech, all reference laboratory testing and consulting services set out in this directory are governed by our standard terms and conditions, available at antechdiagnostics.com/terms-service.

Privacy Policy

As part of Mars Petcare, we are committed to A Better World for Pets.™ Antech respects and is committed to protecting your privacy. Please visit mars.com/privacy to see a copy of our privacy statement and how we handle the privacy of your data.

NOTE

The tests and services listed in this directory have been organized in a manner determined to be either appropriate medically or for ease of use by veterinary practitioners in the estimation of our medical operations teams. Any headings or groupings of tests or services are for reference purposes only. Additionally, all turnaround times and prices set forth in this directory are subject to change at our sole discretion. Please contact Customer Service at 1-800-872-1001 if you wish to confirm any information in this directory.

Table of Contents

14	Support Services
15	Account Information
16	Quick Reference
18	Consulting Services
19	Essential Panels
20	Preventative Care
41	CBC Chemistry Profiles
73	Hematology
79	Chemistry
87	Urine
91	Endocrinology
99	Infectious
109	Microbiology
117	Pathology
123	Molecular
135	General
151	Avian • Exotics
161	Equine • Large Animal
185	Test Index

NOTES | FREQUENTLY USED TEST CODES

SUPPORT SERVICES

Customer Service

1-800-872-1001, Dial 0

M – F	5 a.m. – 9 p.m. PST	8 a.m. – Midnight EST
SAT	5 a.m. – 5 p.m. PST	8 a.m. – 8 p.m. EST
SUN	6 a.m. – 3 p.m. PST	9 a.m. – 6 p.m. EST

Dispatch/Specimen Pickup

1-800-872-1001, Dial 1

AUTOMATED 24/7		
M – F	5 a.m. – 9 p.m. PST	8 a.m. – Midnight EST
SAT	5 a.m. – 5 p.m. PST	8 a.m. – 8 p.m. EST
SUN	6 a.m. – 3 p.m. PST	9 a.m. – 6 p.m. EST

Consultation Services

1-800-872-1001, Dial 2

M – F	5 a.m. – 6 p.m. PST	8 a.m. – 9 p.m. EST
SAT	6 a.m. – 3 p.m. PST	9 a.m. – 6 p.m. EST
SUN	CLOSED	

Billing

1-800-872-1001, Dial 3

M – F	6 a.m. – 5 p.m. PST	9 a.m. – 8 p.m. EST
SAT – SUN	CLOSED	

Equine Support

1-800-872-1001, Dial 4

M – F	5 a.m. – 9 p.m. PST	8 a.m. – Midnight EST
SAT	5 a.m. – 5 p.m. PST	8 a.m. – 8 p.m. EST
SUN	6 a.m. – 3 p.m. PST	9 a.m. – 6 p.m. EST

Antech Imaging Services

1-800-872-1001, Dial 5

AUTOMATED	24/7
------------------	------

Website Support

1-800-872-1001, Dial 6

support@antechdiagnostics.com

M – F	6 a.m. – 4 p.m. PST	9 a.m. – 7 p.m. EST
SAT	6 a.m. – Noon PST	9 a.m. – 3 p.m. EST
SUN	CLOSED	

Sound

1-800-819-5538

M – F	8 a.m. – 5 p.m. PST	11 a.m. – 8 p.m. EST
SAT – SUN	CLOSED	

ACCOUNT INFORMATION

ANTECH ACCOUNT NUMBER

HEALTHTRACKS NUMBER

USERNAME

PASSWORD

HealthTracks®
online.antechdiagnostics.com

HealthTracks provides you access to your diagnostic needs in one innovative digital platform. This digital center simplifies your practice’s diagnostic workflow, ensuring you have the features you need to provide the best care for your patients.

Elevate your diagnostic experience and explore all the HealthTracks features at [antechdiagnostics.com/healthtracks](https://online.antechdiagnostics.com/healthtracks)

TERRITORY SALES MANAGER

E-MAIL

MOBILE

INSIDE ACCOUNT REPRESENTATIVE

E-MAIL

MOBILE

PROFESSIONAL SERVICES VETERINARIAN

E-MAIL

MOBILE

SOUND EQUIPMENT SPECIALIST

E-MAIL

MOBILE

QUICK REFERENCE

Antech's dedicated Customer Service Team is here to help at **1-800-872-1001**

Completion Codes — Submitting a Sample for Profile Completion

Missing one of the sample types? No problem. If you do not have a necessary sample when submitting a profile, you may submit the omitted sample within 30 days of the original profile submission to complete the profile at no additional charge.

1. Create a new Test Requisition Form using the applicable completion code from the list below
2. On the Test Requisition Form, reference the original profile's 12-character accession number in the notes field
3. Print the Test Requisition Form and submit the sample as normal

Completion codes submitted after 30 days, not related to a prior profile submission, or without reference to the original accession number will be invoiced for the applicable test code.

- T991NC — KeyScreen® GI Parasite PCR Completion
- T805NC — Ova and Parasite Completion
- T808NC — Ova and Parasite with *Giardia* Completion
- T760NC — Urinalysis Completion
- T415NC — PT/PTT Completion

Follow-up and Monitoring Tests

We offer you reduced-charge testing for serial patient monitoring and rechecks. Within 30 days of the original submission, order one of the codes below and include the original profile's 12-character accession number on the Test Requisition Form.

- AC100R — Accuplex® Recheck
- RECHECK — Superchem with SDMA, CBC
- T330R — CBC Recheck
- T760R — Urinalysis Recheck
- M130R — Urine Culture Recheck
- M133R — Urinalysis and Urine Culture Recheck

Please note that recheck codes submitted after 30 days, not related to a prior profile submission, or without reference to the original accession number, will be charged the applicable non-recheck test fee.

Add-on Testing

Easily add testing to base profiles at a reduced charge. Add-ons must be ordered during the initial submission to receive discounted add-on pricing. After the original profile is submitted, additional tests will be charged the regular individual test fee.

Test Cancellation

Testing that has not begun can be canceled at no charge by calling Customer Service. Testing that is already in process cannot be canceled.

*Histology Note: all samples received for histology will be processed and evaluated. Histology tests may be canceled only prior to sample processing, and all samples in the accession must be canceled.

See the full reference lab terms of service at antechdiagnostics.com/terms-service

Specimen Retention Times

CLINICAL SUBMISSIONS	
Whole Blood	7 days
Serum/Plasma	7 days
Urine	7 days
EIA (Coggins): AGID and ELISA	30 days
Hematology – Slides (includes Wright’s stained urine smears)	7 days
Microbiology – all specimens	7 days
Microbiology – Slides (gram stain smears)	14 days
PATHOLOGY (CYTOLOGY) SUBMISSIONS	
Biological Specimens	7 days
Slides	90 days
PATHOLOGY (HISTOLOGY) SUBMISSIONS	
Biological Specimens (wet tissue in formalin)	30 days
Blocks	18 months
Slides	180 days

Supplies

Multiple ways to order supplies:

- Online — HealthTracks (preferred): online.antechdiagnostics.com
- Phone — **1-800-872-1001**
- Fax — **1-888-397-8356**

Supplies are typically delivered within 2 – 5 days of ordering.

Test Requisition Forms

Antech offers both electronic and paper test requisition forms (TRFs). Electronic TRFs can be generated directly from your practice management system, or through HealthTracks. Electronic TRFs are the most efficient and convenient way to order diagnostic testing. They have many benefits, including increased data and order integrity and the ability to archive requisitions. Manual TRFs are also available online at antechdiagnostics.com and can be ordered through HealthTracks.

CONSULTING SERVICES

Live, expert, and complimentary
peer-to-peer medical consultations

Looking for a trusted expert to help you interpret diagnostic results or develop a diagnostic and treatment plan? Antech provides current customers with comprehensive consulting services to assist with diagnostic interpretation, clinical application of results, and general case support. With a world-class team of 60+ board-certified specialists across 12 specialties, Antech is committed to partnerships you can count on.

Request A Consult

Please have your Antech account number and the accession number(s) available for the results you wish to discuss. Immediately connect with the first available Internal Medicine specialist consultant or request a call-back.

1-800-872-1001, Dial 2
online.antechdiagnostics.com



Internal Medicine specialists are
available Monday through Saturday

M – F	5 a.m. – 6 p.m. PST	8 a.m. – 9 p.m. EST
SAT	6 a.m. – 3 p.m. PST	9 a.m. – 6 p.m. EST

During business days, all other specialists are
available by call back within 24 – 48 hours



Avian and
Exotics



Behavior



Cardiology



Dentistry



Dermatology



Equine and
Large Animal



Neurology



Nutrition



Oncology



Reproduction

ESSENTIAL PANELS

	SUPERCHEM WITH SDMA SA010	VET SCREEN WITH SDMA SA025	ADULT WELLNESS CHEMISTRY WITH SDMA SA665	PRE-OP SCREEN WITH SDMA SA040
A/G Ratio	●	●		●
Albumin	●	●		●
Alkaline Phosphatase	●	●		●
ALT (SGPT)	●	●		●
Amylase	●			
AST (SGOT)	●	●		
BUN	●	●		●
BUN/Creatinine Ratio	●	●		●
Calcium	●	●		
Chloride	●	●		●
Cholesterol	●	●		
CPK	●	●		
Creatinine	●	●		●
GGT	●			
Globulin	●	●		●
Glucose	●	●		●
Magnesium	●			
NA/K Ratio	●	●		●
Pancreatic Sensitive Lipase (PSL)	●			
Phosphorus	●	●		
Potassium	●	●		●
SDMA	●	●		●
Sodium	●	●		●
Total Bilirubin	●	●		
Total Protein	●	●		●
Triglycerides	●			

PREVENTATIVE CARE

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
AC100	<p>Accuplex®</p> <p>This is a canine vector-borne disease screening for heartworm, Lyme disease (includes screening for antibodies against two C6 peptides), <i>E. canis</i>, and <i>A. phagocytophilum</i>.</p>	<p>0.5 mL serum in red top or serum separator tube</p>	<p>1-2 days</p>
KAC100	<p>Accuplex®, KeyScreen® GI Parasite PCR</p> <p>Accuplex® is a canine vector-borne disease screening for heartworm, Lyme disease (includes screening for antibodies against two C6 peptides), <i>E. canis</i>, and <i>A. phagocytophilum</i>. KeyScreen® is a fecal PCR test that detects 20 individual parasite species. If detected, <i>Toxocara</i> spp. will be further speciated (<i>T. cati</i>, <i>T. canis</i>, and <i>T. leonina</i>), <i>Ancylostoma caninum</i> will be evaluated for benzimidazole resistance, and <i>Giardia</i> evaluated for subtypes A and B.</p>	<p>0.5 mL serum and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, Antech provided fecal container</p>	<p>1-3 days</p>
KA535	<p>Adult Chem with Lytes, SDMA, CBC, Accuplex®, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, Accuplex®, and KeyScreen®</p> <p>Wellness chemistry with electrolytes, complete blood count, SDMA for glomerular filtration rate estimation (see T1035), Accuplex® for vector-borne disease screening (AC100), and KeyScreen® for intestinal parasite detection (see T991).</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	<p>1-3 days</p>
AC535S	<p>Adult Chem with Lytes, SDMA, CBC, O&P, Accuplex®</p> <p>Wellness Chemistry with Electrolytes, CBC, Fecal O&P, Accuplex®, and SDMA</p> <p>Wellness chemistry with electrolytes, SDMA for glomerular filtration rate estimation (see T1035), complete blood count, Accuplex® for vector-borne disease screening (see AC100), and fecal analysis using zinc sulfate with centrifugation/flotation for ova and parasite detection (T805).</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, and 5.0 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	<p>1-2 days</p>

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
KA583	Adult Chem with Lytes, SDMA, CBC, Accuplex®, T4, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Accuplex®, and KeyScreen® Wellness chemistry with electrolytes, SDMA (T1035), complete blood count, Accuplex® for vector-borne disease screening (see AC100), total T4, and KeyScreen® for intestinal parasite detection (see T991).	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
KS519	Adult Chem with Lytes, SDMA, CBC, Feline Heartworm, FIV, FeLV, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, FeLV, FIV, Heartworm Antibody, and KeyScreen® Wellness chemistry with electrolytes, SDMA (T1035), complete blood count, FeLV antigen detection, FIV and heartworm antibody detection, and KeyScreen® for intestinal parasite detection (see T991).	1.25 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
KS673	Adult Chem with Lytes, SDMA, CBC, T4, UA, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Urinalysis, and KeyScreen® A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), total T4, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035).	0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
KS575	Adult Chem with Lytes, SDMA, CBC, FIV, FeLV, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, FeLV, FIV, and KeyScreen® Wellness chemistry with electrolytes, SDMA (T1035), complete blood count, FeLV antigen detection, FIV antibody detection, and KeyScreen® for intestinal parasite detection (see T991).	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
KS534	Adult Chem with Lytes, SDMA, CBC, Heartworm, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, Heartworm Antigen detection, and KeyScreen® Wellness chemistry with electrolytes, SDMA, complete blood count, heartworm antigen detection, and KeyScreen® for intestinal parasite detection (see T991).	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
KS535	<p>Adult Chem with Lytes, SDMA, CBC, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, and KeyScreen®</p> <p>Wellness chemistry with electrolytes, SDMA, complete blood count, and KeyScreen® for intestinal parasite detection (see T991).</p>	<p>0.5 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	1-3 days
KS591	<p>Adult Chem with Lytes, SDMA, CBC, T4, Feline Heartworm, FIV, FeLV, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, T4, FeLV, FIV, Heartworm Antibody, and KeyScreen®</p> <p>Wellness chemistry with electrolytes, SDMA, complete blood count, total T4, FeLV antigen detection, FIV and heartworm antibody detection, and KeyScreen® for intestinal parasite detection (see T991).</p>	<p>1.25 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	1-3 days
KS597	<p>Adult Chem with Lytes, SDMA, CBC, T4, FIV, FeLV, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, T4, FeLV, FIV, and KeyScreen®</p> <p>Wellness chemistry with electrolytes, SDMA, complete blood count, total T4, FeLV antigen detection, FIV antibody detection, and KeyScreen® for intestinal parasite detection (see T991).</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	1-3 days
KS583	<p>Adult Chem with Lytes, SDMA, CBC, T4, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, T4, and KeyScreen®</p> <p>Wellness chemistry with electrolytes, SDMA, complete blood count, total T4, and KeyScreen® for intestinal parasite detection (see T991).</p>	<p>0.5 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	1-3 days
KS619	<p>Adult Chem with Lytes, SDMA, CBC, T4, UA, FIV, FeLV, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Urinalysis, FeLV, FIV, and KeyScreen®</p> <p>A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), total T4, FeLV antigen detection, FIV antibody detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035).</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
KA053	Adult Chem with Lytes, SDMA, CBC, UA, Accuplex®, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, Urinalysis, Accuplex®, and KeyScreen® A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), Accuplex® for vector-borne disease screening (see AC100), KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035).	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
KA673	Adult Chem with Lytes, SDMA, CBC, UA, Accuplex®, T4, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Urinalysis, Accuplex®, and KeyScreen® A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), total T4, Accuplex® for vector-borne disease screening (see AC100), KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035).	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
KS634	Adult Chem with Lytes, SDMA, CBC, UA, FIV, FeLV, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, Urinalysis, FeLV, FIV, and KeyScreen® A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), FeLV antigen detection, FIV antibody detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
KS592	Adult Chem with Lytes, SDMA, CBC, UA, Heartworm, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, Urinalysis, Heartworm, KeyScreen® A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), heartworm antigen detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
KS587	<p>Adult Chem with Lytes, SDMA, CBC, UA, Heartworm, T4, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Urinalysis, Heartworm Antigen test, and KeyScreen®</p> <p>A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), total T4, heartworm antigen detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035).</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days
KS053	<p>Adult Chem with Lytes, SDMA, CBC, UA, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, Urinalysis, and KeyScreen®</p> <p>A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035).</p> <p>Interferences: marked hemolysis and lipemia.</p>	<p>0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days
KS623	<p>Adult Chem with Lytes, SDMA, CBC, UA, Feline Heartworm, FIV, FeLV, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, Urinalysis, FeLV, FIV, Heartworm Antibody, and KeyScreen®</p> <p>A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), FeLV antigen detection, FIV and heartworm antibody detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035).</p>	<p>1.25 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days
KS589	<p>Adult Chem with Lytes, SDMA, CBC, T4, UA, Feline Heartworm, FIV, FeLV, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Urinalysis, FeLV, FIV and Heartworm Antibody detection, and KeyScreen®</p> <p>A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), total T4, FeLV antigen detection, FIV and heartworm antibody detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035).</p>	<p>1.25 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
KA600	Adult Chem with SDMA, Accuplex®, KeyScreen® Wellness Chemistry with SDMA, CBC, Accuplex®, and KeyScreen® Wellness chemistry, complete blood count, Accuplex® for canine vector-borne disease screening (see AC100), KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
AC600	Adult Chem with SDMA, CBC, Accuplex® Wellness Chemistry with SDMA, CBC, and Accuplex® Wellness chemistry, complete blood count, Accuplex® for canine vector-borne disease screening (see AC100), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis and lipemia.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
KA670	Adult Chem with SDMA, CBC, Accuplex®, T4, KeyScreen® Wellness Chemistry with SDMA, CBC, T4, Accuplex®, and KeyScreen® Wellness chemistry, complete blood count, total T4, Accuplex® for canine vector-borne disease screening (see AC100), KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
KS622	Adult Chem with SDMA, CBC, Feline Heartworm, FIV, FeLV, KeyScreen® Wellness Chemistry with SDMA, CBC, FeLV, FIV, Heartworm Antibody, and KeyScreen® Wellness chemistry, complete blood count, FeLV antigen detection, FIV and heartworm antibody detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.25 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
KS675	Adult Chem with SDMA, CBC, FIV, FeLV, KeyScreen® Wellness Chemistry with SDMA, CBC, FeLV, FIV, and KeyScreen® Wellness chemistry, complete blood count, FeLV antigen detection, FIV antibody detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
KS605	Adult Chem with SDMA, CBC, Heartworm, KeyScreen® Wellness Chemistry with SDMA, CBC, Heartworm Antigen, and KeyScreen® Wellness chemistry, complete blood count, heartworm antigen detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
KS685	Adult Chem with SDMA, CBC, Heartworm, T4, KeyScreen® Wellness Chemistry with SDMA, CBC, T4, Heartworm Antigen, and KeyScreen® Wellness chemistry, complete blood count, total T4, heartworm antigen detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
KS600	Adult Chem with SDMA, CBC, KeyScreen® Wellness Chemistry with SDMA, CBC, and KeyScreen® Wellness chemistry, complete blood count, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
AC615	Adult Chem with SDMA, CBC, O&P, Accuplex® Wellness Chemistry with SDMA, CBC, Fecal O&P with Centrifugation, Accuplex® Wellness chemistry, complete blood count, fecal analysis via zinc sulfate centrifugation/floatation for ova and parasite detection (T805), Accuplex® for canine vector disease borne screening (AC100), and SDMA for glomerular filtration rate estimation (see T1035).	1.0 mL serum, 1.0 mL EDTA whole blood, and 5 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days
AC655	Adult Chem with SDMA, CBC, O&P, Giardia, Accuplex® Wellness Chemistry with SDMA, CBC, Fecal O&P with Centrifugation, <i>Giardia</i> , Accuplex® Wellness chemistry, complete blood count, SDMA for glomerular filtration rate estimation (see T1035), fecal analysis via zinc sulfate centrifugation/floatation for ova and parasite detection (T805) and <i>Giardia</i> antigen detection via ELISA. Accuplex® can be reflexed for a nominal fee. Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
KS590	Adult Chem with SDMA, CBC, T4, Feline Heartworm, FIV, FeLV, KeyScreen® Wellness Chemistry with SDMA, CBC, T4, FeLV, FIV, Heartworm Antibody, and KeyScreen® Wellness chemistry, complete blood count, total T4, FeLV antigen detection, FIV and heartworm antibody detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.25 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
KS697	Adult Chem with SDMA, CBC, T4, FIV, FeLV, KeyScreen® Wellness Chemistry with SDMA, CBC, T4, FeLV, FIV, and KeyScreen® Wellness chemistry, complete blood count, total T4, FeLV antigen detection, FIV antibody detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
KS588	Adult Chem with SDMA, CBC, T4, UA, Feline Heartworm, FIV, FeLV, KeyScreen® Wellness Chemistry with SDMA, CBC, T4, Urinalysis, FeLV, FIV, Heartworm Antibody, and KeyScreen® A complete minimum database (chemistry, complete blood count, urinalysis), total T4, FeLV antigen detection, FIV and heartworm antibody detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.25 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
KS621	Adult Chem with SDMA, CBC, T4, UA, FIV, FeLV, KeyScreen® Wellness Chemistry with SDMA, CBC, T4, Urinalysis, FeLV, FIV, and KeyScreen® A complete minimum database (chemistry, complete blood count, urinalysis), total T4, FeLV antigen detection, FIV antibody detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
KS672	Adult Chem with SDMA, CBC, T4, UA, KeyScreen® Wellness Chemistry with SDMA, CBC, T4, Urinalysis, and KeyScreen® A complete minimum database (chemistry, complete blood count, urinalysis), total T4, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
AC607	Adult Chem with SDMA, CBC, UA, Accuplex® Wellness Chemistry with SDMA, CBC, Urinalysis, Accuplex® A complete minimum database (chemistry, complete blood count, urinalysis), Accuplex® for canine vector borne disease screening (AC100), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
KA607	Adult Chem with SDMA, CBC, UA, Accuplex®, KeyScreen® Wellness Chemistry with SDMA, CBC, Urinalysis, Accuplex®, and KeyScreen® A complete minimum database (chemistry, complete blood count, urinalysis), Accuplex® for canine vector borne disease screening (AC100), KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
KA672	Adult Chem with SDMA, CBC, UA, Accuplex®, T4, KeyScreen® Wellness Chemistry with SDMA, CBC, T4, Urinalysis, Accuplex®, and KeyScreen® A complete minimum database (chemistry, complete blood count, urinalysis), total T4, Accuplex® for canine vector-borne disease screening (see AC100), KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
KS624	<p>Adult Chem with SDMA, CBC, UA, Feline Heartworm, FIV, FeLV, KeyScreen®</p> <p>Wellness Chemistry with SDMA, CBC, Urinalysis, FeLV, FIV, Heartworm Antibody, and KeyScreen®</p> <p>A complete minimum database (chemistry, complete blood count, urinalysis), FeLV antigen detection, FIV and heartworm antibody detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035).</p> <p>Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.</p>	<p>1.25 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	1-3 days
KS631	<p>Adult Chem with SDMA, CBC, UA, FIV, FeLV, KeyScreen®</p> <p>Wellness Chemistry with SDMA, CBC, Urinalysis, FeLV, FIV, and KeyScreen®</p> <p>A complete minimum database (chemistry, complete blood count, urinalysis), FeLV antigen detection, FIV antibody detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035).</p> <p>Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days
KS625	<p>Adult Chem with SDMA, CBC, UA, Heartworm, KeyScreen®</p> <p>Wellness Chemistry with SDMA, CBC, Urinalysis, Heartworm Antibody, and KeyScreen®</p> <p>A complete minimum database (chemistry, complete blood count, urinalysis), heartworm antigen detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035).</p> <p>Interferences: marked hemolysis and lipemia.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days
KS687	<p>Adult Chem with SDMA, CBC, UA, Heartworm, T4, KeyScreen®</p> <p>Wellness Chemistry with SDMA, CBC, T4, Urinalysis, Heartworm Antibody, and KeyScreen®</p> <p>A complete minimum database (chemistry, complete blood count, urinalysis), total T4, heartworm antigen detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). KS625 to which T4 has been added.</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
KS607	Adult Chem with SDMA, CBC, UA, KeyScreen® Wellness Chemistry with SDMA, CBC, Urinalysis, and KeyScreen® A complete minimum database (chemistry, complete blood count, urinalysis), KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
AC652	Adult Chem with SDMA, CBC, UA, UPC with Accuplex® Wellness Chemistry with SDMA, CBC, Accuplex®, Urinalysis, Protein/Creatinine Ratio A complete minimum database (chemistry, complete blood count, urinalysis), Accuplex® for canine vector borne disease screening (AC100), SDMA for glomerular filtration rate estimation (see T1035), and a urine protein to creatinine ratio (T775). Interferences: marked hemolysis and lipemia.	1.0 mL EDTA whole blood, 1.0 mL serum, and 6.5 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
KS670	Adult Chem with SDMA, T4, KeyScreen® Wellness Chemistry with SDMA, CBC, T4, and KeyScreen® Wellness chemistry, complete blood count, total T4, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
AC664	Adult Chem with SDMA, Vaccine Titers, Fecal Combo, Accuplex® Wellness Chemistry with SDMA, CBC, Distemper Parvo Vaccinal Titer, O&P with Centrifugation, <i>Giardia</i> , Accuplex® Wellness chemistry, complete blood count (CBC), Accuplex® for canine vector borne disease screening (AC100), Distemper/Parvovirus vaccinal titer (T565), fecal analysis via zinc sulfate centrifugation/ floatation and <i>Giardia</i> Ag ELISA for ova and parasite detection (T808), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL EDTA whole blood, 1.0 mL serum, and 6 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
AC670	Adult Wellness with SDMA, CBC, T4, Accuplex® Wellness Chemistry with SDMA, CBC, Total T4, Accuplex® A comprehensive chemistry, complete blood count, total T4, SDMA for glomerular filtration rate estimation (see T1035), and Accuplex® for canine vector borne disease screening (AC100). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
AC608	Adult Wellness with SDMA, T4, UA, O&P, Giardia, Accuplex® Wellness Chemistry with SDMA, CBC, T4, Urinalysis, O&P with Centrifugation, <i>Giardia</i> , Accuplex® A comprehensive minimum database (wellness chemistry, complete blood count, urinalysis), total T4, Accuplex® for canine vector borne disease screening (AC100), SDMA for glomerular filtration rate estimation (see T1035), and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (T808). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 6 grams of feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
AC710	Basic Wellness Screen with SDMA, Accuplex® Superchem with SDMA, CBC, Total T4, Urinalysis, Accuplex® The most comprehensive chemistry, CBC, Total T4, urinalysis, and Accuplex®. SDMA for glomerular filtration rate estimation (see T1035). Interferences: gross hemolysis and lipemia. Lipemia can falsely decrease T4 results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
AC651	Canine Wellness with SDMA, UA, O&P, Giardia, Accuplex® Wellness Chemistry with SDMA, CBC, Fecal O&P with Centrifugation, <i>Giardia</i> , Urinalysis, Accuplex® A complete minimum database (wellness chemistry, complete blood count, and urinalysis), Accuplex® for canine vector-borne disease screening (see AC100), SDMA for glomerular filtration rate estimation (see T1035), and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (T808). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 6 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
AC770	<p>Superchem with SDMA, CBC, T4, UA, O&P, Giardia, Accuplex®</p> <p>Superchem with SDMA, CBC, Total T4, Urinalysis, Fecal O&P with Centrifugation, <i>Giardia</i>, Accuplex®</p> <p>The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count, and urinalysis), a total T4, Accuplex® for canine vector-borne disease screening (AC100), and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (T808). Accuplex® can be reflexed for a nominal fee.</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 6 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-2 days
AC712	<p>Superchem with SDMA, CBC, T4, UA, UPC, Accuplex®</p> <p>Superchem with SDMA, CBC, Total T4, Urinalysis, urine protein to creatinine ratio, Accuplex®</p> <p>The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), total T4, urine protein to creatinine ratio (T775), and Accuplex® for canine vector borne disease screening (AC100).</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube</p>	1-2 days
AC750	<p>Senior Profile 1 with SDMA, O&P, Accuplex®</p> <p>Superchem with SDMA, CBC, Total T4, Urinalysis, Fecal O&P with Centrifugation, Accuplex®</p> <p>The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, Accuplex® for canine vector-borne disease screening, and fecal analysis using zinc sulfate centrifugation/flotation for ova and parasite detection.</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-2 days
AC705	<p>Senior Profile 1 with SDMA, Accuplex®</p> <p>Superchem with SDMA, CBC, Total T4, Urinalysis, Accuplex®</p> <p>The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, and Accuplex® (canine vector-borne disease screening for heartworm, Lyme disease (includes screening for antibodies against two C6 peptides), <i>E. canis</i> and <i>A. phagocytophilum</i>)).</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube</p>	1-2 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
KA705	Senior Profile 1 with SDMA, Accuplex®, KeyScreen® Superchem with SDMA, CBC, T4, Urinalysis, Accuplex®, and KeyScreen® The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, Accuplex® for canine vector-borne disease screening (see AC100), and KeyScreen® for intestinal parasite detection (see T991). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
KS722	Senior Profile 1 with SDMA, Feline Heartworm, KeyScreen® Superchem with SDMA, CBC, Total T4, Urinalysis, Heartworm Antibody, KeyScreen® The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, heartworm antibody detection, and KeyScreen® for intestinal parasite detection (see T991). Interferences: marked hemolysis and lipemia. Lipemia can falsely decrease T4 results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
KS710	Senior Profile 1 with SDMA, Heartworm, KeyScreen® Superchem with SDMA, CBC, T4, Urinalysis, Heartworm Antigen, and KeyScreen® The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), total T4, heart worm antigen detection, and KeyScreen® for intestinal parasite detection (see T991). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
KS705	Senior Profile 1 with SDMA, KeyScreen® Superchem with SDMA, CBC, T4, Urinalysis, and KeyScreen® The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, and KeyScreen® for intestinal parasite detection (see T991). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
AC720	Senior Profile 2 with SDMA, Accuplex® Superchem with SDMA, CBC, Total T4, Urinalysis, Free T4 By Equilibrium Dialysis, and Accuplex® The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a thyroid panel (T4 and FT4 ED), and Accuplex® (AC100) for canine vector borne disease screening. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-3 days
KA720	Senior Profile 2 with SDMA, Accuplex®, KeyScreen® Superchem with SDMA, CBC, Total T4, Urinalysis, Free T4 by ED, Accuplex®, and KeyScreen® The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a thyroid panel (T4 and FT4 ED), Accuplex® for canine vector-borne disease screening (see AC100), and KeyScreen® for intestinal parasite detection (see T991). Interferences: marked hemolysis and lipemia. Lipemia can falsely decrease T4 results.	1.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
KS780	Senior Profile 2 with SDMA, FeLV, FIV, KeyScreen® Superchem with SDMA, CBC, Total T4, Urinalysis, Free T4 by ED, FeLV, FIV, and KeyScreen® The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a thyroid panel (T4 and FT4 ED), FeLV antigen and FIV antibody detection, and KeyScreen® for intestinal parasite detection (see T991). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
KS720	Senior Profile 2 with SDMA, KeyScreen® Superchem with SDMA, CBC, Total T4, Urinalysis, Free T4 by ED, and KeyScreen® The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a thyroid panel (T4 and FT4 ED), and KeyScreen® for intestinal parasite detection (see T991). Interferences: marked hemolysis and lipemia. Lipemia can falsely decrease T4 results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
AC775	<p>Senior Profile 2 with SDMA, O&P, Giardia, Accuplex®</p> <p>Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, O&P with Centrifugation, <i>Giardia</i>, and Accuplex®</p> <p>The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a thyroid panel (T4 and FT4 ED), Accuplex® (AC100) test for canine vector borne disease screening, and fecal analysis using zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection.</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>1.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 6 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days
AC727	<p>Senior Profile 1 with SDMA, Vaccine Titers, Fecal Combo, Accuplex®</p> <p>Superchem with SDMA, CBC, T4, Urinalysis, Distemper Parvo Vaccinal Titer, O&P with Centrifugation, <i>Giardia</i>, Accuplex®</p> <p>The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, Distemper/Parvovirus vaccinal titer (T565), Accuplex® for canine vector-borne disease screening (see AC100) and fecal analysis using zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA (T808) for ova and parasite detection.</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days
KA114	<p>Superchem with SDMA, CBC, Accuplex®, T4, KeyScreen®</p> <p>Superchem with SDMA, CBC, T4, Accuplex®, and KeyScreen®</p> <p>Most comprehensive chemistry profile (superchemistry with SDMA (SA010)), a complete blood count, total T4, and KeyScreen® for intestinal parasite detection (see T991).</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	1-3 days
KA020	<p>Superchem with SDMA, CBC, Accuplex®, KeyScreen®</p> <p>Superchem with SDMA, CBC, Accuplex®, and KeyScreen®</p> <p>Most comprehensive chemistry profile (superchemistry with SDMA (SA010)), a complete blood count, Accuplex® for canine vector-borne disease screening (see AC100), and KeyScreen® for intestinal parasite detection (see T991).</p> <p>Interferences: marked hemolysis and lipemia.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
KS737	Superchem with SDMA, CBC, Feline Heartworm, FIV, FeLV, KeyScreen® Superchem with SDMA, CBC, FeLV, FIV, Heartworm Antibody, and KeyScreen® Most comprehensive chemistry profile (superchemistry with SDMA (SA010)), a complete blood count, FeLV antigen detection, FIV and heartworm antibody detection, and KeyScreen® for intestinal parasite detection (see T991). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.25 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
KS490	Superchem with SDMA, CBC, FIV, FeLV, KeyScreen® Superchem with SDMA, CBC, FeLV, FIV, and KeyScreen® Most comprehensive chemistry profile (superchemistry with SDMA (SA010)), a complete blood count, FeLV antigen and FIV antibody detection, and KeyScreen® for intestinal parasite detection (see T991). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
KS028	Superchem with SDMA, CBC, Heartworm, KeyScreen® Superchem with SDMA, CBC, Heartworm Antigen, and KeyScreen® Most comprehensive chemistry profile (superchemistry with SDMA (SA010)), a complete blood count, heartworm antigen detection, and KeyScreen® for intestinal parasite detection (see T991). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
KS110	Superchem with SDMA, CBC, Heartworm, T4, KeyScreen® Superchem with SDMA, CBC, T4, Heartworm Antigen, and KeyScreen® Most comprehensive chemistry profile (superchemistry with SDMA (SA010)), a complete blood count, total T4, heartworm antigen detection, and KeyScreen® for intestinal parasite detection (see T991). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
KS020	Superchem with SDMA, CBC, KeyScreen® Superchem with SDMA, CBC, and KeyScreen® Most comprehensive chemistry profile (superchemistry with SDMA (SA010)), a complete blood count, and KeyScreen® for intestinal parasite detection (see T991). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
KS724	<p>Superchem with SDMA, CBC, T4, Feline Heartworm, FIV, FeLV, KeyScreen®</p> <p>Superchem with SDMA, CBC, T4, FeLV, FIV, Heartworm Antibody (Feline) and KeyScreen®</p> <p>Most comprehensive chemistry profile (superchemistry with SDMA (SA010)), a complete blood count, total T4, FeLV antigen detection, FIV and heartworm antibody detection, and KeyScreen® for intestinal parasite detection (see T991).</p> <p>Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.</p>	<p>1.25 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	1-3 days
KS220	<p>Superchem with SDMA, CBC, T4, FIV, FeLV, KeyScreen®</p> <p>Superchem with SDMA, CBC, T4, FeLV, FIV, and KeyScreen®</p> <p>Most comprehensive chemistry profile (superchemistry with SDMA (SA010)), a complete blood count, total T4, FeLV antigen and FIV antibody detection, and KeyScreen® for intestinal parasite detection (see T991).</p> <p>Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	1-3 days
KS120	<p>Superchem with SDMA, CBC, T4, KeyScreen®</p> <p>Superchem with SDMA, CBC, T4, and KeyScreen®</p> <p>Most comprehensive chemistry profile (superchemistry with SDMA (SA010)), a complete blood count, total T4, and KeyScreen® for intestinal parasite detection (see T991).</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>0.5 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	1-3 days
KS715	<p>Superchem with SDMA, CBC, T4, UA, Feline Heartworm, FIV, FeLV, KeyScreen®</p> <p>Superchem with SDMA, CBC, T4, Urinalysis, FeLV, FIV, Heartworm Antibody (Feline), and KeyScreen®</p> <p>The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), total T4, FeLV antigen detection, FIV and heartworm antibody detection, and KeyScreen® for intestinal parasite detection (see T991).</p> <p>Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.</p>	<p>1.25 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
KS700	Superchem with SDMA, CBC, T4, UA, FIV, FeLV, KeyScreen® Superchem with SDMA, CBC, T4, Urinalysis, FeLV, FIV, and KeyScreen® The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), total T4, FeLV antigen and FIV antibody detection, and KeyScreen® for intestinal parasite detection (see T991). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
KA021	Superchem with SDMA, CBC, UA, Accuplex®, KeyScreen® Superchem with SDMA, CBC, Urinalysis, Accuplex®, and KeyScreen® The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), Accuplex® for canine vector-borne disease screening (see AC100) and KeyScreen® for intestinal parasite detection (see T991). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
KS494	Superchem with SDMA, CBC, UA, Feline Heartworm, FIV, FeLV, KeyScreen® Superchem with SDMA, CBC, Urinalysis, FeLV, FIV, Heartworm Antibody, and KeyScreen® The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), FeLV antigen detection, FIV and heartworm antibody detection, and KeyScreen® for intestinal parasite detection (see T991). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.25 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
KS702	Superchem with SDMA, CBC, UA, FIV, FeLV, KeyScreen® Superchem with SDMA, CBC, Urinalysis, FeLV, FIV, and KeyScreen® The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), FeLV antigen and FIV antibody detection, and KeyScreen® for intestinal parasite detection (see T991). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
KS112	<p>Superchem with SDMA, CBC, UA, Heartworm, KeyScreen®</p> <p>Superchem with SDMA, CBC, Urinalysis, Heartworm Antigen, and KeyScreen®</p> <p>The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), heartworm antigen detection, and KeyScreen® for intestinal parasite detection (see T991).</p> <p>Interferences: marked hemolysis and lipemia.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days
KS021	<p>Superchem with SDMA, CBC, UA, KeyScreen®</p> <p>Superchem with SDMA, CBC, Urinalysis, and KeyScreen®</p> <p>The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis) and KeyScreen® for intestinal parasite detection (see T991).</p> <p>Interferences: marked hemolysis and lipemia.</p>	<p>0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days
AC114	<p>Total Body Function with SDMA, O&P, Accuplex®</p> <p>Superchem with SDMA, CBC, Total T4, Fecal O&P with Centrifugation, and Accuplex®</p> <p>The most comprehensive chemistry panel with SDMA (SA010), a complete blood count, total T4, Accuplex® for canine vector borne disease screening (AC100) and fecal analysis using zinc sulfate centrifugation/flotation (T805) for ova and parasite detection.</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, and 5.0 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	1-2 days
AC123	<p>Total Body Function with SDMA, O&P, Giardia, Accuplex®</p> <p>Superchem with SDMA, CBC, Total T4, Fecal O&P with Centrifugation and <i>Giardia</i>, and Accuplex®</p> <p>The most comprehensive chemistry with SDMA (SA010), a complete blood count, total T4, Accuplex® for canine vector borne disease screening (AC100) and fecal analysis using zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA (T808) for ova and parasite detection.</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
KS260	FeLV, FIV, KeyScreen® FeLV, FIV, and KeyScreen® FeLV antigen detection, FIV antibody detection, and KeyScreen® for intestinal parasite detection (see T991). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	0.5 mL serum and 0.3 grams feces Serum in red top or serum separator tube, Antech provided fecal container	1-3 days
KT615	Heartworm, KeyScreen® Heartworm Antigen and KeyScreen® Heartworm antigen detection and KeyScreen® for intestinal parasite detection (see T991). Interferences: marked hemolysis and lipemia.	0.5 mL serum and 0.3 grams feces Serum in red top or serum separator tube, Antech provided fecal container	1-3 days
AC805	O&P with Accuplex® Fecal O&P with Centrifugation, Accuplex® Fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection and an Accuplex® test for canine vector borne disease screening (AC100).	0.5 mL serum and 5 grams feces Serum in red top or serum separator tube, Antech provided fecal container	1-2 days
AC808	O&P, Giardia with Accuplex® Fecal O&P with Centrifugation, <i>Giardia</i> , Accuplex® Fecal analysis via zinc sulfate centrifugation/flotation and a <i>Giardia</i> Ag (ELISA) for ova and parasite detection and an Accuplex® for canine vector borne disease screening (AC100).	0.5 mL serum and 6 grams feces Serum in red top or serum separator, Antech provided fecal container	1-2 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CBC CHEMISTRY PROFILES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
KS630	Adult Chem with SDMA, CBC, Feline Heartworm, UA, KeyScreen® Wellness Chemistry with SDMA, CBC, Heartworm Antibody, Urinalysis, and KeyScreen® A complete minimum database (wellness chemistry, complete blood count, urinalysis), heartworm antibody detection, KeyScreen® for intestinal parasite detection (see T991), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
SA675	Adult Chemistry with SDMA, CBC, FeLV, FIV Wellness Chemistry with SDMA, CBC, FeLV, and FIV Wellness chemistry, complete blood count, FeLV antigen and FIV antibody detection, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis and lipemia.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
SA622	Adult Feline Wellness with SDMA, FeLV, FIV Wellness Chemistry with SDMA, CBC, Feline Heartworm Antibody, FeLV, FIV A comprehensive chemistry panel, a complete blood count, FeLV antigen, FIV and heartworm antibody detection, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
SA686	Adult Wellness Panel with SDMA, Heartworm, O&P, UA Wellness Chemistry with SDMA, CBC, Heartworm Antigen, Fecal O&P with Centrifugation, Urinalysis A comprehensive minimum database (comprehensive chemistry, complete blood count, and urinalysis), heartworm antigen detection, SDMA for glomerular filtration rate estimation (see T1035), and fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection (T805). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA600	Adult Wellness with SDMA Adult Wellness Chemistry with SDMA, CBC A comprehensive chemistry panel, complete blood count, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
SA620	Adult Wellness with SDMA, Ehrlichia, Lyme Wellness Chemistry with SDMA, CBC, <i>Ehrlichia canis</i> , Heartworm Antigen, Lyme Titer IgG A comprehensive chemistry panel, complete blood count, SDMA for glomerular filtration rate estimation (see T1035) and screening for vector borne disease (heartworm antigen detection, detection of antibodies (IFA methodology) suggesting exposure to <i>Ehrlichia canis</i> or <i>Borrelia burgdorferi</i> (natural or vaccinal)). Interferences: marked hemolysis and lipemia.	1.25 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
SA684	Adult Wellness with SDMA, FeLV, FIV, O&P, Giardia Wellness Chemistry with SDMA, CBC, FeLV, FIV, O&P with Centrifugation, <i>Giardia</i> A comprehensive chemistry, complete blood count, FeLV antigen detection, FIV antibody detection, SDMA for glomerular filtration rate estimation (see T1035), and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (T808). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, 6 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days
SA651	Adult Wellness with SDMA, Heartworm, UA, O&P, Giardia Wellness Chemistry with SDMA, CBC, Heartworm Antigen, Fecal O&P with Centrifugation, <i>Giardia</i> , Urinalysis A comprehensive minimum database (wellness chemistry, complete blood count, urinalysis), heartworm antigen detection, SDMA for glomerular filtration rate estimation (see T1035), and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (T808). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 6 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
SA615	Adult Wellness with SDMA, O&P Wellness Chemistry with SDMA, CBC, Fecal O&P with Centrifugation A comprehensive chemistry panel, complete blood count, SDMA for glomerular filtration rate estimation (see T1035), and fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection (T805). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 5 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA655	Adult Wellness with SDMA, O&P, <i>Giardia</i> Wellness Chemistry with SDMA, CBC, Fecal O&P with Centrifugation, <i>Giardia</i> A comprehensive chemistry, complete blood count, SDMA for glomerular filtration rate estimation (see T1035), and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (T808). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, 6 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days
SA670	Adult Wellness with SDMA, T4 Wellness Chemistry with SDMA, CBC, Total T4 A comprehensive chemistry, complete blood count, total T4, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
SA674	Adult Wellness with SDMA, T4, O&P Wellness Chemistry with SDMA, CBC, T4, O&P with Centrifugation A comprehensive chemistry, complete blood count, total T4, SDMA for glomerular filtration rate estimation (see T1035), and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> ELISA for ova and parasite detection (T808). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 5 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days
SA676	Adult Wellness with SDMA, T4, O&P, UA Wellness Chemistry with SDMA, CBC, T4, O&P with Centrifugation, Urinalysis A comprehensive minimum database (wellness chemistry, complete blood count and urinalysis), total T4, SDMA for glomerular filtration rate estimation (see T1035), and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> ELISA for ova and parasite detection (T808). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
SA672	Adult Wellness with SDMA, T4, UA Wellness Chemistry with SDMA, CBC, Total T4, Urinalysis A comprehensive minimum database (wellness chemistry, complete blood count, urinalysis), total T4, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	24 hours Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA625	Canine Adult Wellness with SDMA, UA Wellness Chemistry with SDMA, CBC, Heartworm Antigen, Urinalysis, A comprehensive minimum database (wellness chemistry, complete blood count, urinalysis), heartworm antigen detection, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
SA100	Canine Comprehensive with SDMA Superchem, CBC, Total T4, Free T4, T3, The most comprehensive chemistry, a complete blood count, thyroid evaluation (T3, T4, FT4), SDMA for glomerular filtration rate estimation (see T1035).	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-3 days
SA710	Canine Senior Profile with SDMA, Heartworm Superchem with SDMA, CBC, Total T4, Urinalysis, Heartworm Antigen The most comprehensive minimum database (superchemistry, complete blood count, and urinalysis), total T4, heartworm antigen detection, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
SA765	Canine Senior with SDMA, Fecal Combo Superchem with SDMA, CBC, Total T4, Urinalysis, Heartworm Antigen, Fecal O&P with Centrifugation, <i>Giardia</i> The most comprehensive minimum database (superchemistry, complete blood count, and urinalysis), total T4, heartworm antigen detection, SDMA for glomerular filtration rate estimation (see T1035), and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (T808). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 6 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
SA635	Canine Wellness with SDMA, O&P Wellness Chemistry with SDMA, CBC, Heartworm Antigen, Fecal O&P with Centrifugation A comprehensive chemistry, complete blood count, heartworm antigen detection, SDMA for glomerular filtration rate estimation (see T1035), and fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection (T805). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, and 5 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA220	Cat Scan Plus with SDMA Superchem with SDMA, CBC, Total T4, FeLV, and FIV The most comprehensive chemistry, complete blood count, total T4, FeLV antigen detection, FIV antibody detection and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
SA230	Cat Scan with SDMA Superchem with SDMA, CBC, T4, and FeLV The most comprehensive chemistry, complete blood count, total T4, FeLV antigen detection, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
SA1348	Cat Scan With SDMA, T4 Superchem with SDMA, CBC, T4, and FeLV A comprehensive chemistry (Vet Screen), complete blood count, total T4, FeLV antigen detection, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
SA610	Cat Wellness Profile with SDMA Wellness Chemistry with SDMA, CBC, and Feline Heartworm Antibody A comprehensive chemistry panel, a complete blood count, heartworm antibody detection, and SDMA for glomerular filtration rate estimation (see T1035).	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
SA132	Cocci Profile with SDMA, T4 Superchem with SDMA, CBC, T4, and Coccidioidomycosis, Screen and Titer The most comprehensive chemistry, complete blood count, total T4, SDMA for glomerular filtration rate estimation (see T1035), and <i>Coccidioides</i> antibody assessment by AGID (T535). Interferences: lipemia. T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL of serum and 1.0 mL of whole blood Serum in red top or serum separator tube, lavender top	1-5 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA105	D1 with SDMA, FT4-ED Superchem with SDMA, CBC, T3, Total T4, Free T4 by ED The most comprehensive chemistry panel, a complete blood count, a thyroid panel (T4, T3, FT4 ED), and SDMA for glomerular filtration rate estimation (see T1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-3 days
SA504	Desert Disease Panel with SDMA, T4 Superchem with SDMA, CBC, T4, Coccidioidomycosis, Screen and Titer, and <i>Ehrlichia canis</i> The most comprehensive chemistry, a complete blood count, total T4, SDMA for glomerular filtration rate estimation (see T1035), and screening for exposure to <i>Ehrlichia canis</i> (IFA) and Coccidioidomycosis (IgM and IgG). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-5 days
SA800	Diabetes Monitoring Panel with SDMA CBC, Urinalysis, Fructosamine Assay, Chemistry panel: Total Protein, ALT (SGPT), Alk Phos, Total Bilirubin, BUN, Creatinine, Glucose, PSL, and SDMA A minimum database (chemistry (ZW85710), complete blood count, urinalysis), SDMA for glomerular filtration rate estimation (see T1035), and fructosamine. Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	24 hours
SA605	Dog Wellness Profile with SDMA Wellness Chemistry with SDMA, CBC, Heartworm Antigen A comprehensive chemistry panel, a complete blood count, heartworm antigen detection, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis and lipemia.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
SA650	Dog Wellness Profile with SDMA, O&P, Giardia Wellness Chemistry with SDMA, CBC, Heartworm Antigen, Fecal O&P with Centrifugation, <i>Giardia</i> A comprehensive chemistry, complete blood count, heartworm antigen detection, SDMA for glomerular filtration rate estimation (see T1035), and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (T808). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 grams of feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA685	Dog Wellness with SDMA, T4 Wellness Chemistry with SDMA, CBC, Heartworm Antigen, Total T4 A comprehensive chemistry, complete blood count, total T4, heartworm antigen detection, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum 1.0 mL EDTA whole blood Serum in red top or serum separator, lavender top	1-2 days
SA640	Feline Adult Wellness Immune Profile with SDMA Superchem with SDMA, CBC, O&P with Centrifugation, Panleukopenia Vaccine Titer (Feline Only), Urinalysis A comprehensive minimum database (wellness chemistry, complete blood count, and urinalysis), SDMA for glomerular filtration rate estimation (see T1035), feline panleukopenia vaccinal titer, and fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection (T805). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-5 days
SA660	Feline Adult Wellness with SDMA, Fecal Combo Superchem with SDMA, CBC, Feline Heartworm Antibody, O&P with Centrifugation, <i>Giardia</i> A comprehensive chemistry, complete blood count, heartworm antibody detection, SDMA for glomerular filtration rate estimation (see T1035) and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (T808). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days
SA661	Feline Adult Wellness with SDMA, O&P, <i>Giardia</i>, UA Superchem with SDMA, CBC, Feline Heartworm Antibody, O&P with Centrifugation, <i>Giardia</i> , Urinalysis A comprehensive minimum database (wellness chemistry, complete blood count and urinalysis), heartworm antibody detection, SDMA for glomerular filtration rate estimation (see T1035) and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (T808). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA696	Feline Adult Wellness with SDMA, O&P, UA Wellness Chemistry with SDMA, CBC, Heartworm Antibody, Fecal O&P with Centrifugation, and Urinalysis A comprehensive minimum database (wellness chemistry, complete blood count and urinalysis), heartworm antibody detection, SDMA for glomerular filtration rate estimation (see T1035) and fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection (T805). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
SA629	Feline Adult Wellness with SDMA, O&P, UA, FeLV, FIV Wellness Chemistry with SDMA, CBC, Feline Heartworm Antibody, FeLV, FIV, Fecal O&P with Centrifugation, Urinalysis SA696 to which FeLV antigen and FIV antibody detection has been added. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	2.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
SA630	Feline Adult Wellness with SDMA, UA Wellness Chemistry with SDMA, CBC, Feline Heartworm Antibody, Urinalysis A comprehensive minimum database (wellness chemistry, complete blood count, urinalysis), heartworm antibody detection, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
SA624	Feline Adult Wellness with SDMA, UA, FeLV, FIV Wellness Chemistry with SDMA, CBC, FeLV, FIV, Urinalysis, and Heartworm Antibody (Feline) SA630 to which FeLV antigen and FIV antibody detection has been added. Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
SA718	Feline Comprehensive Wellness with SDMA, O&P Superchem with SDMA, CBC, Total T4, FIV Antibody, Urinalysis-Complete, FeLV Antigen ELISA, Heartworm Antibody, Feline, Fecal O&P The most comprehensive minimum database (superchemistry, complete blood count, urinalysis), a total T4, FeLV antigen detection, FIV and heartworm antibody detection, SDMA for glomerular filtration rate estimation (see T1035), and fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection (T805). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA200	Feline Comprehensive Plus with SDMA Superchem, CBC, Total T4, Free T4, T3, FeLV, FIV, Feline Coronavirus Titer Most comprehensive chemistry panel, a complete blood count, thyroid evaluation (T3, T4, FT4), FeLV antigen detection, FIV antibody detection, SDMA for glomerular filtration rate estimation (see T1035), and feline coronavirus titer. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	2.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-3 days
SA210	Feline Comprehensive with SDMA Superchem with SDMA, CBC, FeLV, FIV, FCV Titer, <i>Toxoplasma</i> Ab - IgG/IgM Most comprehensive chemistry panel, a complete blood count, FeLV antigen detection, FIV antibody detection, SDMA for glomerular filtration rate estimation (see T1035), feline coronavirus titer, and <i>Toxoplasma</i> titer (IgM and IgG). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-3 days
SA682	Feline EDP with SDMA, FeLV, FIV, O&P Wellness Chemistry with SDMA, CBC, FeLV, FIV, O&P with Centrifugation A comprehensive chemistry, complete blood count, FeLV antigen and FIV antibody detection, SDMA for glomerular filtration rate estimation (see T1035), and fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection (T805). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, and 5 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days
SA908	Feline Health Check with SDMA Vet Screen with SDMA, CBC, FeLV, FCV, FIV, T4 A comprehensive chemistry, a complete blood count, total T4, FeLV antigen and FIV antibody detection, SDMA for glomerular filtration rate estimation (see T1035), and feline coronavirus titer. SA914 to which a T4 has been added. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-3 days
SA914	Feline Health Profile with SDMA Vet Screen with SDMA, CBC, FeLV, FCV, FIV A comprehensive chemistry, a complete blood count, FeLV antigen and FIV antibody detection, SDMA for glomerular filtration rate estimation (see T1035), and feline coronavirus titer. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-3 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA460	Feline Panel #101 with SDMA Superchem with SDMA, CBC, FeLV, FIV, and FCV Titer The most comprehensive chemistry panel, a complete blood count, FeLV antigen and FIV antibody detection, SDMA for glomerular filtration rate estimation (see T1035), and feline coronavirus titer. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-3 days
SA722	Feline Senior Profile with SDMA Superchem with SDMA, CBC, Total T4, Urinalysis, Feline Heartworm Antibody The most comprehensive minimum (superchemistry, complete blood count and urinalysis), total T4, heartworm antibody detection test and and SDMA for glomerular filtration rate estimation (see T1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
SA782	Feline Senior 2 with SDMA, Panleukopenia Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, Panleukopenia Vaccinal Titer (Feline Only) The most comprehensive minimum database (superchemistry, complete blood count, and urinalysis), thyroid panel (T4, FT4 by ED), SDMA for glomerular filtration rate estimation (see T1035), and Panleukopenia vaccinal titer. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.25 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-5 days
SA190	Feline Total Health Check with SDMA Superchem with SDMA, CBC, T4, FIV, FeLV, FCV, and <i>Toxoplasma</i> Ab - IgG/IgM Most comprehensive chemistry panel, a complete blood count, total T4, FeLV antigen and FIV antibody detection, SDMA for glomerular filtration rate estimation (see T1035), feline coronavirus and <i>Toxoplasma</i> titer (IgM and IgG). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA786	Feline Total Health Check with SDMA, UA Superchem with SDMA, CBC, T4, Urinalysis, FeLV, FIV, FCV, <i>Toxoplasma</i> Ab - IgG/IgM The most comprehensive minimum database (superchemistry, complete blood count, and urinalysis), total T4, FeLV antigen and FIV antibody detection, SDMA for glomerular filtration rate estimation (see T1035), feline coronavirus and <i>Toxoplasma</i> titer (IgM and IgG). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-3 Days
SA180	Feline Total Health Plus with SDMA Superchem with SDMA, CBC, T4, FIV, Feline Heartworm Antibody, FeLV, FCV, <i>Toxoplasma</i> Ab - IgG/IgM Most comprehensive chemistry panel, a complete blood count, total T4, FeLV antigen detection, FIV and heartworm antibody detection, SDMA for glomerular filtration rate estimation (see T1035), feline coronavirus and <i>Toxoplasma</i> titer (IgM and IgG). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-3 days
SA700	Feline Wellness 3 with SDMA Superchem with SDMA, CBC, T4, Urinalysis, FeLV, FIV The most comprehensive minimum database (superchemistry, complete blood count, and urinalysis), total T4, FeLV antigen and FIV antibody detection, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
SA902	Feline Wellness 4 with SDMA Vet Screen with SDMA, CBC, T4, FeLV, FIV, Feline Heartworm Antibody A comprehensive chemistry panel, a complete blood count, total T4, FeLV antigen detection, FIV and heartworm antibody detection, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
SA683	Feline Wellness with SDMA, FeLV, FIV, O&P Wellness Chemistry with SDMA, CBC, Feline Heartworm Antibody, FeLV, FIV, O&P with Centrifugation A comprehensive chemistry panel, a complete blood count, FeLV antigen detection, FIV and heartworm antibody detection, SDMA for glomerular filtration rate estimation (see T1035), and fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection (T805). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, and 5 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA910	Geriatric Feline with SDMA Vet Screen with SDMA, CBC, FeLV, FIV, T4 A comprehensive chemistry, a complete blood count, total T4, FeLV antigen and FIV antibody detection, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
SA924	Geriatric Panel with SDMA, Feline Vet Screen with SDMA, CBC, T4, FeLV, FIV, and Urinalysis A comprehensive minimum database (vetscreen, complete blood count, and urinalysis), a total T4, FeLV antigen and FIV antibody detection, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
SA037	Geriatric Profile with SDMA Vet Screen with SDMA, CBC, Total T4, Urinalysis A comprehensive minimum database (vetscreen, complete blood count, and urinalysis), a total T4, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	24 hours Performed each shift
SA1204	GHP Chem with Lytes, SDMA, CBC GHP Chem with Lytes (Total Protein, Albumin, Globulin, ALT (SGPT), Alk Phos, T Bilirubin, BUN, Creatinine, Phosphorus, Glucose, Calcium, Sodium, Potassium, Chloride, Cholesterol, Amylase) with SDMA, and CBC A chemistry panel containing electrolytes (without AST, CPK, GGT, triglyceride, PSL), complete blood count, and SDMA for glomerular filtration rate estimation (see T1035).	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
SA1212	GHP Chem with Lytes, SDMA, CBC, T4 GHP Chemistry with Lytes and SDMA, CBC, and Total T4 A chemistry panel containing electrolytes (without AST, CPK, GGT, triglyceride, PSL), a complete blood count, total T4, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA1214	GHP Chem with Lytes, SDMA, CBC, T4, UA GHP Chemistry with Lytes and SDMA, CBC, Total T4, Urinalysis A complete minimum database (SA1202 (chemistry panel with SDMA), complete blood count, and urinalysis), and total T4. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	24 hours Performed each shift
SA1104	GHP Chem with SDMA, CBC GHP Chemistries with SDMA and CBC A chemistry profile (without electrolytes (sodium, potassium, chloride), AST, CPK, GGT, triglyceride, PSL), a complete blood count, and SDMA for glomerular filtration rate estimation (see T1035).	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
SA1108	GHP Chem with SDMA, CBC, T4 GHP Chem (Total Protein, Albumin, Globulin, ALT (SGPT), Alk Phos, T Bilirubin, BUN, Creatinine, Phosphorus, Glucose, Calcium, Cholesterol, Amylase) with SDMA, CBC, and T4 A chemistry profile (without electrolytes (sodium, potassium chloride), AST, CPK, GGT, triglyceride, PSL), a complete blood count, total T4, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
SA1114	GHP Chem with SDMA, CBC, T4, UA GHP Chemistries with SDMA, CBC, T4, Urinalysis A minimum database (chemistry with SDMA (SA1102), complete blood count, and urinalysis), and total T4. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	24 hours Performed each shift
SA1210	GHP with Lytes, SDMA, CBC, UA GHP Chemistries with Lytes and SDMA, CBC, Urinalysis A complete minimum database (SA1202 (chemistry panel with SDMA), complete blood count, and urinalysis).	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	24 hours Performed each shift
SA235	Hyperthyroid Feline with SDMA Superchem with SDMA, CBC, Total T4, and Free T4 by ED Most comprehensive chemistry profile (superchemistry with SDMA (SA010)), a complete blood count, Total T4, and FT4 ED.	1.25 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA805	Hyperthyroid Panel (with GGT), SDMA CBC, T4, Nsaid Chemistries with GGT AST (SGOT), ALT (SGPT), Alk Phos, BUN, Creatinine GGT A useful panel to evaluate the response to methimazole after a full minimum database has been run to make the initial diagnosis of hyperthyroidism. Includes an abbreviated chemistry (crea, BUN, AST, ALT, ALP and GGT), complete blood count, total T4, and SDMA. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
SA805G	Hyperthyroid Panel with Glucose, SDMA NSAID Chemistries with Glucose, CBC, and T4 A useful panel to evaluate the response to methimazole after a full minimum database has been run to make the initial diagnosis of hyperthyroidism. Includes an abbreviated chemistry, complete blood count, total T4, and SDMA. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
SA440	Hyperthyroid Panel with SDMA NSAID Chemistries AST (SGOT), ALT (SGPT), Alk Phos, BUN, Creatinine, CBC, T4 A useful panel to evaluate the response to methimazole after a full minimum database has been run to make the initial diagnosis of hyperthyroidism. Includes an abbreviated chemistry (crea, BUN, AST, ALT, ALP), complete blood count, total T4, and SDMA. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
SA328	Liver Chem, CBC with Phenobarbital Liver Chemistries, CBC, and Phenobarbital A liver chemistry panel (SA324), complete blood count and phenobarbital level.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top tube, lavender top	1-2 days
SA327	Liver Chemistry, CBC A liver chemistry panel (SA326) to which a complete blood count has been added.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA1502	Liver Profile Chemsitry Panel 112 (Total Protein, Albumin, Globulin, A/G Ratio, AST (SGOT), ALT (SGPT), ALK, T. Bilirubin, BUN, Cholesterol) and CBC A comprehensive liver chemistry which includes AST, ALT, ALP and parameters to assess synthetic capacity/function (albumin, globulin, A/G ratio, BUN, cholesterol, glucose, and bilirubin) and a complete blood count.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
SA320	Liver Profile Standard Liver Chemistry, Complete Blood Count, Bile Acids (Pre and Post) A comprehensive liver panel (SA324), complete blood count and pre/post bile acids (T220). Interferences: marked hemolysis or lipemia. Ursodiol administration may cause bile acids to be increased. Consider withholding ursodiol for five days prior to testing.	0.5 mL serum for resting and post-prandial samples (fasted sample labeled Pre and 2-hour post sample labeled Post) and 1.0 mL EDTA whole blood 2 serum in red top or serum separator tubes (pre and post), lavender top	1-2 days
SA1004	Mini Early Detection Chem with SDMA, CBC Mini Early Detection Chemistries (Total Protein, ALT (SGPT), ALK PHOS, BUN, Creatinine, Glucose) with SDMA and CBC A chemistry (TP, ALT, ALP, BUN, creatinine, and glucose), complete blood count and SDMA.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
SA1010	Mini Early Detection Chem with SDMA, CBC, O&P, Heartworm Mini Early Detection Chemistries (Total Protein, ALT (SGPT), ALK PHOS, BUN, Creatinine, Glucose) with SDMA, CBC, Heartworm Antigen, and O&P with Centrifugation A chemistry (TP, ALT, ALP, BUN, creatinine, and glucose), complete blood count, heartworm antigen detection, SDMA, and fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection.	0.5 mL serum, 1.0 mL EDTA whole blood, and 5 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days
SA070	Miniscreen with CBC Miniscreen 4 Chemistries, CBC Miniscreen chemistry (SA060) and a complete blood count.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
SA072	Miniscreen with SDMA, CBC Miniscreen 11 Chemsitries (Total Protein, Albumin, Globulin, A/G Ratio, ALT (SGPT), Alkaline Phosphatase, T. Bilirubin, BUN, Creatinine, BUN/Creat Ratio, Glucose) and CBC A comprehensive chemistry with SDMA (SA071) and a complete blood count.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA1519	Miniscreen with SDMA, T4 Miniscreen 4 Chemistries, CBC, and T4 A comprehensive chemistry with SDMA (D4C), complete blood count and total T4. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
SA810	NSAID 1 with SDMA Chemistry Panel (Total Protein, ALT (SGPT), Alkaline Phosphatase, Total Bilirubin, BUN, Creatinine) with SDMA, CBC, Urinalysis A minimum database (chemistry with SDMA (SA804), complete blood count and urinalysis) limited to renal and liver values used to monitor a patient on NSAID therapy. Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	24 hours Performed each shift
SA815	NSAID 2 with SDMA Chemistry Panel (Total Protein, ALT (SGPT), Alkaline Phosphatase, Total Bilirubin, BUN, Creatinine) with SDMA, CBC, Urinalysis, Bile Acids A minimum database (chemistry with SDMA (SA804), complete blood count and urinalysis) and a single bile acids test. Interferences: marked hemolysis and lipemia. Ursodeoxycholic acid may be detected by bile acid assay causing falsely elevated values.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
SA820	NSAID 3 with SDMA Chemistry Panel (Total Protein, ALT (SGPT), Alkaline Phosphatase, Total Bilirubin, BUN, Creatinine) with SDMA, CBC, Urinalysis, Bile Acid: Creatinine Ratio A minimum database (chemistry with SDMA (SA804), complete blood count and urinalysis) and a urine transport tube bile acid: creatinine ratio (T227). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
SA054	Pre-Op Panel Plus with Lytes, SDMA Pre-Op Chemistries with Electrolytes and SDMA, CBC, PT, and aPTT A chemistry with electrolytes and SDMA (SA043), complete blood count, prothrombin time, and partial thromboplastin time. Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 0.5 mL of citrated plasma Serum in red top or serum separator tube, lavender top, citrated plasma (blue top)	24 hours Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA050	Pre-Op Panel Plus with SDMA Pre-Op Screen with SDMA, CBC, PT and aPTT A smaller chemistry profile with SDMA (SA040), a complete blood count, prothrombin time, and partial thromboplastin time. Interferences: marked hemolysis and lipemia. Ensure that the submission of a blue top tube that has been appropriately filled to the designated mark.	0.5 mL serum, 1.0 mL EDTA whole blood, and 0.5 mL of citrated plasma Serum in red top or serum separator tube, lavender top, citrated plasma (blue top)	24 hours Performed each shift
SA053	Pre-Op Panel with Electrolytes, SDMA and Urinalysis Pre-Op Chem with Electrolytes and SDMA, CBC, Urinalysis A minimum database that includes a smaller chemistry with SDMA (SA043), a complete blood count and urinalysis. Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	24 hours Performed each shift
SA920	Pre-Op Panel with SDMA, T4, UA Pre-Op Screen with SDMA, CBC, T4, Urinalysis A minimum database (smaller chemistry profile with SDMA (SA040), a complete blood count, urinalysis) and total T4. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	24 hours Performed each shift
SA512	Pre-Op Plus with SDMA, CBC, Heartworm Pre-Op Screen with SDMA, CBC, PT, PTT, Heartworm Antigen A smaller chemistry profile with SDMA (SA040), complete blood count, heartworm antigen detection, and prothrombin time and partial thromboplastin time. Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.5 mL of citrated plasma Serum in red top or serum separator tube, lavender top, citrated plasma (blue top)	1-2 days
SA510	Pre-Op Plus with SDMA, CBC, FeLV, FIV Pre-Op Screen with SDMA, CBC, PT, PTT, FIV, FeLV A smaller chemistry profile with SDMA (SA040), complete blood count, FeLV antigen and FIV antibody detection, and prothrombin and partial thromboplastin time. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.5 mL citrated plasma Serum in red top or serum separator tube, lavender top, citrated plasma (or blue top)	1-2 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA508	Pre-Op Screen with SDMA, CBC, T4 Pre-Op Screen with SDMA, CBC, and Total T4 A smaller chemistry profile with SDMA (SA040), a complete blood count, and total T4 (T495). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
SA514	Pre-Op Screen with SDMA, CBC, T4, FeLV, FIV Pre-Op Screen with SDMA, CBC, T4, FIV, and FeLV SA508 to which FeLV antigen and FIV antibody detection has been added. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
SA052	Pre-Op Screen with SDMA, CBC, UA Pre-Op Screen with SDMA, CBC, and Urinalysis A minimum database including smaller chemistry profile with SDMA (SA040), a complete blood count, and urinalysis. Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	24 hours Performed each shift
SA516	Pre-Op Screen with Electrolytes, SDMA, CBC Pre-Op Chemistry with Electrolytes and SDMA, CBC A smaller chemistry profile which includes electrolytes and SDMA (SA043) and a complete blood count. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
SA055	Pre-Op Screen with SDMA, CBC Pre-Op Screen with SDMA and CBC A smaller chemistry profile with SDMA (SA040) and a complete blood count. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
RECHECK	Recheck Profile Standard with SDMA Superchem with SDMA and CBC This is a superchemistry with SDMA (SA010) and complete blood count submitted within 30 days of a previous accession to assess/monitor trends in values. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
RECHECK2	Recheck Profile Vetscreen with SDMA, CBC Vet Screen with SDMA and CBC This is a vet screen chemistry with SDMA (SA025) and complete blood count submitted within 30 days of a previous accession to assess/monitor trends in values. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
SA313	Renal Profile Standard with SDMA Chemistry Renal Profile with SDMA and CBC A chemistry panel (total protein, albumin, globulin, BUN, creatinine, phosphorous, calcium, sodium, and potassium) with SDMA and a complete blood count. Interferences: marked hemolysis or lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
SA310	Renal Profile with SDMA Chemistry Renal Profile with SDMA, CBC, and Urinalysis A minimum database including a chemistry profile with SDMA (T7008), complete blood count and urinalysis. Interferences: marked hemolysis or lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	24 hours Performed each shift
SA724	Superchem with SDMA, CBC, T4, FeLV, FIV, HeartwormAB Superchem with SDMA, CBC, Total T4, FeLV, FIV, and Heartworm Antibody, Feline The most comprehensive chemistry with SDMA (SA010), a complete blood count, total T4, FeLV antigen, and FIV and heartworm antibody detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum and 1.0 mL of EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
SA098	Superchem with SDMA, CBC, T4, FT4ED, Heartworm Superchem with SDMA, CBC, Total T4, Free T4 by ED, and Heartworm Antigen The most comprehensive chemistry with SDMA (SA010), a complete blood count, a thyroid screening panel (T4, FT4 ED) and heartworm antigen detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum and 1.0 mL of EDTA whole blood Serum in red top or serum separator tube, lavender top	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA712	<p>Senior Care Basic Plus Profile with SDMA</p> <p>Superchem with SDMA, CBC, Total T4, Urinalysis, and Protein/Creatinine Ratio</p> <p>The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, and urine protein to creatinine ratio (T775).</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>0.5 mL serum and 1.0 mL of EDTA whole blood</p> <p>Serum in red top or serum separator tube, lavender top</p>	<p>24 hours</p>
SA080	<p>Senior Comprehensive Plus with SDMA</p> <p>Superchem with SDMA, CBC, Total T4, Free T4 by ED, and TSH</p> <p>The most comprehensive chemistry panel with SDMA (SA010), a complete blood count, and thyroid panel (T4, FT4 ED, and TSH).</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>1.5 mL serum and 1.0 mL EDTA whole blood</p> <p>Serum in red top or serum separator tube, lavender top</p>	<p>1-3 days</p>
SA090	<p>Senior Comprehensive with SDMA</p> <p>Superchem with SDMA, CBC, Total T4, and Free T4 by ED</p> <p>The most comprehensive chemistry panel with SDMA (SA010), a complete blood count, and thyroid panel (T4 and FT4 ED).</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>1.0 mL serum and 1.0 mL EDTA whole blood</p> <p>Serum in red top or serum separator tube, lavender top</p>	<p>1-3 days</p>
SA086	<p>Senior Comprehensive with SDMA, D.Bilirubin, I.Bilirubin, Anion, Osmolality</p> <p>Superchem with Direct & Indirect Bilirubin and SDMA, Anion gap, Osmolality, CBC, Total T4, and Free T4 by ED</p> <p>The most comprehensive chemistry with SDMA (SA010) to which indirect bilirubin, direct bilirubin, anion gap, serum osmolality has been added, a complete blood count, and thyroid panel (T4 and FT4 ED).</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>1.0 mL serum and 1.0 mL EDTA whole blood</p> <p>Serum in red top or serum separator tube, lavender top</p>	<p>1-3 days</p>
SA740	<p>Senior Feline Plus with SDMA</p> <p>Superchem with SDMA, CBC, Total T4, FeLV, FIV Antibody, Urinalysis, and O&P with Centrifugation</p> <p>The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, FeLV antigen and FIV antibody detection, and fecal analysis using zinc sulfate centrifugation/flotation for ova and parasite detection.</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	<p>1-2 days</p>

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA755	Senior Panel 2 with SDMA, O&P Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, O&P with Centrifugation The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), thyroid panel (T4 and FT4 ED), and fecal analysis using zinc sulfate centrifugation/flotation for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
SA746	Senior Panel 2 with SDMA, UMIC Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, Culture The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), thyroid panel (T4 and FT4 ED), and urine culture. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 7.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube	1-3 days
SA760	Senior Panel 2 with SDMA, Vaccine Titers Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, Distemper/Parvo Vaccinal Titer The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), thyroid panel (T4 and FT4 ED), and Distemper/Parvovirus vaccine titer (T565). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-3 days
SA954	Senior Profile 1 with SDMA, Feline Heartworm Ab, FeLV, FIV, Fecal Combo Superchem with SDMA, CBC, T4, Urinalysis, Heartworm Antibody, Feline, FeLV, FIV Antibody, O&P with Centrifugation and <i>Giardia</i> The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, FeLV antigen detection, FIV and heartworm antibody detection, and fecal analysis using zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 6 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA750	Senior Profile 1 with SDMA, O&P Superchem with SDMA, CBC, Total T4, Urinalysis, Fecal O&P with Centrifugation The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, and fecal analysis using zinc sulfate centrifugation/flotation for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
SA741	Senior Profile 1 with SDMA, Virals, Fecal Combo Superchem with SDMA, CBC, T4, Urinalysis, FeLV, FIV Antibody, and O&P with Centrifugation and <i>Giardia</i> The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, FeLV antigen and FIV antibody detection, and fecal analysis using zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 6 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
SA705	Senior Profile 1 with SDMA Superchem with SDMA, CBC, Total T4, Urinalysis The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis) and a total T4. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	24 hours Performed each shift
SA566	Senior Profile 1 with SDMA, UMIC Superchem with SDMA, CBC, Total T4, Urinalysis, Culture The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, and urine culture is added. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 7.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-4 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA770	Senior Profile 1 with SDMA, Fecal Combo Superchem with SDMA, CBC, Total T4, Urinalysis, Fecal O&P with Centrifugation, <i>Giardia</i> The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, and fecal analysis using zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 6 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
SA719	Senior Profile 1 with SDMA, Feline Heartworm Ag Superchem with SDMA, CBC, T4, Urinalysis, Feline Heartworm Antigen The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, and heartworm antigen detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
SA735	Senior Profile 1 with SDMA, FT4ED, Heartworm Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, Heartworm Antigen The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a thyroid panel (T4 and FT ED), and heartworm antigen detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-3 days
SA720	Senior Profile 2 with SDMA Superchem with SDMA, CBC, Total T4, Urinalysis, and Free T4 by ED The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis) and a thyroid panel (T4 and FT4 ED). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA775	Senior Profile 2 with SDMA, Fecal Combo Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, and O&P with Centrifugation, <i>Giardia</i> The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a thyroid panel (T4 and FT4 ED), and fecal analysis using zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 6 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
SA766	Senior Profile 2 with SDMA, Fel Heartworm Ab Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, and Feline Heartworm Antibody The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a thyroid panel (T4 and FT4 ED), and heartworm antibody detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-3 days
SA780	Senior Profile 2 with SDMA, FeLV, FIV Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, FIV, and FeLV The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a thyroid panel (T4 and FT4 ED), FeLV antigen and FIV antibody detection. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-3 days
SA790	Senior Profile 2 with SDMA, O&P, Fel Heartworm Ab Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, O&P with Centrifugation, Feline Heartworm Antibody The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a thyroid panel (T4 and FT4 ED), heartworm antibody detection, and fecal analysis using zinc sulfate centrifugation/flotation for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA762	Senior Profile 2 with SDMA, O&P, Heartworm Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, Heartworm Antigen, and O&P with Centrifugation The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a thyroid panel (T4 and FT4 ED), heartworm antigen detection, and fecal analysis using zinc sulfate centrifugation/flotation for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
SA730	Senior Profile with SDMA, Vaccine Titers Superchem with SDMA, CBC, T4, Urinalysis, and Distemper/Parvo Vaccinal Titer The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, and Distemper/Parvovirus vaccine titer. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-3 days
SA768	Senior Profile with SDMA, Feline Heartworm Ab, Fecal Combo Superchem with SDMA, CBC, T4, Urinalysis, Feline Heartworm Antibody, and O&P with Centrifugation, <i>Giardia</i> The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, heartworm antibody detection, and fecal analysis using zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 6 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
SA725	Senior Profile with SDMA, Panleukopenia Superchem with SDMA, CBC, T4, Urinalysis, and Panleukopenia Vaccine Titer (Feline Only) The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, and Panleukopenia vaccine titer. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-5 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA745	Senior Super Plus with SDMA Superchem with SDMA, CBC, Total T4, Urinalysis, Heartworm Antigen, and Fecal O&P with Centrifugation The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, heartworm antigen detection test, and fecal analysis using zinc sulfate centrifugation/flotation for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
SA764	Senior Wellness with SDMA, Feline Heartworm Ab, O&P Superchem with SDMA, CBC, T4, Urinalysis, Feline Heartworm Antibody, O&P with Centrifugation The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a total T4, heartworm antibody detection, and fecal analysis using zinc sulfate centrifugation/flotation for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
SA081	Senior Comprehensive Plus with SDMA, UA Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, and TSH The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), and thyroid panel (T4, FT4 ED, and TSH). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-3 days
SA784	Senior Comprehensive with SDMA, Heartworm, UA, UPCR Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, Heartworm Antigen, and Protein/Creatinine Ratio The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), thyroid panel (T4, FT4 ED), heartworm antigen detection, and urine protein to creatinine ratio. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum, 1.0 mL EDTA whole blood, and 6.5 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA715	Senior Profile 1 with SDMA, Virals, Feline Heartworm Ab Superchem with SDMA, CBC, Total T4, Urinalysis, FeLV, FIV, Feline Heartworm Antibody The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), a total T4, FeLV antigen, FIV and heartworm antibody detection. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
SA1601	Superchem with SDMA, CBC, Free T4 ED Superchem with SDMA, CBC, and Free T4 by ED Most comprehensive chemistry profile (superchemistry with SDMA (SA010)), a complete blood count and FT4 ED. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-3 days
SA020	Superchem with SDMA, CBC Superchem with SDMA and CBC Most comprehensive chemistry profile which includes SDMA for estimation of the glomerular filtration rate (T1035) and a complete blood count. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
SA028	Superchem with SDMA, CBC, Heartworm Antigen Superchem with SDMA, CBC, and Heartworm Antigen The most comprehensive chemistry with SDMA (SA010), a complete blood count and heart worm antigen detection. Interferences: marked hemolysis and lipemia.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
SA450	Superchem with SDMA, CBC, Thyroid Panel (T3, T4, FT4ED, TSH, TGAA) Superchem with SDMA, CBC, T4, T3, Free T4 by ED, TSH, and Thyroglobulin Auto Antibody Most comprehensive chemistry profile (superchemistry with SDMA (SA010)), a complete blood count, and thyroid panel (T3, T4, FT4 ED, TSH and thyroglobulin autoantibody (TGAA)). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top, urine transport tube	1-3 days
SA021	Superchem with SDMA, CBC, UA Superchem with SDMA, CBC, Urinalysis The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis). Interferences: marked hemolysis and lipemia (serum and whole blood).	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	24 hours Performed each shift

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA490	Superchem with SDMA, CBC, Viral Panel Superchem with SDMA, CBC, FeLV Antigen and FIV Antibody Most comprehensive chemistry profile (superchemistry with SDMA (SA010)), a complete blood count, FeLV antigen and FIV antibody detection. Interferences: marked hemolysis or lipemia.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
SA788	Total Body Function Plus with SDMA, Fecal Combo Superchem with SDMA, CBC, T4, Heartworm Antigen, and O&P with Centrifugation and <i>Giardia</i> The most comprehensive chemistry panel with SDMA (SA010), a complete blood count, total T4, heartworm antigen detection and fecal analysis using zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA (T808) for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days
SA114	Total Body Function Plus with SDMA, O&P Superchem with SDMA, CBC, T4, Heartworm Antigen, and O&P with Centrifugation The most comprehensive chemistry panel with SDMA (SA010), a complete blood count, total T4, heartworm antigen detection and fecal analysis using zinc sulfate centrifugation/flotation (T805) for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 5.0 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days
SA120	Total Body Function with SDMA Superchem with SDMA, CBC, and Total T4 The most comprehensive chemistry panel with SDMA (SA010), a complete blood count, and a total T4. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift
SA110	Total Body Function with SDMA, Heartworm Superchem with SDMA, CBC, T4, Heartworm Antigen The most comprehensive chemistry panel with SDMA (SA010), a complete blood count, a total T4, and heartworm antigen detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA225	Total Body Function with SDMA, Feline Virals Superchem with SDMA, CBC, T4, FeLV, FIV Antibody, and FCV Titer The most comprehensive chemistry panel with SDMA (SA010), a complete blood count, a total T4, FeLV antigen and FIV antibody detection, and feline coronavirus titer. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-3 days
SA122	Total Body Function with SDMA, O&P Superchem with SDMA, CBC, T4, and O&P with Centrifugation The most comprehensive chemistry panel with SDMA (SA010), a complete blood count, total T4, and fecal analysis using zinc sulfate centrifugation/flotation (T805) for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, and 5.0 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days
SA1633	Total Body Function with SDMA, TSH Superchem with SDMA, CBC, T4, and TSH The most comprehensive chemistry panel with SDMA (SA010), a complete blood count and a thyroid panel (total T4 and TSH). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
SA035	Vet Screen 2 Plus with SDMA, CBC Vet Screen with SDMA with Amylase and CBC A comprehensive chemistry profile with SDMA (SA025), with amylase evaluation added, and a complete blood count. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours
SA030	Vet Screen with SDMA, CBC Vet Screen with SDMA and CBC A comprehensive chemistry profile with SDMA (SA025) and a complete blood count. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours
SA530	Vet Screen with SDMA, CBC, Heartworm, T4 Vet Screen with SDMA, CBC, Heartworm Antigen, and T4 A comprehensive chemistry panel with SDMA (SA025), complete blood count, heartworm antigen detection, and total T4. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA906	Vet Screen with SDMA, CBC, FeLV, FIV, UA Vet Screen with SDMA, CBC, FeLV, FIV, Urinalysis A comprehensive minimum database (comprehensive chemistry with SDMA (SA025), a complete blood count and urinalysis), FeLV antigen and FIV antibody detection. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
SA1645	Vet Screen with SDMA, CBC, FeLV, T4 Vet Screen with SDMA, CBC, FeLV, and T4 A comprehensive chemistry with SDMA (SA025), complete blood count, total T4 and FeLV antigen detection. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
SA039	Vet Screen with SDMA, CBC, Heartworm Vet Screen with SDMA, CBC, Heartworm Antigen A comprehensive chemistry profile with SDMA (SA025), a complete blood count, and evaluation for heartworm antigen. SA030 to which heartworm antigen detection has been added. Interferences: marked hemolysis and lipemia.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
SA905	Vet Screen with SDMA, CBC, O&P, Giardia Vet Screen with SDMA, CBC, O&P with Centrifugation and <i>Giardia</i> A comprehensive chemistry with SDMA (SA025), a complete blood count, and fecal analysis using zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA (T808) for ova and parasite detection. Interferences: marked hemolysis or lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days
SA031	Vet Screen with SDMA, CBC, PT, PTT Vet Screen with SDMA, CBC, PT and APTT A comprehensive chemistry profile with SDMA (SA025), a complete blood count, Prothrombin Time (PT) and Activated Partial Thromboplastin Time (APTT). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 0.5 mL citratd plasma Serum in red top or serum separator tube, lavender top, citrated plasma (blue top)	24 hours Performed each shift
SA034	Vet Screen with SDMA, CBC, T4 Vet Screen with SDMA, CBC, and Total T4 A comprehensive chemistry profile with SDMA (SA025), complete blood count, and a total T4. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	24 hours Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA036	Vet Screen with SDMA, CBC, Urinalysis Vet Screen with SDMA, CBC, Urinalysis A complete minimum database (comprehensive chemistry and SDMA (SA025), complete blood count, urinalysis). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	24 hours Performed each shift
SA912	Vet Screen with SDMA, Thyroids Vet Screen with SDMA, T4, CBC, Free T4 by ED A comprehensive chemistry profile with SDMA (SA025), a complete blood count, and thyroid profile (T4 and a FT4 ED). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-3 days
SA904	Vet Screen with SDMA, CBC, UA, O&P, Heartworm Vet Screen with SDMA, CBC, Urinalysis, O&P with Centrifugation, Heartworm Antigen A complete minimum database (comprehensive chemistry with SDMA (SA025), complete blood count, and urinalysis), heartworm antigen detection and fecal analysis using zinc sulfate centrifugation/flotation (T805) for ova and parasite detection. Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5.0 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
SA647	Wellness Profile with SDMA, O&P, UA Wellness Chemistry with SDMA, CBC, Fecal O&P with Centrifugation, Urinalysis The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis) and fecal analysis using zinc sulfate centrifugation/flotation (T805) for ova and parasite detection. Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5.0 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

HEMATOLOGY

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
ADD03	Add-on CBC Includes WBC, RBC, Hemoglobin, Hematocrit, MCV, MCH, MCHC, Platelet count and Differential: Percent and Absolute Neutrophils, Bands, Lymphocytes, Monocytes, Eosinophils, Basophils Assessment of the quantity and morphology of the erythrocytes, leukocytes and platelets. Interferences: marked hemolysis.	1.0 mL EDTA whole blood in lavender top tube	Daily Performed each shift
T395	APTT Partial Thromboplastin Time (PTT) measures the integrity of the intrinsic and common components of the coagulation cascade. Interferences: marked hemolysis or lipemia. Partially full blue top tube may falsely increase sample's coagulation time. Note: clotting of the sample may preclude the analysis. The blue top tube needs to be filled to 2/3 or more of its capacity.	0.5 mL citrated plasma in non-additive transport tube.	24 hours Performed each shift
S16100	Blood Type, Canine - Full Panel Blood Type: DEA 1, DEA 4, DEA 5, DEA 7 This panel is used for RBC typing for identifying blood donors. It evaluates the presence of the blood group antigens; DEA 1.1, 1.2, 1.3, 4, 5, and 7. Antibody screen can be performed at no extra charge if clients submit serum with LTT.	1.0 mL EDTA whole blood in lavender top tube	2-4 days
T315	Blood Type, Canine DEA 1.1 Only This test determines the appropriateness of canine blood donors prior to transfusion. Dog Erythrocyte Antigen (DEA) 1.1 is the most antigenic blood group in dogs, and canine blood donors are often screened to see if they are positive or negative for this blood group. Interferences: freezing of sample may preclude testing. Note: keep sample refrigerated prior to transport and send with an ice pack.	1.0 mL EDTA whole blood in lavender top tube	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T330R	CBC Recheck Includes WBC, RBC, Hemoglobin, Hematocrit, MCV, MCH, MCHC, Platelet count and Differential: Percent and Absolute Neutrophils, Bands, Lymphocytes, Monocytes, Eosinophils, Basophils A complete blood count reevaluation. Resubmission must be within 4 weeks of original accession. Interferences: marked hemolysis and lipemia.	1.0 mL EDTA whole blood in lavender top tube 	24 hours Performed each shift
AE275	CBC Small Mammalian Includes WBC, RBC, Hemoglobin, Hematocrit, MCV, MCH, MCHC, Platelet count and Differential: Percent and Absolute Neutrophils, Heterophils, Bands, Lymphocytes, Monocytes, Eosinophils, Basophils Interferences: marked hemolysis.	0.5 mL EDTA whole blood in lavender top tube 	24 Hours Performed each shift
SA290	Coagulation Panel with D-Dimer CBC, PT and aPTT, Fibrinogen (Quantitative), and D-Dimer Panel used to evaluate a patient exhibiting unexplained bleeding and/or thromboembolic disorders. Additionally, it can be used to evaluate risk of bleeding secondary to provocative procedures. Includes a CBC, PT/PTT, fibrinogen and D-dimer. Interferences: marked hemolysis, lipemia or clotting.	1.0 mL EDTA whole blood and 1.0 mL citrated plasma Lavender top, citrated plasma (or filled blue top)	24 hours Performed each shift
SA300	Coagulation Profile 2 PT and aPTT, Fibrinogen (Quantitative), D-Dimer, and Platelet Count Panel used to evaluate a patient exhibiting unexplained bleeding and/or thromboembolic disorders. Additionally, it can be used to evaluate risk of bleeding secondary to provocative procedures. Includes PT/PTT, platelet count, fibrinogen and D- dimer. Interferences: marked hemolysis, lipemia or clotting.	1.0 mL EDTA whole blood and 1.0 mL citrated plasma Lavender top, citrated plasma (or filled blue top)	24 hours Performed each shift

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA305	Coagulation Profile 3 PT and aPTT, and Platelet Count Panel used to evaluate a patient exhibiting unexplained bleeding. Additionally, it can be used to evaluate risk of bleeding secondary to provocative procedures. Includes PT/PTT and platelet count. Interferences: marked hemolysis, marked lipemia or clotting.	1.0 mL EDTA whole blood and 1.0 mL citrated plasma Lavender top, citrated plasma (or filled blue top)	24 hours Performed each shift
T330 Add-on Equivalent ADD03	Complete Blood Count Includes WBC, RBC, Hemoglobin, Hematocrit, MCV, MCH, MCHC, Platelet count and Differential: Percent and Absolute Neutrophils, Bands, Lymphocytes, Monocytes, Eosinophils, Basophils Includes WBC, RBC, HGB, HCT, MCV, MCH, MCHC, platelet count and estimate, WBC differential, RBC and WBC morphology. Interferences: marked hemolysis and lipemia.	1.0 mL EDTA whole blood in lavender top tube	24 hours Performed each shift
T337	Complete Blood Count with Retic Count CBC and Reticulocyte Count This includes WBC, RBC, HGB, HCT, MCV, MCH, MCHC, platelet count and estimate, WBC differential, RBC, WBC morphology and a reticulocyte count. Interferences: marked hemolysis and lipemia.	1.0 mL EDTA whole blood in lavender top tube	24 hours Performed each shift
T540	Direct Coombs' Test Warm Used to investigate hemolytic anemia. Interferences: marked hemolysis and lipemia. Test should be performed within 24 hours of sample collection.	1.0 mL EDTA whole blood in lavender top tube	24 hours

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T320	Feline Blood Typing Evaluation of feline blood type as Type A, Type B, or AB. Note: keep sample refrigerated and send with an ice pack.	1.0 mL EDTA whole blood in lavender top tube	2-10 days
T365 Add-on Equivalent ADD40	Fibrinogen Quantitative This test is used to measure the concentration of functional fibrinogen in the plasma. Interferences: clotted sample precludes analysis. Citrated Whole Blood (blue top tube) or citrated plasma is the only acceptable sample. The tube should be greater than 2/3rds filled. If submitting separated citrated plasma, label it as Citrated Plasma.	0.5 mL citrated plasma collected as whole blood in blue top tube at least 2/3rds full to the fill line	24 hours Performed each shift
RE VW	Path Review CBC This code is used to add-on a clinicopathologist review to a CBC.	0.5 mL EDTA whole blood in lavender top tube	1-2 days
T400 Add-on Equivalent ADD120	Platelet Count Interferences: marked hemolysis.	1.0 mL EDTA whole blood in lavender top tube	24 hours Performed each shift

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T401	Platelet Count Manual Interferences: marked hemolysis.	1.0 mL EDTA whole blood in lavender top tube	24 Hours Performed each shift
T410	Prothrombin Time Prothrombin Time (PT) measures the integrity of the extrinsic and common components of the coagulation cascade. Interferences: marked hemolysis or lipemia. Partially full blue top tube may falsely increase sample's coagulation time. Note: clotting of the sample may preclude the analysis. The blue top tube needs to be filled to 2/3 or more of its capacity.	0.5 mL citrated plasma in non-additive transport tube.	24 hours Performed each shift
T415	PT and aPTT Prothrombin Time and Activated Partial Thromboplastin Time. Prothrombin time measures the integrity of the extrinsic and common components of the coagulation cascade. Partial Thromboplastin Time (PTT) measures the integrity of the intrinsic and common components of the coagulation cascade. Interferences: marked hemolysis or lipemia. Partially full blue top tube may falsely increase sample's coagulation time. Note: clotting of the sample may preclude the analysis. The blue top tube needs to be filled to 2/3 or more of its capacity. Partially full blue top tubes may falsely increase the coagulation times.	0.5 mL citrated plasma in non additive transport tube.	24 hours Performed each shift
ADD290	PT and APTT Add-on Prothrombin Time and Activated Partial Thromboplastin Time. Prothrombin time measures the integrity of the extrinsic and common components of the coagulation cascade. Partial Thromboplastin Time (PTT) measures the integrity of the intrinsic and common components of the coagulation cascade. Interferences: marked hemolysis or lipemia. Partially full blue top tube may falsely increase sample's coagulation time. Note: clotting of the sample may preclude the analysis. The blue top tube needs to be filled to 2/3 or more of its capacity. Partially full blue top tubes may falsely increase the coagulation times.	0.5 mL citrated plasma in non-additive transport tube.	24 hours Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T425	Reticulocyte Count	1.0 mL EDTA whole blood in lavender top tube	24 hours
Add-on Equivalent ADD140	<p>Reticulocytes, Absolute Reticulocytes</p> <p>Includes an unadjusted reticulocyte percent and absolute reticulocyte count. Only aggregate reticulocytes are counted.</p> <p>Interferences: marked hemolysis.</p>		Performed each shift
ADD140	Reticulocyte Count Add-on	1.0 mL EDTA whole blood in lavender top tube	24 hours
	<p>Reticulocytes, Absolute Reticulocytes</p> <p>Includes an unadjusted reticulocyte percent and absolute reticulocyte count. Only aggregate reticulocytes are counted.</p> <p>Interferences: unable to be performed if sample has marked hemolysis.</p>		Performed each shift
T331	Special CBC, Path Review	1.0 mL EDTA whole blood in lavender top tube	1-3 days
	<p>This includes WBC, RBC, HGB, HCT, MCV, MCH, MCHC, platelet count and estimate, WBC differential, RBC, and WBC morphology. A clinical pathologist will then review results and blood smears and provide an interpretative comment of the complete blood count.</p> <p>Interferences: marked hemolysis and lipemia. Clotting or freezing of sample may preclude analysis.</p>		
S17123	Von Willebrand Factor	1.5 mL citrated plasma collected as whole blood in blue top tube, spun and plasma placed in non-additive tube (labeled as Citrated plasma). Freeze citrated plasma.	5-7 days
	<p>To evaluate for Von Willebrands factor deficiency. Results will be expressed as a percentage. The percentage will indicate the risk of bleeding problems.</p> <p>Note: for genetic screening, do not test bitches in season, pregnant or lactating. Do not test unhealthy animals (with no evidence of primary clotting issues), those on medication for a recent illness, or those vaccinated within 14 days. Large breed pups can be tested at seven weeks of age. Small breed pups should be over 12 weeks of age prior to testing.</p>		

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CHEMISTRY

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
FRUCTADD	Add-on Fructosamine Interferences: hemolysis and lipemia may affect results.	0.5 mL serum in red top or serum separator tube	24 hours
SA665	Adult Wellness Chemistry with SDMA Wellness Chemistry (ALT (SGPT), Alkaline Phosphatase, Total Protein, Albumin, Globulin, A/G Ratio, BUN, Creatinine, BUN/Creatinine Ratio, Glucose, Sodium, Potassium, Na/K Ratio, and Chloride) with Electrolytes and SDMA A comprehensive chemistry panel and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	24 hours
T225	Bile Acids Single sample labelled as either pre or post bile acids. Interferences: marked hemolysis or lipemia. Note: ursodeoxycholic acid may be detected by bile acid assay, causing falsely elevated values. 1.0 mL of serum labeled either as fasted (10-12 hours fast) or post prandial (2 hours after feeding a maintenance diet meal).	0.5 mL serum in red top or serum separator tube	1-2 days
T220	Bile Acids Pre and Post Fasted bile acids, 2 hour post feeding bile acids Pre and post prandial bile acids Interferences: marked hemolysis or lipemia. Ursodeoxycholic acid may be detected by bile acid assay, causing falsely elevated values. Note: suggested protocol: 1. Obtain a fasting serum sample (1 mL) and label the tube preprandial. 2. Obtain a second serum sample (1 mL) 2 hours after feeding the animal a maintenance diet and label the tube postprandial.	0.5 mL serum in red top or serum separator tube (fasted sample labeled Pre) and 0.5 serum in red top or serumseparator tube (2-hour post sample labeled Post)	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T110	Calcium Total Calcium Interferences: marked hemolysis or lipemia. Lipemia can result in falsely elevated total calcium values.	0.5 mL serum in red top tube or serum separator tube	24 hours Performed each shift
S18537	Calcium, Ionized Used to further explore the relevance of abnormalities in total calcium by assessing its ionized fraction. Interferences: aerobic exposure, or severe lipemia. Note: sample must be anaerobically transferred from spun collection tube into plain red top.	2.0 mL serum anaerobically transferred from spun red top or serum separator tube into a plain, unopened red top tube.	1-4 days
SA804	Chemistry Panel with SDMA Total Protein, ALT (SGPT), Alkaline Phosphatase, Total Bilirubin, BUN, Creatinine, Sodium, Potassium, Chloride, and SDMA An abbreviated chemistry panel including ALT, ALP, total bilirubin, BUN, creatinine, total protein, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
T7008	Chemistry Renal Profile with SDMA Total Protein, Albumin, Globulin, A/G Ratio, BUN, Creatinine, BUN/Creat Ratio, Phosphorus, Calcium, Corrected Calcium, Sodium, Potassium, Na/K Ratio, Chloride, and SDMA A comprehensive profile panel that can be used to trend response to therapy when treating renal disease (glomerular or tubulointerstitial) which includes SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
CSFPR	CSF Protein This is a component of the CSF fluid analysis and cytology but can be requested alone. Interferences: marked hemolysis and lipemia.	0.5 mL cerebrospinal fluid in lavender or red top tube	24 hours

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T140	Electrolyte Screen	0.5 mL serum in red top or serum separator tube	24 hours
	Sodium, potassium, sodium/potassium ratio, chloride, bicarbonate and anion gap		Performed each shift
	Evaluation of serum sodium, potassium, sodium/potassium ratio, chloride, bicarbonate, and anion gap. Interferences: marked hemolysis and lipemia.		
T145	GGTP	0.5 mL serum in red top or serum separator tube	24 hours
	Gamma-glutamyl transferase		Performed each shift
	Interferences: marked hemolysis or lipemia.		
SA1202	GHP Chem with Lytes, SDMA	0.5 mL serum in red top or serum separator tube	24 hours
	Total Protein, Albumin, Globulin, AL:T (SGPT), Alkaline Phosphatase, Total Bilirbin, BUN, Creatinine, Phosphorus, Glucose, Calcium, Sodium, Potassium, Chloride, Cholesterol, Amylase, and SDMA		Performed each shift
	A chemistry panel containing electrolytes (without AST, CPK, GGT, triglyceride, PSL) and SDMA for glomerular filtration rate estimation (see T1035).		
SA1102	GHP Chem with SDMA	0.5 mL serum in red top or serum separator tube	24 hours
	Total Protein, Albumin, Globulin, AL:T (SGPT), Alkaline Phosphatase, Total Bilirbin, BUN, Creatinine, Phosphorus, Glucose, Calcium, Sodium, Potassium, Chloride, Cholesterol, Amylase, and SDMA		Performed each shift
	A chemistry profile (without electrolytes (sodium, potassium, chloride), AST, CPK, GGT, triglyceride, PSL), and SDMA for glomerular filtration rate estimation (see T1035).		
SA321	Liver Chemistries	0.5 mL serum in red top or serum separator tube	24 hours
	Total Protein, Albumin, Globulin, A/G Ratio, AST (SGOT), ALT (SGPT), Alkaline Phosphatase, GGT, Total Bilirubin, Direct Bilirubin, Sodium, Potassium, Chloride, and Cholesterol		Performed each shift
	A chemistry panel that includes all cholestatic and hepatocellular liver enzymes, bilirubin, albumin, globulin and cholesterol. Interferences: marked hemolysis or lipemia.		

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA324	Liver Chemistry Total Protein, Albumin, Globulin, A/G Ratio, AST (SGOT), ALT (SGPT), Alkaline Phosphatase, GGT, Total Bilirubin, Bun, Sodium, Potassium, Chloride, Glucose A chemistry panel that includes all cholestatic and hepatocellular liver enzymes, bilirubin, albumin, globulin, BUN and glucose. Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
SA326	Liver Chemistry Panel Total Protein, Albumin, Globulin, A/G Ratio, AST (SGOT), ALT (SGPT), Alkaline Phosphatase, GGT, Total Bilirubin, Direct and Indirect Bilirubin, Bun, Sodium, Potassium, Chloride, Glucose, Cholesterol The most comprehensive liver chemistry which includes all cholestatic and hepatocellular liver enzymes and parameters to assess synthetic capacity/function (albumin, BUN, cholesterol, glucose, and bilirubin). Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
SA329	Liver Chemistry with Phenobarbital Total Protein, Albumin, Globulin, A/G Ratio, AST (SGOT), ALT (SGPT), Alkaline Phosphatase, GGT, Total Bilirubin, Direct Bilirubin, Bun, Sodium, Potassium, Chloride, Cholesterol, and Phenobarbital A liver chemistry panel (SA321) to which a phenobarbital level has been added. Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top	24 hours Performed each shift
SA1002	Mini Early Detection Chem with SDMA Total Protein, ALT (SGPT), ALK PHOS, BUN, Creatinine, Sodium, Potassium, Chloride, Glucose, and SDMA A chemistry which includes TP, ALT, ALP, BUN, creatinine, and glucose. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	24 hours
SA071	Miniscreen 11 Chem with SDMA Total Protein, Albumin, Globulin, A/G Ratio, ALT (SGPT), Alkaline Phosphatase, T. Bilirubin, BUN, Creatinine, BUN/Creat Ratio, Sodium, Potassium, Chloride, Glucose, and SDMA A comprehensive chemistry including total protein, albumin, globulin, A/G ratio, ALT, ALP, bilirubin, BUN, creatinine, BUN/creatinine ratio, blood glucose and SDMA. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	24 hours

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA060	Miniscreen 4 Chem	0.5 mL serum in red top or serum separator tube	24 hours
	Total Protein, ALT (SGPT), Bun, Sodium, Potassium, Chloride, Glucose		
	Miniscreen chemistry that allows assessment of four analytes (TP, ALT, BUN, and glucose).		
	Interferences: marked hemolysis and lipemia.		
SA822	NSAID Chemistries with SDMA	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
	AST (SGOT), ALT (SGPT), Alkaline Phosphatase, BUN, Sodium, Potassium, Chloride, Creatinine with SDMA		
	A chemistry panel including AST, ALT, ALP, BUN, creatinine and SDMA.		
	Interferences: marked hemolysis and lipemia.		
T85360	Pre NSAID Use Panel 3 with SDMA	0.5 mL serum in red top or serum separator tube	24 hours
	Pre NSAID Panel 3 with Electrolytes and SDMA		
	A chemistry panel that includes albumin, AST, ALT, ALP, GGT, BUN, creatinine and SDMA.		
SA043	Pre-Op Chem with Electrolytes, SDMA	0.5 mL serum in red top or serum separator tube	24 hours
	Pre-Op Chemistries with Electrolytes and SDMA		
	A smaller chemistry with electrolyte assessment. Includes albumin, globulin, ALT, ALP, BUN, creatinine, glucose, sodium, potassium, chloride and SDMA.		
	Interferences: marked hemolysis and lipemia.		
T165 Add-on Equivalent ADD90	PrecisionPSL	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
	PrecisionPSL (Pancreatic Sensitive Lipase) is used to help diagnose pancreatitis in patients with consistent clinical signs.		
	Interferences: marked hemolysis or lipemia.		

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
ADD90	PrecisionPSL Add-on PrecisionPSL (Pancreatic Sensitive Lipase) is used to help diagnose acute pancreatitis in patients with consistent clinical signs. Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
SA040	Pre-Op Screen with SDMA Includes Total Protein, Albumin, Globulin, A/G Ratio, ALT (SGPT), Alkaline Phosphatase, BUN, Creatinine, BUN/Creat Ratio, Sodium, Potassium, Chloride, Na/K Ratio, Glucose, and SDMA A smaller chemistry which includes total protein, albumin, globulin, ALT, ALP, BUN, creatinine, glucose and SDMA. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
T240 Add-on Equivalent ADD130	Protein Electrophoresis (Serum) Total protein, Albumin, Globulin, Alpha 1, Alpha 2, Beta 1, and Gamma 1 fraction assessment with interpretation. An evaluation of the globulin fraction of the serum (alpha 1, alpha 2, beta, and gamma) to determine if the globulin fraction is monoclonal based on these components. Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	2-4 days
RECC	Recheck Profile with SDMA Total Protein, Albumin, Globulin, Albumin/Globulin Ratio, AST (SGOT), ALT (SGPT), Alkaline Phosphatase, GGT, Total Bilirubin, BUN, Creatinine, BUN/Creatinine Ratio, Phosphorous, Glucose, Calcium, Magnesium, Sodium, Potassium, Sodium/Potassium Ratio, Chloride, Cholesterol, Triglyceride, Amylase, PrecisionPSL, CPK, and SDMA	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
SA311	Renal Profile with SDMA Total Protein, Albumin, Globulin, Albumin/globulin Ratio, BUN, Creatinine, Phosphorous, Glucose, Calcium, Sodium, Potassium, Sodium/potassium Ratio, Chloride, and SDMA A chemistry panel which includes albumin, globulin, BUN, creatinine, phosphorous, glucose, calcium, sodium, potassium, chloride and SDMA. Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T1035	<p>SDMA</p> <p>SDMA (symmetric dimethylarginine) is freely filtered by the kidneys, and elevation in SDMA is indicative of a reduced glomerular filtration rate.</p>	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
SA010	<p>Superchem with SDMA</p> <p>Total Protein, Albumin, Globulin, Albumin/Globulin Ratio, AST (SGOT), ALT (SGPT), Alkaline Phosphatase, GGT, Total Bilirubin, BUN, Creatinine, BUN/Creatinine Ratio, Phosphorous, Glucose, Calcium, Magnesium, Sodium, Potassium, Sodium/Potassium Ratio, Chloride, Cholesterol, Triglyceride, Amylase, PrecisionPSL, CPK, and SDMA</p> <p>Most comprehensive chemistry profile which includes SDMA for estimation of the glomerular filtration rate (T1035).</p> <p>Interferences: marked hemolysis and lipemia.</p>	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
SA120C	<p>Superchem with SDMA, T4</p> <p>Superchem with SDMA and T4</p> <p>The most comprehensive chemistry with SDMA (SA010) and a total T4.</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
S85783	<p>Troponin I</p> <p>Interferences: marked hemolysis.</p>	1.0 mL frozen serum (should be fasting, non-hemolyzed sample)	3-5 days (includes shipping time to Reference Lab)
SA025	<p>Vet Screen with SDMA</p> <p>Total Protein, Albumin, Globulin, Albumin/Globulin Ratio, AST (SGOT), ALT (SGPT), Alkaline Phosphatase, GGT, Total Bilirubin, BUN, Creatinine, BUN/Creatinine Ratio, Phosphorous, Glucose, Calcium, Sodium, Potassium, Sodium/Potassium Ratio, Chloride, Cholesterol, CPK, and SDMA</p> <p>A comprehensive chemistry profile with SDMA.</p> <p>Interferences: marked hemolysis and lipemia.</p>	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA026	Vet Screen with Amylase, SDMA Vet Screen with SDMA and Amylase A comprehensive chemistry profile with SDMA (SA025) to which amylase evaluation has been added. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
SA1656	Vet Screen with SDMA, T4 Vet Screen with SDMA and T4 A comprehensive chemistry with SDMA (SA025) and total T4. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
S16016	Vitamin D	1.0 mL serum in red top or serum separator tube. Overnight fast recommended. Store and transport refrigerated or frozen (on ice packs)	4-10 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

URINE

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S16735	<p>Crystallographic Stone Analysis</p> <p>Nidus Composition, Stone Composition and Shell Composition</p> <p>The report includes a comprehensive stone analysis using optical crystallography with a region-by-region stone fracture. Stone composition is provided for nidus, stone, shell, and surface.</p> <p>Note: Do not send in formalin.</p>	Dry stone in red top tube or plain container	5-7 days
M134	<p>MIC Urine If Indicated</p> <p>Urinalysis with Culture, if indicated</p> <p>Includes a urinalysis followed by a culture if there is any of the following: glucosuria, any bacteria visible, greater than 4 wbcs per hpf, or a urine specific gravity less than 1.020. In the case of any of the above criteria the urine will be plated by standard plate methodology and evaluated.</p>	7.0 mL urine in urine transport tube	1-4 days
T764	<p>Path Review (Urinalysis)</p> <p>This test can be added onto a complete urinalysis and includes sediment evaluation by a clinical pathologist.</p>		1-2 days
ADD220	<p>Urinalysis Add-on</p> <p>Includes chemical analysis (Color, Appearance, Specific Gravity, pH, Protein, Glucose, Ketones, Bilirubin, Occult blood) and microscopic evaluation of urine sediment</p> <p>The complete urinalysis includes a physical (color, appearance, USG), chemical (pH, protein, glucose, ketones, bilirubin, occult blood), and microscopic (WBC, RBC, casts crystals, bacteria, transitional and squamous epithelial cells, fat droplets) exam of the urine.</p> <p>Interferences: visible levels of hemolysis, drugs containing dyes, nitrofurantoin, or riboflavin.</p>	6.0 mL urine in urine transport tube	24 hours Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T925	Urinalysis And UPC Ratio Urinalysis and Urine Protein/Creatinine Ratio A complete urinalysis (T760) and a urine protein to creatinine ratio (T775).	6.5 mL urine in urine transport tube	1-2 days
T760R	Urinalysis Recheck Includes chemical analysis (Color, Appearance, Specific Gravity, pH, Protein, Glucose, Ketones, Bilirubin, Occult blood) and microscopic evaluation of urine sediment A complete urinalysis which includes urine specific gravity, urine chemistry, and full microscopic evaluation as a follow up to a previous assessment within the last month. Interferences: visible levels of hemolysis, drugs containing dyes, nitrofurantoin, or riboflavin.	6.0 mL urine in urine transport tube	24 hours Performed each shift
T760 Add-on Equivalent ADD220	Urinalysis-Complete Includes chemical analysis (Color, Appearance, Specific Gravity, pH, Protein, Glucose, Ketones, Bilirubin, Occult blood) and microscopic evaluation of urine sediment The complete urinalysis includes a physical (color, appearance, USG), chemical (pH, protein, glucose, ketones, bilirubin, occult blood), and microscopic (WBC, RBC, casts crystals, bacteria, transitional and squamous epithelial cells, fat droplets) exam of the urine. Interferences: visible levels of hemolysis, drugs containing dyes, nitrofurantoin, or riboflavin.	6.0 mL urine in urine transport tube	24 hours Performed each shift
T227	Urine Bile Acid: Creatinine Ratio Evaluation of the amount of urinary bile acid present with respect to creatinine. Note: urine containing ascorbic acid may result in false positive results.	1.0 mL urine in urine transport tube	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T765	<p>Urine Clearance Ratios</p> <p>Urinary Clearance Creatinine Ratios (Phosphorous, Calcium, Sodium, Potassium, Chloride)</p> <p>Determination of urine fractional excretion is indicated to assess renal clearance of specific substances.</p> <p>Interferences: marked hemolysis or lipemia in serum.</p>	<p>0.5 mL urine and 0.5 mL serum</p> <p>Urine transport tube, serum in red top or serum separator tube</p>	1-2 days
T770	<p>Urine Cortisol/Creatinine Ratio</p> <p>Urine Cortisol, Urine Creatinine, and Urine Cortisol and Creatinine Ratio</p> <p>The urine cortisol/creatinine ratio is most commonly used as a screening test for canine hyperadrenocorticism.</p>	<p>0.5 mL urine in urine transport tube</p>	1-2 days
T830C	<p>Urine Microalbumin Canine</p> <p>Measures albumin concentration in the urine.</p>	<p>0.5 mL urine in urine transport tube</p>	24 hours
T830F	<p>Urine Microalbumin Feline</p> <p>Measures albumin concentration in the urine.</p>	<p>0.5 mL urine in urine transport tube</p>	24 hours

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
ADD230	Urine Protein/Creatinine Add-on Urine Protein, Urine Creatinine and Urine Protein and Creatinine Ratio Urine Protein/Creatinine Ratio (UPC) to determine the magnitude of the proteinuria.	0.5 mL urine in urine transport tube	1-2 days
T775 Add-on Equivalent ADD230	Urine Protein/Creatinine Ratio Urine Protein, Urine Creatinine and Urine Protein and Creatinine Ratio Urine Protein/Creatinine Ratio (UPC) to determine the magnitude of the proteinuria.	0.5 mL urine in urine transport tube	1-2 days
T310	Urine Uric Acid/Creatinine Urine Creatinine, Urine Uric Acid, and Urine Acid and Creatinine Ratio Used to monitor the success of preventive treatment for urate bladder stones in Dalmatian dogs. This test is not preferred to a 24-hour urine uric acid concentration. urine spot testing (urine urate: creatinine ratios) are too variable and do not correlate well with 24-hour urine urate excretion.	0.5 mL urine in urine transport tube	1-2 days
M133	UTI Cystitis Profile with MIC Urine Urinalysis-Complete with Urine Culture A complete urinalysis (T760), and urine culture (M130). Culture performed using <i>FIRST</i> Tract methodology allowing for rapid growth within 5 hours and subsequent reporting of a positive or negative result in under 24 hours. A positive result will be immediately plated for bacterial identification and antimicrobial sensitivity testing. Interferences: patients should be off antibiotics for seven to 10 days prior to urine culture. Samples obtained via free catch and catheterization need to be interpreted carefully in light of the method of collection. Cystocentesis is the preferred sampling method for culture.	6.5 mL urine in urine transport tube	1-4 days

ENDOCRINOLOGY

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
ADD260	Add-on Fructosamine Serum fructosamine concentrations reflect the mean blood glucose concentrations during the preceding one to two weeks. The measurement can be used to differentiate stress hyperglycemia from diabetes mellitus. Interferences: hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
S14410	Anti-Mullerian Hormone (C/F) A qualitative AMH test for dogs and cats is meant to distinguish between spayed and ovarian intact bitches and queens once they have reached the age of sexual maturity. This test is to be used on canine and feline species only. Interferences: marked hemolysis or lipemia. Note: serum should be separated and placed in a plain, non-additive tube. Sample must be kept cold — refrigerate or freeze after collection and send on ice packs.	0.5 mL serum in non-additive tube (plastic preferred). Avoid use of serum separator tubes due to possible assay interference.	5-10 days
ADD300	Canine Post Pill T4 Add-on Used for monitoring thyroid status in patients previously diagnosed as hypothyroid, being supplemented on synthetic thyroid hormone (4-6 hours post pill). Note: refrigerate and send on ice. Ideal sampling time is typically 4-6 hours post pill administration.	0.5 mL serum in red top or serum separator tube	1-2 days
T445	Cortisol Cortisol Assay, Single sample Single cortisol level Interferences: marked hemolysis may affect results, and marked lipemia (if unable to be cleared by centrifugation) may falsely decrease cortisol levels.	0.5 mL serum in red top or serum separator tube	1-2 days
ACTH2	Cortisol Serial ACTH 2 Cortisol Assay, 2 Samples: Pre (Baseline) and Post ACTH Stimulation This is an ACTH stimulation test that requires two samples: one pre, one post. Results include baseline cortisol and one post-ACTH cortisol. Interferences: marked hemolysis may affect results, and marked lipemia (if unable to be cleared by centrifugation) may falsely decrease cortisol levels.	0.5 mL serum in red top or serum separator tube for baseline sample (labeled Pre) and 0.5 serum in red top or serum separator tube for post sample (labeled Post)	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
ACTH3	<p>Cortisol Serial ACTH 3</p> <p>Cortisol Assay, 3 Samples: Pre (Baseline) and 2 Post ACTH Stimulation</p> <p>This is an ACTH stimulation test that requires three samples. One pre cortrosyn administration and then two post at established intervals.</p> <p>Interferences: marked hemolysis may affect results, and marked lipemia (if unable to be cleared by centrifugation) may falsely decrease cortisol levels.</p>	<p>0.5 mL serum in red top or serum separator tube for baselin esample (labeled Pre) and 0.5 serum in red top or serum separator tube for each post sample (post samples labeled Post)</p>	1-2 days
ACTH4	<p>Cortisol Serial ACTH 4</p> <p>Cortisol Assay, 4 Samples: Pre (Baseline) and 3 Post ACTH Stimulation</p> <p>This is an ACTH stimulation test that requires four samples. On pre cortrosyn administration and then 3 post at set intervals.</p> <p>Interferences: marked hemolysis may affect results, and marked lipemia (if unable to be cleared by centrifugation) may falsely decrease cortisol levels.</p>	<p>0.5 mL serum in red top or serum separator tube for baseline sample (labeled Pre) and 0.5 serum in red top or serum separator tube for each post sample (post samples labeled Post)</p>	1-2 days
DEX2	<p>Cortisol Serial DEX 2</p> <p>Cortisol Assay, 2 Samples: Pre (Baseline) and Post Dexamethasone Suppression</p> <p>This test measures cortisol concentrations before and after administration of exogenous dexamethasone (0 hours followed by an 8-hour sample). This is only a two cortisol assay. The low dose dexamethasone suppression test, a three cortisol sample test (mnemonic CDEX3), is preferred for diagnosing hyperadrenocorticism. The high-dose dexamethasone suppression test, a three cortisol sample test, is used to differentiate pituitary-dependent from adrenal tumor hyperadrenocorticism.</p> <p>Interferences: marked hemolysis may affect results, and marked lipemia (unable to be cleared by centrifugation) may falsely decrease cortisol levels.</p>	<p>0.5 mL serum in red top or serum separator tube for baseline sample (labeled Pre) and 0.5 serum in red top or serum separator tube for post sample (labeled Post)</p>	1-2 days
DEX3	<p>Cortisol Serial DEX 3</p> <p>Cortisol Assay, 3 Samples: Pre (Baseline) and 2 Post Dexamethasone Suppression</p> <p>This test measures cortisol concentrations before and after administration of exogenous dexamethasone. This suppression test requires three samples; baseline cortisol prior to dexamethasone administration and two cortisol levels post dexamethasone administration (4 and 8 hours). The low-dose dexamethasone suppression test is preferred for diagnosing hyperadrenocorticism. The high-dose dexamethasone suppression test is used to differentiate pituitary-dependent from adrenal tumor hyperadrenocorticism.</p> <p>Interferences: marked hemolysis may affect results, and marked lipemia (unable to be cleared by centrifugation) may falsely decrease cortisol levels.</p>	<p>0.5 mL serum in red top or serum separator tube for baseline sample (labeled Pre) and 0.5 serum in red top or serum separator tube for each post sample (post samples labeled Post)</p>	1-2 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
DEX4	Cortisol Serial DEX 4 Cortisol Assay, 4 Samples: Pre (Baseline) and 3 Post Dexamethasone Suppression This test measures cortisol concentrations before and after administration of exogenous dexamethasone. This suppression test includes four samples; baseline cortisol prior to dexamethasone administration and three additional cortisol samples post dexamethasone administration (last occurring at 8 hours). Samples taken at 2 or 6 hours are interpreted similarly to those obtained at 4 hours. Interferences: marked hemolysis may affect results, and marked lipemia (unable to be cleared by centrifugation) may falsely decrease cortisol levels.	0.5 mL serum in red top or serum separator tube for baselin esample (labeled Pre) and 0.5 serum in red top or serum separator tube for each post sample (post samples labeled Post)	1-2 days
T435	Endogenous ACTH Endogenous ACTH used for small animal testing. Note: endogenous ACTH is extremely labile. Immediately separated and frozen EDTA-plasma without aprotinin may be acceptable but is NOT the preferred sample. Serum or heparinized plasma is not acceptable for endogenous ACTH testing. Call lab to obtain an aprotinin-treated lavender top tube and transfer tube.	1.0 mL aprotinin treated EDTA plasma in non-additive transport tube (labeled as AP treated plasma). Transport refrigerated (on ice packs)	2-3 days
T460 Add-on Equivalent ADD50	Free T4 By Equilibrium Dialysis Note: this test should not be performed as an add-on to samples older than 72 hours.	0.5 mL serum in red top or serum separator tube	2-3 days
ADD50	Free T4 Equilibrium Dialysis Add-on	0.5 mL serum in red top or serum separator tube	2-3 days
S16345 Add-on Equivalent ADD260	Fructosamine Assay Fructosamine Interferences: hemolysis and lipemia may affect results.	0.5 mL serum in red top or serum separator tube	24 hours
T497 Add-on Equivalent ADD300	Hypothyroid Post Pill T4 Used for monitoring thyroid status in patients previously diagnosed as hypothyroid, being supplemented on synthetic thyroid hormone (4-6 hours post pill). Note: draw serum 4-6 hours post pill.	0.5 mL serum in red top or serum separator tube	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T470	Insulin with Glucose Interferences: marked hemolysis and lipemia. Note: generally used to evaluate for the presence of an insulinoma. Concurrent blood glucose concentration should be below 60 mg/dL (3.3 mmol/L) when serum insulin levels are measured. Not to be used on patients receiving exogenous insulin without an appropriate withdrawal period.	0.5 mL serum in red top or spun serum separator tube	1-2 days
T470F	Insulin Feline with Glucose Submission of a single serum sample where a blood glucose (ideal BG less than 50 mg/dL or 2.8 mmol/L) is assessed along side an insulin level. Interferences: marked hemolysis. Note: generally used to evaluate for the presence of an insulinoma. Concurrent blood glucose concentration should be below 60 mg/dL (3.3 mmol/L) when serum insulin levels are measured. Not to be used on patients receiving exogenous insulin without an appropriate withdrawal period.	0.5 mL serum in red top or serum separator tube	3-7 days
S16520	Luteinizing Hormone Interferences: marked hemolysis or lipemia. Store refrigerated up to 24 hours. Freeze for longer storage.	0.5 mL serum in red top or serum separator tube	1-3 days
S86698	Malignancy Profile Ionized Calcium, Parathyroid Hormone, and Parathormone Related Protein Evaluation of the ionized fraction of calcium concurrent with parathyroid hormone and parathyroid hormone related protein.	Ionized calcium: 2.0 mL serum anaerobically transferred from spun red top or serum separator tube into a plain, unopened red top tube. PTH: 1 mL frozen serum in plain transport tube. PTHrP: 1.0 mL frozen EDTA plasma in plain transport tube. Include patient history with submission.	7-10 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA1528	OFA Thyroid Panel OFA, OFA Interpretation, Free T4 by ED, TSH, Thyroglobulin Auto Antibody The Orthopedic Foundation for Animals (OFA) thyroid panel for dogs includes a FT4 by ED, thyroid stimulation hormone and thyroglobulin auto antibody. Include OFA form in submission. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	2.0 mL serum in red top tube. Do not use serum separator tube for submission.	OFA will issue report (with results) directly to submitting client
S16596	Parathormone Related Protein Determines the presence of the parathormone-related protein (PTHrp). Interferences: marked hemolysis or lipemia Note: submit frozen plasma and label submisson tube as plasma PTHrP.	1.0 mL frozen EDTA plasma in plain transport tube, collected as whole blood in lavender top tube (spin and separate plasma into plain transport tube). Include patient history. Frozen plasma (Call customer service for specimen handling protocol: 1-800-872-1001, dial 0)	7-14 days
T448	Post Cortisol A single cortisol level, submitted as a post, used to complete an ACTH stimulation or dexamethasone suppression test (8 hours). Interferences: marked lipemia.	0.5 serum in red top or serum separator tube (labeled Post)	1-2 days
T475	Progesterone Progesterone level Interferences: serum separator gel may interfere with test.	0.5 mL serum in red top tube. Submission of serum separator tube is not recommended.	1-2 days
S16595	PTH, Ionized Calcium Evaluates the total calcium, ionized calcium, and parathyroid hormone to further define the cause of abnormalities in calcium homeostasis. Note: serum samples must be anaerobically obtained (exposure to air may artifactually decrease ionized calcium). Samples should not be transported in serum separator tubes. Transfer serum into two plain red top tubes and label as serum ionized calcium and serum PTH. Samples should be frozen and sent on ice.	2.0 mL serum anaerobically transferred from spun red top or serum separator tube into a plain, unopened red top tube, stored and transported frozen along with patient history	7-10 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T480	T3 Tri-iodothyronine (T3)	0.5 mL serum in red top or serum separator tube	1-2 days
SA430	T3 Suppression T4, T3, Post Pill T3 and Post Pill T4 The submission of two separate serum samples (baseline T3/T4 and a post T3/T4). The second sample (post T3/T4) is obtained following a specific protocol of oral liothyronine sodium administration.	0.5 mL serum (labeled Pre) and 0.5 serum (labeled Post) in red top or serum separator tubes	1-2 days
T495 Add-on Equivalent ADD190	T4 Total T4 Interferences: marked hemolysis and moderate to marked lipemia. Lipemia can falsely decrease T4 results.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
ADD190	T4 Add-on Total T4 as an add-on test to a previous profile. Interferences: marked hemolysis and moderate to marked lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
SA401	T4, TSH Thyroid panel including a total T4 and thyroid stimulating hormone (TSH). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum in red top or serum separator tube	1-2 days
S16760	Testosterone Single testosterone level	0.5 mL serum in red top or serum separator tube	3-5 days
T505	Thyroglobulin Auto Antibody Canine Thyroglobulin Auto Antibody Level Thyroglobulin auto antibody is used to confirm the presence of autoimmune thyroiditis.	0.5 mL serum in red top or serum separator tube	1-5 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA360	Thyroid 1 T4 and T3 Used to evaluate T4 and T3 levels. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.75 mL serum in red top or serum separator tube	1-2 days
SA370	Thyroid Profile 2 T4, Free T4 by Equilibrium Dialysis Panel includes a total T4 and free T4 by Equilibrium Dialysis (FT4 ED). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum in red top or serum separator tube	1-3 days
SA380	Thyroid Profile 3 Total T4, Free T4 By Equilibrium Dialysis, and TSH Panel includes a total T4, free T4 by Equilibrium Dialysis (FT4 ED), and thyroid stimulating hormone (TSH). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum in red top or serum separator tube	1-3 days
SA400	Thyroid Profile 5 T4, Free T4 By Equilibrium Dialysis, TSH, and Thyroglobulin Auto Antibody Panel includes a total T4, free T4 by Equilibrium Dialysis (FT4 ED), thyroid stimulating hormone (TSH), and thyroglobulin autoantibody (TGAA). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	2.0 mL serum in red top or serum separator tube	1-5 days
SA410	Thyroid Profile 6 T3, T4, Free T4 By Equilibrium Dialysis, T3 Autoantibodies, T4 Autoantibodies, TSH, and Thyroglobulin Auto Antibody Panel includes total T3, T3 autoantibody, total T4, free T4 by Equilibrium Dialysis (FT4 ED), T4 autoantibody, thyroid stimulating hormone (TSH), and thyroglobulin autoantibody (TGAA). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	2.5 mL serum in red top or serum separator tube	1-5 days
T510 Add-on Equivalent ADD200	TSH Thyroid stimulating hormone (TSH) Interferences: marked lipemia.	0.5 mL serum in red top or serum separator tube	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
ADD200	TSH Add-on Thyroid stimulating hormone (TSH) Interferences: marked lipemia.	0.5 mL serum in red top or serum separator tube	1-2 days
SA390	TSH, Free T4 ED TSH, Free T4 ED by Equilibrium Dialysis Panel that includes a free T4 by Equilibrium Dialysis (FT4 ED) and thyroid stimulating hormone (TSH). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum in red top or serum separator tube	1-3 days
Add-on Equivalent SA390ADD			

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

INFECTIOUS

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
AC100	Accuplex®	0.5 mL serum in red top or serum separator tube	1-2 days
Add-on Equivalent ADD111	Heartworm, <i>Borrelia burgdorferi</i> , <i>E. canis</i> , <i>Anaplasma phagocytophilum</i> This is a canine vector-borne disease screening test for Heartworm (Ag), Lyme disease (includes screening for antibodies against two C6 peptides), <i>E. canis</i> , and <i>A. phagocytophilum</i> .		
AC617	Accuplex® with Microfilaria Knotts Includes a canine vector-borne disease screening test for Heartworm (Ag), Lyme disease (includes screening for antibodies against two C6 peptides), <i>E. canis</i> , and <i>A. phagocytophilum</i> (AC100) to which a Knotts is added to detect Microfilaria.	0.5 mL serum, and EDTA 1.0 mL whole blood Serum in red top or serum separator tube, lavender top	1-2 days
S16872	<i>Anaplasma phagocytophilum</i> Ab Used to detect antibodies against <i>Anaplasma phagocytophilum</i> by IFA methodology. Exposure does not necessarily indicate that clinical signs are caused by infection. Canine only.	1.0 mL serum in red top or serum separator tube	1-4 days
S16070	<i>Babesia canis</i> Detects the presence of antibody to <i>Babesia canis</i> . Performed by IFA.	0.5 mL serum in red top or serum separator tube	7-10 days
S16075	<i>Babesia gibsoni</i> Detects the presence of antibody to <i>Babesia gibsoni</i> . Performed by IFA.	0.5 mL serum in red top or serum separator tube	7-10 days
S16502	<i>Babesia gibsoni</i> IFA, Export This test evaluates the presence of antibodies to <i>Babesia gibsoni</i> , as required for export to specific countries. Ensure export test requirements prior to test submission.	1.0 mL serum in red top or serum separator tube	7-10 days
S85889	<i>Bartonella henselae</i>, ELISA Evaluation of antibody indicating previous exposure to <i>Bartonella henselae</i> or <i>Bartonella clarridgeiae</i> (ELISA methodology). Note: recently infected cats may not have seroconverted. Titer magnitude does not prove current infection or relevant clinical disease. Active infection is suggested by a rising titer (a four-fold increase over two weeks).	0.5 mL serum in red top or serum separator tube	5-7 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T530	Brucella canis - Screen This test is species-specific (dogs only). This is a screening test performed at 1:50 dilution. A titer is not performed. Note: canine only. Test will not detect antibodies against <i>B. suis</i> , <i>B. melitensi</i> , and <i>B. abortus</i> . Use the <i>Brucella</i> AGID test if a dog may have been infected with <i>Brucella</i> from other species (e.g., from cattle, sheep, goats, pigs, sheep). Do not use this test for export.	0.5 mL serum in red top tube or serum separator tube	1-3 days
S16131	Brucella Screen Multiplex <i>Brucella</i> Screen Multiplex. Includes confirmatory testing with <i>Brucella</i> Slide Agglutination and AGID for samples with non-negative screen results An antibody detection for <i>Babesia canis</i> antigen via the Canine <i>Brucella</i> Multiplex Assay. Samples undergo sequential testing dependent upon initial results. If the Multiplex assay is positive, further testing is performed by Canine <i>Brucella</i> Slide Agglutination and AGID testing at no additional cost. Note: canine only. This test can detect antibodies against <i>B. canis</i> , <i>B. suis</i> , <i>B. melitensi</i> and <i>B. abortus</i> .	1.0 mL serum in red top tube or serum separator tube	7-10 days
SA265	C5 Cat Viral FeLV Antigen ELISA, Feline Coronavirus Titer, FIV Antibody Panel that assesses possible infectious etiologies including FeLV (antigen detection), FIV (antibody detection), and feline coronavirus (antibody detection) for the exhibited clinical signs. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	0.75 mL serum in red top tube or serum separator tube	1-3 days
S16135	Calici Virus - IFA Calici Virus A serology test evaluating for IgG antibody levels by IFA suggesting Feline calicivirus exposure. Interferences: sample submitted in tubes with separator gel.	0.5 mL serum in red top tube. Submission in serum separator tube is not recommended.	10-12 days
T555	Canine Distemper Titer IgG, IgM This test provides an IgM and IgG titer for Canine Distemper. Titer pattern needs to be interpreted in light of the patient (duration of illness, clinical signs, vaccine history). Paired serology may be needed to confirm recent exposure. Interferences: marked hemolysis or marked lipemia if unable to clear by centrifugation.	0.5 mL serum in red top or serum separator tube	1-3 days
ADD04	Cocci AGID Add-on Serology by agar gel immunodiffusion (AGID) used most frequently to screen for exposure to <i>Coccidioides immitis</i> . This is an antibody test evaluating both IgM and IgG antibody. Interferences: lipemia.	0.5 mL serum in red top or serum separator tube	3-5 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T535 Add-on Equivalent ADD04	Coccidioidomycosis, Screen and Titer Serology by agar gel immunodiffusion is used most frequently to screen for exposure to <i>Coccidioides</i> . This is a titer test for antibodies directed against <i>Coccidioides immitis</i> (aka Valley Fever). Interferences: lipemia.	0.5 mL serum in red top or serum separator tube	3-5 days
T550	Cryptococcal Antigen Used as a diagnostic test for Cryptococcosis. Detects the presence of <i>Cryptococcal</i> antigen. Interferences: marked hemolysis and lipemia. Note: serum and CSF submitted in a plain tube (white top or red top) is the only acceptable specimens. Plasma (any type) is not an acceptable specimen for this test. Body cavity fluid (e.g., CSF) cannot be submitted in a tube containing anticoagulant (e.g., LT) because the anticoagulant may invalidate the test.	0.5 mL serum or fluid in red top tube or plain (non-additive) transport tube	2-3 days
T570 Add-on Equivalent ADD05	Ehrlichia canis This test detects antibodies (IgG) directed against <i>Ehrlichia canis</i> . Interferences: marked hemolysis and lipemia. Canine only. Test is not valid for export purposes.	0.5 mL serum in red top or serum separator tube	1-2 days
ADD05	Ehrlichia canis Add-on Antibodies for <i>Ehrlichia canis</i> evaluated for by IFA. Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	1-2 days
S16900	Ehrlichiosis Serology Panel Canine <i>Ehrlichia canis</i> , <i>Anaplasma phagocytophilum</i> AB, <i>Neorickettsia risticii</i> AB Evaluates for the presence of antibodies indicating exposure to <i>Ehrlichia canis</i> , <i>Anaplasma phagocytophilum</i> , and <i>Neorickettsia risticii</i> . Interferences: marked hemolysis and lipemia.	2.0 mL serum in red top or serum separator tube	1-7 days
T812	Fecal Combo with Heartworm Heartworm Antigen, Fecal O&P, Centrifugation, <i>Giardia</i> Detection of heartworm antigen and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> ELISA for ova and parasite detection (T808). Note: sample should be evaluated within 24 hours of collection. If a worm has been identified in the sample, separate the worm and place in a container labeled Worm in black marker. Additionally, indicate Worm on TRF type included with submission.	0.5 mL serum and 6 grams feces Serum in red top or serum separator tube, Antech provided fecal container	1-2 days
T593	Feline Coronavirus Exposure Titer This test is run at a dilution of 1:25. The purpose of this test is to determine if a cat has any evidence of exposure to any coronavirus (FIP or enteric coronavirus).	0.5 mL serum in red top or serum separator tube	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T595	Feline Coronavirus Titer This test determines if a cat has antibody titers against Feline Coronavirus (FCV). Titers are determined at 1:400 and 1:1600. A result of less than 1:400 is not the same as a negative titer. A titer of less than 1:400 indicates that the cat does not have high antibody concentrations against FCV. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	1-3 days
SA270	Feline Serology 2 FeLV Antigen ELISA, FIV Antibody, Feline Coronavirus Titer, <i>Cryptococcal</i> Antigen, <i>Toxoplasma</i> Ab - IgG/IgM Panel that assesses possible infectious etiologies including FeLV (T580), FIV (T610), feline coronavirus (T595), <i>Cryptococcus neoformans</i> (T550), and <i>Toxoplasma gondii</i> (T720) for the exhibited clinical signs. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum in red top or serum separator tube	1-3 days
SA262	Feline Serology, FeLV, FIV, Toxo AB IgG/IgM FeLV Antigen ELISA, FIV Antibody, <i>Toxoplasma</i> Ab - IgG/IgM Panel that assesses possible infectious etiologies including FeLV (T580), FIV (T610), and <i>Toxoplasma gondii</i> (T720) for the exhibited clinical signs. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum in red top or serum separator tube	1-2 days
T580 Add-on Equivalent ADD06	FeLV Antigen ELISA Test detects FeLV antigen by ELISA methodology. Interferences: marked hemolysis or lipemia may cause false positive results.	0.5 mL serum in red top or serum separator tube	1-2 days
ADD06	FeLV Antigen ELISA Add-on Test detects FeLV antigen by ELISA methodology. Interferences: marked hemolysis or lipemia. marked hemolysis or lipemia may result in false positive results.	0.5 mL serum in red top or serum separator tube	1-2 days
T585	FeLV IFA This is an immunofluorescence assay run on blood or bone marrow smears to look for evidence that Feline Leukemia Virus (FeLV) has infected the cat's bone marrow. Note: requires unstained blood or bone marrow smears.	2 freshly prepared, unstained blood smears or 2 bone marrow smears.	1-3 days
S6234	FeLV PCR Used to detect FeLV using PCR in cases where ELISA and IFA are inconclusive. Results are reported as positive or negative. Note: samples should be refrigerated	1 to 5 grams bone marrow in leakproof container OR 1.0 mL EDTA whole blood in lavender top tube	7-10 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
ADD07	FeLV, FIV ELISA Add-on FeLV Antigen ELISA, FIV AB Evaluation for the presence of FeLV antigen (ELISA) and FIV antibodies (IFA). Interferences: marked hemolysis and lipemia. marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	0.5 mL serum in red top or serum separator tube	1-2 days
SA260 Add-on Equivalent ADD07	FeLV, FIV Special FeLV Antigen ELISA, FIV Antibody Evaluation for the presence of FeLV antigen (ELISA) and FIV antibodies (IFA). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	0.5 mL serum in red top or serum separator tube	1-2 days
SA266	FeLV, FIV, Feline Heartworm Antigen FeLV Antigen ELISA, FIV Antibody, Heartworm Antigen, Feline Evaluation for the presence of FeLV antigen (ELISA), antibodies to FIV (IFA), and heartworm antigen. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	0.5 mL serum in red top or serum separator tube	1-2 days
SA269	FeLV, FIV Heartworm Ab FeLV Antigen ELISA, FIV Antibody, Heartworm Antibody, Feline Evaluation for the presence of FeLV antigen (ELISA), FIV (IFA) and heartworm antibody. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	0.5 mL serum in red top or serum separator tube	1-2 days
T610 Add-on Equivalent ADD15	FIV Antibody Evaluation for the presence of FIV antibody by IFA.	0.5 mL serum in red top or serum separator tube	1-2 days
ADD15	FIV Antibody Add-on Evaluation for the presence of FIV antibody by IFA.	0.5 mL serum in red top or serum separator tube	1-2 days
T613	Heartworm Ag-Heat-Treated Serum Heartworm Antigen using heat-treated serum This test is to be done on HW negative samples that are thought to be false negative because of antigen/antibody complexes. This test is NOT recommended for samples that are already HW antigen-positive or borderline positive. Interferences: marked hemolysis or lipemia. Marked hemolysis may cause false positive result. Marked lipemia may cause false negative results.	1.0 mL serum in red top or serum separator tube	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T625	Heartworm Antibody, Feline Heartworm Antibody, Feline This test evaluates for antibodies consistent with heartworm exposure. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	1-2 days
T630	Heartworm Antibody, Antigen Feline This test evaluates for antibodies consistent with heartworm exposure and antigens consistent with the presence of female adult heartworm. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	1-2 days
T615 Add-on Equivalent ADD70	Heartworm Antigen This test evaluates for the presence of heartworm antigen. Interferences: marked hemolysis and lipemia.	1.0 mL serum in red top or serum separator tube	1-2 days
T620	Heartworm Antigen This test evaluates for the presence of heartworm antigen. Interferences: marked hemolysis and lipemia.	1.0 mL serum in red top or serum separator tube	1-2 days
T617	Heartworm Antigen, Microfilaria Microfilaria - Knotts, Heartworm Antigen This test evaluates for the presence of heartworm antigen (female adult heartworm) and a Knotts test is performed to evaluate for microfilariae. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube and 1.0 mL EDTA whole blood in lavender top tube	1-2 days
ADD70	Heartworm Antigen Add-on This test evaluates for the presence of heartworm antigen. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	1-2 days
T618	Heartworm, Ova and Parasite Heartworm Antigen, Fecal O&P with Centrifugation Heartworm antigen detection and fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection (T805).	0.5 mL serum and 5 grams feces Serum in red top or serum separator tube, Antech provided fecal container	1-2 days
T380 Add-on Equivalent ADD20	Hemotropic <i>Mycoplasma</i> <i>Mycoplasma haemofelis</i> Microscopic evaluation of a fresh blood smear for hemotropic <i>Mycoplasma</i> . Interferences: marked hemolysis.	1.0 mL EDTA whole blood in lavender top tube	24 hours

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
ADD20	Hemotropic <i>Mycoplasma</i> Add-on Microscopic evaluation of a fresh blood smear for hemotropic <i>Mycoplasma</i> . Interferences: marked hemolysis.	1.0 mL EDTA whole blood in lavender top tube	24 hours
T640	<i>Histoplasma</i> Antibody Evaluates for antibodies consistent with exposure to <i>Histoplasma</i> spp. (AGID). Results are reported as positive or negative. Interferences: marked hemolysis and lipemia. Note: the <i>Histoplasma</i> antigen test is preferred to the antibody test for serologic diagnosis of histoplasmosis.	0.5 mL serum in red top or serum separator tube	2-3 days
S86569	<i>Histoplasma capsulatum</i> Ag, EIA	2.0 mL urine in urine transport tube	5-7 days
S16405	<i>Histoplasma</i> Titer Detection of antibody indicating exposure to <i>Histoplasma</i> spp. Titer provided.	1.0 mL serum in red top or serum separator tube	9-12 days
S86096	Influenza - Canine Acute H3N8 and H3N2 influenza antibody detection Used to detect for the presence of antibody to H3N8 and H3N2 influenza viruses.	1.0 mL serum in red top tube	10-14 days
S16510	Leptospirosis <i>L. pomona</i> , <i>L. icterohemorrhagiae</i> , <i>L. canicola</i> , <i>L. grippotyphossa</i> , <i>L. hardjo</i> , <i>L. autumnalis</i> , <i>L. bratislava</i> The results of this test include a semi-quantitative titer for serovars. Interferences: recent vaccination (within one month) may interfere with testing. Note: test is not species specific.	1.0 mL serum in red top tube or serum separator tube	3-5 days
T670	Lyme Titer IgG Detection of IgG consistent with exposure (natural or vaccinal) to <i>Borrelia burgdorferi</i> . Interferences: marked hemolysis and lipemia. Note: test does not distinguish between natural and vaccinal Lyme exposure.	0.5 mL serum in red top or serum separator tube	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T390	Microfilaria - Knotts This is a concentration test involving whole blood that is more sensitive than the examination of a peripheral blood smear to determine the presence of microfilariae. Interferences: marked hemolysis.	1.0 mL EDTA whole blood in lavender top tube	1-2 days
S16560	Neospora Caninum - IFA Results are expressed as an antibody titer consistent with exposure to <i>Neospora caninum</i> .	1.0 mL serum in red top or serum separator tube	7-10 days
T695	Parvovirus Antigen Detects parvovirus antigen in fecal samples. Note: modified live virus vaccination for parvovirus may give false-positive results for about two weeks after vaccination. A negative test does not rule out parvovirus infection.	5 grams feces in Antech provided fecal container	1-2 days
S1204	Rabies Diagnostic Non-Export Rabies Antibody Titer, Non-export Results are reported as an antibody titer. Note: not to be used for export. In the case of export, the individual country's guidelines for testing need to be reviewed. With regards to Rabies vaccination guidelines and the utility of titers, the clinician is recommended to refer to state/provincial guidelines for vaccination requirements.	1.5 mL serum in red top or serum separator tube	3-4 weeks
S17108	Rabies Export - FAVN Rabies Antibody Titer, Export by FAVN Evaluation of rabies titer by fluorescent antibody virus neutralization. Interferences: marked hemolysis or lipemia. Note: complete FAVN form. For the post rabies vaccination, wait at least 3-4 weeks prior to drawing the sample.	1.5 mL serum in red top or serum separator tube (with FAVN form)	3-4 weeks
T715	Rocky Mt Spotted Fever IFA Evaluates for IgG antibody consistent with exposure to <i>Rickettsia rickettsii</i> . Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	1-2 days
T720	Toxoplasma Ab - IgG/IgM (Feline) <i>Toxoplasma</i> Ab - IgG/IgM (Feline) Detection of IgG and IgM antibodies to <i>Toxoplasma gondii</i> . Interferences: marked hemolysis and lipemia. This code is for domestic cats only. Alternative codes are available for canines and exotic animals (including felines).	0.5 mL serum in red top or serum separator tube	1-2 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S18708	<i>Toxoplasma gondii</i> PCR Used to evaluate for the presence of <i>Toxoplasma gondii</i> in cases where toxoplasmosis is suspected.	0.5 mL EDTA whole blood OR transtracheal wash, OR BAL, OR lymph node aspirate	7-10 days
S85030	Toxoplasmosis IgG/IgM <i>Toxoplasma</i> IgG, <i>Toxoplasma</i> IgM Detection of IgG and IgM antibodies to <i>Toxoplasma gondii</i> . Note: when submitting CSF, it is recommended that serum be submitted as well.	1.0 mL serum (canine) in red top or serum separator tube	5-7 days
S85819	<i>Tritrichomonas</i> PCR Used to evaluate for the presence of <i>Tritrichomonas foetus</i> in the feces of a patient with the appropriate clinical signs. Note: feces should be diarrheic and not be contaminated with cat litter. Fecal samples collected using a fecal loop are preferred.	0.3 grams feces in Antech provided container	1-2 days
S16581	Vaccine Panel Panleukopenia, Rhino, Calici Feline Panleukopenia, Rhinotracheitis, and Calicivirus Vaccinal Titers This is a semi-quantitative titer that is correlated to the amount of protective humoral antibody present in relation to previous feline panleukopenia (S16053), rhinotracheitis (S16702), and calicivirus vaccination (S16112).	2.5 mL serum in red top or serum separator tube	3-14 days
T820	<i>Giardia</i> ELISA The <i>Giardia</i> antigen capture Enzyme- linked Immunosorbent Assay (ELISA) is a fecal procedure designed to detect a <i>Giardia</i> specific antigen.	5 grams feces in Antech provided fecal container	1-2 days
Add-on Equivalent ADD250			
ADD250	<i>Giardia</i> ELISA Add-on The <i>Giardia</i> antigen capture Enzyme- linked Immunosorbent Assay (ELISA) is a fecal procedure designed to detect a <i>Giardia</i> specific antigen. Interferences: <i>Giardia</i> ELISA should be interpreted in conjunction with fecal flotation/centrifugation and clinical signs.	5 grams feces in Antech provided fecal container	1-2 days
T806	O&P Centrifugation, Smear Fecal ova and parasites and direct smear evaluation. O&P samples are appropriately mixed with zinc sulfate solution, centrifuged, followed by flotation and slide evaluation. Fecal direct smear evaluation is used to check for adult protozoal organisms. Note: sample should be evaluated within 24 hours of collection for O&P. For protozoal organisms, smears should be fresh and are best evaluated in-house immediately after sample collection. If a worm has been identified in the sample, separate the worm and place it in a container labeled Worm in black. Additionally, mark Worm on the TRF being submitted with the sample.	5 grams feces in Antech provided fecal container	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T805	<p>Ova and Parasite</p> <p>Fecal ova and parasite detection. Samples are appropriately mixed with zinc sulfate solution, centrifuged, followed by flotation and slide evaluation.</p> <p>Note: sample should be evaluated within 24 hours of collection. If a worm has been identified in the sample, separate the worm and place it in a container labeled Worm in black, additionally mark Worm on the TRF being submitted with the sample.</p>	<p>5 grams feces in Antech provided fecal container</p>	<p>1-2 days</p>
T808	<p>Ova and Parasite with <i>Giardia</i></p> <p>Fecal ova and parasites and <i>Giardia</i> Ag detection by ELISA. O&P samples are appropriately mixed with zinc sulfate solution, centrifuged, followed by flotation and slide evaluation. <i>Giardia</i> ELISA is considered more sensitive than a single floatation test.</p> <p>Note: sample should be evaluated within 24 hours of collection. If a worm has been identified in the sample, separate the worm and place in a container labeled Worm in black. Additionally, mark Worm on the TRF being submitted with the sample.</p>	<p>6 grams feces in Antech provided fecal container</p>	<p>1-2 days</p>

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

MICROBIOLOGY

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
M010	Acid Fast Stain - Micro	1 air-dried smear	1-3 days
	Acid-fast stain applied to an air-dried smear to determine if mycobacteria present in sample provided.		
	Note: air dried smear is the preferred sample submission.		
M040	Culture Aerobic and Anaerobic	Culturette (aerobic and anaerobic) from fluid (body cavity fluids), TTW, BAL, wound, lesion, or skin. Tissue sample in saline.	3-5 days
	This test is used when a bacterial infection is suspected in a tissue or fluid, but it is uncertain whether an aerobic or an anaerobic organism is the cause.		
	Note: fluids should be submitted in a red top tube. Tissue should be submitted in an air-tight container for anaerobic culture. Samples should be refrigerated prior to submission. Sensitivities are not performed, but drugs of choice for anaerobes are provided as a guide.		
M050	Culture Aerobic and Fungal	Aerobic culture: Culturette from body fluid, TTW, BAL, wound, lesion, or skin; tissue in saline. Fungal culture: Dry hair and skin scraping collected in a sterile red top or other sterile container without additive.	3-21 days
	This test is used when an aerobic bacterial infection or fungal infection is suspected in a tissue or fluid.		
	Note: DTM plates are acceptable for fungal ID. However, please indicate the source of the specimen plated.		

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
M060	<p>Culture, Blood Aerobic Only</p> <p>This test is utilized when bacteremia with an aerobic organism is suspected. Submit a single sample in a BACTEC bottle.</p> <p>Interferences: antibiotics.</p> <p>Note: submit in BACTEC aerobic culture bottle. Clip fur and scrub venipuncture site for aseptic collection. Aim to collect 1.0 mL for cats and small dogs, 2-3 mL for larger dogs to be inoculated into 20 mL bottle. Don't unscrew caps on bottles. Remove the protective top and wipe visible parts of the rubber stopper with 70% ethanol. Allow stopper to dry or wipe with sterile gauze. Replace the drawing needle with a sterile needle before puncturing the rubber stopper, fill until vacuum stops, then gently invert the bottle to mix. Anticoagulants in the media will prevent blood from clotting.</p>	<p>Whole blood collected in BD Bactec Blood Culture Bottle (BCB)</p>	<p>3-5 days</p>
M061	<p>Culture Blood Aerobic/Anaerobic</p> <p>Aerobic and Anaerobic Blood Culture</p> <p>This test is utilized when bacteremia is suspected. Submit a single sample in two separate BACTEC bottles as required for aerobic and anaerobic culture.</p> <p>Interferences: prior to culture submission, patient should be off antibiotics for seven to ten days.</p> <p>Note: submit two BACTEC culture bottles (pink label and gold label). Clip fur and scrub venipuncture site for aseptic collection. Aim to collect 1.0 mL for cats and small dogs, 2-3 mL for larger dogs to be inoculated into BacTec Peds Plus pink cap/strip bottle. Anaerobic blood culture bottles with orange cap/gold strip require 3-10 mL of blood. Don't unscrew caps on bottles. Remove the protective top and wipe visible parts of the rubber stopper with 70% ethanol. Allow stopper to dry or wipe with sterile gauze. Replace the drawing needle with a sterile needle before puncturing the rubber stopper, fill until the vacuum stops, then gently invert the bottle to mix. Anticoagulants in the media will prevent blood from clotting.</p>	<p>Whole blood collected in two separate BD Bactec Blood Culture Bottle (BCB): 1 aerobic BD Bactec BCB and 1 anaerobic BD Bactec BCB</p>	<p>5 days</p>

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
M062	<p>Culture Blood Aerobic/Anaerobic (Serial Blood Culture x2)</p> <p>Bacterial culture and MIC sensitivity testing</p> <p>Utilize this test when bacteremia is suspected. Aerobic/ anaerobic cultures are obtained at two separate time points. Each single time point requires two separate BACTEC bottles as indicated for aerobic and anaerobic culture.</p> <p>Interferences: antibiotics.</p> <p>Note: submit two BACTEC culture bottles (pink label and gold label) for each time point. For each collection, clip fur and scrub venipuncture site for aseptic collection. Aim to collect 1.0 mL for cats and small dogs, 2-3 mL for larger dogs to be inoculated into BacTec Peds Plus pink cap/strip bottle. Anaerobic blood culture bottles with orange cap/gold strip require 3-10 mL of blood. Don't unscrew caps on bottles. Remove the protective top and wipe visible parts of the rubber stopper with 70% ethanol. Allow stopper to dry or wipe with sterile gauze. Replace the drawing needle with a sterile needle before puncturing the rubber stopper, fill until the vacuum stops, then gently invert the bottle to mix. Anticoagulants in the media will prevent blood from clotting.</p>	<p>2 sets of aerobic and anaerobic BD Bactec blood culture bottles (one set for each single time point)</p>	<p>5 days</p>
M063	<p>Culture Blood Aerobic/Anaerobic (Serial Blood Culture x3)</p> <p>Bacterial culture and MIC sensitivity testing</p> <p>Utilize this test when bacteremia is suspected. Aerobic/ anaerobic cultures are obtained at three separate time points. Each single time point requires two separate BACTEC bottles as indicated for aerobic and anaerobic culture.</p> <p>Interferences: antibiotics.</p> <p>Note: submit two BACTEC culture bottles (pink label and gold label) for each time point. For each collection, clip fur and scrub venipuncture site for aseptic collection. Aim to collect 1.0 mL for cats and small dogs, 2-3 mL for larger dogs to be inoculated into BacTec Peds Plus pink cap/strip bottle. Anaerobic blood culture bottles with orange cap/gold strip require 3-10 mL of blood. Don't unscrew caps on bottles. Remove the protective top and wipe visible parts of the rubber stopper with 70% ethanol. Allow stopper to dry or wipe with sterile gauze. Replace the drawing needle with a sterile needle before puncturing the rubber stopper, fill until the vacuum stops, then gently invert the bottle to mix. Anticoagulants in the media will prevent blood from clotting.</p>	<p>3 sets of aerobic and anaerobic BD Bactec blood culture bottles (one set for each single time point)</p>	<p>5 days</p>

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
M070	<p>Culture ID</p> <p>This test is for bacterial culture and identification only. No sensitivity is performed.</p>	<p>Culturette from fluid (body cavity fluids), TTW, BAL, wound, lesion, skin, or environmental sample. Tissue sample in saline.</p> <p>Culturette, red top or white top tube with fluid, or other sterile container. Samples collected in EDTA tube are not acceptable.</p>	<p>2-4 days</p> <p>If a fastidious culture organism is observed, the listed turnaround time may be extended.</p>
M020	<p>Culture, Aerobic</p> <p>This test is used when an aerobic bacterial infection is suspected in a tissue or fluid. In the case of a positive culture of relevant bacteria, sensitivities are performed and reported.</p> <p>Interferences:</p> <ol style="list-style-type: none"> 1. Patient should be off antibiotics for at least 7-10 days. 2. Fluid in lavender top tube is unacceptable (EDTA inhibits growth). <p>Note: tissue should be submitted in an RTT with a few drops of saline (to keep moist). Samples should be refrigerated prior to transportation to the laboratory. A separate sample should be submitted for each type of culture needed.</p>	<p>Culturette from fluid (body cavity fluids), TTW, BAL, wound, lesion, or skin. Tissue sample in saline.</p> <p>Culturette, red top or white top tube with fluid, or other sterile container. Samples collected in EDTA tube are not acceptable.</p>	<p>3-4 days</p> <p>Preliminary report available every 24 hours. Final culture result available in 72 hours. If a fastidious organism is observed, the listed turnaround time may be extended.</p>
M110	<p>Culture, Mycoplasma</p> <p><i>Mycoplasma</i> culture</p> <p>Use this test when <i>Mycoplasma</i> infection is suspected.</p> <p>Interferences: antibiotics.</p>	<p>Culturette from body fluid, TTW, semen, fresh urine. Uterine, cervical, vaginal or lung tissue.</p> <p>Culturette, red top or white top tube with fluid, or other sterile container. Samples collected in EDTA tube are not acceptable.</p>	<p>7-10 days</p>

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
M130 Add-on Equivalent ADD210	Culture, Urine All urine cultures deemed positive by <i>FIRST</i> Tract, a highly accurate and rapid urine culture technique, will then undergo traditional culture for identification and susceptibility. Interferences: patient should be off antibiotics for seven to ten days prior to urine culture.	0.5 mL cystocentesis, clean catch or catheterized urine in sterile red top tube or urine transport tube.	1-3 days Preliminary reports available every 24 hours. Canine urine culture final result available in 48 hours. Feline urine culture final result available in 72 hours. If a fastidious organism is observed, the listed turnaround time may be extended.
ADD210	Culture, Urine Add-on <i>FIRST</i> Tract Urine Culture followed by plating in the case of a positive result for bacterial identification and sensitivity All urine cultures deemed positive by <i>FIRST</i> Tract, a highly accurate and rapid urine culture technique, will then undergo traditional culture for identification and susceptibility. Interferences: patient should be off antibiotics for seven to ten days prior to urine culture.	0.5 mL cystocentesis, clean catch or catheterized urine in sterile red top tube or urine transport tube.	1-3 days Preliminary reports available every 24 hours. Canine urine culture final result available in 48 hours. Feline urine culture final result available in 72 hours. If a fastidious organism is observed, the listed turnaround time may be extended.

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
M130R	<p>Culture, Urine Recheck</p> <p>All urine cultures deemed positive by <i>FIRST</i>Tract, a highly accurate and rapid urine culture technique, will then undergo traditional culture for identification and susceptibility. Used as a follow up code in a patient, where a culture has been previously submitted within the last 30 days.</p> <p>Interferences: patient should be off antibiotics for seven to ten days prior to urine culture.</p>	<p>0.5 mL cystocentesis, clean catch or catheterized urine in sterile red top tube or urine transport tube.</p>	<p>1-3 days</p> <p>Preliminary reports available every 24 hours. Canine urine culture final result available in 48 hours. Feline urine culture final result available in 72 hours. If a fastidious organism is observed, the listed turnaround time may be extended.</p>
M100	<p>Culture, Acid Fast Bacilli</p> <p>Mycobacterial culture and acid fast staining only</p> <p>Use this test when a mycobacterial infection of fluid/tissue is suspected.</p> <p>Interferences: antimicrobials.</p>	<p>Urine, CSF, sputum, or tissue in sterile container or Copan swab (without contact with gel)</p> <p>Call customer service for specimen handling protocol: 1-800-872-1001, dial 0</p>	<p>8 weeks</p> <p>Culture requires 8 weeks to be finalized. No preliminary report is provided, only final report.</p>
M030	<p>Culture, Anaerobic</p> <p>To be used when an anaerobic bacterial infection is suspected in tissue/fluid.</p> <p>Interferences: sample with exposure to air or if dried will preclude accurate testing.</p> <p>Note: fluids should be submitted in an RTT. Tissue should be submitted in an air-tight, sterile container (e.g., urine cup). No saline is necessary. Sensitivities are not performed, but drugs of choice for anaerobes are provided as a guide.</p>	<p>Anaerobic culturette of fluid (body cavity fluids), TTW or BAL, wound, lesion. Tissue sample (at least 2 cm x 2 cm) in sterile, air tight container.</p> <p>Culturette, red top or white top tube with fluid, or other sterile container. Samples collected in EDTA tube are not acceptable.</p>	<p>3-5 days</p> <p>Preliminary report issued 3rd day. Final culture result issued 4th day.</p>
CAMP	<p>Culture, <i>Campylobacter</i></p> <p>Use this test when a <i>Campylobacter</i> infection is suspected to be the cause of gastroenteritis that is not responding to conservative management.</p> <p>Interferences: dried swab.</p>	<p>2 grams of feces in sterile container or in <i>Campylobacter</i>-Thioglycollate broth (Campy-Thio broth), or fecal swab.</p>	<p>3-4 days</p>

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
M240	Culture, Dermatophytes Use this test to investigate ringworm (dermatophyte fungal infections) as a cause for skin lesions.	Hair, nails, or skin scraping in plain container DTM bottle	3 weeks Preliminary report every 7 days. Final report available at 3 weeks (or upon identification).
M160	Culture, Feces Culture specifically evaluates for <i>Salmonella</i> , <i>Shigella</i> and <i>Campylobacter</i> spp. Fecal PCR testing is more sensitive and tests for a broader array of potential pathogens.	Fecal culturette or 5 grams feces in Antech provided fecal container. Submission in in lavender top tube (with EDTA) is not acceptable.	3-4 days Preliminary report every 24 hours. Final report available in 72 hours.
M080	Culture, Fungal This test is for fungal culture and identification. Not for suspected ringworm infections.	Dry hair, nails, skin scraping, body fluid, or lesion material collected in a sterile red top or other sterile container without additive, or on a culturette. DTM bottle, culturette, container	21 days Preliminary report every 7 days. Final report available at 21 days (3 weeks), or upon identification.
M255	Ear Culture Ear Aerobic Culture and Sensitivity An aerobic culture of material obtained from the ear canal. Interferences: antibiotics.	Culturette from ear	3-5 days Preliminary report available every 24 hours. Final culture result available in 72 hours.
M076	Environmental Culture	Culturette of environmental sample	3-4 days Preliminary report available every 24 hours. Final culture result available in 72 hours.

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
M082	Fungal Identification Used after fungal growth has occurred in order to identify the fungi as dermatophytes versus saprophytic fungi.	DTM fungal culture bottle (with fungal growth), or fungal culture plate (with fungal growth)	2 days If further culturing is required, testing may be extended 7-10 days
M090	Gram Stain This test is used to investigate the presence and morphology of bacteria, yeast, and fungi. Note: place specimen on a slide or submit material to prepare a slide. Store at room temperature.	Prepared smear, culturette from fluid (body cavity fluids, TTW or BAL), feces, urine, wound or lesion material. Culturette, slide holder	1-2 days
M280	Ophthalmic Aerobic Culture An aerobic culture and eye-specific antibiotic susceptibility.	Submission of 2 Copan swabs is preferred.	3-5 days Preliminary report available every 24 hours. Final culture result available in 72 hours.
M120	Salmonella, Shigella, Campylobacter <i>Campylobacter</i> Culture, Fecal Culture (<i>Salmonella</i> and <i>Shigella</i>) Culture specifically evaluates for <i>Salmonella</i> , <i>Shigella</i> and <i>Campylobacter</i> . Fecal PCR testing is more sensitive and tests for a broader array of potential pathogens. Note: submit feces in a sterile container, a culturette, or in a Campy thio broth.	Fecal culturette or 5 grams feces in Antech provided fecal container. Submission in in lavender top tube (with EDTA) is not acceptable.	3-4 days Preliminary report available every 24 hours. Final culture result available in 72 hours.
M135	Urine Culture Aerobic and Anaerobic Aerobic culture performed by <i>FIRST</i> tract methodology followed by standard plating in cases of positive growth. Anaerobic culture performed by standard plating but sensitivities are not performed.	0.5 mL cystocentesis, clean catch or catheterized urine in sterile red top tube or urine transport tube.	3-5 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

PATHOLOGY

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
BONECBC	Bone Marrow and CBC <p>This test includes a microscopic evaluation of bone marrow. The report includes descriptions of cellularity, assessment of cell lineage, and pathologist comments. The CBC submitted at the time of the bone marrow sampling is utilized in the interpretation of the bone marrow evaluation to allow for a more complete assessment of the hematopoietic system. A detailed history and time course of events will further aid the pathologist in the interpretation of the findings.</p> <p>Note: ensure the appropriate labeling of lavender top tubes (bone marrow and peripheral blood). Provide a detailed history and any relevant clinicopathologic data for the pathologist's interpretation.</p>	1.0 mL EDTA whole blood, AND bone marrow sample in lavender top tube or 4 to 5 air-dried bone marrow slides.	1-3 days
BMBC	Bone Marrow Core Biopsy <p>The histopathological evaluation of a bone marrow core biopsy. A CBC submitted at the time of the bone marrow sampling is utilized in the interpretation of the bone marrow evaluation to allow for a more complete assessment of the hematopoietic system. A detailed history and time course of events will further aid the pathologist in the interpretation of the findings.</p>	Bone marrow tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid <p>To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin</p>	3-5 business days
BONE	Bone Marrow Cytology <p>This test includes a microscopic evaluation of bone marrow. The report includes description of cellularity, assessment of cell lineage, and pathologist comments.</p> <p>Note: for maximal diagnostic utility, a CBC should be submitted at the time of bone marrow sampling. Additionally, submission of a core biopsy along with an aspirate is necessary for a complete bone marrow evaluation. If submitting bone marrow aspirate slides and a core biopsy sample (formalin), ensure that the formalin sample is in a separate zip lock bag.</p>	4 to 5 air dried bone marrow slides or bone marrow sample in lavender top tube	1-3 days
T325	Buffy Coat For Mast Cells <p>Used to assess for the presence of mast cells, evaluate for hemoparasites/infectious agents, or to do a differential on a complete blood count with low white cells.</p>	1.0 mL EDTA whole blood in lavender top tube	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S86793	Canine Melanoma Diagnostic Panel	Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin	2 Weeks
S86792	Canine Melanoma Prognostic Panel	Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin	2 Weeks
CYTO	Cytology This test includes a microscopic evaluation of cells. The report includes cytologic interpretation, diagnosis, and comments regarding etiology and biological behavior where applicable. Note: single source.	2 or more air-dried, unstained slides with clinical history Lavender top with fluid, red top with fluid, air-dried unstained slides	1-3 days
DERM	Dermatopathology This specialty biopsy is supported by a team of board-certified anatomic pathologists with a specialization in dermatopathology.	Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin	7 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
FLUA	<p>Fluid Analysis</p> <p>Includes preparation of submitted sample, Cell Count (WBC & RBC), Specific Gravity, Protein measurement, and microscopic interpretation by clinical pathologist</p> <p>This is the appropriate test for the evaluation of fluids from a cavity or tissue space. The fluid will be characterized based on protein, specific gravity, white and red cells, and a pathologist will evaluate the fluid and provide a microscopic description and, when possible, diagnosis.</p> <p>Note: a comprehensive history should be included with the fluid analysis submission.</p>	<p>1.0 mL body cavity fluid in lavender or red top tube with 2 unstained smears prepared from fluid</p> <p>Red top with fluid, lavender top with fluid, air dried unstained slides</p>	<p>1-3 days</p>
FL	<p>Fluid Part Only</p> <p>This portion of a fluid analysis includes characterization of the submitted body cavity or tissue spave fluid based on color, clarity, specific gravity, protein and white and red cell counts.</p>	<p>0.5 mL body cavity fluid in lavender or red top tube</p>	<p>1-3 days</p>
FBX	<p>Histopathology</p> <p>Preparation of submitted sample and microscopic interpretation by anatomic pathologist</p> <p>A boarded pathologist will evaluate the tissue. A full written biopsy report will be provided, including source, history, submitted tissue type, microscopic findings, pathologist's comments, and microscopic description. Where applicable, it will also include margin evaluation, grading, interpretation of special stains, and recommendations for further testing if needed.</p> <p>Note: firmly and evenly tighten the lid of the formalin jar and check for any leaks prior to placing the sample in a ziplock bag with the test requisition form. A comprehensive history should be included with each biopsy submission.</p>	<p>Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid</p> <p>To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin</p>	<p>3-5 business days</p>
IHC1	<p>Immunohistochemistry 1 Stain</p> <p>This test includes a single immunohistochemistry stain. It's an antibody-based method to detect a specific protein.</p> <p>Note: typically an add-on test after histopathological evaluation of submitted tissue based on pathologist recommendations.</p>	<p>Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid</p> <p>To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin</p>	<p>10-14 days</p>

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
IHC2	Immunohistochemistry 2 Stains This test includes two immunohistochemistry stains. It's an antibody-based method to detect specific proteins. Note: typically an add-on test after histopathological evaluation of submitted tissue based on pathologist recommendations.	Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin	10-14 days
IHC3	Immunohistochemistry 3 Stains This test includes three immunohistochemistry stains. It's an antibody-based method to detect specific proteins. Note: typically an add-on test after histopathological evaluation of submitted tissue based on pathologist recommendations.	Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin	10-14 days
IHC4	Immunohistochemistry 4 Stains This test includes four immunohistochemistry stains. It's an antibody-based method to detect specific proteins. Note: typically an add-on test after histopathological evaluation of submitted tissue based on pathologist recommendations.	Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin	10-14 days
IHC5	Immunohistochemistry 5 Stains This test includes five immunohistochemistry stains. It's an antibody-based method to detect specific proteins. Note: typically an add-on test after histopathological evaluation of submitted tissue based on pathologist recommendations.	Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin	10-14 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S86601	Mast Cell Panel	Formalin Fixed Paraffin Blocks accompanied by an Antech biopsy report. All submissions need to have an FBX completed by Antech (including any FBX type test, including 2nd opinion by Antech pathologist using block made by some other lab), prior to sending to MSU for an additional fee.	7-10 days
FBXORL	Oral Path Biopsy Preparation of each submitted oral tissue sample and microscopic interpretation by anatomic pathologist A boarded pathologist with a special interest in oral cavity pathology will evaluate the tissue. A full written biopsy report will be provided, including source, history, description of submitted tissue, microscopic description, and diagnosis where possible. Additionally, a dentist will provide a treatment plan based on the histopathology or recommendations for further diagnostic evaluation. Note: firmly and evenly tighten the lid of the formalin jar and check for any leaks prior to placing the sample in a ziplock/ sealable bag with the test requisition form. A comprehensive history should be included with each biopsy submission.	Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin	3-5 business days
FBXNEO	Pet Cancer Specialty Biopsy This specialty biopsy service is comprised of a team of six board-certified anatomic pathologists with specialization in cancer. The biopsy is reviewed by three of these pathologists including one report written with a consensus opinion. Note: firmly and evenly tighten the lid of the formalin jar and check for any leaks prior to placing the sample in a ziplock/ sealable bag with the test requisition form. A comprehensive history should be included with each biopsy submission.	Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin	3-5 business days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S11068 - SearchLight DNA Histology Add-on S10058 - SearchLight DNA Cytology Add-on	<p>SearchLight DNA®</p> <p>SearchLight DNA is an add-on test to any biopsy or cytology where cancer is suspected or diagnosed. This test provides precision diagnostic, prognostic, and therapeutic guidance for multiple types of cancers.</p> <p>This canine cancer genomic panel identifies mutations in 120 relevant cancer genes and provides insights into cancer's origin, behavior, and the optimal approach for treatment.</p> <p>Report includes:</p> <ul style="list-style-type: none"> • Information about mutations in 120 cancer genes including diagnostic, prognostic, and therapeutic biomarker associations with supporting evidence from peer-reviewed literature, clinical consensus, and inference from human FDA guidance. • Pharmacogenomic marker (MDR1) mutation status. • A list of targeted therapeutic drug(s) if the patient's tumor mutations are a match for an existing drug. • Information on clinical trials by tumor type curated from public databases and individual academic centers by Vidium Animal Health. <p>To order SearchLight DNA, please contact customer service.</p> <p>Note: SearchLight DNA can only be added to a biopsy or cytology performed at Antech where neoplasia has been diagnosed or suspected.</p>	<p>Fine-Needle Aspirate (FNA) Samples: FNA Slides - ≥2 unstained or stained slides + 1 Diffquick FNA Liquid - ≥0.5 mL in sterile vial/tube with no additives</p> <p>Biopsy Samples: Formalin-Fixed, Paraffin-Embedded (FFPE) Scrolls – 10 x 10- micron sections with 1 adjacent hematoxylin and eosin (H&E) slide Unstained FFPE Slides – 10 x 10-micron sections with 1 adjacent H&E slide</p> <p>Bone Tumor – Non-decalcified sample</p>	12-17 days
RCYTO	<p>Recheck Cytology</p> <p>Used for a second cytology submission, from the same site as the original submission, in cases where the initial submission was of insufficient cellularity to allow for a diagnosis/description.</p>	Fluid in red top or lavender top tube, or air-dried,unstained slides with clinical history	1-3 days
STAT	<p>Stat Charge Biopsy</p> <p>Additional fee for histopathology results that client would like to get an expedited return on.</p>	Not applicable	2-3 days
CYSTAT	<p>Stat Charge-Cyto</p> <p>An additional charge applied when a STAT assessment of cytological sample is required.</p>	Not applicable	24 hours
FBXLSP	<p>Superliver (Liver Biopsy with Liver Staining Panel)</p> <p>The tissue will be evaluated by a boarded pathologist with a special interest in hepatic disease. A full written biopsy report will be provided, including source, history, description of submitted tissue, microscopic description, and diagnosis where possible. Trichrome, rhodanine, and reticulum staining and interpretation are included. The report will include the staging of fibrosis, evaluation of the degree of parenchymal collapse, and a qualitative copper grade. Internal pathologist to pathologist opinions are sought when needed. There is an option for further discussion with an internal medicine consultant with a special interest in hepatic disease.</p> <p>Note: firmly and evenly tighten the lid of the formalin jar and check for any leaks prior to placing the sample in a ziplock/ sealable bag with the test requisition form. A comprehensive history should be included with each biopsy submission, including the breed, age, and time when liver values were first noted to be elevated. Additional diagnostics performed as well as responses to any specific treatments should be included in the history.</p>	<p>Liver tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid</p> <p>To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin</p>	3-5 business days

MOLECULAR

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T1025	<p>CADET BRAF</p> <p>CADET <i>BRAF</i>. Includes CADET <i>BRAF</i>-PLUS when reflexed based on CADET <i>BRAF</i> result</p> <p>A non-invasive, highly sensitive assay that evaluates free-catch urine samples from dogs for the presence of cells harboring the b-raf mutation or specific copy number variations associated with Transitional Cell Carcinoma/ Urothelial Carcinoma/ Prostatic carcinoma. The assay identifies 95% of TCC/UC cases. It can be used to diagnose TCC/UC and monitor dogs undergoing treatment for remission and relapse.</p> <p>Interferences: cystocentesis drawn urine.</p> <p>Radiation and Chemotherapy are not interferences for the Cadet <i>BRAF</i>/<i>BRAF</i> Plus test assay; however, they may affect the tumor itself which may influence results.</p> <p>Note: urine sample should be free catch. Do not freeze.</p>	<p>40 mL urine collected in CADET BRAF urine container</p> <p>CADET <i>BRAF</i> container (urine must be put in <i>BRAF</i> container within 15 minutes of collection and can be collected over multiple days)</p>	3-6 days
T980	<p>Canine Ehrlichia/Anaplasma PCR</p> <p><i>Anaplasma phagocytophilum</i>, <i>Anaplasma platys</i>, <i>Ehrlichia canis</i>, <i>Ehrlichia chaffeensis</i>, <i>Ehrlichia ewingii</i></p> <p>A comprehensive PCR panel used when the clinical signs are suggestive of Anaplasmosis or Erhlichiosis.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results.</p>	<p>1.0 mL whole blood, synovial fluid or CSF</p>	2-4 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T950	<p>Canine GI PCR Panel</p> <p><i>Campylobacter jejuni/coli</i>, canine enteric coronavirus, canine parvovirus, <i>Clostridium difficile</i> toxins A/B, <i>Clostridium perfringens</i> enterotoxin, <i>Cryptosporidium</i> spp., <i>Giardia</i> spp., <i>Salmonella</i> spp.</p> <p>Highly sensitive and specific panel used to evaluate for viral (canine enteric coronavirus, canine parvovirus, and canine rotavirus), bacterial (<i>Campylobacter</i> spp., <i>Clostridium</i> spp. and <i>Salmonella</i> spp.) and protozoal (<i>Giardia</i> spp. and <i>Cryptosporidium</i> spp.) causes for the gastrointestinal signs exhibited by the patient.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results for protozoal or bacterial organisms.</p> <p>Note: follow-up cultures are immediately performed on submissions found to be PCR-positive for <i>Salmonella</i> (when a Copan fecal swab is received by the laboratory). Additionally, samples found to be PCR-positive for DNA of <i>C. perfringens</i> enterotoxin or <i>C. difficile</i> toxins are immediately tested for these toxins by ELISA.</p>	<p>0.5 grams of feces and 1 copan fecal swab</p> <p>Culturette, Antech provided fecal container</p>	2-4 days
T953	<p>Canine GI PCR with SARS CoV-2 PCR Panel</p> <p>Canine GI PCR Panel, SARS-CoV-2 PCR add-on Panel</p> <p>Canine GI PCR panel to which evaluation for SARS-CoV-2 PCR has been added to determine whether the gastrointestinal signs being exhibited could be related to the typical etiologic agents or potentially SARS-CoV-2.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results for protozoal or bacterial organisms.</p> <p>Note: follow-up cultures are immediately performed on submissions found to be PCR-positive for <i>Salmonella</i> (when a Copan fecal swab is received by the laboratory). Additionally, samples found to be PCR-positive for DNA of <i>C. perfringens</i> enterotoxin or <i>C. difficile</i> toxins are immediately tested for these toxins by ELISA.</p>	<p>0.5 grams of feces and 1 copan fecal swab, conjunctival or nasal swab</p> <p>Culturette, Antech provided fecal container</p>	2-11 days
S8710	<p>Canine Parvovirus PCR</p> <p>Detection of canine parvovirus in whole blood, feces or tissues.</p> <p>Note: vaccination for canine parvovirus within the previous 2-3 weeks may result in false positive results.</p>	<p>2.0-3.0 mL EDTA whole blood in lavender top tube OR 1-2 grams feces OR small tissue in red top</p>	7-10 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T995	<p>Canine Respiratory PCR Panel</p> <p>Canine influenza virus (H3N8), Pan Influenza (detects all Influenza Type A, including H3N2; when positive, confirmation by H3N2 specific qPCR), H5N1 influenza virus, canine adenovirus type 2, canine distemper virus, canine herpes virus, canine parainfluenza virus, canine resp. coronavirus, <i>Bordetella bronchiseptica</i>, <i>Mycoplasma cynos</i>, <i>Streptococcus equi</i> subsp. <i>zooepidemicus</i></p> <p>A highly sensitive and specific panel covering the more common etiological agents causing acute respiratory signs in dogs. Includes Distemper, Parainfluenza, Respiratory coronavirus, Canine influenza (H3N8), Influenza (H5N1), PanInfluenza, Herpes, Adenovirus Type 2, <i>Bordetella bronchiseptica</i>, <i>Mycoplasma cynos</i>, and <i>Streptococcus zooepidemicus</i>.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results.</p>	<p>2 Sterile Swabs (nasal, deep pharyngeal or conjunctival) for PCR (submitted dry without transport media) and 1 Copan Swab in Gel for <i>Bordatella</i> culture. All swabs sampled from same site.</p> <p>Culturette and 2 swabs in white tube</p>	2-4 days
T998	<p>Canine Respiratory PCR with SARS CoV-2 PCR Panel</p> <p>Canine Respiratory PCR Panel(Canine influenza virus (H3N8), Pan Influenza (detects all Influenza Type A, including H3N2; when positive, confirmation by H3N2 specific qPCR), H5N1 influenza virus, canine adenovirus type 2, canine distemper virus, canine herpes virus, canine parainfluenza virus, canine resp. coronavirus, <i>Bordetella bronchiseptica</i>, <i>Mycoplasma cynos</i>, <i>Streptococcus equi</i> subsp. <i>zooepidemicus</i>), SARS-CoV-2 PCR Add-on Panel</p> <p>Canine Respiratory PCR panel to which evaluation for SARS-CoV-2 PCR has been added to determine whether the respiratory signs being exhibited could be related to the typical etiologic agents (see T995) or potentially SARS-CoV-2.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results.</p>	<p>2 Sterile Swabs (nasal, deep pharyngeal or conjunctival) for PCR (submitted dry without transport media) and 1 Copan Swab in Gel for <i>Bordatella</i> culture. All swabs sampled from same site.</p> <p>Culturette and 2 swabs in white tube</p>	2-11 days
T960	<p>Canine Tick Borne PCR Panel</p> <p><i>Anaplasma phagocytophilum</i>, <i>Anaplasma platys</i>, <i>Babesia canis</i>, <i>Babesia</i> sp. (Coco), <i>Babesia conradae</i>, <i>Babesia gibsoni</i>, <i>Bartonella henselae</i>, <i>Bartonella vinsonii</i>, <i>Ehrlichia canis</i>, <i>Ehrlichia chaffeensis</i>, <i>Ehrlichia ewingii</i>, <i>Mycoplasma haemocanis</i>, <i>Cand. M. haematoparuvum</i>, <i>Neorickettsia risticii</i></p> <p>A highly sensitive and specific panel used to evaluate for the evidence of infection in cases suspected of having a vector-borne disease. Includes <i>Anaplasma phagocytophilum</i>, <i>Anaplasma platys</i>, <i>Babesia canis</i>, <i>Babesia</i> spp. (non <i>canis</i>), <i>Bartonella henselae</i>, <i>Bartonella vinsonii</i>, <i>Ehrlichia canis</i>, <i>Ehrlichia</i> spp. (non-<i>canis</i>), <i>Mycoplasma heamocanis/</i> <i>heamatoparvum</i>, <i>Neorickettsia risticii</i>, and <i>Rickettsia rickettsii</i>.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results.</p>	<p>1.0 mL of EDTA whole blood in lavender top tube</p>	2-4 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T955	<p>Feline GI PCR Panel</p> <p><i>Campylobacter coli</i>, <i>Campylobacter jejuni</i>, <i>Clostridium difficile</i> toxins A/B, <i>Clostridium perfringens</i> enterotoxin, <i>Cryptosporidium</i> spp. and <i>C. felis</i>, feline parvovirus, <i>Giardia</i> spp. <i>Salmonella</i> spp. <i>Tritrichomonas foetus</i></p> <p>Highly sensitive and specific panel used to evaluate for viral (feline panleukopenia), bacterial (<i>Campylobacter</i> spp, <i>Clostridium</i> spp toxins, and <i>Salmonella</i> spp.) and protozoal (<i>Giardia</i> spp., <i>Tritrichomonas foetus</i> and <i>Cryptosporidium</i> spp.) causes for the gastrointestinal signs exhibited by the patient by PCR.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results for protozoal or bacterial organisms.</p> <p>Note: follow-up cultures are immediately performed on submissions found to be PCR-positive for <i>Salmonella</i> (when a Copan fecal swab is received by the laboratory). Additionally, samples found to be PCR-positive for DNA of <i>C. perfringens</i> enterotoxin or <i>C. difficile</i> toxins are immediately tested for these toxins by ELISA.</p>	<p>0.5 grams of feces and 1 copan fecal swab</p> <p>Culturette, Antech provided fecal container</p>	<p>2-4 days</p>
T958	<p>Feline GI PCR with SARS CoV-2 PCR Panel</p> <p>Feline GI PCR Panel, SARS-CoV-2 PCR Add-on Panel</p> <p>Feline GI PCR panel (T955) to which evaluation for SARS-CoV-2 PCR has been added to determine whether the gastrointestinal signs being exhibited could be related to the typical etiologic agents or potentially SARS-CoV-2.</p> <p>Note: follow-up cultures are immediately performed on submissions found to be PCR-positive for <i>Salmonella</i> (when a Copan fecal swab is received by the laboratory). Additionally, samples found to be PCR-positive for DNA of <i>C. perfringens</i> enterotoxin or <i>C. difficile</i> toxins are immediately tested for these toxins by ELISA.</p>	<p>0.5 grams of feces and 1 copan fecal swab</p> <p>Culturette, Antech provided fecal container</p>	<p>2-11 days</p>
T985	<p>Feline Hemoplasma PCR Panel</p> <p><i>Mycoplasma haemofelis</i>, <i>Candidatus M. haemominutum</i>, <i>Candidatus M. turicensis</i></p> <p>Used in the evaluation of a regenerative anemia where hemotropic <i>Mycoplasma</i> is a differential. PCR testing for the presence of <i>Mycoplasma haemofelis</i>, <i>Mycoplasma haemominutum</i>, and <i>Mycoplasma turicensis</i>.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results.</p>	<p>1.0 mL EDTA whole blood in lavender top tube</p>	<p>2-4 days</p>

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T990	<p>Feline Respiratory PCR Panel</p> <p>Feline calicivirus (FCV), Feline herpes virus-1 (FHV-1), <i>Bordetella bronchiseptica</i>, <i>Chlamydomphila felis</i>, <i>Mycoplasma felis</i>, H1N1 influenza virus.</p> <p>A high sensitive and specific panel covering the more common etiological agents associated with upper respiratory signs in cats. Includes feline calicivirus, feline herpes virus 1, <i>Chlamydia felis</i>, <i>Mycoplasma felis</i>, <i>Bordetella bronchiseptica</i>, Influenza H1N1.</p> <p>Interferences: previous or current antibiotic/antiviral usage may interfere with PCR test results.</p> <p>Note: complimentary culture/susceptibility testing on all submissions PCR-positive for <i>Bordetella bronchiseptica</i> (when the laboratory receives a Copan swab).</p>	<p>2 Sterile Swabs (nasal, deep pharyngeal or conjunctival) for PCR (submitted dry without transport media) and 1 Copan Swab in Gel for <i>Bordetella</i> culture. All swabs sampled from same site.</p> <p>Culturette and 2 swabs in white tube (swab the location of the predominant clinical signs: nasal exudates, deep pharyngeal swab and/or conjunctival swab)</p>	2-4 days
T993	<p>Feline Respiratory PCR with SARS CoV-2 PCR Panel</p> <p>Feline Respiratory PCR Panel, SARS-CoV-2 PCR Add-on Panel</p> <p>Feline Respiratory PCR panel (T990) to which evaluation for SARS-CoV-2 PCR has been added to determine whether the respiratory signs being exhibited could be related to the typical etiologic agents or potentially SARS-CoV-2.</p> <p>Interferences: previous or current antibiotic/ antiviral usage may interfere with PCR test results.</p>	<p>2 Sterile Swabs (nasal, deep pharyngeal or conjunctival) for PCR (submitted dry without transport media) and 1 Copan Swab in Gel for <i>Bordetella</i> culture. All swabs sampled from same site.</p> <p>Culturette and 2 swabs in white tube</p>	2-11 days
T600	<p>FIP mRNA PCR</p> <p>This PCR test detects mRNA of the M gene of all known feline coronavirus strains in any sample. However, for diagnosis of FIP, only the detection of mRNA outside of the intestinal tract is indicative since active replication of the virus in circulating mononuclear cells is typical for FIP.</p> <p>Note: most common sample submission is peritoneal effusion that on fluid analysis is consistent with FIP. Other sample options include CSF, tissue, and tissue aspirates (lymph node, etc.).</p>	<p>1.0 mL fresh fluid (body cavity, CSF, joint, milk, sputum, or urine), EDTA whole blood, tissue samples in plain, non-additive tube (+/- saline), or unstained slides from aspirated tissue (i.e., lymph node, liver or kidney)</p>	5-7 days
ADD350	<p>KeyScreen® GI Parasite PCR Add-on</p> <p>Detects the following parasites by PCR: <i>Ancylostoma</i> spp. (if positive evaluated for benzimidazole resistance by PCR), <i>Uncinaria stenocephala</i>, <i>Toxocara/Toxascaris</i> spp. (if <i>Toxocara</i> spp. positive further speciated into <i>T. canis</i>, <i>T. cati</i>, <i>T. leonina</i>), <i>Baylisascaris procyonis</i>, <i>Trichuris vulpis</i>, <i>Dipylidium caninum</i>, <i>Echinococcus granulosus</i>, <i>Echinococcus multilocularis</i>, <i>Taenia</i> spp., <i>Giardia duodenalis</i>, <i>Giardia zoonotic</i> (if positive evaluated for potentially zoonotic strains A and B), <i>Cystoisospora</i> spp., <i>Eimeria</i> spp., <i>Cryptosporidium canis</i>, <i>Cryptosporidium felis</i>, <i>Toxoplasma gondii</i> (feline specific), <i>Neospora caninum</i>, <i>Tritrichomonas blagburni</i> (feline specific).</p>	<p>0.3 grams of fresh feces in Antech provided fecal container (minimum 0.15 grams)</p>	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T991 Add-on Equivalent ADD350	KeyScreen® GI Parasite PCR Panel Detects the following parasites by PCR: <i>Ancylostoma</i> spp. (if positive evaluated for benzimidazole resistance by PCR), <i>Uncinaria stenocephala</i> , <i>Toxocara/Toxascaris</i> spp. (if <i>Toxocara</i> spp. positive further speciated into <i>T. canis</i> , <i>T. cati</i> , <i>T. leonina</i>), <i>Baylisascaris procyonis</i> , <i>Trichuris vulpis</i> , <i>Dipylidium caninum</i> , <i>Echinococcus granulosus</i> , <i>Echinococcus multilocularis</i> , <i>Taenia</i> spp., <i>Giardia duodenalis</i> , <i>Giardia</i> zoonotic (if positive evaluated for potentially zoonotic strains A and B), <i>Cystoisospora</i> spp., <i>Eimeria</i> spp., <i>Cryptosporidium canis</i> , <i>Cryptosporidium felis</i> , <i>Toxoplasma gondii</i> (feline specific), <i>Neospora caninum</i> , <i>Tritrichomonas blagburni</i> (feline specific). Detects 20 individual parasite species. Hookworms: <ul style="list-style-type: none">• <i>Ancylostoma caninum</i>• <i>Uncinaria stenocephala</i>• <i>Ancylostoma</i> benzimidazole resistance Roundworms: <ul style="list-style-type: none">• <i>Toxocara</i> spp.• <i>Toxocara canis</i>• <i>Toxocara cati</i>• <i>Toxascaris leonina</i>• <i>Baylisascaris procyonis</i> Whipworms: <ul style="list-style-type: none">• <i>Trichuris vulpis</i> Tapeworms: <ul style="list-style-type: none">• <i>Dipylidium caninum</i>• <i>Echinococcus granulosus</i>• <i>Echinococcus multilocularis</i>• <i>Taenia</i> spp. <i>Giardia</i> <ul style="list-style-type: none">• <i>Giardia duodenalis</i>• <i>Giardia</i> zoonotic (strain A or B) Coccidia <ul style="list-style-type: none">• <i>Cystoisospora</i> spp.• <i>Eimeria</i> spp. Additional Protozoa <ul style="list-style-type: none">• <i>Cryptosporidium canis</i>• <i>Cryptosporidium felis</i>• <i>Toxoplasma gondii</i> (feline specific)• <i>Neospora caninum</i>• <i>Tritrichomonas blagburni</i> (feline specific) If detected, <i>Toxocara</i> spp. will be identified at the species level as <i>T. cati</i> , <i>T. canis</i> , and <i>T. leonina</i> . <i>Ancylostoma caninum</i> will be evaluated for benzimidazole resistance, and <i>Giardia</i> evaluated for the potentially zoonotic strains A or B.	0.3 grams of fresh feces in Antech provided fecal container (minimum 0.15 grams)	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T974	Leptospira PCR Blood Canine <i>Leptospira</i> PCR - Blood: <i>Leptospira interrogans</i> serovars <i>icterohaemorrhagiae</i> , <i>canicola</i> , <i>pomona</i> , <i>australis</i> , <i>bratislava</i> , <i>autumnalis</i> , <i>ballum</i> and <i>pyrogenes</i> , <i>L. kirschneri</i> serovar <i>grippotyphosa</i> , <i>L. interrogans/borgpetersenii</i> serovar <i>sejroe</i> A highly sensitive and specific test used to evaluate patients suspected of having acute leptospirosis. <i>Leptospira</i> bacteremia is brief and is followed by bacteriuria. In cases where a patient has been ill anywhere from five days to two weeks, submission of blood and urine PCR will increase diagnostic sensitivity. Interferences: previous or current antibiotic usage may interfere with PCR test results.	0.5 mL EDTA whole blood	2-4 days
T978	Leptospira PCR Blood/Urine Canine <i>Leptospira</i> PCR Blood/Urine A highly sensitive and specific test used to evaluate patients suspected of having leptospirosis. <i>Leptospira</i> bacteremia is followed by bacteriuria. In cases where a patient has been ill anywhere from five days to two weeks, submission of blood and urine PCR will increase diagnostic sensitivity. Interferences: previous or current antibiotic usage may interfere with PCR test results.	0.5 mL EDTA whole blood and 2.0 mL urine	2-4 days
T976	Leptospira PCR Urine Canine <i>Leptospira</i> PCR - urine: <i>Leptospira interrogans</i> serovars <i>icterohaemorrhagiae</i> , <i>canicola</i> , <i>pomona</i> , <i>australis</i> , <i>bratislava</i> , <i>autumnalis</i> , <i>ballum</i> and <i>pyrogenes</i> , <i>L. kirschneri</i> serovar <i>grippotyphosa</i> , <i>L. interrogans/borgpetersenii</i> serovar <i>sejroe</i> A highly sensitive and specific test used to evaluate patients suspected of having leptospirosis. <i>Leptospira</i> bacteremia is followed by bacteriuria. urine usually becomes positive day 7-10 post-exposure in clinically ill patients. Interferences: previous or current antibiotic usage may interfere with PCR test results.	2.0 mL urine	2-4 days
S14493	OncoK9® OncoK9® - Cancel Signal OncoK9® is a non-invasive multi-cancer early detection (MCED) liquid biopsy test for the detection and characterization of cancer-associated genomic alterations in DNA isolated from canine blood. Interferences: refrigeration, freezing; pregnancy; clots; hemolysis; trauma/surgery.	Specific collection kit for this test is required and can be ordered on Antech Online. Two cell-free DNA collection tubes, each with a minimum fill of 7 mL of whole blood, total of 14-17 mL whole blood.	7-15 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S85503	<p>PARR</p> <p>PARR (PCR for Antigen Receptor Rearrangement) is an assay that amplifies DNA. This test helps differentiate between monoclonal cells (most consistent with neoplasia) or polyclonal cells (most consistent with a reactive process).</p> <p>Note: PARR cannot be run on formalin-fixed, paraffin-embedded samples/slides or slides with coverslips attached. Non-diagnostic samples are obtained when there is too little DNA or an in the presence of an inhibitor of the PCR reaction.</p>	<ol style="list-style-type: none">1. Blood and bone marrow: 0.5 mL EDTA whole blood or bone marrow sample in lavendertop tube. 4-5 very cellular bone marrow slides can also besubmitted. Peripheral blood smears are not acceptable. If available, please provide a complete blood count withpathology review or cytology/histology report of the bonemarrow.2. Lymph node and other organ aspirates: 4-5 cellular slides (smears or cytopspins) or fluid in EDTA(lavender top tube). No additional fluid (i.e.,saline/serum) needs to be added to aspirated fluid. Gluedcover-slipped slides cannot be used. Stained cytology slidesare acceptable. If available, please provide a copy of relevant cytologyand/or fluid analysis with submission.3. Body cavity fluid: 4-5 cellular slides (smears or cytopspins) or fluid in EDTA (lavender top tube). Stained cytology slides are acceptable. If available, please provide a copy of relevant cytologyand/or fluid analysis with submission.4. CSF: Multiple cytopspin preparations are preferred (stained orunstained) or fluid in EDTA (lavender top tube). For thesample to be diagnostic, it is estimated that 50,000 lymphoid cells are needed. If the CSF has a lymphocyte countof 100/uL, at least 0.5 mL of fluid is needed, or the cellsfrom 0.5 mL spun onto slides. If available, please provide a copy of relevant cytologyand/or fluid analysis with submission.5. Formalin-Fixed Paraffin Embedded (FFPE) Tissue: For this sample type, use code S86965 (Feline IntestinalLymphoma Panel performed at MSU) or S85562 (Lymphoma PanelPARR, performed at MSU). If specifically requesting that PARR testing be performed on FFPE tissue at CSU, please use test code S85503 and specify "FFPE tissue for CSU PARR" <p>Tissue or fluid aspirated into EDTA Lavender top, slides in slide holder</p>	5-7 business days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T984	Rabbit/Small Rodent Ringworm PCR Panel <i>Microsporum</i> spp., <i>M. equinum</i> / <i>canis</i> , <i>M. gypseum</i> (<i>Arthroderma gypseum</i> , <i>A. fulvum</i> , <i>A. incurvatum</i>), <i>Trichophyton</i> spp., <i>T. benhamiae</i> / <i>mentagrophytes</i> (<i>Arthroderma benhamiae</i> , <i>A. vanbreuseghemii</i>) A highly sensitive and specific test used when clinical signs are suggestive of ringworm infection. Note: plucked hair with roots (10-20 hairs), a skin scraping, or a toothbrush sample submitted in a sterile, dry container, free of liquids and preservative or culture media.	Minimum of 10 plucked hair with roots, skin scraping, or tooth brush sample in sterile, dry container free of liquids or preservative 	1-5 days
T988	Rabbit/Small Rodent Ringworm PCR with Dermatophyte Culture <i>Microsporum</i> spp., <i>M. equinum</i> / <i>canis</i> , <i>M. gypseum</i> (<i>Arthroderma gypseum</i> , <i>A. fulvum</i> , <i>A. incurvatum</i>), <i>Trichophyton</i> spp., <i>T. benhamiae</i> / <i>mentagrophytes</i> (<i>Arthroderma benhamiae</i> , <i>A. vanbreuseghemii</i>), Culture, Dermatophytes Evaluation for ringworm infection via both PCR and standard culture.	Minimum of 12 plucked hair with roots, skin scraping, or tooth brush sample. For dermatophyte culture can also submit preinoculated DTM bottle. Hair, toothbrush, DTM	1-21 days
T996	SARS-CoV-2 PCR Add-on Panel This test may be added to either a GI PCR or Respiratory PCR panel, where history may suggest that SARS-CoV-2 may be playing an etiological role. Note: not a standalone test. Must be an add-on to Canine/Feline Respiratory or GI PCR.	Conjunctival and deep pharyngeal swabs, nasal swab when combined with 2nd swab, nasopharyngeal or oropharyngeal aspirates or washes, bronchoalveolar lavage, tracheal aspirates, or fecal sample. Culturette, fecal sample	2-11 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S11068 - SearchLight DNA Histology Add-on S10058 - SearchLight DNA Cytology Add-on	<p>SearchLight DNA®</p> <p>SearchLight DNA® is an add-on test to any biopsy or cytology where cancer is suspected or diagnosed. This test provides precision diagnostic, prognostic, and therapeutic guidance for multiple types of cancers.</p> <p>This canine cancer genomic panel identifies mutations in 120 relevant cancer genes and provides insights into cancer’s origin, behavior, and the optimal approach for treatment.</p> <p>Report includes:</p> <ul style="list-style-type: none">• Information about mutations in 120 cancer genes including diagnostic, prognostic, and therapeutic biomarker associations with supporting evidence from peer-reviewed literature, clinical consensus, and inference from human FDA guidance.• Pharmacogenomic marker (MDR1) mutation status.• A list of targeted therapeutic drug(s) if the patient’s tumor mutations are a match for an existing drug.• Information on clinical trials by tumor type curated from public databases and individual academic centers by Vidium Animal Health. <p>To order SearchLight DNA, please contact customer service.</p> <p>Note: SearchLight DNA can only be added to a biopsy or cytology performed at Antech where neoplasia has been diagnosed or suspected.</p>	<p>Fine-Needle Aspirate (FNA) Samples:</p> <p>FNA Slides - ≥2 unstained or stained slides + 1 Diffquick</p> <p>FNA Liquid - ≥0.5 mL in sterile vial/tube with no additives</p> <p>Biopsy Samples:</p> <p>Formalin-Fixed, Paraffin-Embedded (FFPE) Scrolls – 10 x 10- micron sections with 1 adjacent hematoxylin and eosin (H&E) slide</p> <p>Unstained FFPE Slides – 10 x 10-micron sections with 1 adjacent H&E slide</p> <p>Bone Tumor – Non-decalcified sample</p>	12-17 days
S14497 Add-on Equivalent S14515	<p>Canine Wisdom Panel™ Premium</p> <p>Canine Wisdom Panel™ Premium DNA panel including:</p> <ul style="list-style-type: none">• 267 genetic health-associated variants (inc. MDR1 mutation, vWD)• 50+ physical trait variants• Breed background detection (350+ breeds)• Genetic diversity (heterozygosity) scoring• Genetic relatives matching with other dogs in the database• Support with interpretation available from Wisdom Panel™	1 Canine Wisdom swab kit (2 swabs per kit)	14-21 days
S14498 Add-on Equivalent S14516	<p>Feline Wisdom Panel™ Complete</p> <p>Feline Wisdom Panel™ Complete DNA panel including:</p> <ul style="list-style-type: none">• 45 genetic health condition-related variants (inc. MDR1 mutation, Polycystic Kidney Disease)• 25 physical traits variants• Identifies genetic blood type• Provides breed (mix) identification (70+ breeds)• Genetic diversity (heterozygosity) scoring, and• Genetic relatives matching with other cats in the database.• Support with interpretation available from Wisdom veterinarians.	1 Feline Wisdom swab kit (2 swabs per kit)	14-21 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S14515	Canine Wisdom Panel™ Premium, Add-on Canine Wisdom Panel™ Premium DNA panel including: <ul style="list-style-type: none"> • 267 genetic health-associated variants (inc. MDR1 mutation, vWD) • 50+ physical trait variants • Breed background detection (350+ breeds) • Genetic diversity (heterozygosity) scoring • Genetic relatives matching with other dogs in the database • Support with interpretation available from Wisdom Panel™ 	1 Canine Wisdom swab kit (2 swabs per kit)	14-21 days
S14516	Feline Wisdom Panel™ Complete, Add-on Feline Wisdom Panel™ Complete DNA panel including: <ul style="list-style-type: none"> • 45 genetic health condition-related variants (inc. MDR1 mutation, Polycystic Kidney Disease) • 25 physical traits variants • Identifies genetic blood type • Provides breed (mix) identification (70+ breeds) • Genetic diversity (heterozygosity) scoring, and • Genetic relatives matching with other cats in the database. • Support with interpretation available from Wisdom veterinarians. 	1 Feline Wisdom swab kit (2 swabs per kit)	14-21 days
KS14497	Canine Wisdom™ Panel with Keyscreen® Canine Wisdom Panel™ Premium and 2x KeyScreen® GI Parasite PCR Panel For Canine only. Wisdom Panel™ and first Keyscreen® must be submitted at the same time. Follow up submissions must be for the same dog and include matching patient information. The follow up Keyscreen® can be submitted within 365 days of the original. Call customer service for submission details: 1-800-872-1001, dial 0.	1 Canine Wisdom swab kit (2 swabs per kit), 0.3 grams of fresh feces in Antech provided fecal container	1-21 days 1-3 days for KeyScreen, 14-21 days for Wisdom Panel

GENERAL

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S16005	Acetylcholine Receptor Antibody Assessment for the presence of antibodies against the acetylcholine receptor. A positive result supports a diagnosis of generalized acquired Myasthenia Gravis. Note: immunosuppressive therapy with corticosteroids instituted prior to submission of results may lower antibody levels and give a false negative result.	1.0 mL serum in red top or serum separator tube	7-9 days
T010	Albumin Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	24 hours
T020	Alkaline Phosphatase Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
T215	Alkaline Phosphatase Isoenzymes Alkaline phosphatase, steroid induced alkaline phosphatase, percent steroid induced alkaline phosphatase This evaluates the percentage of steroid-induced alkaline phosphatase in the serum relative to the total serum alkaline phosphatase. Elevation of steroid-induced alkaline phosphatase is not specific for hyperadrenocorticism and is of limited diagnostic utility in the diagnosis of this disease.	0.5 mL serum in red top or serum separator tube	1-2 days
T030	ALT SGPT Alanine aminotransferase Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
T040	Amylase Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T050	Amylase, Lipase Amylase and Lipase (PSL) Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
S16040	Anticoagulant Screen Use test when ingestion of warfarin, bromadiolone, coumachlor, brodifacoum, diphacinone, chlorophacinone, or difethialone is suspected. Note: send out test requiring serum, stomach contents, bait, or frozen liver tissue.	5.0 mL serum, 20 grams stomach contents or bait, or 10 grams of frozen unfixed liver	10-14 days
T515	Antinuclear Antibodies ANA Antinuclear Antibody test is a titer test used as part of the Systemic Lupus Erythematosus (SLE) evaluation in patients with appropriate history, clinical signs, and other laboratory findings. Note: canine and feline only.	0.5 mL serum in red top or serum separator tube	1-3 days
T060	AST SGOT Aspartate aminotransferase Interferences: marked hemolysis and lipemia. Serum or plasma should be separated within 1 hour of draw and centrifuged	0.5 mL serum in red top or spun serum separator tube	24 hours Performed each shift
SA280	Autoimmune Profile CBC, Antinuclear Antibodies ANA, and Direct Coombs' Test Warm An autoimmune panel evaluating a complete blood count, Direct Coombs test warm, and Antinuclear antibodies. Canine and feline only. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-3 days
SA170	Autoimmune Profile Standard CBC, Antinuclear Antibodies ANA, and Direct Coombs' Test Warm An autoimmune panel evaluating a complete blood count, Direct Coombs test warm, and Antinuclear antibodies. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-3 days
T115	Bicarbonate Interferences: marked hemolysis or lipemia. Serum stable for 1 hour in an ice bath prior. Due to the instability of CO2 in specimens, low levels may not be accurate. Most accurately done within minutes after drawing and on blood gas analyzer.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
DBIL	Bilirubin, Direct Direct bilirubin Interferences: marked hemolysis and lipemia. Lipemia can falsely increase results.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T090	Bilirubin, Total Total bilirubin Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
T520	Bladder Tumor Analytes This test may be used as an adjunct in the diagnosis of transitional cell carcinomas (TCC). The test is a reasonably sensitive (~90%) but not a specific test for TCC. Consider CADET BRAF (T1025). Interferences: false-positive results may occur in the presence of proteinuria, aciduria, urinary tract infection, and hematuria. An alternative test to consider is the CADET BRAF if looking for a non-invasive diagnostic test with a far higher sensitivity (95%) and specificity (99%).	2.0 mL urine in urine transport tube	1-2 days
T525	Blastomyces Antibody To investigate the possibility of blastomycosis as the cause of clinical signs. Note: test is reported as positive or negative. No titer is given.	0.5 mL serum in red top or serum separator tube	2-3 days
S86293	Blastomyces Quantitative Ag Assay <i>Blastomyces dermatitidis</i> Ag, EIA An antigen test used to diagnose blastomycosis.	2.0 mL urine in urine transport tube	5-7 days
T340	Blood Cross Match Includes Donor ID, Major and/or Minor Cross Match Screen for pre-existing antibodies directed against red cell antigens to determine serological compatibility prior to transfusion. If serum, as well as whole blood, is submitted both minor and major cross-match will be performed. This test is not intended for use in determining pre-breeding compatibility testing, neonatal isoerythrolysis, or the prevention of neonatal isoerythrolysis. Interferences: marked hemolysis and lipemia.	1.0 mL of whole blood in lavender top tube and 1.0 mL of serum in red top or serum separator tube for patient and each donor	1-3 days
T345	Cross Match, Additional Donor 2 Donor ID, Major and Minor Cross Match Screens for pre-existing antibodies directed against red cell antigens to determine serological compatibility prior to transfusion. If serum, as well as whole blood, is submitted both minor and major cross-match will be performed. This test is not intended for use in determining pre-breeding compatibility testing, neonatal isoerythrolysis, or the prevention of neonatal isoerythrolysis. Interferences: marked hemolysis and lipemia, clotting or freezing of sample may preclude analysis. Note: this test must be used in addition to Cross Match - T340.	1.0 mL EDTA whole blood in lavender top tube and 1.0 mL serum in red top tube for recipient and each donor	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S16003	Brucella Titer For Export (KSU) Evaluation for antibody suggestive of <i>Brucella canis</i> exposure. <i>Brucella canis</i> antibody evaluation is performed by serum agglutination. The test is typically used when required for export. Canine only. Note: clearly write microchip number on test requisition form (TRF) and on submitted sample tube.	1.0 mL serum in red top tube or serum separator tube	7-10 days
T100	BUN Blood Urea Nitrogen Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top tube or serum separator tube	24 hours Performed each shift
T107	Bun, Creatinine, SDMA, T4 BUN and Creatinine with SDMA, and T4 BUN, creatinine, SDMA for glomerular filtration rate estimation (see T1035) and total T4. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum in red top tube or serum separator tube	24 hours Performed each shift
T105	BUN, Creatinine, SDMA BUN, Creatinine with SDMA and BUN/Creat Ratio BUN, creatinine, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top tube or serum separator tube	24 hours Performed each shift
SA140	Canine Heartworm Program Miniscreen 4 Chem with Electrolytes, Heartworm Antigen Miniscreen chemistry and heartworm antigen detection. Interferences: marked hemolysis and lipemia.	1.0 mL serum in red top or serum separator tube	1-2 days
SA130	Canine Heartworm Program Plus Miniscreen 4 Chem with Electrolytes, CBC, Heartworm Antigen Miniscreen chemistry, a complete blood count, and heartworm antigen detection. Interferences: marked hemolysis and lipemia.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
T690	Canine Parvovirus Antibody Titer Parvovirus IgG, Parvovirus IgM Evaluation of IgM and IgG titers consistent with Canine Parvovirus exposure. Titers must be interpreted in light of previous vaccination and clinical signs. Interferences: marked hemolysis or marked lipemia if unable to clear by centrifugation.	0.5 mL serum in red top or serum separator tube	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T1010	Cardio BNP-Canine Antech Cardio BNP - Canine This test provides a quantitative assessment of canine heart health by measuring the concentration of the carboxy terminal of the BNP peptide (cBNP), which is released by cardiac myocytes in response to stretch and stress. Note: draw 2.0 mL of whole blood into an LT tube, centrifuge for plasma, draw 0.5 mL plasma from the tube, and inject into closed BNP tube. Keep sample cold. It requires a specific BNP tube for submission.	0.5 mL EDTA plasma injected into BNP tube	2-3 days
T120	Chloride Interferences: marked lipemia. Serum should be separated within 1 hour after draw.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
T125	Cholesterol Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
T235	Cholinesterase - Serum Test used to diagnose organophosphate toxicity.	0.5 mL serum in red top or serum separator tube	5-7 days
T840	Cobalamin B12 Generally used as part of a larger panel when evaluating for malassimilation/maldigestion. Note: ideally freeze and ship on ice. Otherwise submit fresh on ice.	0.5 mL serum in red top or serum separator tube	1-2 days
S16195	Cobalamin B12/Folate Used to further evaluate gastrointestinal signs where malassimilation and maldigestion are suspected. Additional measurement of TLI is recommended to ensure the most accurate interpretation of results. Note: pets should fast overnight prior to sampling. Hemolysis can affect folate results. Previous vitamin supplementation will cause elevated concentrations. Send/ship test on ice packs.	0.5 mL serum in red top or serum separator tube	1-2 days
SA160	Cobalamine, Folate, TLI Canine Used to further evaluate gastrointestinal signs where malassimilation and maldigestion are suspected. Note: pet should fast overnight. Hemolysis of sample may affect folate results. Previous vitamin supplementation will cause elevated concentrations.	1.0 mL serum in red top or serum separator tube	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA275	Cobalamine, Folate, TLI Feline Used to further evaluate gastrointestinal signs where malassimilation and maldigestion are suspected. Note: pet should fast overnight. Hemolysis of sample may affect folate results. Previous vitamin supplementation will cause elevated concentrations.	1.5 mL serum in red top or serum separator tube	1-9 days
S16210	Copper Level This is used to detect exposure to excessive copper. Note: submit sample on ice. In the case of a primary copper, hepatopathy only a quantitative copper on liver tissue combined with histopathology can confirm the diagnosis.	1.0 mL serum or plasma, 5 grams fresh liver, 50 mg fresh liver biopsy (approx. 3 Tru-Cut samples), 3.0 mL urine, or 500 g Feed. For birds: 0.3 mL serum or plasma is adequate. For ruminants and camelids: 5 grams fresh kidney is required.	10-14 days
S16215	Copper Storage Disease Tissue Copper Quantification To be used when a quantitative copper level is required on hepatic tissue. Note: a minimum of 5 grams of fresh liver in saline is the ideal sample.	5 grams fresh liver tissue placed in saline	10-14 days
S16225	Corona Virus IgG and IgM Evaluation for IgM and IgG antibody to canine coronavirus.	0.5 mL serum in red top or serum separator tube	3-5 days
T130	CPK Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
T135	Creatinine Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
SA350	Crypto, Giardia, Clostridium Enterotoxin <i>Clostridium perfringens</i> Enterotoxin, <i>Giardia</i> (FA), <i>Cryptosporidium</i> (FA), <i>Giardia</i> ELISA A panel that can be considered when suspecting <i>Clostridium perfringens</i> , <i>Cryptosporidium</i> spp., or <i>Giardia</i> spp. as the cause for gastroenteritis. Submission includes <i>Cryptosporidium</i> spp. <i>Giardia</i> spp. evaluated by FA, <i>Clostridium perfringens</i> enterotoxin assessment and a <i>Giardia</i> ELISA Antigen detection. Note: the specimen must arrive cold (send on ice) and must be tested within 24 hours of collection for accuracy with regards to Clostridial enterotoxin detection.	10 grams feces in Antech provided fecal container	1-2 days
T350	D-Dimer This test is used to measure the concentration of D dimers. Interferences: marked hemolysis or lipemia. Note: moderately to markedly lipemic samples can cause D dimer results to be spuriously decreased. Citrated Whole Blood (Blue Top Tube) or citrated plasma is the only acceptable sample. The tube should be >2/3rds filled. If submitting separated citrated plasma, label it as Citrated Plasma.	0.5 mL citrated plasma or citrated whole blood in blue top tube at least 2/3rds full to the fill line	24 hours Performed each shift
ADD270	Distemper Parvo Titer Add-on This test detects IgG antibodies against canine parvovirus and canine distemper virus using ELISA. Results are reported as positive (protective) or negative (insufficient neutralizing antibody detected). This test is not to be used in patient suspected of having distemper or parvovirus infection. Interferences: marked hemolysis or lipemia. This test does not give an endpoint titer result. This test is only meant to be used to determine if the amount of antibody present would be sufficient to protect the patient in the face of exposure. This test should not be used on dogs with current or recent clinical signs of distemper or parvovirus infection.	0.5 mL serum in red top or serum separator tube	1-3 days
T565 Add-on Equivalent ADD270	Distemper Parvo Vaccinal Titer This test detects IgG antibodies against canine parvovirus and canine distemper virus using ELISA. Results are reported as positive (protective) or negative (insufficient neutralizing antibody detected). This test is not to be used in patient suspected of having distemper or parvovirus infection. Interferences: marked hemolysis and lipemia. This test should not be used on dogs with current or recent clinical signs of distemper or parvovirus infection as a positive may indicate a response to infection rather than protection.	0.5 mL serum in red top or serum separator tube	1-3 days
S16250	Distemper Smear, FA This test is helpful in making a diagnosis of distemper, in the absence of recent vaccination, when used in conjunction with clinical signs.	EDTA whole blood in lavender top tube, or unstained peripheral blood, buffy coat or urine sediment smears	1-5 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T560	Distemper Vaccinal Titer This test detects IgG antibodies against canine distemper virus. The result is reported as positive or negative. A titer is not reported. This test is not to be used in patients suspected to be ill from distemper virus infection. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	1-2 days
T675	Distemper/Parvo IgG End Point Titer Distemper IgG End PT Titer, Parvo IgG End Point This test provides an endpoint IgG titer for Canine Distemper and Canine Parvovirus exposure. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	1-2 days
S16865	FIV Antibody, Western Blot This is a western blot test to detect antibodies directed against the Feline Immunodeficiency Virus. It may be used as a confirmatory test for the presence of antibodies against FIV. Note: does not differentiate cats infected with FIV from cats vaccinated against or kittens with passive transfer of antibodies against FIV from the queen.	0.5 mL serum in red top or serum separator tube	1-2 days
SFUN	Fungal Serology Fungal Serology (<i>Histoplasma</i> , <i>Blastomyces</i> , <i>Aspergillus</i>) Detects the presence of antibodies to <i>Histoplasma</i> , <i>Blastomyces</i> , and <i>Aspergillus</i> . Detection of antigen is more sensitive than antibody for the diagnosis of histoplasmosis and blastomycosis.	0.5 mL serum in red top or serum separator tube	2-3 days
SA340	Fungal Serology with Cocci Fungal Serology (<i>Histoplasma</i> , <i>Blastomyces</i> , <i>Aspergillus</i>), Coccidioidomycosis Screen and Titer Detects the presence of antibodies to <i>Histoplasma</i> , <i>Blastomyces</i> , <i>Aspergillus</i> , and <i>Coccidioides</i> spp. Detection of antigen is more sensitive than antibody for the diagnosis of histoplasmosis and blastomycosis. Interferences: marked lipemia. A titer is only reported for <i>Coccidioides</i> . All other results are reported as positive or negative for antibody detection. Both the <i>Blastomyces</i> and <i>Histoplasma</i> antigen tests are preferred over antibody testing as a diagnostic test.	1.0 mL serum in red top or serum separator tube	2-5 days
T790	Giardia (FA) and Cryptosporidium (FA) Giardia IFA and Cryptosporidium IFA This assay uses fluorescence-labeled monoclonal antibodies to detect <i>Giardia</i> and <i>Cryptosporidium</i> cysts in stool samples. Note: <i>Giardia</i> FA is less sensitive than the <i>Giardia</i> ELISA and O&P tests. <i>Cryptosporidium</i> FA is likely more sensitive than an O&P test. The GI PCR panel would be preferred to either FA test with regards to sensitivity and specificity.	6 grams fresh fecal specimen in Antech provided fecal container	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T150	Glucose Interferences: hemolysis or lipemia. Hemolysis may cause false decreases in glucose. Lipemia may cause false increases in glucose.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
S16400	Herpes - Conjunctival Smear Conjunctival smears are used to evaluate for the presence of herpes virus antigen in epithelial cells via immunofluorescence. Note: ensure slides have adequate cellularity.	2 conjunctival smears	9-12 days
S86022	Herpes Antibody, IFA Detection of Herpes antibody by IFA.	0.3 mL serum in red top or serum separator tube	5-7 days
S17029	Heska ALLERCEPT Environmental & Food Panel	3.0 mL serum in red top or serum separator tube	4-5 business days
S17026	Heska ALLERCEPT Environmental Panel	2.0 mL serum in red top or serum separator tube	4-5 business days
S17028	Heska ALLERCEPT Food Panel	1.0 mL serum in red top or serum separator tube	4-5 business days
S17027	Heska Equine ALLERCEPT Panel	3.0 mL serum in red top or serum separator tube	4-5 business days
T660	Immunoglobulins A, G, M This test quantifies IgG, IgM and IgA levels by radial immunodiffusion in dogs. Interferences: marked hemolysis or lipemia. Note: reference intervals for adult dogs and cats have been established. However, if submitting serum for a young animal, age-matched controls should be submitted concurrently. References intervals for other species are provided by the manufacturer and are age-dependent. It is imperative to provide the entire signalment, including age, when submitting samples.	0.5 mL serum in red top tube	7-9 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T155	Iron, Serum Interferences: marked hemolysis.	0.5 mL serum in red top tube or serum separator tube	1-2 days
T160	LDH Lactic Dehydrogenase Lactate dehydrogenase (LDH) Interferences: hemolysis. Hemolysis or delayed separation can falsely elevate results.	0.5 mL serum in red top tube or serum separator tube	24 hours
T745	Lead Level Used to detect exposure to lead.	0.5 mL of EDTA whole blood in lavender top tube	3-5 days
T170	Magnesium Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
S16535	Masticatory Muscle Myositis Masticatory Muscle Myositis - 2M Antibody This test is used to diagnose masticatory myositis by detecting autoantibodies specifically directed against masticatory muscle proteins. Note: this assay may be negative if immunosuppressive dosages of cortisocosteroids have been given for longer that 7-10 days prior to test submission.	0.5 mL serum in red top or serum separator tube	12-15 days
T810	Occult Blood, Feces This is a test for the detection of occult blood in fecal samples. The test is not highly sensitive and is dependent on the amount of heme present and how evenly the heme is spread throughout the sample. Note: diets of red meat, chicken, or fresh/uncooked vegetables may cause false-positive results and should not be eaten for 48 hours prior to testing.	3 grams feces in Antech provided fecal container	1-2 days
S16575	Osmolality - Serum This may be of utility in the evaluation of the causes of PU/PD.	1.0 mL serum in red top or serum separator tube	7-10 days
S85364	Pancreatic Lipase Immunoreactivity Measures concentrations of pancreatic lipase in the serum. Interferences: hemolysis. Note: a 12 hour fast is recommended prior to sample submission.	0.5 mL serum in red top or serum separator tube	5-7 business days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S86468	<p>Pancreatitis Profile - Canine</p> <p>Cobalamin B12/Folate, Canine cTLI, Pancreatic Lipase Immunoreactivity</p> <p>Includes cobalamin (B12), folate, trypsin like immunoreactivity (TLI), and pancreatic lipase immunoreactivity (cPLI).</p> <p>Note: pet should fast overnight. Hemolysis of sample may affect folate results. Previous vitamin supplementation will cause elevated concentrations.</p>	1.5 mL serum in red top or serum separator tube	2-10 days
S86288	<p>Pancreatitis Profile - Feline</p> <p>Cobalamin B12/Folate, Feline fTLI, Pancreatic Lipase Immunoreactivity</p> <p>Includes cobalamin (B12), folate, trypsin like immunoreactivity (TLI), and pancreatic lipase immunoreactivity (fPLI).</p>	1.5 mL serum in red top or serum separator tube	2-10 days
S16580	<p>Panleukopenia Titer IgG, IgM</p> <p>This test evaluates for the presence of IgM and IgG antibodies to the feline panleukopenia virus. It will be positive either due to infection or vaccination against panleukopenia. Kittens may test positive due to transmammary transfer of antibodies. This test will also detect antibodies directed against canine parvovirus in cats infected with CPV-2a and CPV-2b.</p> <p>Interferences: marked hemolysis or lipemia (if unable to clear with centrifugation).</p>	0.5 mL serum in red top or serum separator tube	1-5 days
S16053	<p>Panleukopenia Vaccinal Titer</p> <p>The test detects IgG antibodies against Panleukopenia via Immunofluorescent antibody assessment (IFA). A titer of 1:5 or greater, in the absence of clinical signs of disease, indicates an immunological response to vaccination. A vaccine titer of less than 1:5 indicates a low level of circulating antibodies.</p> <p>Interferences: marked hemolysis or lipemia (if unable to clear with centrifugation).</p>	0.5 mL serum in red top or serum separator tube	3-5 days
Add-on Equivalent ADD280			
ADD280	<p>Panleukopenia Vaccinal Titer Add-on</p> <p>The test detects IgG antibodies against Panleukopenia via Immunofluorescent antibody assessment (IFA). A titer of 1:5 or greater, in the absence of clinical signs of disease, indicates an immunological response to vaccination. A vaccine titer of less than 1:5 indicates a low level of circulating antibodies.</p> <p>Interferences: marked hemolysis or lipemia (if unable to clear with centrifugation).</p>	0.5 mL serum in red top or serum separator tube	3-5 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T825	Parasite Identification Note: document where on the animal the parasite was noted.	Parasite sent in 70% ethanol. Include host species and body location where parasite was collected. If submitting slides with worms, place slides in a leak proof container with 70% ethanol. Other representative sample types include feather, intestine, liver, fresh tissue, urine, wash/lavage, feces or vomitus Jar with formalin or alcohol, slides in 70% alcohol	10-12 days
T700	Parvo Antibody, Antigen Parvovirus Antigen, Canine Parvovirus Antibody Titer Detects parvovirus antigen in fecal samples and evaluates a concurrent serum sample for IgM and IgG titers consistent with Canine Parvovirus exposure. Note: modified live virus vaccination for parvovirus may give false-positive results for about two weeks after vaccination. A negative test does not rule out parvovirus infection.	0.5 mL serum and 5 grams feces Serum in red top or serum separator tube, Antech provided fecal container	1-3 days
T705	Parvovirus Vaccinal Titer Used to evaluate current IgG levels to Parvovirus in light of previous vaccination. Interferences: marked hemolysis or lipemia. Note: not to be used in dogs with current or recent clinical signs of parvovirus infection as a positive may indicate response to infection rather than protection.	0.5 mL serum in red top or serum separator tube	1-3 days
T180	Phosphorus Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
T185	Potassium Interferences: hemolysis and marked lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T245	<p>Protein Electrophoresis (Urine)</p> <p>Urine Protein Electrophoresis (Total Protein, Albumin, Globulin, Alpha 1, Alpha 2, Beta 1, Gamma 1 with Interpretation)</p> <p>Evaluation for Bence Jones proteins in Urine.</p>	1.0 mL urine in urine transport tube	2-4 days
ADD130	<p>Protein Electrophoresis Add-on</p> <p>Total protein, Albumin, Globulin, Alpha 1, Alpha 2, Beta 1, and Gamma 1 fraction assessment with interpretation.</p> <p>An evaluation of the globulin fraction of the serum (alpha 1, alpha 2, beta, and gamma) to determine if the globulin fraction is monoclonal based on these components.</p> <p>Interferences: hemolysis may interfere with results. Serum must be less than 48 hours of age to add-on this test.</p>	0.5 mL serum in red top or serum separator tube	2-4 days
T9810	<p>Relaxin</p> <p>Relaxin is a hormone produced by the placenta, and detection of relaxin in serum or plasma is a sensitive and specific test for pregnancy in dogs and cats.</p> <p>Interferences: marked hemolysis.</p>	0.2 mL serum in red top or serum separator tube	2-3 days
S16702	<p>Rhinotrachitis Feline</p> <p>Viral Neutralization Antibody Titer</p> <p>This is a semi-quantitative titer that is correlated to the amount of protective humoral antibody present in relation to previous feline rhinotracheitis vaccination (feline herpes virus).</p>	1.0 mL serum in red top or serum separator tube	10-12 days
S16730	<p>Selenium Level</p> <p>Interferences: marked hemolysis or lipemia.</p>	1.0 mL serum in red top or serum separator tube (spun), 1.0mL EDTA whole blood in lavender top tube, or 10.0 grams of feed (representative sample)	7-10 days
T195	<p>Sodium</p> <p>Interferences: marked hemolysis or lipemia.</p>	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
T200	<p>Sodium And Potassium Combo</p> <p>Sodium, Potassium, and NA/K Ratio</p> <p>Sodium and potassium evaluation</p> <p>Interferences: hemolysis and marked lipemia. Hemolysis will falsely elevate results.</p>	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T485	T3 Autoantibodies Detection of T3 autoantibody is used to assess the presence of autoimmune thyroid disease. This test is less sensitive to autoimmune thyroid disease than the thyroglobulin autoantibody (TGAA) test.	0.5 mL serum in red top or serum separator tube	2-3 days
T500	T4 Autoantibodies Detection of T4 autoantibody is used to assess the presence of autoimmune thyroid disease. This test is less sensitive to autoimmune thyroid disease than the thyroglobulin autoantibody (TGAA) test.	0.5 mL serum in red top or serum separator tube	2-4 days
S16755	Taurine Measurement of taurine concentration in whole blood via HPLC. Note: do not fast prior to sample submission. Whole blood submitted in a green-top refrigerated is preferred. EDTA may also be used but taurine concentrations in EDTA anticoagulated blood will be slightly lower than in a GTT.	1.0 mL heparinized whole blood in green top tube	10-14 days
SA330	Tick Borne Disease Panel <i>Ehrlichia canis</i> , Lyme Titer IgG, and Rocky Mountain Spotted Fever IFA Detects the presence of antibodies to <i>Ehrlichia canis</i> , <i>Borrelia burgdorferi</i> (Lyme), and Rocky Mountain Spotted fever. Interferences: marked hemolysis or lipemia.	1.0 mL serum in red top or serum separator tube	1-2 days
S16800	TLI, Feline Trypsin-Like Immunoreactivity Use to confirm feline exocrine pancreatic insufficiency (EPI). Note: recommend fasting for 8-12 hours prior to drawing sample. Spin and separate serum from clot.	1.0 mL serum in red top or serum separator tube	7-10 days
T190	Total Protein Interferences: marked hemolysis or lipemia	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
T127	Triglyceride, Chol. Panel Triglyceride and cholesterol Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
T205	Triglycerides Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T230	Trypsin-Like Immunoreac Canine cTLI Canine Trypsin- Like Immunoreactivity (cTLI is) used to diagnose exocrine pancreatic insufficiency. Note: a 12-hour fast is required. Canine only.	0.5 mL serum in red top or serum separator tube	1-2 days
S16870	Zinc Evaluation of blood zinc level Note: plasma or serum needs to be separated from cells and transferred to a plastic vial. Do not use tubes with rubber stoppers.	0.3 mL of serum in royal blue top tube or plastic transfer tube	7-10 days
T730 Add-on Equivalent ADD320	Bromide This test is used to monitor bromide therapy. Therapeutic serum bromide concentrations should be measured at three weeks after initiating therapy, but steady-state concentrations may fluctuate among dogs due to differences in drug clearance and bioavailability. Interferences: gel may interfere with test. Do not use a serum separator tube for sample draw or submission.	0.5 mL serum in red top tube. Submission in serum separator tube is not recommended.	1-2 days
ADD320	Bromide Add-on This test is used to monitor bromide therapy. Therapeutic serum bromide concentrations should be measured at three weeks after initiating therapy, but steady-state concentrations may fluctuate among dogs due to differences in drug clearance and bioavailability.	0.5 mL serum in red top tube. Submission in serum separator tube is not recommended.	1-2 days
S18702	Cyclosporine Used to assess the cyclosporine drug level being achieved, at the current drug dosing, in a given patient. Note: peak drug level occurs 2 hours post medication. Trough drug level will occur just prior to the next 12-hour dose.	2.0 mL EDTA whole blood in lavender top tube	7-10 days
T735	Digoxin Digoxin level Interferences: hemolysis, gel may interfere with test. Do not use a serum separator tube for sample draw or submission.	0.5 mL serum in red top tube collected 8 hours post pill	1-2 days
S86541	Keppra Level Interferences: hemolysis and gel in serum separator tube may interfere with test.	2.0 mL serum in red top or non-additive tube. Submission in serum separator tube is not recommended	5-7 days
T750 Add-on Equivalent ADD315	Phenobarbital Phenobarbital level Interferences: hemolysis and gel in serum separator tube may interfere with test. Do not use a serum separator tube for sample draw or submission.	0.5 mL serum in red top tube. Submission in serum separator tube is not recommended.	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
ADD315	Phenobarbital Add-on Phenobarbital level Interferences: hemolysis and gel in serum separator tube may interfere with test.	0.5 mL serum in red top tube. Submission in serum separator tube is not recommended.	1-2 days
SA830	Phenobarbital Panel Plus Liver Chemistry with Electrolytes, CBC, Bile Acids, Phenobarbital Includes a chemistry (SA324), complete blood count, a single bile acids, and a phenobarbital level. Interferences: marked hemolysis and lipemia. All therapeutic drug monitoring should be submitted as serum in a red top tube. A serum separator tube should not be used in sample collection as the drug concentration being measured may be erroneously decreased by as much as 30% due to binding with gel. Ursodeoxycholic acid may be detected by bile acid assay, causing falsely elevated values.	1.5 mL serum and 1.0 mL EDTA whole blood Serum in red top tube, lavender top	1-2 days
S86480	Zonisamide Zonegran Zonisamide concentration Interferences: hemolysis, gel may interfere with test.	0.5 mL serum in red top tube. Submission in serum separator tube is not recommended.	3-5 days

*SEE [ANTECHDIAGNOSTICS.COM/TERMS-SERVICE](https://www.antechediagnostics.com/terms-service) FOR MORE DETAILS ON TURNAROUND TIMES

AVIAN • EXOTICS

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S17116	Adrenal Androgen Panel For Ferrets Estradiol, 17-OH Progesterone, Androstenedione	0.5 mL heparinized plasma spun in green microvial with gel separator	10-14 days
S16011	Aspergillus AB Avian (AFMP1P ELISA) Avian <i>Aspergillus</i> Antibody	0.5 mL serum in microtainer or regular serum separator tube (spun)	5-7 days
S85358	Aspergillus Ag, Galactomannan Assay <i>Aspergillus</i> Antigen, Galactomannan Assay (Avian) Antigen detection test (EIA immunoenzymatic sandwich microplate assay) for galactomannan.	Canine or Feline: 0.8 mL serum in red top or serum separator tube Avian: 0.1 mL serum in red top or serum separator tube	7-10 days
S85359	Aspergillus Profile Avian <i>Aspergillus</i> Antibody Avian, <i>Aspergillus</i> Ag. Galactomannan Assay, Avian Protein Electrophoresis	1.0 mL serum or 1.0 mL heparinized plasma 2 green top micro vials (with gel) OR 2 serum in red top or serum separator tubes	2-10 days
AE290	Avian Blood Lead	0.5 mL of heparinized whole blood collected in green top tube	3-5 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
AE021	Avian Comprehensive Profile with Bile Acid	0.75 mL of heparinized plasma, 0.5 mL of heparinized whole blood and 2 freshly prepared blood smears 2 green top micro vials (with gel), green top micro vial (no gel) or 2 hematocrit tubes, slide holders	24 hours Performed each shift
AE025	Avian Comprehensive Profile with Bile Acid, EPH Comprehensive Avian Profile with Bile Acids (Avian/ Exotic), and Avian Protein Electrophoresis	0.75 mL of heparinized plasma, 0.5 mL of heparinized whole blood and 2 freshly prepared blood smears 2 green top micro vials (gel), green top micro vial (no gel) or 2 hematocrit tubes, slide holders	1-3 days
AE010	Avian Comprehensive Chemistries Total Protein, Albumin, Globulin, SGOT (AST), Phosphorus, Glucose, Calcium, Sodium, Potassium, Chloride, Cholesterol, CPK, Uric Acid	0.5 mL heparinized plasma in spun green top micro vial with gel	24 hours Performed each shift
AE300	Avian Protein Electrophoresis Includes Total Protein, Pre-Albumin, Albumin, Alpha 1, Alpha 2, Beta and Gamma	0.5 mL heparinized plasma in green micro vial (with gel)	2-3 days
AE040	Avian Standard Chemistries Includes Total Protein, SGOT (AST), Phosphorus, Glucose, Calcium, CPK, and Uric Acid	0.5 mL heparinized plasma in green micro vial (with gel)	24 hours Performed each shift

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
AE051	Avian Standard Profile with Bile Acid Standard Avian Profile with Electrolytes, Avian/ Exotic CBC, and Avian/Exotic Bile Acids	0.75 mL of heparinized plasma, 0.5 mL of heparinized whole blood and 2 freshly prepared blood smears 2 green top micro vials (with gel), green top micro vial (no gel) or 2 hematocrit tubes, slide holders	1-2 days
AE055	Avian Standard Profile with Bile Acid, EPH Standard Avian Profile with Electrolytes, Bile Acids Avian/ Exotic, and Avian Protein Electrophoresis	0.75 mL of heparinized plasma, 0.5 mL of heparinized whole blood and 2 freshly prepared blood smears 2 green top micro vials (with gel), green top micro vial (no gel) or 2 hematocrit tubes, slide holders	1-3 days
S16012	Avian Zinc Assay	0.5 mL serum in yellow micro vial with gel	7-10 days
AE270	Avian/Exotic CBC Includes WBC Estimate, Thrombocyte Estimate, Hematocrit, Blood Parasites, and Differential: Percent and Absolute Heterophils, Bands, Lymphocytes, Monocytes, Eosinophils, Basophils and Azurophilic Monocytes Interferences: marked hemolysis. EDTA whole blood is not preferred because it tends to cause hemolysis in reptiles and in some avian species	0.5 mL of heparinized whole blood in green top tube with 2 freshly prepared blood smears	24 hours Performed each shift
AE260	Bile Acids Avian/Exotic Single Bile Acids, Avian/Exotic Interferences: marked hemolysis or lipemia. Ursodeoxycholic acid may be detected by bile acid assay, causing falsely elevated values.	0.5 mL serum in red top or serum separator tube	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S16671	Chlamydia Antibody Titer IFA	0.5 mL serum in red top or serum separator tube	7-10 days
S16788	Chlamydia PCR Blood	0.5 mL heparinized whole blood in green top tube	10-12 days
S16672	Chlamydia PCR Swab	Avian: Swab in red top tube (combined choanal & cloacal) Feline & other species: conjunctival swab in red top	7-10 days
S85206	Chlamydia Profile <i>Chlamydophila</i> PCR Blood, <i>Chlamydophila</i> Antibody Titer IFA, <i>Chlamydophila</i> Titer EBA, <i>Chlamydophila</i> PCR Swab, <i>Chlamydophila</i> Titer EBA, <i>Chlamydophila</i> PCR Swab	0.4 mL heparinied whole blood, 0.4 mL serum, swab in red top tube (combined choanal & cloacal) Green top micro vial (with gel), green top micro vial (no gel), swab (Combined Choanal & Cloacal)	7-14 days
S16670	Chlamydia Titer EBA Psittacosis	0.5 mL serum in serum separator tube	7-14 days
AE030	Comprehensive Avian Post Purchase Avian Comprehensive Chemistries, <i>Chlamydophila</i> Titer EBA, <i>Giardia</i> ELISA, Gram Stain, Psittacine Beak and Feather Disease, Polyoma PCR Swab, Avian Protein Electrophoresis, Avian/Exotic CBC	1.0 mL heparinized plasma, 1.0 mL heparinized whole blood, 0.5 grams feces, and 2 fecal swabs Antech provided fecal container, 3 green top micro vials (with gel), green top micro vial (whole blood), 2 hematocrit tubes, 2 slides, swabs in white tube (2 fecals)	1-10 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
AE020	Comprehensive Avian Profile Avian/Exotic CBC, Avian Comprehensive Chemistries	0.5 mL of heparinized plasma, 0.5 mL of heparinized whole blood and 2 freshly prepared blood smears Green top micro vial (with gel), green top micro vial (no gel) or 2 hematocrit tubes, slide holder	24 hours Performed each shift
AE200	Comprehensive Mammalian Profile CBC Small Mammalian, Mammalian Comprehensive Chemistries	0.5 mL serum and 1.0 mL whole blood Serum in red top or serum separator OR green top micro vial (with gel), lavender or green top vial (no gel)	24 hours Performed each shift
RREP	Comprehensive Reptilian (Recheck) Avian/Exotic CBC, Reptilian Comprehensive Chemistries	0.5 mL of heparinized plasma, 0.5 mL of heparinized whole blood, and 2 freshly prepared blood smears Green top micro vial (with gel), green top micro vial (no gel) or 2 hematocrit tubes, slide holder	24 hours Performed each shift
AE160	Comprehensive Reptilian Profile Avian/Exotic CBC, Reptilian Comprehensive Chemistries	0.5 mL of heparinized plasma, 0.5 mL of heparinized whole blood, and 2 freshly prepared blood smears Green top micro vial (with gel), green top micro vial (no gel) or 2 hematocrit tubes, slide holder	24 hours Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
AE070	Diarrhea Profile Avian Comprehensive Avian Profile, Gram Stain, Culture, Aerobic Culture, <i>Giardia</i> ELISA Interferences: marked hemolysis and lipemia.	0.5 mL heparinized plasma, 0.5 mL whole blood, 2 fecal slides, 0.5 grams of feces, and culturette Green top micro vial (with gel), green top micro vial (no gel) or 2 hematocrit tubes, slide holder, culturette	1-2 days
S16501	Distemper PCR	2.0 mL urine in urine transport tube, EDTA whole blood in lavender top tube, CSF OR other body cavity fluid in red top or lavender top tube	5-9 days
S16107	Distemper Titer-Virus Neutraliz.	0.5 mL serum in red top or serum separator tube	14-21 days
S16877	<i>Encephalitozoon cuniculi</i> IgG Ab Interferences: marked hemolysis.	0.10 mL spun serum or spun heparinized plasma in plain, non-additive transport tube	3-5 days
AE080	Feather Picker Profile Comprehensive Avian Profile, Culture (Aerobic), <i>Giardia</i> ELISA, Gram Stain, Psittacine Beak and Feather Disease, Avian Protein Electrophoresis	1.0 mL heparinized plasma, 1.0 mL heparinized whole blood and 2 freshly prepared blood smear, culturette, 2 fecal slides and 0.5 grams of feces 2 green top micro vial (with gel), green top micro vial (whole blood) or 2 hematocrit tubes, slide holders, Antech provided fecal container, culturette	1-7 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
AE230	Geriatric Weak Ferret Profile Mammalian Comprehensive Chemistries, Insulin Ferret, Urinalysis, CBC Small Mammalian	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator OR 2 green top micro vial (gel), lavender or green top vial (no gel), urine transport tube	1-3 days
AE090	Hepatic Profile Comprehensive Avian Profile, Bile Acids Avian/Exotic, <i>Chlamydophila</i> Titer EBA, Avian Protein Electrophoresis	1.0 mL heparinized plasma (in green microtainer with gel separator) and 0.5 mL heparinized whole blood	1-14 days
AE190	Mammalian Comprehensive Chemistries Total Protein, Albumin, Globulin, SGOT (AST), SGPT (ALT), Alk Phos, T. Bilirubin, Bun, Creatinine, Phosphorus, Glucose, Calcium, Sodium, Potassium, Chloride, Cholesterol, CPK	0.5 mL serum in red top or serum separator OR green top micro vial (gel)	24 hours Performed each shift
AE210	Mammalian Standard Chemistries Total Protein, SGPT (ALT), Alk Phos, T. Bilirubin, Bun, Creatinine, Phosphorus, Glucose, Calcium	0.5 mL serum in red top or serum separator OR green top micro vial (gel)	24 hours Performed each shift
AE060	Mini Avian Post Purchase Comprehensive Avian Profile, Gram Stain, Avian Protein Electrophoresis	0.75 mL heparinized plasma, 0.5 mL heparinized whole blood and 2 freshly blood smears, and fecal slides Green top micro vial (with gel), green top micro vial (no gel) or 2 hematocrit tubes, slide holder, culturette	1-3 days
S16789	Mycoplasma PCR	Tissue (fresh) OR fluid (ocular, conjunctival, BAL) OR nasal swab	10-12 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S16085	Psittacine Beak and Feather Disease	0.5 mL EDTA or heparinized whole blood in lavender or green top tube	7-10 days
AE110	PU/PD Profile Avian Comprehensive Chemistries, Avian/Exotic CBC, Avian Protein Electrophoresis, Urinalysis	0.5 mL heparinized plasma, 0.5 mL heparinized whole blood and 2 freshly preapred blood smear, and 1 mL urine Green top micro vial (with gel), green top micro vial (no gel) or 2 hematocrit tubes, slide holder, urine transport tube	1-3 days
AE240	Rabbit Neurologic Profile CBC Small Mammalian, Mammalian Comprehensive Chemistries, <i>Encephalitozoon cuniculi</i> IgG Ab, Pasteurella Antibody Titer	1.0 mL of serum or heparinized plasma and 1.0 mL EDTA whole blood 3 green top micro vials (with gel), lavender top or green top vial (no gel)	1-10 days
RECHECKAV	Recheck Comprehensive Avian Avian/Exotic CBC, Avian Comprehensive Chemistries	0.5 mL heparinized plasma and 0.5 mL heparinized whole blood Green top micro top vial (gel), green top micro vial (no gel) or 2 hematocrit tubes, slide holders	24 hours Performed each shift
AE150	Reptilian Comprehensive Chemistries Total Protein, Albumin, Globulin, SGOT (AST), Bun, Phosphorus, Glucose, Calcium, Sodium, Potassium, Chloride, CPK, and Uric Acid Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator OR green top micro vial (gel)	24 hours Performed each shift
AE170	Reptilian Standard Chemistries Total Protein, SGOT (AST), Phosphorus, Glucose, Calcium, CPK, and Uric Acid Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator OR green top micro vial (gel)	24 hours Performed each shift

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
AE050	Standard Avian Profile Avian/Exotic CBC, Avian Standard Chemistries with Electrolytes Interferences: marked hemolysis or lipemia.	0.5 mL heparinized plasma and 0.5 mL heparinized whole blood Green top micro vial (gel), 1 green top vial (no gel) or 2 hematocrit tubes with 2 slides	24 hours Performed each shift
AE220	Standard Mammalian Profile CBC Small Mammalian, Mammalian Standard Chemistries with Electrolytes Interferences: marked hemolysis or lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator OR green micro vial (gel), lavender or green vial (no gel)	24 hours Performed each shift
AE180	Standard Reptilian Profile Avian/Exotic CBC, Reptilian Standard Chemistries with Electrolytes Interferences: marked hemolysis or lipemia.	0.5 mL heparinized plasma and 0.5 mL heparinized whole blood with 2 freshly prepared blood smears Green top micro vial (with gel), 1 green top vial (no gel) or 2 hematocrit tubes, slide holder	24 hours Performed each shift
S16792	Toxoplasma antibody - Exotics	1.0 mL serum in red top or serum separator tube	7-14 days
S85448	West Nile Titer (PRNT & IgM ELISA)	1.0 mL serum in red top or serum separator tube (mammals) OR 0.1 mL heparinized plasma (avians)	7-14 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S85449	West Nile Virus PCR	Equine: 1.0 mL EDTA whole blood in lavender top tube, CSF, brain or spinal cord. Avian: Kidney in sterile container, EDTA or heparinized whole blood, or tissue (brain or heart) Other mammals: EDTA whole blood, CSF, brain or spinal cord	3-7 days
S16095	Zoogen Avian DNA Sexing	Blood spot on ASEX TRF 0.3 mL EDTA whole blood in lavender top tube or heparinized whole blood in green top tube or green top micro vial, heparinized hematocrit tube (1/3rd full), or egg shell (with blood membranes attached). Feathers are acceptable, but the DNA extraction rate is low.	3-6 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

EQUINE • LARGE ANIMAL ESSENTIAL PANELS

	EQUINE CHEMISTRY L010	EQUINE PERFORMANCE HORSE CHEMISTRY L601	EQUINE HEPATIC SCREEN L225	EQUINE RENAL SCREEN L240	RUMINANT CHEMISTRY L1001
A/G Ratio	●	●			●
Albumin	●	●	●	●	●
Alkaline Phosphatase	●	●	●		●
AST (SGOT)	●	●	●		●
BUN	●	●	●	●	●
BUN/Creatinine Ratio	●	●			●
Calcium	●			●	●
Chloride	●	●	●	●	●
Cholesterol	●				●
CPK	●	●			●
Creatinine	●	●		●	●
Direct Bilirubin	●		●		●
GGT	●	●	●		●
Globulin	●	●	●		●
Glucose	●			●	●
LDH	●		●		●
Magnesium					●
NA/K Ratio	●				●
Phosphorus	●			●	●
Potassium	●	●	●	●	●
Sodium	●	●	●	●	●
Total Bilirubin	●		●		●
Total Protein	●	●	●	●	●
Triglyceride	●				●

EQUINE

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
EQUINE HEALTH PROGRAM (EHP)			
L510	Equine Health Program Equine Chemistry, CBC, Fibrinogen, EIA (AGID), FEC (MST)	2.0 mL serum, 1.0 mL EDTA whole blood, 0.5 mL of citrated plasma, and 10 grams feces	1-4 days
L510E	Equine Health Program GVL Equine Chemistry, CBC, Fibrinogen, EIA (AGID) by GVL, FEC (MST)	2.0 mL serum, 1.0 mL EDTA whole blood, 0.5 mL of citrated plasma, and 10 grams feces	1-4 days
L510NE	Equine Health Program (no EIA) Equine Chemistry, CBC, Fibrinogen, FEC (MST)	1.0 mL serum, 1.0 mL EDTA whole blood, 0.5 mL of citrated plasma, and 10 grams feces	1-2 days
L511	Equine Health Program (ELISA) Equine Chemistry, CBC, Fibrinogen, EIA (ELISA), FEC (MST)	2.0 mL serum, 1.0 mL EDTA whole blood, 0.5 mL of citrated plasma, and 10 grams feces	1-3 days
L511E	Equine Health Program GVL (ELISA) Equine Chemistry, CBC, Fibrinogen, EIA (ELISA) by GVL, FEC (MST)	2.0 mL serum, 1.0 mL EDTA whole blood, 0.5 mL of citrated plasma, and 10 grams feces	1-3 days
L610	Equine Performance Horse Health Program Equine Performance Horse Chemistry, CBC, Fibrinogen, EIA (AGID), FEC (MST)	2.0 mL serum, 1.0 mL EDTA whole blood, 0.5 mL of citrated plasma, and 10 grams feces	1-4 days
L610E	Equine Performance Horse Health Program GVL Equine Performance Horse Chemistry, CBC, Fibrinogen, EIA (AGID) by GVL, FEC (MST)	2.0 mL serum, 1.0 mL EDTA whole blood, 0.5 mL of citrated plasma, and 10 grams feces	1-4 days
L610NE	Equine Performance Horse Health Program (no EIA) Equine Performance Horse Chemistry, CBC, Fibrinogen, FEC (MST)	1.0 mL serum, 1.0 mL EDTA whole blood, 0.5 mL citrated plasma, and 10 grams of feces	1-2 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
L615	Equine Senior Health Program Equine Chemistry, CBC, Fibrinogen, Endogenous ACTH, Insulin, EIA (AGID), FEC (MST) See L525 and T470E for additional drawing instructions.	3.0 mL of serum, 1.0 mL EDTA whole blood, 0.5 mL of citrated plasma, 10 grams feces, and 1.0 mL EDTA plasma	1-4 days
L615E	Equine Senior Health Program AGID GVL Equine Chemistry, CBC, Fibrinogen, Endogenous ACTH, Insulin, EIA (AGID) by GVL, FEC (MST) See L525 and T470E for additional drawing instructions.	3.0 mL of serum, 1.0 mL EDTA whole blood, 0.5 mL of citrated plasma, 10 grams feces, and 1.0 mL EDTA plasma	1-4 days
L615NE	Equine Senior Health Program (no EIA) Equine Chemistry, CBC, Fibrinogen, Endogenous ACTH, Insulin, FEC (MST) See L525 and T470E for additional drawing instructions.	2.0 mL of serum, 1.0 mL EDTA whole blood, 0.5 mL of citrated plasma, 10 grams feces, and 1.0 mL EDTA plasma	1-3 days
L615T	Equine Senior Health Program with TRH Equine Chemistry, CBC, Fibrinogen, TRH Stimulation Test for PPID, Insulin, EIA (AGID), FEC (MST) See L535 and T470E for additional drawing instructions.	3.0 mL of serum, 1.0 mL EDTA whole blood, 0.5 mL of citrated plasma, 10 grams feces, 1.0 mL Pre EDTA plasma and 1.0 mL Post EDTA plasma	1-4 days
L110	Equine Infectious Anemia (EIA) AGID Note: <ul style="list-style-type: none"> Equine Infectious Anemia (EIA) testing submission must be accompanied by the currently accepted Federal VS 10-11 FEB 2018 official form. All blocks of the official form must be completed. Signed by the submitting veterinarian with a current Federal Category II Accreditation number. Sample tube labeling must match official form by either horse's name or tube number EXACTLY. 	1.0 mL serum in red top or serum separator tube	2-3 days
L111	Equine Infectious Anemia (EIA) AGID by GVL Note: sample tube labeling must match official form by either horse's name or tube number EXACTLY.	1.0 mL serum in red top or serum separator tube	2-3 days
L120	Equine Infectious Anemia (EIA) ELISA Note: <ul style="list-style-type: none"> Equine Infectious Anemia (EIA) testing submission must be accompanied by the currently accepted Federal VS 10-11 FEB 2018 official form. All blocks of the official form must be completed. Signed by the submitting veterinarian with a current Federal Category II Accreditation number. Sample tube labeling must match official form by either horse's name or tube number EXACTLY. 	1.0 mL serum in red top or serum separator tube	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
L121	Equine Infectious Anemia (EIA) ELISA by GVL Note: sample tube labeling must match official form by either horse's name or tube number EXACTLY.	1.0 mL serum in red top or serum separator tube	1-3 days
CBC AND CHEMISTRY			
T332	CBC And Fibrinogen Interferences: marked hemolysis or lipemia.	1.0 mL EDTA whole blood and 0.5 mL citrated plasma (preferred for Fibrinogen)	24 hours Performed each shift
L365	Fibrinogen Quantitative	0.5 mL citrated plasma collected as whole blood in blue top tube at least 2/3rds full to the fill line	24 Hours
L070	Equine Inflammatory Profile CBC, Fibrinogen, Plasma Protein Interferences: marked hemolysis or lipemia.	1.0 mL EDTA whole blood and 0.5 mL citrated plasma	24 hours Performed each shift
T415	PT and aPTT Prothrombin Time. Activated Partial Thromboplastin Time Prothrombin Time and Activated Partial Thromboplastin Time. Prothrombin time measures the integrity of the extrinsic and common components of the coagulation cascade. Partial Thromboplastin Time (PTT) measures the integrity of the intrinsic and common components of the coagulation cascade. Interferences: marked hemolysis or lipemia. Partially full blue top tube may falsely increase sample's coagulation time. Note: clotting of the sample may preclude the analysis. The blue top tube needs to be filled to 2/3 or more of its capacity. Partially full blue top tubes may falsely increase the coagulation times.	0.5 mL citrated plasma or citrated whole blood in blue top tube at least 2/3rds full to the fill line	24 hours Performed each shift
L010	Equine Chemistry Total Protein, Albumin, Globulin, A/G Ratio, AST (SGOT), Alk Phos, GGT, T. Bilirubin, D. Bilirubin, BUN, Creatinine, BUN/Creatinine Ratio, Phosphorus, Glucose, Calcium, Sodium, Potassium, NA/K Ratio, Chloride, Cholesterol, Triglyceride, CPK, LDH Panel is recommended for all large animal patients except for Bovine (use code L1001). Sample Handling: serum separator samples should be spun prior to submission. Red top tubes should be spun and serum transferred to a non-additive tube and marked as SERUM. If the serum is not separated from the red blood cells, chemistry values can be affected.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
L601	Equine Performance Horse Chemistry Total Protein, Albumin, Globulin, A/G Ratio, AST (SGOT), Alk Phos, GGT, BUN, Creatinine, BUN/Creatinine Ratio, Sodium, Potassium, Chloride, CPK Sample Handling: serum separator samples should be spun prior to submission. Red top tubes should be spun and serum transferred to a non-additive tube and marked as SERUM. If the serum is not separated from the red blood cells, chemistry values can be affected.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
L050	Equine Chemistry, CBC	0.5 mL serum and 1.0 mL EDTA whole blood	24 hours Performed each shift
L040	Equine Chemistry, CBC, Fib Panel is recommended for all large animal patients except for Bovine (use code L080).	0.5 mL serum, 1.0 mL EDTA whole blood, and 0.5 mL citrated plasma	24 hours Performed each shift
L040R	Equine Chemistry, CBC, Fib (RECHECK) Note: resubmission must be within 30 days of original accession and previous accession number provided on the new request form.	0.5 mL serum, 1.0 mL EDTA whole blood, and 0.5 mL citrated plasma	24 hours Performed each shift
L640	Equine Performance Horse Panel Equine Performance Horse Chemistry, CBC, Fibrinogen	1.0 mL of serum, 1.0 mL of whole blood, and 0.5 mL of citrated plasma	24 hours Performed each shift
L565	Equine Chemistry, CBC, Fib, T4, fT4 (ED) Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results.	2.0 mL serum, 1.0 mL EDTA whole blood, and 0.5 mL citrated plasma	1-4 days
L035	Equine Chemistry, CBC, Fib, fT4 (ED)	1.5 mL serum, 1.0 mL EDTA whole blood, and 0.5 mL citrated plasma	1-4 days
L030	Equine Chemistry, CBC, Fib, T4 Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results.	0.5 mL serum, 1.0 mL EDTA whole blood, and 0.5 mL citrated plasma	24 hours Performed each shift
L290	Neonatal Foal Panel Equine Chemistry, CBC, Fibrinogen, Equine IgG Total	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.5 mL citrated plasma	1-3 days
L090	Equine IgG Total IgG is generally used to assess immunoglobulin concentration and passive transfer of immunity in foals or immunoglobulin status in adult horses. This test is species specific for horses. Note: IgG testing is species specific (see Large Animal Section for all non-equine species) Methodology: Immunoturbidimetric method	0.5 mL serum in red top or serum separator tube, 1.0 mL EDTA whole blood, or plasma from a EDTA whole blood sample *If submitting milk transfer to a non-additive tube and clearly marked as MILK	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
L225	Equine Hepatic Screen Total Protein, Albumin, Globulin, AST (SGOT), Alk Phos, GGT, T. Bilirubin, D. Bilirubin, I. Bilirubin, BUN, Sodium, Potassium, Chloride, LDH	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
L230	Equine Hepatic Panel Equine Hepatic Screen with Electrolytes, Bile Acids	1.0 mL serum in red top or serum separator tube	1-2 days
T225	Bile Acids Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	1-2 days
L1025	Equine Hepatic Profile Equine Hepatic Screen, Bile Acids, SDH Note: for SDH must separate serum and place in a non-additive tube marked SERUM. Sample must be kept cold or frozen and submitted within 24 hrs.	1.5 mL serum and 0.5 mL cold or frozen serum (for SDH testing) in red top or serum separator tube	1-2 days
T250	Sorbital Dehydrogenase SDH) Note: separate serum and place in a non-additive tube marked SERUM. Sample must be kept cold or frozen and submitted within 24 hrs.	0.5 mL serum in red top or serum separator tube (cold or frozen)	1-2 days
L240	Equine Renal Screen Total Protein, Albumin, BUN, Creatinine, Phosphorus, Glucose, Calcium, Sodium, Potassium, Chloride Sample Handling: serum separator samples should be spun prior to submission. Red top tubes should be spun and serum transferred to a non-additive tube and marked as SERUM. If the serum is not separated from the red blood cells, chemistry values can be affected.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
L275	Equine Muscle Enzyme Screen AST (SGOT), CPK Sample Handling: serum separator samples should be spun prior to submission. Red top tubes should be spun and serum transferred to a non-additive tube and marked as SERUM. If the serum is not separated from the red blood cells, chemistry values can be affected.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
L280	Equine Muscle Enzyme Screen, CBC, FIB AST (SGOT), CPK, CBC, Fibrinogen	1.0 mL of serum, 1.0 mL EDTA whole blood, and 0.5 mL of citrated plasma	24 hours Performed each shift
L190	Equine Rhabdomyolysis Screen AST (SGOT), BUN, Creatinine, Phosphorus, Calcium, Sodium, Potassium, Chloride, CPK, LDH Sample Handling: serum separator samples should be spun prior to submission. Red top tubes should be spun and serum transferred to a non-additive tube and marked as SERUM. If the serum is not separated from the red blood cells, chemistry values can be affected.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T240 Add-on Equivalent ADD130	Protein Electrophoresis (Serum) Total Protein, Albumin, Globulin, Alpha 1, Alpha 2, Beta 1, Gamma 1, fractional assessment with interpretation. An evaluation of the globulin fraction of the serum (alpha 1, alpha 2, beta, and gamma) to determine if the globulin fraction is monoclonal based on these components. Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	2-4 days
S1680	Equine Pre-Purchase Drug Screen Drug Screen (performed at TVMDL) includes the following: Corticosteroids (LC/MS): Betamethasone, dexamethasone, methylprednisolone, prednisolone, prednisone and triamcinolone acetonide NSAIDs: Acetaminophen, acetylsalicylic acid, carprofen, celcoxib, deracoxib, diclofenac, diflunisal, eltenac, ethacrynic acid, etodolac, fenbufen, fenoprofen, firocoxib, flufenamic acid, flunixin, flurbiprofen, ibuprofen, indomethacin, indoprofen, ketoprofen, ketorolac, meclofenamic acid, mefanamic acid, meloxicam, nabumetone, naproxen, oxyphenbutazone, phenylbutazone, piroxicam, salicylic acid, tenoxicam, tolfenamic acid and tolmetin Tranquilizers: trazodone, acepromazine, fluphenazine, fluoxetine, guanabenz, reserpine, romifidine, xylazine, and detomidine Muscle relaxant: Methocarbamol Note: urine is the preferred specimen but serum is acceptable.	10 mL urine in urine transport tube OR 5.0 mL serum in red top or serum separator (Mark accordingly)	7-10 days
URINE			
T760 Add-on Equivalent ADD220	Urinalysis-Complete The complete urinalysis includes a physical (color, appearance, USG), chemical (pH, protein, glucose, ketones, bilirubin, occult blood), and microscopic (WBC, RBC, casts crystals, bacteria, transitional and squamous epithelial cells, fat droplets) exam of the urine. Interferences: visible levels of hemolysis, drugs containing dyes, nitrofurantoin, or riboflavin.	6.0 mL urine in urine transport tube	24 hours Performed each shift
L340	Fractional Excretion of Electrolytes Serum: Calcium, Chloride, Creatinine, Phosphorus, Sodium, Potassium Urine: Calcium, Chloride, Creatinine, Phosphorus, Sodium, Potassium Note: serum and urine must be drawn and submitted at the same time.	1.0 mL serum in red top or serum separator tube AND 1.0 mL urine in urine transport tube (mark each sample as URINE or SERUM)	2-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
FECAL			
T826	FEC: Modified Stoll's Technique <i>Strongyle</i> sp. and <i>Parascaris</i> sp. ova (EPG) performed by Modified Stoll's Technique (MST) and qualitative ID of all other parasites Fecal Egg Count by Modified Stoll's Technique test provides quantitative eggs counts for <i>Strongyle</i> & <i>Parascaris</i> ova reported as eggs per gram (EPG). The lowest detection limit is 2 EPG. Other eggs or parasites seen will be reported. Note: fresh fecal samples should be kept refrigerated until submission. Specimens should be tested within 72 hours of collection.	10 grams feces in Antech provided fecal container	1-2 days
T828	FEC: McMasters Method (Large Animal) <i>Strongyle</i> sp. and <i>Parascaris</i> sp. ova (EPG) performed by McMaster's Method (MM) and qualitative ID of all other parasites Fecal Egg Count by McMaster's Method test provides an estimated egg count for <i>Strongyle</i> & <i>Parascaris</i> ova reported as Eggs per gram (EPG). The lowest detection limit is 100 EPG. Other eggs or parasites seen will be reported. Note: fresh fecal samples should be kept refrigerated until submission. Specimens should be tested within 72 hours of collection.	10 grams feces in Antech provided fecal container	1-2 days
L86181	<i>Clostridium Difficile</i> Toxins A/B *See L950 Equine PCR Gastrointestinal (Diarrhea) Panel also.	5 grams feces in Antech provided fecal container	1-2 days
T16007	<i>Clostridium Perfringens</i> Enterotoxin *See L950 Equine PCR Gastrointestinal (Diarrhea) Panel also.	5 grams feces in Antech provided fecal container Send on ice to keep cold	1-2 days
M160	Culture, Feces Culture specifically evaluates for <i>Salmonella</i> , <i>Shigella</i> and <i>Campylobacter</i> spp. Fecal PCR testing is more sensitive and tests for a broader array of potential pathogens.	5 grams feces in Antech provided fecal container. Submission in lavender top tube (with EDTA) is not acceptable.	3-4 days Preliminary report every 24 hours. Final report available in 72 hours.
M121	Culture, <i>Salmonella</i> Note: a negative <i>Salmonella</i> spp. culture result does not rule out intermittent or lower levels of fecal shedding of <i>Salmonella</i> spp., five-day serial submissions are recommended.	5 grams of feces in Antech provided fecal container	3-4 days Preliminary report available every 24 hours. Final report available in 72 hours.
L496	Acute Diarrhea Panel Fecal Culture, <i>Clostridium perfringens</i> enterotoxin, <i>Clostridium difficile</i> toxins A/B	10 grams feces in Antech provided fecal container	2-4 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
L492	Foal Diarrhea Panel FEC (MST), Fecal Culture, <i>Clostridium perfringens</i> enterotoxin, <i>Clostridium difficile</i> toxins A/B	10 grams feces in Antech provided fecal container	2-4 days
L420	Chronic Diarrhea Panel Equine Chemistry, CBC, Fibrinogen, FEC (MST), Fecal Culture, <i>Clostridium perfringens</i> enterotoxin, <i>Clostridium difficile</i> toxins A/B	1.0 mL serum in red top or serum separator tube, 1.0 mL EDTA whole blood, 0.5 mL citrated plasma, and 10 grams of feces	1-4 days
ENDOCRINOLOGY			
T495	T4 Interferences: marked hemolysis and moderate to marked lipemia. Lipemia can falsely decrease T4 results.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
T460	Free T4 By Equilibrium Dialysis Note: this test should not be performed as an add-on to samples older than 5 days.	0.5 mL serum in red top or serum separator tube	2-3 days
SA370	Thyroid Profile 2 Total T4, Free T4 By Equilibrium Dialysis Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results.	1.0 mL serum in red top or serum separator tube	1-3 days
L590	Equine TRH Stim (Thyroid Function) Baseline T3, Baseline T4, 2-hr Post T3, 4-hr Post T4	1.0 mL serum labeled Pre, 0.5 mL serum labeled Post T3, and 0.5 mL serum labeled Post T4. All samples should be in red top or serum separator tube.	3-5 days
L500	Pituitary Pars Intermedia Dysfunction Monitoring Panel Equine Chemistry, CBC, Fibrinogen, Endogenous ACTH (Equine), Insulin (Equine), T4 Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Note: see L525 and T470E for further sample handling instructions.	1.0 mL serum in red top or serum separator tube, 1.0 mL EDTA whole blood, 0.5 mL citrated plasma, and 1.0 mL EDTA plasma	1-2 days
L500TRH	Pituitary Pars Intermedia Dysfunction (PPID) with TRH Stim Equine Chemistry, CBC, Fibrinogen, TRH Stimulation for PPID, Insulin (Equine), T4 See L535 and T470E for further sample handling instructions.	1.0 mL serum in red top or serum separator tube, 1.0 mL EDTA whole blood, 0.5 mL citrated plasma, 1.0 mL Pre EDTA plasma and 1.0 mL Post EDTA plasma	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
L560	Pituitary Pars Intermedia Dysfunction (PPID) with ft4ED Equine Chemistry, CBC, Fibrinogen, Endogenous ACTH (Equine), Insulin (Equine), T4, Free T4 by Equilibrium Dialysis Sample Handling: see additional information under L525 and T470E.	2.0 mL serum in red top or serum separator tube, 1.0 mL EDTA whole blood, 0.5 mL citrated plasma, and 1.0 mL EDTA plasma	1-3 days
L525	Endogenous ACTH (Baseline) Equine Plasma ACTH levels show a seasonal variation for horses in the northern hemisphere. Seasonally elevated ACTH levels are usually observed between Mid-July to Mid-November. Horses with PPID may have exaggerated levels during this time. Interpretation range values will vary by season. Sample Handling: <ul style="list-style-type: none"> • 12 hour fast with one flake of hay left in the stall • Draw sample using EDTA whole blood tube • Spin sample with 8 hours of draw (lab will not run the test on unspun samples) • Transfer the plasma to new non-additive tube and mark as PLASMA • Send sample on ice pack 	1.0 mL EDTA plasma in non-additive transport tube (labeled as EDTA plasma)	1-2 days
L540	Endogenous ACTH, Insulin Panel Sample Handling: see additional information under L525 and T470E.	1.0 mL serum (for insulin test) AND 1.0 mL EDTA plasma (for Endogenous ACTH test)	1-2 days
L575	Endogenous ACTH, Insulin/Glucose, Leptin (Equine) See additional handling instruction under L525 and T470E.	2.0 mL serum and 1.0 mL EDTA plasma	1-10 days
S14402	Leptin Equine Note: transport refrigerated (on ice packs)	2.0 mL serum in red top or non-additive tube	7-9 Days
L535	TRH Stimulation Test for PPID Endogeneous ACTH (Pre and 10 minute Post) Tests for Cushings Disease (PPID). If Endogenous ACTH results (L525) are normal but disease is strongly suspected, this panel is recommended. Interferences: serum or whole blood submissions will preclude analysis and will NOT be accepted. Sample Handling: see additional information under L525.	1.0 mL EDTA plasma in plain, non-additive tube labeled as Pre EDTA plasma, and 1.0 mL EDTA plasma in plain, non-additive tube labeled as Post EDTA plasma. All samples should be kept cold.	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
T470E	<p>Equine Insulin and Glucose</p> <p>In horses, a high baseline insulin concentration is indicative of Insulin Dysregulation (ID). However, many horses will have a normal baseline insulin with an abnormal response to sugar ingestion. If baseline results are normal but disease is strongly suspected, an oral sugar test (OST) is recommended see code L545 or L550.</p> <p>Sample Handling:</p> <ul style="list-style-type: none">• 12 hour fast with one flake of hay left in the stall• Draw sample in a red top or serum separator tube• Spin sample within 8 hours of draw• Transfer serum into new non-additive tube and mark as SERUM	0.5 mL serum spun and separated in red top or spun serum separator tube	1-2 days
L580	<p>TRH Stimulation with Insulin Glucose</p> <p>See L535 and T470E for additional drawing instructions.</p>	1.0 mL spun and separated serum in red top or spun serum separator tube, 1.0 mL Pre EDTA plasma and 1.0 mL Post EDTA plasma	1-2 days
L545	<p>Oral Sugar Test (2 samples)</p> <p>Insulin (Equine)/Glucose (60 and 90 minute Post)</p> <p>Oral sugar test assesses the insulin response to ingested sugars for a diagnosis of Insulin Dysregulation (ID).</p> <p>See T470E for additional drawing instructions.</p>	1.0 mL serum in red top or serum separator tube (labeled 60 min Post) AND 1.0 mL serum in red top or serum separator tube (labeled 90 min Post)	1-2 days
L550	<p>Oral Sugar Test (3 samples)</p> <p>Insulin (Equine)/Glucose (Pre, 60 and 90 minute Post)</p> <p>Oral sugar test assesses the insulin response to ingested sugars for a diagnosis of Insulin Dysregulation (ID).</p> <p>See T470E for additional drawing instructions.</p>	1.0 mL serum in red top or serum separator tube (labeled Baseline), 1.0 mL serum in red top or serum separator tube (labeled Serum 60 min Post), AND 1.0 mL serum in red top or serum separator tube (labeled Serum 90 min Post)	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
NEUROLOGY			
S14388	EPM SAG 2, 4/3 ELISA The ELISA format contains surface antigen proteins (SAG 2,4/3) of the <i>S. neurona</i> parasite. A negative serum result is a good rule out for EPM caused by <i>S. neurona</i> . Note: if submitting BOTH serum and CSF, should be submitted separately as two different requests in two different bags.	1.0 mL serum in red top or serum separator OR 1.0 mL CSF in red top tube	3-5 days
S14390	EPM SAG 2, 4/3 ELISA Serum/CSF Ratio The specific antibody index is performed on paired serum and CSF samples in conjunction with the <i>S. neurona</i> SAG 2,4/3 ELISA. This test can help discern an intrathecal IgG response from contamination by utilizing a ratio of serum and CSF albumin and IgG levels. Note: both serum and CSF MUST be submitted together in the same bag.	1.0 mL serum in red top or serum separator AND 1.0 mL CSF in red top tube	3-5 days
S14392	EPM SAG 2, 4/3 ELISA, N hughesi ELISA The <i>Neospora hughesi</i> ELISA utilizes a specific protein, SAG 1, found on the <i>N. hughesi</i> parasite. Serum and/or CSF samples can be tested for an IgG response to the recombinant NhSAG 1 protein. The ELISA generates an endpoint titer. Note: if submitting BOTH serum and CSF, should be submitted separately as two different requests in two different bags.	1.0 mL serum OR 1.0 mL CSF in red top tube	4-8 days
S16275	Encephalitis EEE, WEE, VEE (PRNT) and EEE (IgM Capture ELISA)	1.0 mL serum in red top or serum separator tube	10-14 days
S17500	Equine Viral Encephalitis Panel Encephalitis, West Nile Virus Ab (PRNT/IgM), EHV-I Ab	3.0 mL serum in red top or serum separator tube	10-14 days
S14477	Equine Comprehensive Neurological Panel Equine Viral Encephalitis Panel, EPM SAG 2,4/3 ELISA, West Nile titer (PRNT/IgM), EHV-I PCR	3.0 mL serum in red top or serum separator tube. For EHV-1 PCR component: 5.0 - 7.0 mL EDTA whole blood in lavender top tube, 5.0 mL nasal wash in red top tube, and/or nasal swab in non-additive tube with no media.	5-14 days
S85448	West Nile Titer (PRNT & IgM ELISA)	1.0 mL serum in red top or serum separator tube (mammals) OR 0.1 mL heparinized plasma (avians)	7-14 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
PCR			
L950	<p>Equine Gastrointestinal PCR Panel</p> <p><i>Clostridium difficile</i> A & B, <i>Clostridium perfringens</i>, <i>C. perfringens</i> toxin A, <i>C. perfringens</i> toxin B, <i>C. perfringens</i> toxin E, <i>C. perfringens</i> toxin NetF, Equine Coronavirus, <i>N. risticii</i> (Potomac Horse Fever), Equine Rotavirus, and <i>Salmonella</i> spp.</p> <p>Note: A single negative <i>Salmonella</i> spp. PCR result does NOT rule out intermittent or lower levels of fecal shedding of <i>Salmonella</i> spp. Three (3) serial submissions are recommended (use code S14416 for sample 2 and/or sample 3).</p> <p>Interferences: wooden swabs can interfere with testing and plastic swabs are preferred.</p>	5 grams of feces in Antech provided fecal container	4-5 days
L955	<p>Equine PCR Reproductive Panel</p> <p>EVA, EHV-1, <i>Strep equi</i> subsp. <i>zooepidemicus</i>, <i>Leptospira</i> spp., <i>Klebsiella pneumoniae</i>, <i>E. coli</i></p> <p>Interferences: wooden swabs can interfere with testing and plastic swabs are preferred.</p>	Uterine lavage fluid, semen, endometrial or cervical swab in red top tube (without media), OR placenta	4-5 days
L960	<p>Equine PCR Respiratory Panel</p> <p><i>Streptococcus equi</i>, <i>Streptococcus zooepidemicus</i>, <i>Rhodococcus equi</i>, Equine Herpes Virus 1 (EHV1), Equine Herpes Virus 4 (EHV4), Equine Influenza Virus (EIV), Equine Rhinitis Virus A and B (ERAV & ERBV)</p> <p>Interferences: wooden swabs can interfere with testing and plastic swabs are preferred.</p>	Trans tracheal wash, bronchoalveolar lavage, guttural pouch wash, nasal swab (in non-additive tube without media) or respiratory tract tissue	4-5 days
L965	<p>Equine PCR FUO Panel</p> <p>Swab: <i>Streptococcus equi</i>, Equine Herpes Virus 1, Equine Herpes Virus 4, Equine Influenza Virus, Equine Rhinitis Virus A & B</p> <p>EDTA whole blood: Equine Herpes Virus 1, <i>Anaplasma phagocytophilum</i>, <i>Neorickettsia risticii</i> (PHF)</p> <p>Fecal: <i>Neorickettsia risticii</i> (PHF), Equine Coronavirus (ECoV)</p> <p>Interferences: wooden swabs can interfere with testing and plastic swabs are preferred.</p>	1.0 mL EDTA whole blood, nasal swab in non-additive tube (without media), and 5 grams of feces in Antech provided fecal container	4-5 days
L970	<p>Equine PCR FUO Mini</p> <p>EDTA whole blood: <i>Anaplasma phagocytophilum</i>, <i>Neorickettsia risticii</i> (PHF)</p> <p>Fecal: <i>Neorickettsia risticii</i> (PHF), Equine Coronavirus (ECoV)</p>	1.0 mL EDTA whole blood and 5 grams of feces in Antech provided fecal container	4-5 days
S14421	<p><i>Anaplasma phagocytophilum</i> PCR</p>	1.0 mL EDTA whole blood	3-5 days
S14414	<p>Equine Coronavirus (ECoV)</p>	5.0 grams of feces in Antech provided fecal container	4-5 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S14394	Equine Herpes Virus 1 (EHV1) PCR Viremia can be detected in the blood as soon as four days post-infection. Viral shedding in nasal secretions can be detected as soon as six days post-infection and for a longer time than in blood. Sample Handling: submission of BOTH nasal swab AND EDTA whole blood is recommended.	5.0 mL EDTA whole blood, 5.0 mL nasal wash in red top tube, or nasal swab in non-additive tube without media	3-5 days
L974	Equine <i>Leptospira</i> PCR Blood	0.5 mL EDTA whole blood	3-5 days
L976	Equine <i>Leptospira</i> PCR Urine	2.0 mL urine in urine transport tube	3-5 days
L978	Equine <i>Leptospira</i> PCR Blood/Urine	0.5 mL EDTA whole blood and 2.0 mL of urine in urine transport tube	3-5 days
L973	Equine Lepto PCR Aq/Vitreous Humor	1.5 mL aqueous or vitreous humor fluid in lavender top tube	2-4 days
S14479	<i>Neorickettsia risticii</i> PCR (PHF) PCR testing will be done on both EDTA whole blood and fecal samples for no additional charge when both samples are submitted.	1.0 mL EDTA whole blood and/or 5 grams of feces in Antech provided fecal container	4-5 days
S14396	<i>Rhodococcus equi</i> PCR Identifies strains of <i>R. equi</i> carrying virulence plasmid gene vapA. Transtracheal wash sample is the recommended specimen. Note: a fecal sample can be tested but negative results will be unreliable. Interferences: wooden swabs can interfere with testing and plastic swabs are preferred.	Transtracheal wash, bronchoalveolar lavage, and/or nasal swab in non-additive tube without media.	4-5 days
T983	Equine Ringworm PCR Panel <i>Microsporum</i> spp., <i>M. equinum</i> / <i>canis</i> , <i>M. gypseum</i> (<i>Arthroderma gypseum</i> , <i>A. fulvum</i> , <i>A. incurvatum</i>), <i>Trichophyton</i> spp., <i>T. equinum</i> / <i>mentagrophytes</i> (<i>Arthroderma benhamiae</i> , <i>A. vanbreuseghemii</i>)	Minimum of 10 plucked hair with roots, skin scraping, or tooth brush sample in sterile, dry container free of liquids or preservatives	1-5 days
T987	Equine Ringworm PCR with Dermatophyte Culture <i>Microsporum</i> spp., <i>M. equinum</i> / <i>canis</i> , <i>M. gypseum</i> (<i>Arthroderma gypseum</i> , <i>A. fulvum</i> , <i>A. incurvatum</i>), <i>Trichophyton</i> spp., <i>T. equinum</i> / <i>mentagrophytes</i> (<i>Arthroderma benhamiae</i> , <i>A. vanbreuseghemii</i>), and Dermatophyte Culture	Minimum of 15 plucked hair with roots, skin scraping, or tooth brush sample in sterile, dry container free of liquids or preservatives	1-21 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S14416	Salmonella Spp. PCR Detects DNA from pathogenic species of <i>Salmonella</i> . <i>Salmonella</i> PCR testing can also be used as a bio-surveillance tool to identify asymptomatic/shedding horses and monitor in-patients in a hospital setting. Note: a single negative <i>Salmonella</i> spp. PCR result does not rule out intermittent or lower levels of fecal shedding of <i>Salmonella</i> spp. Three (3) serial submissions are recommended.	5 grams of feces in Antech provided fecal container	4-5 days
S86308	Streptococcus equi equi PCR Detects <i>S. equi</i> bacterial DNA and is used primarily to identify asymptomatic carriers. Nasal washes, nasal swabs, or guttural pouch washes are appropriate samples for the <i>S. equi</i> PCR. Back-up culture is also available. Interferences: wooden swabs can interfere with testing and plastic swabs are preferred.	5.0 mL nasal wash, OR 5.0 mL guttural pouch wash, OR 1 pharyngeal swab in non-additive tube without media	4-5 days
MARE REPRODUCTION			
S16295	Estradiol Recommended for testing if greater than 120 days post-breeding in conjunction with progesterone (L140) using code L470.	1.0 mL serum in red top or serum separator tube	7-10 days
S16300	Estrone Sulfate Recommended for testing if greater than 100 days post breeding.	1.0 mL serum in red top or serum separator tube	7-10 days
S16635	Pregnant Mare Serum Gonadotropin (PMSG) Recommended for testing at day 45 to 120 days post-breeding in conjunction with progesterone (L140) using test code L460.	2.0 mL serum in red top or serum separator tube	7-9 days
L140	Progesterone Interferences: serum separator gel may interfere with test.	0.5 mL serum in red top tube. Submission of serum separator tube is not recommended.	1-2 days
L460	Progesterone/PMSG	2.5 mL serum in red top tube. Submission of serum separator tube is not recommended.	1-9 days
L470	Progesterone/Estradiol Panel (Equine)	1.5 mL serum in red top tube. Submission of serum separator tube is not recommended.	1-10 days
L480	Progesterone/PMSG/Estradiol	3.0 mL serum in red top tube. Submission of serum separator tube is not recommended.	1-10 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
GRANULOSA CELL TUMOR TESTING			
S85857	Equine Granulosa Cell Tumor Progesterone, Testosterone, Inhibin	3.0 mL serum in red top tube. Submission of serum separator tube is not recommended.	2-21 days
S4131	Inhibin Inhibin is a heat-sensitive protein. Interferences: marked hemolysis or lipemia.	1.0 mL serum in red top tube or serum separator tube	2-3 weeks
S14320	Anti-Mullerian Hormone (Equine)	1.0 mL serum in red top or serum separator tube	5-11 days
STALLION/GELDING REPRODUCTION			
S16760	Testosterone Single testosterone level	0.5 mL serum in red top or serum separator tube	3-5 days
SEMEN	Semen Analysis Sperm Count (including Total Volume and pH), with Cytologic Evaluation	2.0 mL semen (fresh or extended) in lavender top tube	1-2 days
CRYPTORCHID TESTING			
S16300	Estrone Sulfate >3 yrs of age	1.0 mL serum in red top or serum separator tube	7-10 days
S14320	Anti-Mullerian Hormone (Equine) <3 yrs of age	1.0 mL serum in red top or serum separator tube	5-11 days
S85530	Testosterone Panel Pre, 20 min post, 1 hr post, 2 hr post HCG Note: spin serum and separate in non-additive tube marked SERUM. Equine HCG Stimulation Testing Protocol: 1. Take baseline serum sample for testosterone. 2. Administer 6000 IU (total dose) of HCG IV. 3. Take post-HCG testosterone at 20 minutes, 1 hour and 2 hour post administration. 4. If testicular tissue is present, the post values should be greater than 0.4 nmol/L and/or double the baseline level.	0.5 mL serum for resting and each post sample (baseline sample labeled Pre and each post sample labeled accordingly)	3-5 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
MICROBIOLOGY			
M020	<p>Culture, Aerobic</p> <p>This test is used when an aerobic bacterial infection is suspected in a tissue or fluid. In the case of a positive culture of relevant bacteria, sensitivities are performed and reported.</p> <p>Interferences:</p> <ol style="list-style-type: none">1. Patient should be off antibiotics for at least 7-10 days.2. Fluid in lavender top tube is unacceptable (EDTA inhibits growth). <p>Note: tissue should be submitted in an RTT with a few drops of saline (to keep moist). Samples should be refrigerated prior to transportation to the laboratory. A separate sample should be submitted for each type of culture needed.</p>	<p>Culturette from fluid (body cavity fluids), TTW, BAL, wound, lesion, or skin. Tissue sample in saline.</p> <p>Culturette, red top or white top tube with fluid, or other sterile container. Samples collected in EDTA tube are not acceptable.</p>	<p>3-4 days</p> <p>Preliminary report available every 24 hours. Final culture result available in 72 hours. If a fastidious organism is observed, the listed turnaround time may be extended.</p>
M030	<p>Culture, Anaerobic</p> <p>To be used when an anaerobic bacterial infection is suspected in tissue/fluid.</p> <p>Interferences: sample with exposure to air or if dried will preclude accurate testing.</p> <p>Note: fluids should be submitted in an RTT. Tissue should be submitted in an air-tight, sterile container (e.g., urine cup). No saline is necessary. Sensitivities are not performed, but drugs of choice for anaerobes are provided as a guide.</p>	<p>Anaerobic culturette of fluid (body cavity fluids), TTW or BAL, wound, lesion. Tissue sample (at least 2 cm x 2 cm) in sterile, air tight container.</p> <p>Culturette, red top or white top tube with fluid, or other sterile container. Samples collected in EDTA tube are not acceptable.</p>	<p>3-5 days</p> <p>Preliminary report issued 3rd day. Final culture result issued 4th day.</p>
M040	<p>Culture Aerobic and Anaerobic</p> <p>This test is used when a bacterial infection is suspected in a tissue or fluid, but it is uncertain whether an aerobic or an anaerobic organism is the cause.</p> <p>Note: fluids should be submitted in a red top tube. Tissue should be submitted in an air-tight container for anaerobic culture. Samples should be refrigerated prior to submission. Sensitivities are not performed, but drugs of choice for anaerobes are provided as a guide.</p>	<p>Culturette (aerobic and anaerobic) from fluid (body cavity fluids), TTW, BAL, wound, lesion, or skin. Tissue sample in saline.</p>	<p>3-5 days</p>

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
M060	Culture, Blood Aerobic Only This test is utilized when bacteremia with an aerobic organism is suspected. Submit a single sample in a BACTEC bottle. Interferences: antibiotics. Note: submit in BACTEC aerobic culture bottle. Clip fur and scrub venipuncture site for aseptic collection. Aim to collect 1.0 mL for cats and small dogs, 2-3 mL for larger dogs to be inoculated into 20 mL bottle. Don't unscrew caps on bottles. Remove the protective top and wipe visible parts of the rubber stopper with 70% ethanol. Allow stopper to dry or wipe with sterile gauze. Replace the drawing needle with a sterile needle before puncturing the rubber stopper, fill until vacuum stops, then gently invert the bottle to mix. Anticoagulants in the media will prevent blood from clotting.	Whole blood collected in BD Bactec Blood Culture Bottle (BCB)	3-5 days
M220	Culture, Enhanced Fluid Aerobic Note: if submitting multiple joints, must be submitted separately in different bags and forms.	Fluid in BD Bactec Blood Culture Bottle (BCB)	3-4 days
M080	Culture, Fungal This test is for fungal culture and identification. Not for suspected ringworm infections.	Dry hair, nails, skin scraping, body fluid, or lesion material collected in a sterile red top or other sterile container without additive, or on a culturette. DTM bottle, culturette, container	21 days Preliminary report every 7 days. Final report available at 21 days (3 weeks), or upon identification.
M130 Add-on Equivalent ADD210	Culture, Urine All urine cultures deemed positive by <i>FIRST</i> tract, a highly accurate and rapid urine culture technique, will then undergo traditional culture for identification and susceptibility. Interferences: patient should be off antibiotics for seven to ten days prior to urine culture.	0.5 mL cystocentesis, clean catch or catheterized urine in sterile red top tube or urine transport tube.	1-3 days Preliminary reports available every 24 hours. If a fastidious organism is observed, the listed turnaround time may be extended.

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CYTOLOGY			
CYTO	<p>Cytology</p> <p>Includes preparation of submitted sample and microscopic interpretation by clinical pathologist</p> <p>This test includes a microscopic evaluation of cells. The report includes cytologic interpretation, diagnosis, and comments regarding etiology and biological behavior where applicable.</p> <p>Note: single source.</p>	<p>Fluid in red top or lavender top tube, or air-dried, unstained slides with clinical history</p>	<p>1-3 days</p>
FLUA	<p>Fluid Analysis</p> <p>Includes preparation of submitted sample, Cell Count (WBC & RBC), Specific Gravity, Protein measurement, and microscopic interpretation by clinical pathologist</p> <p>This is the appropriate test for the evaluation of fluids from a cavity or tissue space. The fluid will be characterized based on protein, specific gravity, white and red cells, and a pathologist will evaluate the fluid and provide a microscopic description and, when possible, diagnosis.</p> <p>Note: a comprehensive history should be included with the fluid analysis submission.</p>	<p>1.0 mL body fluid in lavender top or red top with 2 unstained smears prepared from fluid</p>	<p>1-3 days</p>
CSF	<p>Fluid Analysis CSF</p> <p>Includes preparation of submitted sample, CSF Cell Count (WBC & RBC), Specific Gravity, Protein & Glucose measurement, and microscopic interpretation by clinical pathologist</p>	<p>1.0 mL cerebrospinal fluid in lavender or red top tube</p>	<p>1-3 days</p>
SYFLUA	<p>Synovial Fluid Analysis</p> <p>Synovial Fluid Analysis (Color, Clarity, Specific Gravity, WBC, RBC), Total Protein, and Cytologic Evaluation</p> <p>Note: single source.</p>	<p>1.0 mL synovial fluid in lavender top tube</p>	<p>1-3 days</p>

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
BAL	Airway Wash Includes preparation of submitted sample and microscopic interpretation by clinical pathologist	1.0 mL airway wash fluid in red top or lavender top tube, or 2 unstained air-dried smears prepared from fluid	2-4 days
TTW	Airway Wash Includes preparation of submitted sample and microscopic interpretation by clinical pathologist	1.0 mL airway wash fluid in red top or lavender top tube, or 2 unstained air-dried smears prepared from fluid	2-4 days
FBX	Histopathology Preparation of submitted sample and microscopic interpretation by anatomic pathologist A boarded pathologist will evaluate the tissue. A full written biopsy report will be provided, including source, history, submitted tissue type, microscopic findings, pathologist's comments, and microscopic description. Where applicable, it will also include margin evaluation, grading, interpretation of special stains, and recommendations for further testing if needed. Note: firmly and evenly tighten the lid of the formalin jar and check for any leaks prior to placing the sample in a ziplock bag with the test requisition form. A comprehensive history should be included with each biopsy submission.	Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin	3-5 business days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
INDIVIDUAL TESTS			
S14466	<p>Equine 11 Metal & Mineral Panel</p> <p>Serum: Cobalt, Copper, Iron, Magnesium, Molybdenum, Selenium, Zinc</p> <p>EDTA Whole Blood: Arsenic, Cadmium, Lead, Thallium</p> <p>Serum selenium provides current-day levels, while whole blood selenium indicates selenium status over a longer period of time.</p>	<p>One or more of the following:</p> <ol style="list-style-type: none">1.0 mL serum in non-additive tube AND EDTA whole blood in lavender top tube5 grams fresh liver sample in sealed, leak-proof container50 mg fresh liver biopsy (3 Tru-Cut samples, no fluid added) in sealed, leak-proof container500 grams of feed (representative sample) in sealed container	5-10 days
S16285	<p>Equine Viral Arteritis</p>	<p>1.0 mL serum in red top or serum separator tube</p>	7-10 days
S16510	<p>Leptospirosis</p> <p><i>L. pomona</i>, <i>L. icterohemorrhagiae</i>, <i>L. canicola</i>, <i>L. grippotyphossa</i>, <i>L. hardjo</i>, <i>L. autumnalis</i>, <i>L. bratislava</i></p> <p>The results of this test include a semi-quantitative titer for serovars.</p> <p>Interferences: recent vaccination (within one month) may interfere with testing.</p> <p>Note: test is not species specific.</p>	<p>1.0 mL serum in red top tube or serum separator tube</p>	3-5 days
T672	<p>Lyme Equine Multiplex</p>	<p>0.5 mL serum in red top or serum separator tube</p>	5-10 days
S16848	<p>Piroplasmosis, Equine</p> <p><i>Theileria (Babesia) equi</i>, <i>Babesia caballi</i>.</p>	<p>2.0 mL serum in red top or serum separator tube</p>	10-12 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S16740	Streptococcus Equi ELISA The SeM ELISA test assesses the IgG titer against the <i>S. equi</i> specific M protein. It is helpful for making vaccination decisions on horses with existing titers, identifying horses at risk of complications due to elevated titers or for horses with aberrant abscesses (bastard strangles). The SeM ELISA should not be used to determine the infection status of horses with clinical signs of strangles.	5.0 mL serum in red top or serum separator tube	7-10 days
S16850	Vitamin E Note: cover sample to prevent exposure to light.	1.5 mL serum in red top or serum separator tube (spun)	7-10 days
S17505	Vitamin E, Selenium	1.5 mL serum in red top or serum separator tube (spun)	7-10 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE FOR MORE DETAILS ON TURNAROUND TIMES

LARGE ANIMAL

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CBC AND CHEMISTRY			
T332	CBC And Fibrinogen Interferences: marked hemolysis or lipemia.	1.0 mL EDTA whole blood and 0.5 mL citrated plasma (preferred for Fibrinogen)	24 hours Performed each shift
L1001	Ruminant Chemistry with Magnesium Total Protein, Albumin, Globulin, A/G Ratio, AST (SGOT), Alk Phos, GGT, T. Bilirubin, D. Bilirubin, BUN, Creatinine, BUN/Creat Ratio, Phosphorus, Glucose, Calcium, Magnesium, Sodium, Potassium, NA/K Ratio, Chloride, Cholesterol, Triglyceride, CPK, LDH Sample Handling: serum separator samples should be spun prior to submission. Red top tubes should be spun and serum transferred to a non-additive tube and marked as SERUM. If the serum is not separated from the red blood cells, chemistry values can be affected.	0.5 mL serum in red top or serum separator tube	24 hours Performed each shift
L080	Ruminant Chem with Mag, CBC	0.5 mL serum, 1.0 mL EDTA whole blood, and 0.5 mL citrated plasma	24 hours Performed each shift
INDIVIDUAL TESTS			
S16035	Anaplasma CF, Bovine Interferences: severe hemolysis or lipemia will preclude testing.	2.0 mL serum in red top or serum separator tube	5-10 days
S16124	Babesia Bovis	1.0 mL serum in red top or serum separator tube	7-10 days
T785	Baermann Baermann test is to detect lungworm larvae in fecal samples. Sample Handling: specimen should be tested within 24 hours of collection.	10 grams feces in Antech provided fecal container	2-4 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
S14456	Bovine Herpes Virus Note: transport refrigerated (on ice packs).	1.0 mL serum in red top or serum separator tube	5-7 days
S16844	Bovine Viral Diarrhea (BVD) - ELISA	1.0 mL of serum from precolostral newborn calves or calves older than three months are suitable for testing. Ear notches from animals of all ages may also be tested	7-10 days
S86551	Bovine Viral Diarrhea (BVD) PCR	2.0 mL EDTA whole blood in lavender top tube, fresh tissue, semen, or milk	7-10 days
S16145	Caprine Arthritis Encephalitis (CAE)	1.0 mL serum in red top or serum separator tube	7-10 days
S16425	IgG Bovine Test is species specific for bovine. Methodology is the RID Method.	1.0 mL serum in red top or serum separator tube	7-10 days
S16430	IgG Llama Test is species specific for llama. Methodology is the RID Method.	1.0 mL serum in red top or serum separator tube	2-3 days
S16302	Johne's Disease Antibody Whole herd screening.	1.0 mL serum in red top or serum separator tube	5-7 days
T805	Ova and Parasite Fecal ova and parasite detection. Samples are appropriately mixed with zinc sulfate solution, centrifuged, followed by flotation and slide evaluation. Note: sample should be evaluated within 24 hours of collection. If a worm has been identified in the sample, separate the worm and place it in a container labeled Worm in black, additionally mark Worm on the TRF being submitted with the sample.	5 grams feces in Antech provided fecal container	1-2 days

TEST INDEX

TEST NAME	CODE	PAGE
Accuplex®	AC100	20, 99
Accuplex® with Microfilaria Knotts	AC617	99
Accuplex®, KeyScreen® GI Parasite PCR	KAC100	20
Acetylcholine Receptor Antibody	S16005	135
Acid Fast Stain - Micro	M010	109
Acute Diarrhea Panel	L496	168
Add-on CBC	ADD03	73, 75
Add-on Fructosamine	ADD260	91, 93
Add-on Fructosamine	FRUCTADD	79
Adrenal Androgen Panel For Ferrets	S17116	151
Adult Chem with Lytes, SDMA, CBC, Accuplex®, KeyScreen®	KA535	20
Adult Chem with Lytes, SDMA, CBC, Accuplex®, T4, KeyScreen®	KA583	21
Adult Chem with Lytes, SDMA, CBC, Feline Heartworm, FIV, FeLV, KeyScreen®	KS519	21
Adult Chem with Lytes, SDMA, CBC, FIV, FeLV, KeyScreen®	KS575	21
Adult Chem with Lytes, SDMA, CBC, Heartworm, KeyScreen®	KS534	21
Adult Chem with Lytes, SDMA, CBC, KeyScreen®	KS535	22
Adult Chem with Lytes, SDMA, CBC, O&P, Accuplex®	AC535S	20
Adult Chem with Lytes, SDMA, CBC, T4, Feline Heartworm, FIV, FeLV, KeyScreen®	KS591	22
Adult Chem with Lytes, SDMA, CBC, T4, FIV, FeLV, KeyScreen®	KS597	22
Adult Chem with Lytes, SDMA, CBC, T4, KeyScreen®	KS583	22
Adult Chem with Lytes, SDMA, CBC, T4, UA, Feline Heartworm, FIV, FeLV, KeyScreen®	KS589	24
Adult Chem with Lytes, SDMA, CBC, T4, UA, FIV, FeLV, KeyScreen®	KS619	22
Adult Chem with Lytes, SDMA, CBC, T4, UA, KeyScreen®	KS673	21
Adult Chem with Lytes, SDMA, CBC, UA, Accuplex®, KeyScreen®	KA053	23
Adult Chem with Lytes, SDMA, CBC, UA, Accuplex®, T4, KeyScreen®	KA673	23
Adult Chem with Lytes, SDMA, CBC, UA, Feline Heartworm, FIV, FeLV, KeyScreen®	KS623	24
Adult Chem with Lytes, SDMA, CBC, UA, FIV, FeLV, KeyScreen®	KS634	23
Adult Chem with Lytes, SDMA, CBC, UA, Heartworm, KeyScreen®	KS592	23
Adult Chem with Lytes, SDMA, CBC, UA, Heartworm, T4, KeyScreen®	KS587	24
Adult Chem with Lytes, SDMA, CBC, UA, KeyScreen®	KS053	24
Adult Chem with SDMA, Accuplex®, KeyScreen®	KA600	25
Adult Chem with SDMA, CBC, Accuplex®	AC600	25

TEST NAME	CODE	PAGE
Adult Chem with SDMA, CBC, Accuplex®, T4, KeyScreen®	KA670	25
Adult Chem with SDMA, CBC, Feline Heartworm, FIV, FeLV, KeyScreen®	KS622	25
Adult Chem with SDMA, CBC, Feline Heartworm, UA, KeyScreen®	KS630	41
Adult Chem with SDMA, CBC, FIV, FeLV, KeyScreen®	KS675	25
Adult Chem with SDMA, CBC, Heartworm, KeyScreen®	KS605	26
Adult Chem with SDMA, CBC, Heartworm, T4, KeyScreen®	KS685	26
Adult Chem with SDMA, CBC, KeyScreen®	KS600	26
Adult Chem with SDMA, CBC, O&P, Accuplex®	AC615	26
Adult Chem with SDMA, CBC, O&P, <i>Giardia</i> , Accuplex®	AC655	26
Adult Chem with SDMA, CBC, T4, Feline Heartworm, FIV, FeLV, KeyScreen®	KS590	27
Adult Chem with SDMA, CBC, T4, FIV, FeLV, KeyScreen®	KS697	27
Adult Chem with SDMA, CBC, T4, UA, Feline Heartworm, FIV, FeLV, KeyScreen®	KS588	27
Adult Chem with SDMA, CBC, T4, UA, FIV, FeLV, KeyScreen®	KS621	27
Adult Chem with SDMA, CBC, T4, UA, KeyScreen®	KS672	28
Adult Chem with SDMA, CBC, UA, Accuplex®	AC607	28
Adult Chem with SDMA, CBC, UA, Accuplex®, KeyScreen®	KA607	28
Adult Chem with SDMA, CBC, UA, Accuplex®, T4, KeyScreen®	KA672	28
Adult Chem with SDMA, CBC, UA, Feline Heartworm, FIV, FeLV, KeyScreen®	KS624	29
Adult Chem with SDMA, CBC, UA, FIV, FeLV, KeyScreen®	KS631	29
Adult Chem with SDMA, CBC, UA, Heartworm, KeyScreen®	KS625	29
Adult Chem with SDMA, CBC, UA, Heartworm, T4, KeyScreen®	KS687	29
Adult Chem with SDMA, CBC, UA, KeyScreen®	KS607	30
Adult Chem with SDMA, CBC, UA, UPC with Accuplex®	AC652	30
Adult Chem with SDMA, T4, KeyScreen®	KS670	30
Adult Chem with SDMA, Vaccine Titers, Fecal Combo, Accuplex®	AC664	30
Adult Chemistry with SDMA, CBC, FeLV, FIV	SA675	41
Adult Feline Wellness with SDMA, FeLV, FIV	SA622	41
Adult Wellness Chemistry with SDMA	SA665	79
Adult Wellness Panel with SDMA, Heartworm, O&P, UA	SA686	41
Adult Wellness with SDMA	SA600	42
Adult Wellness with SDMA, CBC, T4, Accuplex®	AC670	31
Adult Wellness with SDMA, <i>Ehrlichia</i> , Lyme	SA620	42
Adult Wellness with SDMA, FeLV, FIV, O&P, <i>Giardia</i>	SA684	42
Adult Wellness with SDMA, Heartworm, UA, O&P, <i>Giardia</i>	SA651	42
Adult Wellness with SDMA, O&P	SA615	42
Adult Wellness with SDMA, O&P, <i>Giardia</i>	SA655	43
Adult Wellness with SDMA, T4	SA670	43
Adult Wellness with SDMA, T4, O&P	SA674	43

TEST NAME	CODE	PAGE
Adult Wellness with SDMA, T4, O&P, UA	SA676	43
Adult Wellness with SDMA, T4, UA	SA672	43
Adult Wellness with SDMA, T4, UA, O&P, <i>Giardia</i> , Accuplex®	AC608	31
Adult Wellness with SDMA, UA, Fecal Combo	SA649	44
Adult Wellness with SDMA, Urinalysis	SA607	44
Airway Wash	BAL	180
Airway Wash	TTW	180
Albumin	T010	135
Alkaline Phosphatase	T020	135
Alkaline Phosphatase Isoenzymes	T215	135
ALT SGPT	T030	135
Amylase	T040	135
Amylase, Lipase	T050	136
<i>Anaplasma</i> CF, Bovine	S16035	183
<i>Anaplasma phagocytophilum</i> Ab	S16872	99
<i>Anaplasma phagocytophilum</i> PCR	S14421	173
Anti-Mullerian Hormone (C/F)	S14410	91
Anti-Mullerian Hormone (Equine)	S14320	176
Anticoagulant Screen	S16040	136
Antinuclear Antibodies ANA	T515	136
APTT	T395	73
<i>Aspergillus</i> AB Avian (AFMP1P ELISA)	S16011	151
<i>Aspergillus</i> Ag, Galactomannan Assay	S85358	151
<i>Aspergillus</i> Profile Avian	S85359	151
AST SGOT	T060	136
Autoimmune Profile	SA280	136
Autoimmune Profile Standard	SA170	136
Avian Blood Lead	AE290	151
Avian Comprehensive Chemistries	AE010	152
Avian Comprehensive Profile with Bile Acid	AE021	152
Avian Comprehensive Profile with Bile Acid, EPH	AE025	152
Avian Protein Electrophoresis	AE300	152
Avian Standard Chemistries	AE040	152
Avian Standard Profile with Bile Acid	AE051	153
Avian Standard Profile with Bile Acid, EPH	AE055	153
Avian Zinc Assay	S16012	153
Avian/Exotic CBC	AE270	153
<i>Babesia</i> Bovis	S16124	183

TEST NAME	CODE	PAGE
<i>Babesia canis</i>	S16070	99
<i>Babesia gibsoni</i>	S16075	99
<i>Babesia gibsoni</i> IFA, Export	S16502	99
Baermann	T785	183
<i>Bartonella henselae</i> , ELISA	S85889	99
Basic Wellness Screen with SDMA, Accuplex®	AC710	31
Bicarbonate	T115	136
Bile Acids	T225	79, 166
Bile Acids Avian/Exotic	AE260	153
Bile Acids Pre and Post	T220	79
Bilirubin, Direct	DBIL	136
Bilirubin, Total	T090	137
Bladder Tumor Analytes	T520	137
<i>Blastomyces</i> Antibody	T525	137
<i>Blastomyces</i> Quantitative Ag Assay	S86293	137
Blood Cross Match	T340	137
Blood Type, Canine - Full Panel	S16100	73
Blood Type, Canine DEA 1.1 Only	T315	73
Bone Marrow and CBC	BONECBC	117
Bone Marrow Core Biopsy	BMBC	117
Bone Marrow Cytology	BONE	117
Bovine Herpes Virus	S14456	184
Bovine Viral Diarrhea (BVD) - ELISA	S16844	184
Bovine Viral Diarrhea (BVD) PCR	S86551	184
Bromide	T730	149
Bromide Add-on	ADD320	149
<i>Brucella canis</i> - Screen	T530	100
<i>Brucella</i> Screen Multiplex	S16131	100
<i>Brucella</i> Titer For Export (KSU)	S16003	138
Buffy Coat For Mast Cells	T325	117
BUN Blood Urea Nitrogen	T100	138
BUN, Creatinine, SDMA	T105	138
Bun, Creatinine, SDMA, T4	T107	138
C1 with SDMA, FT4-ED	SA205	44
C5 Cat Viral	SA265	100
CADET® <i>BRAF</i>	T1025	123
Calcium	T110	80
Calcium, Ionized	S18537	80

TEST NAME	CODE	PAGE
Calici Virus - IFA	S16135	100
Canine Adult Wellness Immune Profile with SDMA	SA645	44
Canine Adult Wellness with SDMA, UA	SA625	45
Canine Comprehensive with SDMA	SA100	45
Canine Distemper Titer IgG, IgM	T555	100
Canine <i>Ehrlichia</i> / <i>Anaplasma</i> PCR	T980	123
Canine GI PCR Panel	T950	124
Canine GI PCR with SARS CoV-2 PCR Panel	T953	124
Canine Heartworm Program	SA140	138
Canine Heartworm Program Plus	SA130	138
Canine Melanoma Diagnostic Panel	S86793	118
Canine Melanoma Prognostic Panel	S86792	118
Canine Parvovirus Antibody Titer	T690	138
Canine Parvovirus PCR	S8710	124
Canine Post Pill T4 Add-on	ADD300	91
Canine Respiratory PCR Panel	T995	125
Canine Respiratory PCR with SARS CoV-2 PCR Panel	T998	125
Canine Senior Profile with SDMA, Heartworm	SA710	45
Canine Senior with SDMA, Fecal Combo	SA765	45
Canine Tick Borne PCR Panel	T960	125
Canine Tick Borne PCR with Lyme	T961	126
Canine Wellness with SDMA, O&P	SA635	45
Canine Wellness with SDMA, UA, O&P, <i>Giardia</i> , Accuplex®	AC651	31
Canine Wisdom Panel™ Premium	S14497	133
Canine Wisdom Panel™ Premium, Add-on	S14515	134
Canine Wisdom Panel™ with Keyscreen®	KS14497	134
Canine/Feline Ringworm PCR Panel	T982	126
Canine/Feline Ringworm PCR with Dermatophyte Culture	T986	126
Caprine Arthritis Encephalitis (CAE)	S16145	184
Cardio BNP-Canine	T1010	139
Cat Scan Plus with SDMA	SA220	46
Cat Scan with SDMA	SA230	46
Cat Scan With SDMA, T4	SA1348	46
Cat Wellness Profile with SDMA	SA610	46
CBC And Fibrinogen	T332	164, 183
CBC Recheck	T330R	74
CBC Small Mammalian	AE275	74
Chemistry Panel with SDMA	SA804	80

TEST NAME	CODE	PAGE
Chemistry Renal Profile with SDMA	T7008	80
<i>Chlamydia</i> Antibody Titer IFA	S16671	154
<i>Chlamydia</i> PCR Blood	S16788	154
<i>Chlamydia</i> PCR Swab	S16672	154
<i>Chlamydia</i> Profile	S85206	154
<i>Chlamydia</i> Titer EBA	S16670	154
Chloride	T120	139
Cholesterol	T125	139
Cholinesterase - Serum	T235	139
Chronic Diarrhea Panel	L420	169
<i>Clostridium</i> Difficile Toxins A/B	L86181	168
<i>Clostridium</i> Perfringens Enterotoxin	T16007	168
Coagulation Panel with D-Dimer	SA290	74
Coagulation Profile 2	SA300	74
Coagulation Profile 3	SA305	75
Cobalamin B12	T840	139
Cobalamin B12/Folate	S16195	139
Cobalamine, Folate, TLI Canine	SA160	139
Cobalamine, Folate, TLI Feline	SA275	140
Cocci AGID Add-on	ADD04	100
Cocci Profile with SDMA, T4	SA132	46
Coccidioidomycosis, Screen and Titer	T535	101
Complete Blood Count	T330	75
Complete Blood Count with Retic Count	T337	75
Comprehensive Avian Post Purchase	AE030	154
Comprehensive Avian Profile	AE020	155
Comprehensive Mammalian Profile	AE200	155
Comprehensive Reptilian (Recheck)	REREP	155
Comprehensive Reptilian Profile	AE160	155
Copper Level	S16210	140
Copper Storage Disease	S16215	140
Corona Virus IgG and IgM	S16225	140
Cortisol	T445	91
Cortisol Serial ACTH 2	ACTH2	91
Cortisol Serial ACTH 3	ACTH3	92
Cortisol Serial ACTH 4	ACTH4	92
Cortisol Serial DEX 2	DEX2	92
Cortisol Serial DEX 3	DEX3	92

TEST NAME	CODE	PAGE
Cortisol Serial DEX 4	DEX4	93
CPK	T130	140
Creatinine	T135	140
Cross Match, Additional Donor 2	T345	137
Crypto, <i>Giardia</i> , <i>Clostridium</i> Enterotoxin	SA350	141
<i>Cryptococcal</i> Antigen	T550	101
Crystallographic Stone Analysis	S16735	87
CSF Protein	CSFPR	80
Culture Aerobic and Anaerobic	MO40	109, 177
Culture Aerobic and Fungal	MO50	109
Culture Blood Aerobic/Anaerobic	MO61	110
Culture Blood Aerobic/Anaerobic (Serial Blood Culture x2)	MO62	111
Culture Blood Aerobic/Anaerobic (Serial Blood Culture x3)	MO63	111
Culture ID	MO70	112
Culture, Acid Fast Bacilli	M100	114
Culture, Aerobic	MO20	112, 177
Culture, Anaerobic	MO30	114, 177
Culture, Blood Aerobic Only	MO60	110, 178
Culture, <i>Campylobacter</i>	CAMP	114
Culture, Dermatophytes	M240	115
Culture, Enhanced Fluid Aerobic	M220	178
Culture, Feces	M160	115, 168
Culture, Fungal	MO80	115, 178
Culture, <i>Mycoplasma</i>	M110	112
Culture, <i>Salmonella</i>	M121	168
Culture, Urine	M130	113, 178
Culture, Urine Add-on	ADD210	113
Culture, Urine Recheck	M130R	114
Cyclosporine	S18702	149
Cytology	CYTO	118, 179
D-Dimer	T350	141
D1 with SDMA, FT4-ED	SA105	47
Dermatopathology	DERM	118
Desert Disease Panel with SDMA, T4	SA504	47
Diabetes Monitoring Panel with SDMA	SA800	47
Diarrhea Profile Avian	AE070	156
Digoxin	T735	149
Direct Coombs' Test Warm	T540	75

TEST NAME	CODE	PAGE
Distemper Parvo Titer Add-on	ADD270	141
Distemper Parvo Vaccinal Titer	T565	141
Distemper PCR	S16501	156
Distemper Smear, FA	S16250	141
Distemper Titer-Virus Neutraliz.	S16107	156
Distemper Vaccinal Titer	T560	142
Distemper/Parvo IgG End Point Titer	T675	142
Dog Wellness Profile with SDMA	SA605	47
Dog Wellness Profile with SDMA, O&P, <i>Giardia</i>	SA650	47
Dog Wellness with SDMA, T4	SA685	48
Ear Culture	M255	115
<i>Ehrlichia canis</i>	T570	101
<i>Ehrlichia canis</i> Add-on	ADD05	101
Ehrlichiosis Serology Panel Canine	S16900	101
Electrolyte Screen	T140	81
Encephalitis	S16275	172
<i>Encephalitozoon cuniculi</i> IgG Ab	S16877	156
Endogenous ACTH	T435	93
Endogenous ACTH (Baseline) Equine	L525	170
Endogenous ACTH, Insulin Panel	L540	170
Endogenous ACTH, Insulin/Glucose, Leptin (Equine)	L575	170
Environmental Culture	M076	115
EPM SAG 2, 4/3 ELISA	S14388	172
EPM SAG 2, 4/3 ELISA Serum/CSF Ratio	S14390	172
EPM SAG 2, 4/3 ELISA, <i>N. hughesi</i> ELISA	S14392	172
Equine 11 Metal & Mineral Panel	S14466	181
Equine Chemistry, CBC	L050	165
Equine Chemistry, CBC, Fib	L040	165
Equine Chemistry, CBC, Fib (RECHECK)	L040R	165
Equine Chemistry, CBC, Fib, fT4 (ED)	L035	165
Equine Chemistry, CBC, Fib, T4	L030	165
Equine Chemistry, CBC, Fib, T4, fT4 (ED)	L565	165
Equine Comprehensive Neurological Panel	S14477	172
Equine Coronavirus (ECoV)	S14414	173
Equine Gastrointestinal PCR Panel	L950	173
Equine Granulosa Cell Tumor	S85857	176
Equine Health Program	L510	162
Equine Health Program (ELISA)	L511	162

TEST NAME	CODE	PAGE
Equine Health Program (no EIA)	L510NE	162
Equine Health Program GVL	L510E	162
Equine Health Program GVL (ELISA)	L511E	162
Equine Hepatic Panel	L230	166
Equine Hepatic Profile	L1025	166
Equine Hepatic Screen	L225	166
Equine Herpes Virus 1 (EHV1) PCR	S14394	174
Equine IgG Total	L090	165
Equine Infectious Anemia (EIA) AGID	L110	163
Equine Infectious Anemia (EIA) AGID by GVL	L111	163
Equine Infectious Anemia (EIA) ELISA	L120	163
Equine Infectious Anemia (EIA) ELISA by GVL	L121	164
Equine Inflammatory Profile	L070	164
Equine Insulin and Glucose	T470E	171
Equine Lepto PCR Aq/Vitreous Humor	L973	174
Equine <i>Leptospira</i> PCR Blood	L974	174
Equine <i>Leptospira</i> PCR Blood/Urine	L978	174
Equine <i>Leptospira</i> PCR Urine	L976	174
Equine Muscle Enzyme Screen	L275	166
Equine Muscle Enzyme Screen, CBC, FIB	L280	166
Equine PCR FUO Mini	L970	173
Equine PCR FUO Panel	L965	173
Equine PCR Reproductive Panel	L955	173
Equine PCR Respiratory Panel	L960	173
Equine Performance Horse Chemistry	L601	165
Equine Performance Horse Health Program	L610	162
Equine Performance Horse Health Program (no EIA)	L610NE	162
Equine Performance Horse Health Program GVL	L610E	162
Equine Performance Horse Panel	L640	165
Equine Pre-Purchase Drug Screen	S1680	167
Equine Renal Screen	L240	166
Equine Rhabdomyolysis Screen	L190	166
Equine Ringworm PCR Panel	T983	174
Equine Ringworm PCR with Dermatophyte Culture	T987	174
Equine Senior Health Program	L615	163
Equine Senior Health Program (no EIA)	L615NE	163
Equine Senior Health Program AGID GVL	L615E	163
Equine Senior Health Program with TRH	L615T	163

TEST NAME	CODE	PAGE
Equine TRH Stim (Thyroid Function)	L590	169
Equine Viral Arteritis	S16285	181
Equine Viral Encephalitis Panel	S17500	172
Estradiol	S16295	175
Estrone Sulfate	S16300	175, 176
Feather Picker Profile	AE080	156
FEC: Modified Stoll's Technique	T826	168
FEC: McMasters Method (Large Animal)	T828	168
Fecal Combo with Heartworm	T812	101
Feline Adult Wellness Immune Profile with SDMA	SA640	48
Feline Adult Wellness with SDMA, Fecal Combo	SA660	48
Feline Adult Wellness with SDMA, O&P, <i>Giardia</i> , UA	SA661	48
Feline Adult Wellness with SDMA, O&P, UA	SA696	49
Feline Adult Wellness with SDMA, O&P, UA, FeLV, FIV	SA629	49
Feline Adult Wellness with SDMA, UA	SA630	49
Feline Adult Wellness with SDMA, UA, FeLV, FIV	SA624	49
Feline Blood Typing	T320	76
Feline Comprehensive Plus with SDMA	SA200	50
Feline Comprehensive Wellness with SDMA, O&P	SA718	49
Feline Comprehensive with SDMA	SA210	50
Feline Coronavirus Exposure Titer	T593	101
Feline Coronavirus Titer	T595	102
Feline EDP with SDMA, FeLV, FIV, O&P	SA682	50
Feline Flea and Tick Borne PCR	T965	126
Feline GI PCR Panel	T955	127
Feline GI PCR with SARS CoV-2 PCR Panel	T958	127
Feline Health Check with SDMA	SA908	50
Feline Health Profile with SDMA	SA914	50
Feline Hemoplasma PCR Panel	T985	127
Feline Panel #101 with SDMA	SA460	51
Feline Respiratory PCR Panel	T990	128
Feline Respiratory PCR with SARS CoV-2 PCR Panel	T993	128
Feline Senior 2 with SDMA, Panleukopenia	SA782	51
Feline Senior Profile with SDMA	SA722	51
Feline Serology 2	SA270	102
Feline Serology, FeLV, FIV, Toxo AB IgG/IgM	SA262	102
Feline Total Health Check with SDMA	SA190	51
Feline Total Health Check with SDMA, UA	SA786	52

TEST NAME	CODE	PAGE
Feline Total Health Plus with SDMA	SA180	52
Feline Wellness 3 with SDMA	SA700	52
Feline Wellness 4 with SDMA	SA902	52
Feline Wellness with SDMA, FeLV, FIV, O&P	SA683	52
Feline Wisdom Panel™ Complete	S14498	133
Feline Wisdom Panel™ Complete, Add-on	S14516	133
FeLV Antigen ELISA	T580	102
FeLV Antigen ELISA Add-on	ADD06	102
FeLV IFA	T585	102
FeLV PCR	S6234	102
FeLV, FIV ELISA Add-on	ADD07	103
FeLV, FIV Heartworm Ab	SA269	103
FeLV, FIV Special	SA260	103
FeLV, FIV, Feline Heartworm Antigen	SA266	103
FeLV, FIV, KeyScreen®	KS260	40
Fibrinogen Quantitative	L365	164
Fibrinogen Quantitative	T365	76
FIP mRNA PCR	T600	128
FIV Antibody	T610	103
FIV Antibody Add-on	ADD15	103
FIV Antibody, Western Blot	S16865	142
Fluid Analysis	FLUA	119, 179
Fluid Analysis CSF	CSF	179
Fluid Part Only	FL	119
Foal Diarrhea Panel	L492	169
Fractional Excretion of Electrolytes	L340	167
Free T4 By Equilibrium Dialysis	T460	93, 169
Free T4 Equilibrium Dialysis Add-on	ADD50	93
Fructosamine Assay	S16345	93
Fungal Identification	MO82	116
Fungal Serology	SFUN	142
Fungal Serology with Cocci	SA340	142
Geriatric Feline with SDMA	SA910	53
Geriatric Panel with SDMA, Feline	SA924	53
Geriatric Profile with SDMA	SA037	53
Geriatric Weak Ferret Profile	AE230	157
GGTP	T145	81
GHP Chem with Lytes, SDMA	SA1202	81

TEST NAME	CODE	PAGE
GHP Chem with Lytes, SDMA, CBC	SA1204	53
GHP Chem with Lytes, SDMA, CBC, T4	SA1212	53
GHP Chem with Lytes, SDMA, CBC, T4, UA	SA1214	54
GHP Chem with SDMA	SA1102	81
GHP Chem with SDMA, CBC	SA1104	54
GHP Chem with SDMA, CBC, T4	SA1108	54
GHP Chem with SDMA, CBC, T4, UA	SA1114	54
GHP with Lytes, SDMA, CBC, UA	SA1210	54
<i>Giardia</i> ELISA	T820	107
<i>Giardia</i> ELISA Add-on	ADD250	107
<i>Giardia</i> (FA) and <i>Cryptosporidium</i> (FA)	T790	142
Glucose	T150	143
Gram Stain	MO90	116
Heartworm Ag-Heat-Treated Serum	T613	103
Heartworm Antibody, Antigen Feline	T630	104
Heartworm Antibody, Feline	T625	104
Heartworm Antigen	T615	104
Heartworm Antigen	T620	104
Heartworm Antigen Add-on	ADD70	104
Heartworm Antigen, Microfilaria	T617	104
Heartworm, KeyScreen®	KT615	40
Heartworm, Ova and Parasite	T618	104
Hemotrophic <i>Mycoplasma</i>	T380	104
Hemotropic <i>Mycoplasma</i> Add-on	ADD20	105
Hepatic Profile	AE090	157
Herpes - Conjunctival Smear	S16400	143
Herpes Antibody, IFA	S86022	143
Heska ALLERCEPT Environmental & Food Panel	S17029	143
Heska ALLERCEPT Environmental Panel	S17026	143
Heska ALLERCEPT Food Panel	S17028	143
Heska Equine ALLERCEPT Panel	S17027	143
Histopathology	FBX	119
<i>Histoplasma</i> Antibody	T640	105
<i>Histoplasma capsulatum</i> Ag, EIA	S86569	105
<i>Histoplasma</i> Titer	S16405	105
Hyperthyroid Feline with SDMA	SA235	54
Hyperthyroid Panel (with GGT), SDMA	SA805	55
Hyperthyroid Panel with Glucose, SDMA	SA805G	55

TEST NAME	CODE	PAGE
Hyperthyroid Panel with SDMA	SA440	55
Hypothyroid Post Pill T4	T497	93
IgG Bovine	S16425	184
IgG Llama	S16430	184
Immunoglobulins A, G, M	T660	143
Immunohistochemistry 1 Stain	IHC1	119
Immunohistochemistry 2 Stains	IHC2	120
Immunohistochemistry 3 Stains	IHC3	120
Immunohistochemistry 4 Stains	IHC4	120
Immunohistochemistry 5 Stains	IHC5	120
Influenza - Canine Acute	S86096	105
Inhibin	S4131	176
Insulin Feline with Glucose	T470F	94
Insulin with Glucose	T470	94
Iron, Serum	T155	144
Johne's Disease Antibody	S16302	184
Keppra Level	S86541	149
KeyScreen® GI Parasite PCR Add-on	ADD350	128
KeyScreen® GI Parasite PCR Panel	T991	129
Large Animal Chemistry	L010	164
LDH Lactic Dehydrogenase	T160	144
Lead Level	T745	144
Leptin Equine	S14402	170
<i>Leptospira</i> PCR Blood	T974	130
<i>Leptospira</i> PCR Blood/Urine	T978	130
<i>Leptospira</i> PCR Urine	T976	130
Leptospirosis	S16510	105, 181
Liver Chem, CBC with Phenobarbital	SA328	55
Liver Chemistries	SA321	81
Liver Chemistry	SA324	82
Liver Chemistry Panel	SA326	82
Liver Chemistry with Phenobarbital	SA329	82
Liver Chemistry, CBC	SA327	55
Liver Profile	SA1502	56
Liver Profile Standard	SA320	56
Luteinizing Hormone	S16520	94
Lyme Equine Multiplex	T672	181
Lyme Titer IgG	T670	105

TEST NAME	CODE	PAGE
Magnesium	T170	144
Malignancy Profile	S86698	94
Mammalian Comprehensive Chemistries	AE190	157
Mammalian Standard Chemistries	AE210	157
Mast Cell Panel	S86601	121
Masticatory Muscle Myositis	S16535	144
MIC Urine If Indicated	M134	87
Microfilaria - Knotts	T390	106
Mini Avian Post Purchase	AE060	157
Mini Early Detection Chem with SDMA	SA1002	82
Mini Early Detection Chem with SDMA, CBC	SA1004	56
Mini Early Detection Chem with SDMA, CBC, O&P, Heartworm	SA1010	56
Miniscreen 11 Chem with SDMA	SA071	82
Miniscreen 4 Chem	SA060	83
Miniscreen with CBC	SA070	56
Miniscreen with SDMA, CBC	SA072	56
Miniscreen with SDMA, T4	SA1519	57
<i>Mycoplasma</i> PCR	S16789	157
Neonatal Foal Panel	L290	165
<i>Neorickettsia risticii</i> PCR (PHF)	S14479	174
<i>Neospora caninum</i> - IFA	S16560	106
NSAID 1 with SDMA	SA810	57
NSAID 2 with SDMA	SA815	57
NSAID 3 with SDMA	SA820	57
NSAID Chemistries with SDMA	SA822	83
O&P Centrifugation, Smear	T806	107
O&P with Accuplex®	AC805	40
O&P, <i>Giardia</i> with Accuplex®	AC808	40
Occult Blood, Feces	T810	144
OFA Thyroid Panel	SA1528	95
OncoK9®	S14493	130
Ophthalmic Aerobic Culture	M280	116
Oral Path Biopsy	FBXORL	121
Oral Sugar Test (2 samples)	L545	171
Oral Sugar Test (3 samples)	L550	171
Osmolality - Serum	S16575	144
Ova and Parasite	T805	108, 184
Ova and Parasite with <i>Giardia</i>	T808	108

TEST NAME	CODE	PAGE
Pancreatic Lipase Immunoreactivity	S85364	144
Pancreatitis Profile - Canine	S86468	145
Pancreatitis Profile - Feline	S86288	145
Panleukopenia Titer IgG, IgM	S16580	145
Panleukopenia Vaccinal Titer	S16053	145
Panleukopenia Vaccinal Titer Add-on	ADD280	145
Parasite Identification	T825	146
Parathormone Related Protein	S16596	95
PARR	S85503	131
Parvo Antibody, Antigen	T700	146
Parvovirus Antigen	T695	106
Parvovirus Vaccinal Titer	T705	146
Path Review (Urinalysis)	T764	87
Path Review CBC	REVW	76
Pet Cancer Specialty Biopsy	FBXNEO	121
Phenobarbital	T750	149
Phenobarbital Add-on	ADD315	149
Phenobarbital Panel Plus	SA830	150
Phosphorus	T180	146
Piroplasmosis, Equine	S16848	181
Pituitary Pars Intermedia Dysfunction (PPID) with fT4ED	L560	170
Pituitary Pars Intermedia Dysfunction (PPID) with TRH Stim	L500TRH	169
Pituitary Pars Intermedia Dysfunction Monitoring Panel	L500	169
Platelet Count	T400	76
Platelet Count Manual	T401	77
Post Cortisol	T448	95
Potassium	T185	146
Pre NSAID Use Panel 3 with SDMA	T85360	83
Pre-Op Chem with Electrolytes, SDMA	SA043	83
Pre-Op Panel Plus with Lytes, SDMA	SA054	57
Pre-Op Panel Plus with SDMA	SA050	58
Pre-Op Panel with Electrolytes, SDMA and Urinalysis	SA053	58
Pre-Op Panel with SDMA, T4, UA	SA920	58
Pre-Op Plus with SDMA, CBC, FeLV, FIV	SA510	58
Pre-Op Plus with SDMA, CBC, Heartworm	SA512	58
Pre-Op Screen with Electrolytes, SDMA, CBC	SA516	59
Pre-Op Screen with SDMA	SA040	84
Pre-Op Screen with SDMA, CBC	SA055	59

TEST NAME	CODE	PAGE
Pre-Op Screen with SDMA, CBC, T4	SA508	59
Pre-Op Screen with SDMA, CBC, T4, FeLV, FIV	SA514	59
Pre-Op Screen with SDMA, CBC, UA	SA052	59
PrecisionPSL	T165	83
PrecisionPSL Add-on	ADD90	84
Pregnant Mare Serum Gonadotropin (PMSG)	S16635	175
Progesterone	L140	175
Progesterone	T475	95
Progesterone/Estradiol Panel (Equine)	L470	175
Progesterone/PMSG	L460	175
Progesterone/PMSG/Estradiol	L480	175
Protein Electrophoresis (Serum)	T240	84, 167
Protein Electrophoresis (Urine)	T245	147
Protein Electrophoresis Add-on	ADD130	147
Prothrombin Time	T410	77
Psittacine Beak and Feather Disease	S16085	158
PT and aPTT	T415	77, 164
PT and APTT Add-on	ADD290	77
PTH, Ionized Calcium	S16595	95
PU/PD Profile	AE110	158
Rabbit Neurologic Profile	AE240	158
Rabbit/Small Rodent Ringworm PCR Panel	T984	132
Rabbit/Small Rodent Ringworm PCR with Dermatophyte Culture	T988	132
Rabies Diagnostic Non-Export	S1204	106
Rabies Export - FAVN	S17108	106
Recheck Comprehensive Avian	RECHECKAV	158
Recheck Cytology	RCYTO	122
Recheck Profile Standard with SDMA	RECHECK	59
Recheck Profile Vetscreen with SDMA, CBC	RECHECK2	60
Recheck Profile with SDMA	RECC	84
Relaxin	T9810	147
Renal Profile Standard with SDMA	SA313	60
Renal Profile with SDMA	SA310	60
Renal Profile with SDMA	SA311	84
Reptilian Comprehensive Chemistries	AE150	158
Reptilian Standard Chemistries	AE170	158
Reticulocyte Count	T425	78
Reticulocyte Count Add-on	ADD140	78

TEST NAME	CODE	PAGE
Rhinotrachitis Feline	S16702	147
<i>Rhodococcus equi</i> PCR	S14396	174
Rocky Mt Spotted Fever IFA	T715	106
Ruminant Chem with Mag, CBC	L080	183
Ruminant Chemistry with Magnesium	L1001	183
<i>Salmonella</i> Spp. PCR	S14416	175
<i>Salmonella</i> , <i>Shigella</i> , <i>Campylobacter</i>	M120	116
SARS-CoV-2 PCR Add-on Panel	T996	132
SDMA	T1035	85
SearchLight DNA Cytology Add-on	S10058	122, 133
SearchLight DNA Histology Add-on	S11068	122, 133
Selenium Level	S16730	147
Semen Analysis	SEMEN	176
Senior Care Basic Plus Profile with SDMA	SA712	61
Senior Comprehensive Plus with SDMA	SA080	61
Senior Comprehensive Plus with SDMA, UA	SA081	67
Senior Comprehensive with SDMA	SA090	61
Senior Comprehensive with SDMA, D.Bilirubin, I.Bilirubin, Anion, Osmolality	SA086	61
Senior Comprehensive with SDMA, Heartworm, UA, UPCR	SA784	67
Senior Feline Plus with SDMA	SA740	61
Senior Panel 2 with SDMA, O&P	SA755	62
Senior Panel 2 with SDMA, UMIC	SA746	62
Senior Panel 2 with SDMA, Vaccine Titers	SA760	62
Senior Profile 1 with SDMA	SA705	63
Senior Profile 1 with SDMA, Accuplex®	AC705	32
Senior Profile 1 with SDMA, Accuplex®, KeyScreen®	KA705	33
Senior Profile 1 with SDMA, Fecal Combo	SA770	64
Senior Profile 1 with SDMA, Feline Heartworm Ab, FeLV, FIV, Fecal Combo	SA954	62
Senior Profile 1 with SDMA, Feline Heartworm Ag	SA719	64
Senior Profile 1 with SDMA, Feline Heartworm, KeyScreen®	KS722	33
Senior Profile 1 with SDMA, FT4ED, Heartworm	SA735	64
Senior Profile 1 with SDMA, Heartworm, KeyScreen®	KS710	33
Senior Profile 1 with SDMA, KeyScreen®	KS705	33
Senior Profile 1 with SDMA, O&P	SA750	63
Senior Profile 1 with SDMA, O&P, Accuplex®	AC750	32
Senior Profile 1 with SDMA, UMIC	SA566	63
Senior Profile 1 with SDMA, Vaccine Titers, Fecal Combo, Accuplex®	AC727	35
Senior Profile 1 with SDMA, Virals, Fecal Combo	SA741	63

TEST NAME	CODE	PAGE
Senior Profile 1 with SDMA, Virals, Feline Heartworm Ab	SA715	68
Senior Profile 2 with SDMA	SA720	64
Senior Profile 2 with SDMA, Accuplex®	AC720	34
Senior Profile 2 with SDMA, Accuplex®, KeyScreen®	KA720	34
Senior Profile 2 with SDMA, Fecal Combo	SA775	65
Senior Profile 2 with SDMA, Fel Heartworm Ab	SA766	65
Senior Profile 2 with SDMA, FeLV, FIV	SA780	65
Senior Profile 2 with SDMA, FeLV, FIV, KeyScreen®	KS780	34
Senior Profile 2 with SDMA, KeyScreen®	KS720	34
Senior Profile 2 with SDMA, O&P, Fel Heartworm Ab	SA790	65
Senior Profile 2 with SDMA, O&P, <i>Giardia</i> , Accuplex®	AC775	35
Senior Profile 2 with SDMA, O&P, Heartworm	SA762	66
Senior Profile with SDMA, Feline Heartworm Ab, Fecal Combo	SA768	66
Senior Profile with SDMA, Panleukopenia	SA725	66
Senior Profile with SDMA, Vaccine Titers	SA730	66
Senior Super Plus with SDMA	SA745	67
Senior Wellness with SDMA, Feline Heartworm Ab, O&P	SA764	67
Sodium	T195	147
Sodium And Potassium Combo	T200	147
Sorbital Dehydrogenase (SDH)	T250	166
Special CBC, Path Review	T331	78
Standard Avian Profile	AE050	159
Standard Mammalian Profile	AE220	159
Standard Reptilian Profile	AE180	159
Stat Charge Biopsy	STAT	122
Stat Charge-Cyto	CYSTAT	122
<i>Streptococcus</i> Equi ELISA	S16740	182
<i>Streptococcus equi equi</i> PCR	S86308	175
Superchem with SDMA	SA010	85
Superchem with SDMA, CBC	SA020	68
Superchem with SDMA, CBC, Accuplex®, KeyScreen®	KA020	35
Superchem with SDMA, CBC, Accuplex®, T4, KeyScreen®	KA114	35
Superchem with SDMA, CBC, Feline Heartworm, FIV, FeLV, KeyScreen®	KS737	36
Superchem with SDMA, CBC, FIV, FeLV, KeyScreen®	KS490	36
Superchem with SDMA, CBC, Free T4 ED	SA1601	68
Superchem with SDMA, CBC, Heartworm Antigen	SA028	68
Superchem with SDMA, CBC, Heartworm, KeyScreen®	KS028	36
Superchem with SDMA, CBC, Heartworm, T4, KeyScreen®	KS110	36

TEST NAME	CODE	PAGE
Superchem with SDMA, CBC, KeyScreen®	KS020	36
Superchem with SDMA, CBC, T4, Feline Heartworm, FIV, FeLV, KeyScreen®	KS724	37
Superchem with SDMA, CBC, T4, FeLV, FIV, HeartwormAB	SA724	60
Superchem with SDMA, CBC, T4, FIV, FeLV, KeyScreen®	KS220	37
Superchem with SDMA, CBC, T4, FT4ED, Heartworm	SA098	60
Superchem with SDMA, CBC, T4, KeyScreen®	KS120	37
Superchem with SDMA, CBC, T4, UA, Feline Heartworm, FIV, FeLV, KeyScreen®	KS715	37
Superchem with SDMA, CBC, T4, UA, FIV, FeLV, KeyScreen®	KS700	38
Superchem with SDMA, CBC, T4, UA, O&P, <i>Giardia</i> , Accuplex®	AC770	32
Superchem with SDMA, CBC, T4, UA, UPC, Accuplex®	AC712	32
Superchem with SDMA, CBC, Thyroid Panel (T3, T4, FT4ED, TSH, TGAA)	SA450	68
Superchem with SDMA, CBC, UA	SA021	68
Superchem with SDMA, CBC, UA, Accuplex®, KeyScreen®	KA021	38
Superchem with SDMA, CBC, UA, Feline Heartworm, FIV, FeLV, KeyScreen®	KS494	38
Superchem with SDMA, CBC, UA, FIV, FeLV, KeyScreen®	KS702	38
Superchem with SDMA, CBC, UA, Heartworm, KeyScreen®	KS112	39
Superchem with SDMA, CBC, UA, KeyScreen®	KS021	39
Superchem with SDMA, CBC, Viral Panel	SA490	69
Superchem with SDMA, T4	SA120C	85
Superliver (Liver Biopsy with Liver Staining Panel)	FBXLSP	122
Synovial Fluid Analysis	SYFLUA	179
T3	T480	96
T3 Autoantibodies	T485	148
T3 Suppression	SA430	96
T4	T495	96, 169
T4 Add-on	ADD190	96, 169
T4 Autoantibodies	T500	148
T4, TSH	SA401	96
Taurine	S16755	148
Testosterone	S16760	96, 176
Testosterone Panel	S85530	176
Thyroglobulin Auto Antibody	T505	96
Thyroid 1	SA360	97
Thyroid Profile 2	SA370	97, 169
Thyroid Profile 3	SA380	97
Thyroid Profile 5	SA400	97
Thyroid Profile 6	SA410	97
Tick Borne Disease Panel	SA330	148

TEST NAME	CODE	PAGE
TLI, Feline	S16800	148
Total Body Function Plus with SDMA, Fecal Combo	SA788	69
Total Body Function Plus with SDMA, O&P	SA114	69
Total Body Function with SDMA	SA120	69
Total Body Function with SDMA, Feline Virals	SA225	70
Total Body Function with SDMA, Heartworm	SA110	69
Total Body Function with SDMA, O&P	SA122	70
Total Body Function with SDMA, O&P, Accuplex®	AC114	39
Total Body Function with SDMA, O&P, <i>Giardia</i> , Accuplex®	AC123	39
Total Body Function with SDMA, TSH	SA1633	70
Total Protein	T190	148
<i>Toxoplasma gondii</i> PCR	S18708	107
<i>Toxoplasma</i> Ab - IgG/IgM (Feline)	T720	106
<i>Toxoplasma</i> antibody - Exotics	S16792	159
Toxoplasmosis IgG/IgM	S85030	107
TRH Stimulation Test for PPID	L535	170
TRH Stimulation with Insulin Glucose	L580	171
Triglyceride, Chol. Panel	T127	148
Triglyceride	T205	148
<i>Tritrichomonas</i> PCR	S85819	107
Troponin I	S85783	85
Trypsin-Like Immunoreac Canine cTLI	T230	149
TSH	T510	97
TSH Add-on	ADD200	98
TSH, Free T4 ED	SA390	98
Urinalysis Add-on	ADD220	87
Urinalysis And UPC Ratio	T925	88
Urinalysis Recheck	T760R	88
Urinalysis-Complete	T760	88, 167
Urine Bile Acid: Creatinine Ratio	T227	88
Urine Clearance Ratios	T765	89
Urine Cortisol/Creatinine Ratio	T770	89
Urine Culture Aerobic and Anaerobic	M135	116
Urine Microalbumin Canine	T830C	89
Urine Microalbumin Feline	T830F	89
Urine Protein/Creatinine Add-on	ADD230	90
Urine Protein/Creatinine Ratio	T775	90
Urine Uric Acid/Creatinine	T310	90

TEST NAME	CODE	PAGE
UTI Cystitis Profile with MIC Urine	M133	90
Vaccine Panel Panleukopenia, Rhino, Calici	S16581	107
Vet Screen 2 Plus with SDMA, CBC	SA035	70
Vet Screen with Amylase, SDMA	SA026	86
Vet Screen with SDMA	SA025	85
Vet Screen with SDMA, CBC	SA030	70
Vet Screen with SDMA, CBC, FeLV, FIV, UA	SA906	71
Vet Screen with SDMA, CBC, FeLV, T4	SA1645	71
Vet Screen with SDMA, CBC, Heartworm	SA039	71
Vet Screen with SDMA, CBC, Heartworm, T4	SA530	70
Vet Screen with SDMA, CBC, O&P, <i>Giardia</i>	SA905	71
Vet Screen with SDMA, CBC, PT, PTT	SA031	71
Vet Screen with SDMA, CBC, T4	SA034	71
Vet Screen with SDMA, CBC, UA, O&P, Heartworm	SA904	72
Vet Screen with SDMA, CBC, Urinalysis	SA036	72
Vet Screen with SDMA, T4	SA1656	86
Vet Screen with SDMA, Thyroids	SA912	72
Vitamin D	S16016	86
Vitamin E	S16850	182
Vitamin E, Selenium	S17505	182
Von Willebrand Factor	S17123	78
Wellness Profile with SDMA, O&P, UA	SA647	72
West Nile Titer (PRNT & IgM ELISA)	S85448	159, 172
West Nile Virus PCR	S85449	160
Zinc	S16870	149
Zonisamide Zonegran	S86480	150
Zoogen Avian DNA Sexing	S16095	160

**BETTER DIAGNOSTICS.
BETTER MEDICINE.**



ANTECH[®]

THE

RIGHT DIAGNOSTICS

AT THE RIGHT TIME

We are proud to deliver the most comprehensive portfolio of veterinary products and services to you and your patients, now all offered by Antech.



ANTECH

1-800-872-1001



HESKA

1-800-464-3752



sound

1-800-268-5354



AIS

1-877-727-6800

EXPERIENCE A COMPLETE DIAGNOSTIC SOLUTION

Antech partners with veterinary care teams by providing a flexible, complete diagnostic portfolio that empowers you to practice medicine your way — because better diagnostics, means better medicine.

LEARN HOW TODAY



[ANTECHDIAGNOSTICS.COM/LEARN](https://AntechDiagnostics.com/Learn)

© 2023 Antech Diagnostics, Inc. All rights reserved. Antech, the Antech logo, and all other trademarks used herein are the registered trademarks of Antech Diagnostics, Inc. or its affiliates. 101823