



ANTECH




2024

DIRECTORY OF SERVICES

[ANTECHDIAGNOSTICS.COM](https://antechdiagnostics.com)

CANADA



BETTER DIAGNOSTICS. BETTER MEDICINE.

You deserve a choice when making diagnostic decisions for your patients.

Antech is your diagnostics partner with North America's largest reference laboratory network, best-in-class in-house diagnostics from Heska, the industry's most-trusted imaging equipment from Sound, and breakthrough telemedicine from Antech Imaging Services.

**The right diagnostics at the right time
— all from Antech.**



NORTH AMERICA'S LARGEST REFERENCE LAB NETWORK

- Core diagnostics, including hematology, chemistry, parasitology, microbiology, and clinical and anatomic pathology
- Novel innovations like KeyScreen® GI Parasite PCR — the most advanced GI parasite screening solution in veterinary medicine
- Medical consulting services provided by board-certified specialists across all key specialties



HESKA®

IN-HOUSE DIAGNOSTICS BY HESKA, AN ANTECH COMPANY

- Most comprehensive in-clinic diagnostics portfolio on the market
- Innovative, accurate, and cost-effective core diagnostics including hematology, chemistry, immunodiagnostics, and ova and parasites
- Emergency and critical care diagnostics including blood, gases, and acid-base
- Best-in-class serological allergy assessment
- Rapid infectious disease testing with expanding menu (heartworm, FIV/FelV)



HESKA'S IN-HOUSE DIAGNOSTIC SUITE, NOW OFFERED BY ANTECH



DIAGNOSTIC IMAGING AND EDUCATION FROM SOUND

- Digital radiology for superior image quality and diagnostic capabilities
- A comprehensive, accessible ultrasound portfolio
- Advanced therapeutic laser solutions for impactful treatment
- Efficient mobile and fixed CT solutions
- Innovative PACS to view your images confidently from any modality
- Gold standard for training and education



TELEMEDICINE AND AI-POWERED SOLUTIONS

- Expert interpretations of digital radiographs
- On-the-spot interpretations of ultrasounds
- Consultations by the world's largest team of board-certified veterinary radiologists
- AI-powered interpretations of routine digital radiographs in just minutes in many circumstances



HESKA®

IN-HOUSE DIAGNOSTICS BY HESKA, AN ANTECH COMPANY



The world's only fully automated
fecal and urine point-of-care lab

Element AIM™ works to perfectly prepare, automate, and take images of feces and urine. It finds important cells, bacteria, crystals, casts, and parasites (including roundworm, hookworm, tapeworm, whipworm, *Giardia*, and more) using a fast and sanitary load-and-go workflow. The Element AIM provides quality results and images to you and your clients — faster than ever.

Learn more

heska.com/element-aim

- Sealed test cassettes eliminate the mess and smell of handling slides and cover slips
- Superior lab-quality optics for crystal clear images in just minutes
- Digital images are easily stored in medical records and shared with clients and specialists
- Automation removes user error and wasted time at the microscope and during sample prep



DIAGNOSTIC IMAGING AND EDUCATION FROM SOUND



Sound SmartDR®

An unparalleled digital radiography system



Sound HD Dental

Superior diagnostic intraoral and periodontal images



GE® Ultrasound

Helping veterinarians incorporate ultrasound into practice since 1996



doctorVet Plus Therapy Laser

Everything you need in a therapy laser and nothing that you don't



CereTom® CT

A portable, wireless, self-shielded CT scanner to roll into your surgical suite



OEC C-Arm

High-quality images from a proven, dependable, classic C-Arm



Sound Assurance®

Protecting your investment and your practice

ACADEMY of VETERINARY IMAGING

Education

Possess the training to perform specific ultrasound exams

Sound imaging equipment and services are not only high-tech, but they are also intuitive. We are here for you and your patients, from customer service and tech support to customer education.

Speak with a Sound Representative today

1-800-268-5354

soundvet.com/request-demo



TELEMEDICINE AND AI-POWERED SOLUTIONS

Antech Imaging Services (AIS®) is here to provide your clients and their pets with the best care possible by offering board-certified specialists in radiology, cardiology, internal medicine, avian/exotics, oncology, neurology, and surgery — all by using cutting-edge digital technology partnered with the best practices of traditional veterinary medicine.

AIS PennHIP — the most effective hip screening tool available for dogs.

All dogs can benefit from PennHIP testing. For dogs at risk of developing hip osteoarthritis (OA), early intervention can help prevent or lessen the severity of canine hip dysplasia (CHD).

It gives an estimate of the risk for painful osteoarthritis (OA) of canine hip dysplasia (CHD) later in life. With this information, preventive and palliative strategies can be recommended by the PennHIP-trained veterinarian.

AIS RapidRead™ — accurate radiology reports at the speed of light.

AIS Radiology RapidRead combines an ever-expanding, continuously trained neural network and a world-class team of radiologists and data scientists that delivers:

- 63+ findings and growing
- Contextual assessments for deeper insight so you can diagnose and treat your patient with confidence, sooner
- We continue to add new findings and will soon add ECG, Dental, CT, and MRI to our AIS RapidRead platform
- AIS online picture archiving and communication system (PACS) service is available 24 hours per day, 7 days a week, 365 days a year to securely store and digitally transmit electronic images and clinically relevant reports

Find a specialist that's right for you

1-877-727-6800

**[antechimagingervices.com/
antechweb/contact](https://antechimagingervices.com/antechweb/contact)**

Our goal is to build close working relationships between practitioners and specialists. AIS personalization allows you to pick a preferred specialist, or a group of specialists, based on your preferences.






COMPREHENSIVE GENETIC SCREENING TO IDENTIFY POTENTIAL HEALTH CONDITIONS EARLY — FOR DOGS AND CATS

Wisdom Panel™ Premium enables better risk management in pre-anesthetic and preventative care:

- Detect before clinical signs develop — and proactively counsel pet owners on risk, diet, exercise, and breeding
- Detects MDR1 Medication Sensitivity (cat and dog), von Willebrand's Disease Type 1 (dog), Intervertebral Disc Disease Risk (dog), Polycystic Kidney Disease (cat), Hypertrophic Cardiomyopathy (cat), and many more
- Canine panel includes 265+ disorder-associated genetic health variants
- Feline panel includes 45+ health conditions
- Backed by 35+ published papers
- Includes pet owner-friendly breed and trait identifications



Wisdom's veterinary consultants support you in interpreting genetic results so you can build proactive, tailored healthcare plans for your patients

Canine
Wisdom
Panel™
Premium


CS14497
page 109



Feline
Wisdom
Panel™
Complete

CS14498
page 109





KEYSCREEN® GI PARASITE PCR

The new standard in
parasite diagnostics

CT991
page 105

KeyScreen GI Parasite PCR brings the power of PCR to routine parasite screening. Use KeyScreen to find more GI parasites and treat them with greater precision, speed, and confidence — all at an affordable price.

- Screens for 20 intestinal parasites
- Detects benzimidazole treatment' resistance in hookworms
- Determines the zoonotic potential of *Giardia*

START FINDING MORE GI PARASITES, FASTER, TODAY


antechdiagnostics.com/parasites



Better pet health in one test



Detects
benzimidazole-resistant canine
hookworm infections¹⁻⁵



KeyScreen GI Parasite
PCR found parasites in

1 in 4

samples⁴

Reference lab O&P
only found parasites in

1 in 10

samples⁷



**Treatment
resistant
Hookworms**
found in more dog breeds and
across the US and Canada¹⁻⁵



Giardia is NOT
typically zoonotic
less than **4%** are⁴



Parasites
are evolving
so is KeyScreen GI Parasite PCR^{2, 3, 5, 6}

1. Comparative Study of KeyScreen and traditional fecal flotation methods: <https://parasitesandvectors.biomedcentral.com/articles/10.1186/s13071-023-05904-z>
2. Leutenegger CM, et al. Emergence of Ancylostoma caninum parasites with the benzimidazole resistance F167Y polymorphism in the US dog population. Int. J. Parasitol. Drugs Drug Resist. 2023;14:131-140. <https://www.sciencedirect.com/science/article/pii/S2211320723000015?via%3Dihub>
3. Evason, MD, et al. Emergence of canine hookworm treatment resistance: Novel detection of Ancylostoma caninum anthelmintic resistance markers by fecal PCR in 11 dogs from Canada, Am J Vet Res. 2023 July; <https://doi.org/10.2460/ajvr.23.05.0116>
4. Leutenegger CM, et al. Frequency of intestinal parasites in dogs and cats identified by molecular diagnostics. ACVIM, Philadelphia, June 2023.
5. Leutenegger CM, et al. Association of the novel benzimidazole resistance marker Q134H with F167Y in dogs with Ancylostoma caninum. ACVIM, Philadelphia June 2023.
6. Venkatesan A, et al. Molecular evidence of widespread benzimidazole drug resistance in Ancylostoma caninum from domestic dogs throughout the USA and discovery of a novel β -tubulin benzimidazole resistance mutation. PLoS Pathog. Mar 2023;19:e1011146. <https://journals.plos.org/plospathogens/article?id=10.1371/journal.ppat.1011146>
7. KeyScreen Whitepaper <https://www.antechediagnostics.com/keyscreens/#keyscreens-whitepaper>

ACCUPLEX® WITH C6

Screen for canine vector-borne diseases with confidence

CAC100

page 20

Accuplex, now with C6, tests for Lyme disease, heartworm, *Ehrlichia*, and *Anaplasma*.

- Lyme disease detection with Accuplex includes screening for antibodies to C6, and is in alignment with the 2018 ACVIM consensus statement for Lyme borreliosis.*
- The C6 peptide is only expressed following natural infection or exposure of a host. Therefore, a positive C6 result will not be reported from cross-reactivity with vaccinal antibodies.

EXPERIENCE IT NOW AND MOVE YOUR CANINE VECTOR-BORNE DISEASE SCREENING FROM IN-HOUSE TO OUR HOUSE

antechdiagnostics.com/accuplex

*Meryl P Littman, Bernhard Gerber, Richard E Goldstein, Mary Anna Labato, Michael R Lappin, George E Moore. ACVIM consensus update on Lyme borreliosis in dogs and cats. J Vet Intern Med. 2018 May;32(3):887-903. <https://onlinelibrary.wiley.com/doi/full/10.1111/jvim.15085>



WELCOME TO

ANTECH'S
DIRECTORY
OF SERVICES



Quickly find the test you need, when you need it

- See the Table of Contents to explore all that is available at your fingertips.
- Head to our Test Index ([page 157](#)) to find the exact location of the test or profile you need.
- Find our frequently used phone numbers in the Support Services section ([page 14](#)).
- The Quick Reference section ([page 16](#)) contains all the information you need for scheduling pickups and ordering supplies.

Reference Laboratories Terms of Service

Unless otherwise agreed in an agreement signed by an authorized representative of Antech, all reference laboratory testing and consulting services set out in this directory are governed by our standard terms and conditions, available at antechdiagnostics.com/terms-service.

Privacy Policy

As part of Mars Petcare, we are committed to A Better World for Pets™. Antech respects and is committed to protecting your privacy. Please visit mars.com/privacy to see a copy of our privacy statement and how we handle the privacy of your data.

NOTE

The tests and services listed in this directory have been organized in a manner determined to be either appropriate medically or for ease of use by veterinary practitioners in the estimation of our medical operations teams. Any headings or groupings of tests or services are for reference purposes only. Additionally, all turnaround times and prices set forth in this directory are subject to change at our sole discretion. Please contact Customer Service at 1-800-341-3440 if you wish to confirm any information in this directory.

Table of Contents

14	Support Services
15	Account Information
16	Quick Reference
18	Consulting Services
19	Essential Panels
20	Preventative Care
35	CBC Chemistry Profiles
59	Hematology
65	Chemistry
69	Urine
71	Endocrinology
79	Infectious
87	Microbiology
93	Pathology
99	Molecular
111	General
127	Avian • Exotics
135	Equine • Large Animal
157	Test Index

NOTES | FREQUENTLY USED TEST CODES

SUPPORT SERVICES

Customer Service

1-800-341-3440, Dial 0

M – F	5 a.m. – 9 p.m. PST	8 a.m. – Midnight EST
SAT	5 a.m. – 5 p.m. PST	8 a.m. – 8 p.m. EST
SUN	6 a.m. – 3 p.m. PST	9 a.m. – 6 p.m. EST

Dispatch/Specimen Pickup

1-800-341-3440, Dial 1

AUTOMATED 24/7		
M – F	5 a.m. – 9 p.m. PST	8 a.m. – Midnight EST
SAT	5 a.m. – 5 p.m. PST	8 a.m. – 8 p.m. EST
SUN	6 a.m. – 3 p.m. PST	9 a.m. – 6 p.m. EST

Consultation Services

1-800-341-3440, Dial 2

M – F	5 a.m. – 6 p.m. PST	8 a.m. – 9 p.m. EST
SAT	6 a.m. – 3 p.m. PST	9 a.m. – 6 p.m. EST
SUN	CLOSED	

Billing

1-800-341-3440, Dial 3

M – F	6 a.m. – 5 p.m. PST	9 a.m. – 8 p.m. EST
SAT – SUN	CLOSED	

Equine Support

1-800-341-3440, Dial 4

M – F	5 a.m. – 9 p.m. PST	8 a.m. – Midnight EST
SAT	5 a.m. – 5 p.m. PST	8 a.m. – 8 p.m. EST
SUN	6 a.m. – 3 p.m. PST	9 a.m. – 6 p.m. EST

Antech Imaging Services

1-800-341-3440, Dial 5

AUTOMATED 24/7		
-----------------------	--	--

Website Support

1-800-341-3440, Dial 6

support@antechdiagnostics.com

M – F	6 a.m. – 4 p.m. PST	9 a.m. – 7 p.m. EST
SAT	6 a.m. – Noon PST	9 a.m. – 3 p.m. EST
SUN	CLOSED	

Sound

1-800-819-5538

M – F	8 a.m. – 5 p.m. PST	11 a.m. – 8 p.m. EST
SAT – SUN	CLOSED	

ACCOUNT INFORMATION

ANTECH ACCOUNT NUMBER

HEALTHTRACKS NUMBER

USERNAME

PASSWORD

HealthTracks®
online.antechdiagnostics.com

HealthTracks provides you access to your diagnostic needs in one innovative digital platform. This digital center simplifies your practice’s diagnostic workflow, ensuring you have the features you need to provide the best care for your patients.

Elevate your diagnostic experience and explore all the HealthTracks features at [antechdiagnostics.com/healthtracks](https://online.antechdiagnostics.com/healthtracks)

TERRITORY SALES MANAGER

E-MAIL

MOBILE

PROFESSIONAL SERVICES VETERINARIAN

E-MAIL

MOBILE

SOUND EQUIPMENT SPECIALIST

E-MAIL

MOBILE

QUICK REFERENCE

Antech's dedicated Customer Service Team is here to help at **1-800-341-3440**

Completion Codes — Submitting a Sample for Profile Completion

Missing one of the sample types? No problem. If you do not have a necessary sample when submitting a profile, you may submit the omitted sample within 30 days of the original profile submission to complete the profile at no additional charge.

1. Create a new Test Requisition Form using the applicable completion code from the list below
2. On the Test Requisition Form, reference the original profile's 12-character accession number in the notes field
3. Print the Test Requisition Form and submit the sample as normal

Completion codes submitted after 30 days, not related to a prior profile submission, or without reference to the original accession number will be invoiced for the applicable test code.

- CT991NC — KeyScreen® GI Parasite PCR Completion
- CT805NC — Ova and Parasite Completion
- CT808NC — Ova and Parasite with Giardia Completion
- 40510 CUADELAY — Urinalysis Completion

Follow-up and Monitoring Tests

We offer you reduced-charge testing for serial patient monitoring and rechecks. Within 30 days of the original submission, order one of the codes below and include the original profile's 12-character accession number on the Test Requisition Form.

- CRECHECK — Superchem with SDMA, CBC
- CT330R — Complete Blood Count Recheck
- CT760R — Urinalysis Recheck
- CM130R — Urine Culture Recheck

Please note that recheck codes submitted after 30 days, not related to a prior profile submission, or without reference to the original accession number, will be charged the applicable non-recheck test fee.

Add-on Testing

Add-on codes beginning with CADD are for additional diagnostics only and cannot be used at the time of original submission.

Additional add-on guidelines:

- Add-on codes can only be used for samples with existing results
- Add-on codes can only be added to tests already accessioned
- New samples sent to the lab cannot be added to previously resulted accessions
- Add-on codes do not apply to samples with a hold request
- Samples over seven days old cannot use add-on codes
- You cannot upgrade completed profiles with add-on codes
- The sample must be in the lab for the add-on code to be applied

Test Cancellation

Testing that has not begun can be canceled at no charge by calling Customer Service. Testing that is already in process cannot be canceled.

*Histology Note: all samples received for histology will be processed and evaluated. Histology tests may be canceled only prior to sample processing, and all samples in the accession must be canceled.

See the full reference lab terms of service at antechdiagnostics.com/terms-service-ca

Specimen Retention Times

CLINICAL SUBMISSIONS	
Whole Blood	7 days
Serum/Plasma	7 days
Urine	7 days
EIA (Coggins): AGID and ELISA	30 days
Hematology – Slides (includes Wright’s stained urine smears)	7 days
Hematology – Slides (rechecks and pathology review)	1 year
Microbiology – all specimens	7 days
Microbiology – Slides (gram stain smears)	14 days
PATHOLOGY (CYTOLOGY) SUBMISSIONS	
Biological Specimens	7 days
Slides	90 days
PATHOLOGY (HISTOLOGY) SUBMISSIONS	
Biological Specimens (wet tissue in formalin)	30 days
Blocks	18 months
Slides	180 days

Supplies

Multiple ways to order supplies:

- Online — HealthTracks (preferred): online.antechdiagnostics.com
- Phone — **1-800-341-3440**
- Email — canada.supplies@antechmail.com

Supplies are typically delivered within 2 – 5 days of ordering.

Test Requisition Forms

Antech offers both electronic and paper test requisition forms (TRFs). Electronic TRFs can be generated directly from your practice management system, or through HealthTracks. Electronic TRFs are the most efficient and convenient way to order diagnostic testing. They have many benefits, including increased data and order integrity and the ability to archive requisitions. Manual TRFs are also available online at antechdiagnostics.com and can be ordered through HealthTracks.

CONSULTING SERVICES

Live, expert, and complimentary
peer-to-peer medical consultations

Looking for a trusted expert to help you interpret diagnostic results or develop a diagnostic and treatment plan? Antech provides current customers with comprehensive consulting services to assist with diagnostic interpretation, clinical application of results, and general case support. With a world-class team of 60+ board-certified specialists across 12 specialties, Antech is committed to partnerships you can count on.

Request A Consult

Please have your Antech account number and the accession number(s) available for the results you wish to discuss. Immediately connect with the first available Internal Medicine specialist consultant or request a call-back.

1-800-341-3440, Dial 2
online.antechdiagnostics.com



Internal Medicine specialists are
available Monday through Saturday

M – F	5 a.m. – 6 p.m. PST	8 a.m. – 9 p.m. EST
SAT	6 a.m. – 3 p.m. PST	9 a.m. – 6 p.m. EST

During business days, all other specialists are
available by call back within 24 – 48 hours



Avian and
Exotics



Behavior



Cardiology



Dentistry



Dermatology



Equine and
Large Animal



Neurology



Nutrition



Oncology



Reproduction

ESSENTIAL PANELS

	SUPERCHEM WITH SDMA CSA010	VET SCREEN WITH SDMA CSA025	ADULT WELLNESS CHEMISTRY WITH SDMA CSA665	PRE-OP SCREEN WITH SDMA CSA040
A/G Ratio	●	●		●
Albumin	●	●		●
Alkaline Phosphatase	●	●		●
ALT (SGPT)	●	●		●
Amylase	●			
AST (SGOT)	●	●		
BUN	●	●		●
BUN/Creatinine Ratio	●	●		●
Calcium	●	●		
Chloride	●	●		●
Cholesterol	●	●		
CPK	●	●		
Creatinine	●	●		●
GGT	●			
Globulin	●	●		●
Glucose	●	●		●
Magnesium	●			
NA/K Ratio	●	●		●
Pancreatic Sensitive Lipase (PSL)	●			
Phosphorus	●	●		
Potassium	●	●		●
SDMA	●	●		●
Sodium	●	●		●
Total Bilirubin	●	●		
Total Protein	●	●		●
Triglycerides	●			

PREVENTATIVE CARE

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CAC100	<p>Accuplex®</p> <p>This is a canine vector-borne disease screening for heartworm, Lyme disease (includes screening for antibodies against two C6 peptides), <i>E. canis</i>, and <i>A. phagocytophilum</i>.</p>	0.5 mL serum in red top or serum separator tube	1-2 days
CKAC100	<p>Accuplex®, KeyScreen®</p> <p>Accuplex® is a canine vector-borne disease screening for Heartworm (Ag), Lyme disease (includes screening for antibodies against two C6 peptides), <i>E. canis</i>, and <i>A. phagocytophilum</i>. KeyScreen® is a fecal PCR test that detects 20 individual parasite species. If detected, <i>Toxocara</i> spp. will be further speciated (<i>T. cati</i>, <i>T. canis</i>, and <i>T. leonina</i>), <i>Ancylostoma caninum</i> will be evaluated for benzimidazole resistance, and <i>Giardia</i> evaluated for subtypes A and B.</p> <p>Accuplex® (CAC100) can be reflexed for an additional charge.</p>	<p>0.5 mL serum and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, Antech provided fecal container</p>	1-3 days
CKA535	<p>Adult Chem with Lytes, SDMA, CBC, Accuplex®, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, Accuplex®, and KeyScreen®</p> <p>Wellness chemistry with electrolytes, SDMA (see CT1035), complete blood count, Accuplex® for vector-borne disease screening (CAC100), and KeyScreen® for intestinal parasite detection (see CT991).</p> <p>Accuplex® (CAC100) can be reflexed for an additional charge.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	1-3 days
CKA583	<p>Adult Chem with Lytes, SDMA, CBC, Accuplex®, T4, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Accuplex®, and KeyScreen®</p> <p>Wellness chemistry with electrolytes, SDMA (CT1035), complete blood count, total T4, Accuplex® for vector-borne disease screening (see CAC100), and KeyScreen® for intestinal parasite detection (see CT991).</p> <p>Accuplex® (CAC100) can be reflexed for an additional charge.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	1-3 days
CKS519	<p>Adult Chem with Lytes, SDMA, CBC, Feline Heartworm, FIV, FeLV, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, FeLV, FIV, Heartworm Antibody, and KeyScreen®</p> <p>Wellness chemistry with electrolytes, SDMA (CT1035), complete blood count, FeLV antigen detection, FIV and heartworm antibody detection, and KeyScreen® for intestinal parasite detection (see CT991).</p>	<p>1.25 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	1-3 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CKS575	Adult Chem with Lytes, SDMA, CBC, FIV, FeLV, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, FeLV, FIV, and KeyScreen® Wellness chemistry with electrolytes, SDMA, complete blood count, FeLV antigen detection, FIV antibody detection, and KeyScreen® for intestinal parasite detection (see CT991).	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS534	Adult Chem with Lytes, SDMA, CBC, Heartworm, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, Heartworm Antigen detection, and KeyScreen® Wellness chemistry with electrolytes, SDMA, complete blood count, heartworm antigen detection, and KeyScreen® for intestinal parasite detection (see CT991).	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS586	Adult Chem with Lytes, SDMA, CBC, Heartworm, T4, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, Heartworm Antigen, Total T4, and KeyScreen® Wellness chemistry with electrolytes, SDMA, complete blood count, total T4, heartworm antigen detection, and KeyScreen® for intestinal parasite detection (see CT991).	1.0 mL serum and 0.3 grams feces Serum in red top or serum separator tube, Antech provided fecal container	1-3 days
CKS535	Adult Chem with Lytes, SDMA, CBC, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, and KeyScreen® Wellness chemistry with electrolytes, SDMA, complete blood count, and KeyScreen® for intestinal parasite detection (see CT991).	0.5 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS591	Adult Chem with Lytes, SDMA, CBC, T4, Feline Heartworm, FIV, FeLV, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, T4, FeLV, FIV, Heartworm Antibody, and KeyScreen® Wellness chemistry with electrolytes, SDMA, complete blood count, total T4, FeLV antigen detection, FIV and heartworm antibody detection, and KeyScreen® for intestinal parasite detection (see CT991).	1.25 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CKS597	Adult Chem with Lytes, SDMA, CBC, T4, FIV, FeLV, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, T4, FeLV, FIV, and KeyScreen® Wellness chemistry with electrolytes, SDMA, complete blood count, total T4, FeLV antigen detection, FIV antibody detection, and KeyScreen® for intestinal parasite detection (see CT991).	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS583	Adult Chem with Lytes, SDMA, CBC, T4, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, T4, and KeyScreen® Wellness chemistry with electrolytes, SDMA, complete blood count, total T4, and KeyScreen® for intestinal parasite detection (see CT991).	0.5 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS589	Adult Chem with Lytes, SDMA, CBC, T4, UA, Feline Heartworm, FIV, FeLV, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Urinalysis, FeLV, FIV and Heartworm Antibody detection, and KeyScreen® A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), total T4, FeLV antigen detection, FIV and heartworm antibody detection, , KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035).	1.25 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
CKS619	Adult Chem with Lytes, SDMA, CBC, T4, UA, FIV, FeLV, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Urinalysis, FeLV, FIV, and KeyScreen® A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), total T4, FeLV antigen detection, FIV antibody detection, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035).	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CKA673	Adult Chem with Lytes, SDMA, CBC, UA, Accuplex®, T4, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Urinalysis, Accuplex®, and KeyScreen® A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), total T4, Accuplex® for vector-borne disease screening (see CAC100), KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035). Accuplex® (CAC100) can be reflexed for an additional charge.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
CKS623	Adult Chem with Lytes, SDMA, CBC, UA, Feline Heartworm, FIV, FeLV, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, Urinalysis, FeLV, FIV, Heartworm Antibody, and KeyScreen® A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), FeLV antigen detection, FIV and heartworm antibody detection, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035).	1.25 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
CKS634	Adult Chem with Lytes, SDMA, CBC, UA, FIV, FeLV, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, Urinalysis, FeLV, FIV, and KeyScreen® A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), FeLV antigen detection, FIV antibody detection, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035).	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
CKS592	Adult Chem with Lytes, SDMA, CBC, UA, Heartworm, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, Urinalysis, Heartworm Antigen, KeyScreen® A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), heartworm antigen detection, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis and lipemia. Lipemia can falsely decrease T4 results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CKS587	<p>Adult Chem with Lytes, SDMA, CBC, UA, Heartworm, T4, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Urinalysis, Heartworm Antigen test, and KeyScreen®</p> <p>A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), total T4, heartworm antigen detection, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035).</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days
CKS053	<p>Adult Chem with Lytes, SDMA, CBC, UA, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, Urinalysis, and KeyScreen®</p> <p>A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035).</p> <p>Interferences: marked hemolysis and lipemia.</p>	<p>0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days
CKA600	<p>Adult Chem with SDMA, CBC, Accuplex®, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, Accuplex®, and KeyScreen®</p> <p>Wellness chemistry, complete blood count, Accuplex® for canine vector-borne disease screening (see CAC100), KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035).</p> <p>Interferences: marked hemolysis and lipemia.</p> <p>Accuplex® (CAC100) can be reflexed for an additional charge.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	1-3 days
CKA670	<p>Adult Chem with SDMA, CBC, Accuplex®, T4, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Accuplex®, and KeyScreen®</p> <p>Wellness chemistry, complete blood count, total T4, Accuplex® for canine vector-borne disease screening (see CAC100), KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035).</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p> <p>Accuplex® (CAC100) can be reflexed for an additional charge.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, Antech provided fecal container</p>	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CKS622	Adult Chem with SDMA, CBC, Feline Heartworm, FIV, FeLV, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, FeLV, FIV, Heartworm Antibody, and KeyScreen® Wellness chemistry, complete blood count, FeLV antigen detection, FIV and heartworm antibody detection, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.25 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS675	Adult Chem with SDMA, CBC, FIV, FeLV, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, FeLV, FIV, and KeyScreen® Wellness chemistry, complete blood count, FeLV antigen detection, FIV antibody detection, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS605	Adult Chem with SDMA, CBC, Heartworm, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, Heartworm Antigen, and KeyScreen® Wellness chemistry, complete blood count, heartworm antigen detection, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS685	Adult Chem with SDMA, CBC, Heartworm, T4, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Heartworm Antigen, and KeyScreen® Wellness chemistry, complete blood count, total T4, heartworm antigen detection, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS600	Adult Chem with SDMA, CBC, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, and KeyScreen® Wellness chemistry, complete blood count, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CKS590	Adult Chem with SDMA, CBC, T4, Feline Heartworm, FIV, FeLV, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, T4, FeLV, FIV, Heartworm Antibody, and KeyScreen® Wellness chemistry, complete blood count, total T4, FeLV antigen detection, FIV and heartworm antibody detection, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.25 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS697	Adult chem with SDMA, CBC, T4, FIV, FeLV, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, T4, FeLV, FIV, and KeyScreen® Wellness chemistry, complete blood count, total T4, FeLV antigen detection, FIV antibody detection, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS670	Adult Chem with SDMA, CBC, T4, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, T4, and KeyScreen® Wellness chemistry, complete blood count, total T4, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS588	Adult Chem with SDMA, CBC, T4, UA, Feline Heartworm, FIV, FeLV, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Urinalysis, FeLV, FIV, Heartworm Antibody, and KeyScreen® A complete minimum database (chemistry, complete blood count, urinalysis), total T4, FeLV antigen detection, FIV and heartworm antibody detection, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.25 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CKS621	<p>Adult Chem with SDMA, CBC, T4, UA, FIV, FeLV, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Urinalysis, FeLV, FIV, and KeyScreen®</p> <p>A complete minimum database (chemistry, complete blood count, urinalysis), total T4, FeLV antigen detection, FIV antibody detection, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035).</p> <p>Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days
CKS672	<p>Adult Chem with SDMA, CBC, T4, UA, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Urinalysis, and KeyScreen®</p> <p>A complete minimum database (chemistry, complete blood count, urinalysis), total T4, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035).</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p>	<p>0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days
CKA607	<p>Adult Chem with SDMA, CBC, UA, Accuplex®, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, Urinalysis, Accuplex®, and KeyScreen®</p> <p>A complete minimum database (chemistry, complete blood count, urinalysis), Accuplex® for canine vector borne disease screening (CAC100), KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035).</p> <p>Interferences: marked hemolysis and lipemia.</p> <p>Accuplex® (CAC100) can be reflexed for an additional charge.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days
CKA672	<p>Adult Chem with SDMA, CBC, UA, Accuplex®, T4, KeyScreen®</p> <p>Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Urinalysis, Accuplex®, and KeyScreen®</p> <p>A complete minimum database (chemistry, complete blood count, urinalysis), total T4, Accuplex® for canine vector-borne disease screening (see CAC100), KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035).</p> <p>Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.</p> <p>Accuplex® (CAC100) can be reflexed for an additional charge.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CKS624	Adult Chem with SDMA, CBC, UA, Feline Heartworm, FIV, FeLV, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, Urinalysis, FeLV, FIV, Heartworm Antibody, and KeyScreen® A complete minimum database (chemistry, complete blood count, urinalysis), FeLV antigen detection, FIV and heartworm antibody detection, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.25 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
CKS631	Adult Chem with SDMA, CBC, UA, FIV, FeLV, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, Urinalysis, FeLV, FIV, and KeyScreen® A complete minimum database (chemistry, complete blood count, urinalysis), FeLV antigen detection, FIV antibody detection, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
CKS625	Adult Chem with SDMA, CBC, UA, Heartworm, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, Urinalysis, Heartworm Antibody, and KeyScreen® A complete minimum database (chemistry, complete blood count, urinalysis), heartworm antigen detection, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
CKS687	Adult Chem with SDMA, CBC, UA, Heartworm, T4, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Urinalysis, Heartworm Antibody, and KeyScreen® A complete minimum database (chemistry, complete blood count, urinalysis), total T4, heartworm antigen detection, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035). CKS625 to which T4 has been added. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CKS607	Adult Chem with SDMA, CBC, UA, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, Urinalysis, and KeyScreen® A complete minimum database (chemistry, complete blood count, urinalysis), KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
CKS673	Adult Chem with SDMA, Lytes, CBC, T4, UA, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, T4, Urinalysis, and KeyScreen® A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), total T4, KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035).	0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
CKA053	Adult Chem with SDMA, Lytes, CBC, UA, Accuplex®, KeyScreen® Wellness Chemistry with Electrolytes and SDMA, CBC, Urinalysis, Accuplex®, and KeyScreen® A complete minimum database (chemistry with electrolytes, complete blood count, urinalysis), Accuplex® for vector-borne disease screening (see CAC100), KeyScreen® for intestinal parasite detection (see CT991), and SDMA for glomerular filtration rate estimation (see CT1035). Accuplex® (CAC100) can be reflexed for an additional charge.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
CKA705	Senior Profile 1 with SDMA, Accuplex®, KeyScreen® Superchem with SDMA, CBC, T4, Urinalysis, Accuplex®, and KeyScreen® The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), a total T4, Accuplex® for canine vector-borne disease screening (see CAC100), and KeyScreen® for intestinal parasite detection (see CT991). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results. Accuplex® (CAC100) can be reflexed for an additional charge.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CKS710	Senior Profile 1 with SDMA, Heartworm, KeyScreen® Superchem with SDMA, CBC, T4, Urinalysis, Heartworm Antigen, and KeyScreen® The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), total T4, heart worm antigen detection, and KeyScreen® for intestinal parasite detection (see CT991). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
CKS705	Senior Profile 1 with SDMA, KeyScreen® Superchem with SDMA, CBC, T4, Urinalysis, and KeyScreen® The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), a total T4, and KeyScreen® for intestinal parasite detection (see CT991). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
CKA020	Superchem with SDMA, CBC, Accuplex®, KeyScreen® Most comprehensive chemistry profile (superchemistry with SDMA (CSA010)), a complete blood count, Accuplex® for canine vector-borne disease screening (see CAC100), and KeyScreen® for intestinal parasite detection (see CT991). Interferences: marked hemolysis and lipemia. Accuplex® (CAC100) can be reflexed for an additional charge.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKA114	Superchem with SDMA, CBC, Accuplex®, T4, KeyScreen® Most comprehensive chemistry profile (superchemistry with SDMA (CSA010)), a complete blood count, total T4, and KeyScreen® for intestinal parasite detection (see CT991). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results. Accuplex® (CAC100) can be reflexed for an additional charge.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CKS737	Superchem with SDMA, CBC, Feline Heartworm, FIV, FeLV, KeyScreen® Superchem with SDMA, CBC, FeLV, FIV, Heartworm Antibody, and KeyScreen® Most comprehensive chemistry profile (superchemistry with SDMA (CSA010)), a complete blood count, FeLV antigen detection, FIV and heartworm antibody detection, and KeyScreen® for intestinal parasite detection (see CT991). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.25 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS490	Superchem with SDMA, CBC, FIV, FeLV, KeyScreen® Superchem with SDMA, CBC, FeLV, FIV, and KeyScreen® Most comprehensive chemistry profile (superchemistry with SDMA (CSA010)), a complete blood count, FeLV antigen and FIV antibody detection, and KeyScreen® for intestinal parasite detection (see CT991). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS028	Superchem with SDMA, CBC, Heartworm, KeyScreen® Most comprehensive chemistry profile (superchemistry with SDMA (CSA010)), a complete blood count, heartworm antigen detection, and KeyScreen® for intestinal parasite detection (see CT991). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS110	Superchem with SDMA, CBC, Heartworm, T4, KeyScreen® Most comprehensive chemistry profile (superchemistry with SDMA (CSA010)), a complete blood count, total T4, heartworm antigen detection, and KeyScreen® for intestinal parasite detection (see CT991). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS020	Superchem with SDMA, CBC, KeyScreen® Most comprehensive chemistry profile (superchemistry with SDMA (CSA010)), a complete blood count, and KeyScreen® for intestinal parasite detection (see CT991). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CKS724	Superchem with SDMA, CBC, T4, Feline Heartworm, FIV, FeLV, KeyScreen® Superchem with SDMA, CBC, T4, FeLV, FIV, Heartworm Antibody (Feline) and KeyScreen® Most comprehensive chemistry profile (superchemistry with SDMA (CSA010)), a complete blood count, total T4, FeLV antigen detection, FIV and heartworm antibody detection, and KeyScreen® for intestinal parasite detection (see CT991). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.25 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS220	Superchem with SDMA, CBC, T4, FIV, FeLV, KeyScreen® Most comprehensive chemistry profile (superchemistry with SDMA (CSA010)), a complete blood count, total T4, FeLV antigen and FIV antibody detection, and KeyScreen® for intestinal parasite detection (see CT991). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS120	Superchem with SDMA, CBC, T4, KeyScreen® Most comprehensive chemistry profile (superchemistry with SDMA (CSA010)), a complete blood count, total T4, and KeyScreen® for intestinal parasite detection (see CT991). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-3 days
CKS715	Superchem with SDMA, CBC, T4, UA, Feline Heartworm, FIV, FeLV, KeyScreen® Superchem with SDMA, CBC, T4, Urinalysis, FeLV, FIV, Heartworm Antibody (Feline), and KeyScreen® The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), total T4, FeLV antigen detection, FIV and heartworm antibody detection, and KeyScreen® for intestinal parasite detection (see CT991). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.25 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
CKS700	Superchem with SDMA, CBC, T4, UA, FIV, FeLV, KeyScreen® Superchem with SDMA, CBC, T4, Urinalysis, FeLV, FIV, and KeyScreen® The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), total T4, FeLV antigen and FIV antibody detection, and KeyScreen® for intestinal parasite detection (see CT991). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CKA021	<p>Superchem with SDMA, CBC, UA, Accuplex®, KeyScreen®</p> <p>Superchem with SDMA, CBC, Urinalysis, Accuplex®, and KeyScreen®</p> <p>The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), and KeyScreen® for intestinal parasite detection (see CT991). Accuplex® for canine vector-borne disease screening (see CAC100) can be reflexed on request.</p> <p>Interferences: marked hemolysis and lipemia.</p> <p>Accuplex® (CAC100) can be reflexed for an additional charge.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days
CKS494	<p>Superchem with SDMA, CBC, UA, Feline Heartworm, FIV, FeLV, KeyScreen®</p> <p>Superchem with SDMA, CBC, Urinalysis, FeLV, FIV, Heartworm Antibody, and KeyScreen®</p> <p>The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), FeLV antigen detection, FIV and heartworm antibody detection, and KeyScreen® for intestinal parasite detection (see CT991).</p> <p>Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.</p>	<p>1.25 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days
CKS702	<p>Superchem with SDMA, CBC, UA, FIV, FeLV, KeyScreen®</p> <p>Superchem with SDMA, CBC, Urinalysis, FeLV, FIV, and KeyScreen®</p> <p>The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), FeLV antigen and FIV antibody detection, and KeyScreen® for intestinal parasite detection (see CT991).</p> <p>Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.</p>	<p>1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces</p> <p>Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container</p>	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CKS112	Superchem with SDMA, CBC, UA, Heartworm, KeyScreen® Superchem with SDMA, CBC, Urinalysis, Heartworm Antigen, and KeyScreen® The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), heartworm antigen detection, and KeyScreen® for intestinal parasite detection (see CT991). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
CKS021	Superchem with SDMA, CBC, UA, KeyScreen® Superchem with SDMA, CBC, Urinalysis, and KeyScreen® The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis) and KeyScreen® for intestinal parasite detection (see CT991). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 0.3 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-3 days
CKS260	FelV, FIV, KeyScreen® FelV antigen detection, FIV antibody detection, and KeyScreen® for intestinal parasite detection (see CT991). Interferences: marked hemolysis or lipemia may result in false positive results on the FelV Antigen ELISA.	0.5 mL serum and 0.3 grams feces Serum in red top or serum separator tube, Antech provided fecal container	1-3 days
CKT615	Heartworm, KeyScreen® Heartworm Antigen and KeyScreen® Heartworm antigen detection and KeyScreen® for intestinal parasite detection (see CT991). Interferences: marked hemolysis and lipemia.	0.5 mL serum and 0.3 grams feces Serum in red top or serum separator tube, Antech provided fecal container	1-3 days

CBC CHEMISTRY PROFILES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA600	Adult Wellness with SDMA (No Heartworm) Wellness Chemistry with Electrolytes and SDMA, CBC A comprehensive chemistry panel, complete blood count, and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	Daily
CAC655	Adult Wellness with SDMA, Fecal Combo with Accuplex® Wellness Chemistry with Electrolytes and SDMA, CBC, Fecal O&P with Centrifugation, <i>Giardia</i> , Accuplex® A comprehensive chemistry, complete blood count, SDMA for glomerular filtration rate estimation (see CT1035), and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (CT808). Accuplex® (CAC100). Interferences: marked hemolysis and lipemia. Accuplex® (CAC100) can be reflexed for an additional charge.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days
CSA615	Adult Wellness with SDMA, Fecal O&P Wellness Chemistry with Electrolytes and SDMA, CBC, Fecal O&P with Centrifugation A comprehensive chemistry panel, complete blood count, SDMA for glomerular filtration rate estimation (see CT1035), and fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection (CT805). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 5 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days
CSA655	Adult Wellness with SDMA, Fecal Combo Wellness Chemistry with Electrolytes and SDMA, CBC, Fecal O&P with Centrifugation, <i>Giardia</i> A comprehensive chemistry, complete blood count, SDMA for glomerular filtration rate estimation (see CT1035), and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (CT808). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA682	Adult Wellness with SDMA, FeLV, FIV, O&P Wellness Chemistry with Electrolytes and SDMA, CBC, FeLV, FIV, O&P with Centrifugation A comprehensive chemistry, complete blood count, FeLV antigen detection, FIV antibody detection, SDMA for glomerular filtration rate estimation (see CT1035), and fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection (CT805). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, and 5 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days
CSA684	Adult Wellness with SDMA, FeLV, FIV, O&P, Giardia Wellness Chemistry with Electrolytes and SDMA, CBC, FeLV, FIV, O&P with Centrifugation, <i>Giardia</i> A comprehensive chemistry, complete blood count, FeLV antigen detection, FIV antibody detection, SDMA for glomerular filtration rate estimation (see CT1035), and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (CT808). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days
CSA651	Adult Wellness with SDMA, Heartworm, UA, O&P, Giardia Wellness Chemistry with Electrolytes and SDMA, CBC, Heartworm Antigen, Fecal O&P with Centrifugation, <i>Giardia</i> , Urinalysis A comprehensive minimum database (wellness chemistry, complete blood count, urinalysis), heartworm antigen detection, SDMA for glomerular filtration rate estimation (see CT1035), and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (CT808). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
CSA647	Adult Wellness with SDMA, O&P, UA Wellness Chemistry with Electrolytes and SDMA, CBC, Fecal O&P with Centrifugation, Urinalysis A comprehensive minimum database (wellness chemistry, complete blood count, urinalysis), SDMA for glomerular filtration rate estimation (see CT1035), and fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection (CT805). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
CSA670	Adult Wellness with SDMA, T4 Wellness Chemistry with Electrolytes and SDMA, CBC, Total T4 A comprehensive chemistry, complete blood count, total T4, and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	Daily Performed every shift

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA672	Adult Wellness with SDMA, T4, UA Wellness Chemistry with Electrolytes and SDMA, CBC, Total T4, Urinalysis A comprehensive minimum database (wellness chemistry, complete blood count, urinalysis), total T4, and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	Daily Performed every shift
CAC607	Adult Wellness with SDMA, UA, Accuplex® Wellness Chemistry with Electrolytes and SDMA, CBC, Urinalysis, Accuplex® A comprehensive minimum database (wellness chemistry, complete blood count, and urinalysis), SDMA for glomerular filtration rate estimation (see CT1035), and Accuplex® for canine vector borne disease screening (CAC100). Interferences: marked hemolysis and lipemia. Accuplex® (CAC100) can be reflexed for an additional charge.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
CSA649	Adult Wellness with SDMA, UA, O&P, Giardia Wellness Chemistry with Electrolytes and SDMA, CBC, O&P with Centrifugation, <i>Giardia</i> , Urinalysis A comprehensive minimum database (wellness chemistry, complete blood count, urinalysis), SDMA for glomerular filtration rate estimation (see CT1035), and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (CT808). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, urine, and 6 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
CAC651	Adult Wellness with SDMA, UA, O&P, Giardia, Accuplex® Wellness Chemistry with Electrolytes and SDMA, CBC, Fecal O&P with Centrifugation, <i>Giardia</i> , Urinalysis, Accuplex® A comprehensive minimum database (wellness chemistry, complete blood count, urinalysis), Accuplex® for canine vector borne disease screening (CAC100), SDMA for glomerular filtration rate estimation (see CT1035), and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (CT808). Interferences: marked hemolysis and lipemia. Accuplex® (CAC100) can be reflexed for an additional charge.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 6 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
CSA607	Adult Wellness with SDMA, Urinalysis Wellness Chemistry with Electrolytes and SDMA, CBC, Urinalysis A comprehensive minimum database (wellness chemistry, complete blood count, urinalysis) and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	Daily Performed every shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA710	Basic Wellness Screen with SDMA Superchem with SDMA, CBC, Total T4, Urinalysis, and Heartworm Antigen The most comprehensive minimum database, total T4, heartworm antigen detection, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
CSA625	Canine Adult Wellness with SDMA, Urinalysis Wellness Chemistry with Electrolytes and SDMA, CBC, Heartworm Antigen, Urinalysis A comprehensive minimum database (wellness chemistry, complete blood count, urinalysis), heartworm antigen detection, and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
CSA650	Canine Adult Wellness with SDMA, Heartworm, Fecal Combo Wellness Chemistry with Electrolytes and SDMA, CBC, Heartworm Antigen, Fecal O&P with Centrifugation, <i>Giardia</i> A comprehensive chemistry, complete blood count, heartworm antigen detection, SDMA for glomerular filtration rate estimation (see CT1035), and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (CT808). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 6 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days
CSA620	Canine Adult Wellness with SDMA, Heartworm, Ehrlichia & Lyme IgG Wellness Chemistry with Electrolytes and SDMA, CBC, <i>Ehrlichia canis</i> , Heartworm Antigen, Lyme Titer IgG A comprehensive chemistry panel, a complete blood count, SDMA for glomerular filtration rate estimation (see CT1035) and evaluation for vector borne disease exposure (<i>Ehrlichia canis</i> , <i>Dirofilaria immitis</i> , and <i>Borrelia burgdorferi</i> antibody detection). Interferences: marked hemolysis and lipemia.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-4 days
CSA635	Canine Adult Wellness with SDMA, Heartworm, Ova & Parasite Wellness Chemistry with Electrolytes and SDMA, CBC, Heartworm Antigen, Fecal O&P with Centrifugation A comprehensive chemistry, complete blood count, heartworm antigen detection, SDMA for glomerular filtration rate estimation (see CT1035), and fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection (CT805). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 5 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA685	Canine Adult Wellness with SDMA, T4 Wellness Chemistry with Electrolytes and SDMA, CBC, Heartworm Antigen, Total T4 A comprehensive chemistry, complete blood count, total T4, heartworm antigen detection, and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
CSA680	Canine Adult Wellness with SDMA, UA, O&P Wellness Chemistry with Electrolytes and SDMA, CBC, Heartworm Antigen, Fecal O&P with Centrifugation, Urinalysis A comprehensive minimum database (wellness chemistry, complete blood count, urinalysis), heartworm antigen detection, SDMA for glomerular filtration rate estimation (see CT1035), and fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection (CT805). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
CSA645	Canine Adult Wellness with SDMA, UA, O&P, Vaccine Titer Wellness Chemistry with Electrolytes and SDMA, CBC, O&P with Centrifugation, Distemper/Parvo Vaccine Titer, Urinalysis A comprehensive minimum database (wellness chemistry, complete blood count and urinalysis), SDMA for glomerular filtration rate estimation (see CT1035), fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection (CT805) and Distemper/ Parvovirus vaccine titer assessment. Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-4 days
CSA765	Canine Senior Profile with SDMA, Fecal Combo Superchem with SDMA, CBC, Total T4, Urinalysis, Heartworm Antigen, Fecal O&P with Centrifugation, <i>Giardia</i> The most comprehensive minimum database (superchemistry, complete blood count, and urinalysis), total T4, heartworm antigen detection, SDMA for glomerular filtration rate estimation (see CT1035), and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (CT808). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 6 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA605	Canine Wellness Profile with SDMA Wellness Chemistry with Electrolytes and SDMA, CBC, Heartworm Antigen A comprehensive chemistry panel, a complete blood count, heartworm antigen detection, and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis and lipemia.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
CSA220	Cat Scan Plus with SDMA Superchem with SDMA, CBC, Total T4, FeLV, and FIV The most comprehensive chemistry, complete blood count, total T4, FeLV antigen detection, FIV antibody detection and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
CSA230	Cat Scan with SDMA Superchem with SDMA, CBC, T4, and FeLV The most comprehensive chemistry, complete blood count, total T4, FeLV antigen detection, and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
CSA800	Diabetes Monitoring Panel with SDMA CBC, Urinalysis, Fructosamine Assay, Chemistry panel: Total Protein, ALT (SGPT), Alk Phos, Total Bilirubin, BUN, Creatinine, Glucose, Sodium, Potassium, Chloride, PSL, and SDMA A minimum database (chemistry, complete blood count, urinalysis), SDMA for glomerular filtration rate estimation (see CT1035), and fructosamine. Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
CSA610	Feline Adult Wellness Profile with SDMA Wellness Chemistry with Electrolytes and SDMA, CBC, Feline Heartworm Antibody A comprehensive chemistry panel, a complete blood count, heartworm antibody detection, and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis and lipemia.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA630	Feline Adult Wellness with SDMA, Urinalysis Wellness Chemistry with Electrolytes and SDMA, CBC, Feline Heartworm Antibody, Urinalysis A comprehensive minimum database (wellness chemistry, complete blood count, urinalysis), heartworm antibody detection, and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
CSA622	Feline Adult Wellness with SDMA, FeLV, FIV Wellness Chemistry with Electrolytes and SDMA, CBC, Feline Heartworm Antibody, FeLV, FIV A comprehensive chemistry panel, a complete blood count, FeLV antigen detection, FIV and heartworm antibody detection, and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
CSA660	Feline Adult Wellness with SDMA, Heartworm, Fecal Combo Wellness Chemistry with Electrolytes and SDMA, CBC, Feline Heartworm Antibody, O&P with Centrifugation, <i>Giardia</i> A comprehensive chemistry, complete blood count, heartworm antibody detection, SDMA for glomerular filtration rate estimation (see CT1035) and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (CT808). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-4 days
CSA640	Feline Adult Wellness with SDMA, UA, O&P, Panleukopenia Wellness Chemistry with Electrolytes and SDMA, CBC, O&P with Centrifugation, Panleukopenia Vaccine Titer (Feline Only), Urinalysis A comprehensive minimum database (wellness chemistry, complete blood count, and urinalysis), SDMA for glomerular filtration rate estimation (see CT1035), feline panleukopenia vaccinal titer, and fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection (CT805). Interferences: marked hemolysis and lipemia.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-4 days
CSA205	Feline Comprehensive Plus with SDMA, FT4ED Superchem with SDMA, CBC, T3, T4, Free T4 by ED, FeLV, FIV, FCV Most comprehensive chemistry panel, a complete blood count, thyroid panel (T3, T4, FT4 ED), FeLV antigen detection, FIV antibody detection, SDMA for glomerular filtration rate estimation (see CT1035), and feline coronavirus titer. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-4 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA210	Feline Comprehensive with SDMA Superchem with SDMA, CBC, FeLV, FIV, FCV, <i>Toxoplasma</i> Ab - IgG/IgM Most comprehensive chemistry panel, a complete blood count, FeLV antigen detection, FIV antibody detection, SDMA for glomerular filtration rate estimation (see CT1035), feline coronavirus titer, and <i>Toxoplasma</i> titer (IgM and IgG). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
CSA715	Feline Comprehensive Wellness Screen with SDMA Superchem with SDMA, CBC, Total T4, Urinalysis, FeLV, FIV, Feline Heartworm Antibody The most comprehensive minimum database (superchemistry, complete blood count, urinalysis), total T4, FeLV antigen detection, FIV and heartworm antibody detection, and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-4 days
CSA240	Feline Heartworm Program Plus Miniscreen 4 Chemistry with Electrolytes, CBC, Heartworm Antigen, Feline Heartworm Antibody A small chemistry panel (TP, ALT, BUN, and glucose), a complete blood count, heartworm antibody and heartworm antigen detection. Interferences: marked hemolysis and lipemia.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
CSA190	Feline Total Health Check with SDMA Superchem with SDMA, CBC, T4, FIV, FeLV, FCV, <i>Toxoplasma</i> Ab - IgG/IgM Most comprehensive chemistry panel, a complete blood count, total T4, FeLV antigen and FIV antibody detection, SDMA for glomerular filtration rate estimation (see CT1035), feline coronavirus and <i>Toxoplasma</i> titer (IgM and IgG). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-4 days
CSA786	Feline Total Health Check with SDMA, UA Superchem with SDMA, CBC, T4, Urinalysis, FeLV, FIV, FCV, <i>Toxoplasma</i> Ab - IgG/IgM The most comprehensive minimum database (superchemistry, complete blood count, and urinalysis), total T4, FeLV antigen and FIV antibody detection, SDMA for glomerular filtration rate estimation (see CT1035), feline coronavirus and <i>Toxoplasma</i> titer (IgM and IgG). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-4 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA180	Feline Total Health Plus with SDMA Superchem with SDMA, CBC, T4, FIV, Feline Heartworm Antibody, FeLV, FCV, <i>Toxoplasma</i> Ab - IgG/IgM Most comprehensive chemistry panel, a complete blood count, total T4, FeLV antigen detection, FIV and heartworm antibody detection, SDMA for glomerular filtration rate estimation (see CT1035), feline coronavirus and <i>Toxoplasma</i> titer (IgM and IgG). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-4 days
CSA700	Feline Wellness Profile with SDMA Superchem with SDMA, CBC, T4, Urinalysis, FeLV, FIV The most comprehensive minimum database (superchemistry, complete blood count, and urinalysis), total T4, FeLV antigen and FIV antibody detection, and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
CSA683	Feline Wellness with SDMA, FeLV, FIV, O&P Wellness Chemistry with Electrolytes and SDMA, CBC, Feline Heartworm Antibody, FeLV, FIV, O&P with Centrifugation A comprehensive chemistry panel, a complete blood count, FeLV antigen detection, FIV and heartworm antibody detection, SDMA for glomerular filtration rate estimation (see CT1035), and fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection (CT805). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, and 5 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-4 days
CSA705	General Senior Profile with SDMA Superchem with SDMA, CBC, Total T4, Urinalysis The most comprehensive minimum database , total T4, and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	Daily Performed each shift
CSA440	Hyperthyroid Monitoring Profile with SDMA NSAID Chemistries AST (SGOT), ALT (SGPT), Alk Phos, BUN, Creatinine, Sodium, Potassium, Chloride, CBC, T4 A useful panel to evaluate the response to methimazole after a full minimum database has been run to make the initial diagnosis of hyperthyroidism. Includes an abbreviated chemistry (crea, BUN, AST, ALT, ALP), complete blood count, total T4, and SDMA. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	Daily Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA805	Hyperthyroid Panel with SDMA CBC, T4, NSAID Chemistries with GGT AST (SGOT), ALT (SGPT), Alk Phos, BUN, Creatinine, Sodium, Potassium, Chloride, GGT A useful panel to evaluate the response to methimazole after a full minimum database has been run to make the initial diagnosis of hyperthyroidism. Includes an abbreviated chemistry (crea, BUN, AST, ALT, ALP and GGT), complete blood count, total T4, and SDMA. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	Daily Performed each shift
CSA320	Liver Profile Liver Chemistry with Electrolytes, Complete Blood Count, Bile Acids (Pre and Post) A comprehensive liver panel (CSA324), complete blood count and pre/post bile acids (CT220).	0.5 mL serum for resting and post-prandial samples (fasted sample labeled Pre and 2-hour post sample labeled Post) and 1.0 mL EDTA whole blood 2 serum in red top or serum separator (pre and post), lavender top	1-2 days
CSA070	Miniscreen 4 Chem, CBC Miniscreen 4 Chemistries with Electrolytes (Total Protein, ALT (SGPT), BUN, Sodium, Potassium, Chloride, Glucose) and CBC Miniscreen chemistry (TP, ALT, BUN and glucose) and a complete blood count.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	Daily Performed each shift
CSA810	NSAID 1 with SDMA Chemistry Panel with Electrolytes (Total Protein, ALT (SGPT), Alkaline Phosphatase, Total Bilirubin, BUN, Creatinine, Sodium, Potassium, Chloride) with SDMA, CBC, Urinalysis A minimum database (chemistry with SDMA (CSA804), complete blood count and urinalysis) limited to renal and liver values used to monitor a patient on NSAID therapy. Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	Daily Performed each shift
CSA815	NSAID 2 with SDMA Chemistry Panel with Electrolytes (Total Protein, ALT (SGPT), Alkaline Phosphatase, Total Bilirubin, BUN, Creatinine, Sodium, Potassium, Chloride) with SDMA, CBC, Urinalysis, Bile Acids A minimum database (chemistry with SDMA (CSA804), complete blood count and urinalysis) and a single bile acids test. Interferences: marked hemolysis and lipemia. Ursodeoxycholic acid may be detected by bile acid assay causing falsely elevated values.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA820	NSAID 3 with SDMA Chemistry Panel with Electrolytes(Total Protein, ALT (SGPT), Alkaline Phosphatase, Total Bilirubin, BUN, Creatinine Sodium, Potassium, Chloride) with SDMA, CBC, Urinalysis, Urine Bile Acid: Creatinine Ratio A minimum database (chemistry with SDMA (CSA804), complete blood count and urinalysis) and a urine bile acid: creatinine ratio (CT227). Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
CSA516	Pre-Op Chem with Electrolytes, SDMA, CBC Pre-Op Chemistries with Electrolytes and SDMA, CBC A smaller chemistry profile which includes electrolytes and SDMA (CSA043) and a complete blood count. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	Daily Performed each shift
CSA053	Pre-Op Chem with Electrolytes, SDMA, CBC, Urinalysis Pre-Op Chemistries with Electrolytes and SDMA, CBC, Urinalysis A minimum database that includes a smaller chemistry with SDMA (CSA043), a complete blood count and urinalysis. Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	Daily Performed each shift
CSA050	Pre-Op Panel Plus with SDMA Pre-Op Screen with Electrolytes and SDMA, CBC, PT, and aPTT A smaller chemistry profile with SDMA (CSA040), a complete blood count, prothrombin time, and partial thromboplastin time. Interferences: marked hemolysis and lipemia. Ensure that the submission of a blue top tube that has been appropriately filled to the designated mark.	0.5 mL serum, 1.0 mL EDTA whole blood and 0.5 mL of citrated plasma Serum in red top or serum separator tube, lavender top, citrated plasma (or blue top)	Daily Performed each shift
CSA512	Pre-Op Panel Plus with SDMA, Heartworm Pre-Op Screen with Electrolytes and SDMA, CBC, PT, aPTT, Heartworm Antigen A smaller chemistry profile with SDMA (CSA040), complete blood count, heartworm antigen detection, and prothrombin time and partial thromboplastin time. Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood and 0.5 mL of citrated plasma Serum in red top or serum separator tube, lavender top, citrated plasma (or blue top)	1-2 days
CSA920	Pre-Op Panel with SDMA, T4, UA Pre-Op Screen with Electrolytes and SDMA, CBC, T4, Urinalysis A minimum database including a smaller chemistry profile (CSA040), a complete blood count, urinalysis and a T4. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	Daily Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA510	Pre-Op Screen with SDMA, CBC, PT, PTT, FeLV, FIV Pre-Op Screen with Electrolytes and SDMA, CBC, PT, PTT, FIV, and FeLV A smaller chemistry profile with SDMA (CSA040), complete blood count, FeLV antigen and FIV antibody detection, prothrombin and partial thromboplastin time. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood and 0.5 mL of citrated plasma Serum in red top or serum separator tube, lavender top, citrated plasma (or blue top)	1-2 days
CSA514	Pre-Op Screen with SDMA, CBC, T4, FeLV, FIV Pre-Op Screen with Electrolytes and SDMA, CBC, T4, FIV, and FeLV CSA508 to which FeLV antigen and FIV antibody detection has been added. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
CSA055	Pre-Op Screen with SDMA, CBC Pre-Op Screen with Electrolytes and SDMA and CBC A smaller chemistry profile with SDMA (CSA040) and a complete blood count. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	Daily Performed each shift
CSA056	Pre-Op Screen with SDMA, CBC, Heartworm Pre-Op Screen with Electrolytes and SDMA, CBC, and Heartworm A smaller chemistry panel with SDMA (CSA040), a complete blood count and heartworm antigen detection. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
CSA052	Pre-Op Screen with SDMA, CBC, Urinalysis Pre-Op Screen with Electrolytes and SDMA, CBC, and Urinalysis A minimum database including smaller chemistry profile with SDMA (CSA040), a complete blood count, and urinalysis. Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	Daily Performed each shift
CAC055	Pre-Op with SDMA, CBC with Accuplex® Pre-Op Screen with Electrolytes and SDMA, CBC, Accuplex® A smaller chemistry panel with SDMA (CSA040), a complete blood count, and Accuplex®. Interferences: marked hemolysis and lipemia. Accuplex® (CAC100) can be reflexed for an additional charge.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA051	Pre-Op with SDMA, CBC, FeLV, FIV Pre-Op Screen with Electrolytes and SDMA, CBC, FIV, FeLV A smaller chemistry panel with SDMA (CSA040), a complete blood count, FeLV antigen and FIV antibody detection. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
CSA508	Pre-Op with SDMA, CBC, T4 Pre-Op Screen with Electrolytes and SDMA, CBC, Total T4 A smaller chemistry profile with SDMA (CSA040), a complete blood count, and a total T4. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	Daily Performed each shift
CSA310	Renal Profile with SDMA Chemistry Renal Profile with Electrolytes and SDMA, CBC, and Urinalysis A minimum database including a chemistry profile with SDMA (CT7008), complete blood count and urinalysis. Interferences: marked hemolysis or lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	Daily Performed each shift
CSA111	Superchem with SDMA, CBC, Heartworm, T4 Free by Immulite Superchem with SDMA, CBC, T4 (Free) by Immulite, Heartworm Antigen Most comprehensive chemistry panel with SDMA (CSA010), a complete blood count, T4 (Free) run by Immulite, and heartworm antigen detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
CSA080	Senior Comprehensive Plus with SDMA Superchem with SDMA, CBC, Total T4, Free T4 by ED, and TSH The most comprehensive chemistry panel with SDMA (CSA010), a complete blood count, and thyroid panel (T4, FT4 ED, and TSH). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-6 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA090	Senior Comprehensive Profile with SDMA Superchem with SDMA, CBC, Total T4, and Free T4 by ED The most comprehensive chemistry panel with SDMA (SA010), a complete blood count, and thyroid panel (T4 and FT4 ED). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-6 days
CSA105	Senior Comprehensive with SDMA, T3 Superchem with SDMA, CBC, T3, Total T4, Free T4 by ED The most comprehensive chemistry panel with SDMA (SA010), a complete blood count, and thyroid panel (T3, T4, and FT4 ED). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-6 days
CSA740	Senior Feline Plus with SDMA Superchem with SDMA, CBC, T4, Urinalysis, FeLV, FIV, O&P with Centrifugation The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), a total T4, FeLV antigen and FIV antibody detection, and fecal analysis using zinc sulfate centrifugation/flotation for ova and parasite detection. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
CSA566	Senior Profile 1 with SDMA, UMIC Superchem with SDMA, CBC, Total T4, Urinalysis, Urine Culture The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), a total T4, and urine culture. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 7.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-4 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CAC705	Senior Profile 1 with SDMA, Accuplex® Superchem with SDMA, CBC, Total T4, Urinalysis, Accuplex® The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), a total T4, and Accuplex® (canine vector-borne disease screening for heartworm, Lyme disease (includes screening for antibodies against two C6 peptides), <i>E. canis</i> and <i>A. phagocytophilum</i>). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results. Accuplex® (CAC100) can be reflexed for an additional charge.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
CAC770	Senior Profile 1 with SDMA, Fecal Combo, Accuplex® Superchem with SDMA, CBC, Total T4, Urinalysis, Fecal O&P with Centrifugation, <i>Giardia</i> , Accuplex® The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), a total T4, and fecal analysis using zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection. This profile allows a reflex Accuplex® for canine vector borne disease screening at a nominal fee. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results. Accuplex® (CAC100) can be reflexed for an additional charge.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine, and 6 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
CSA770	Senior Profile 1 with SDMA, Fecal Combo Superchem with SDMA, CBC, Total T4, Urinalysis, Fecal O&P with Centrifugation, <i>Giardia</i> The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), a total T4, and fecal analysis using zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine, and 6 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
CSA745	Senior Profile 1 with SDMA, Heartworm, Ova & Parasite Superchem with SDMA, CBC, Total T4, Urinalysis, Heartworm Antigen, Fecal O&P with Centrifugation The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), a total T4, heartworm antigen detection and fecal analysis using zinc sulfate centrifugation/flotation for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA750	Senior Profile 1 with SDMA, Ova & Parasite Superchem with SDMA, CBC, Total T4, Urinalysis, Fecal O&P with Centrifugation The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), a total T4, and fecal analysis using zinc sulfate centrifugation/flotation for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
CSA730	Senior Profile 1 with SDMA, Vaccine Titers Superchem with SDMA, CBC, T4, Urinalysis, Distemper/Parvo Vaccinal Titer The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), a total T4, and Distemper/Parvovirus vaccine titer. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-4 days
CSA780	Senior Profile 2 with SDMA, FeLV, FIV Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, FIV, and FeLV The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), a thyroid panel (T4 and FT4 ED), FeLV antigen and FIV antibody detection. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-6 days
CSA720	Senior Profile 2 with SDMA Superchem with SDMA, CBC, Total T4, Urinalysis, and Free T4 by ED The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis) and a thyroid panel (T4 and FT4 ED). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-6 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA775	Senior Profile 2 with SDMA, Fecal Combo Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, and O&P with Centrifugation, <i>Giardia</i> The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), a thyroid panel (T4 and FT4 ED), and fecal analysis using zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 6 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-6 days
CSA790	Senior Profile 2 with SDMA, Feline Heartworm Ab, O&P Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, O&P with Centrifugation, and Feline Heartworm Antibody The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), a thyroid panel (T4 and FT4 ED), heartworm antibody detection, and fecal analysis using zinc sulfate centrifugation/flotation for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Antech provided fecal container, lavender top, red top, urine transport tube	1-6 days
CSA755	Senior Profile 2 with SDMA, Ova & Parasite Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, and O&P with Centrifugation The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), a thyroid panel (T4 and FT4 ED), and fecal analysis using zinc sulfate centrifugation/flotation for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-6 days
CSA746	Senior Profile 2 with SDMA, UMIC Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, and Urine Culture The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a thyroid panel (T4 and FT4 ED), and urine culture. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 7.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-6 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA760	Senior Profile 2 with SDMA, Vaccine Titers Superchem with SDMA, CBC, T4, Urinalysis, Free T4 by ED, and Distemper/Parvo Vaccinal Titer The most comprehensive minimum database (superchemistry with SDMA (SA010), complete blood count and urinalysis), a thyroid panel (T4 and FT4 ED), and Distemper/Parvovirus vaccine titer. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-6 days
CSA768	Senior Profile with SDMA, Feline Heartworm Ab, Fecal Combo Superchem with SDMA, CBC, T4, Urinalysis, Feline Heartworm Antibody, and O&P with Centrifugation, <i>Giardia</i> The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis), a total T4, heartworm antibody detection, and fecal analysis using zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 6 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-4 days
CSA955	Superchem with SDMA, CBC, O&P, <i>Giardia</i> Superchem with SDMA, CBC, O&P with Centrifugation, <i>Giardia</i> Most comprehensive chemistry profile (superchemistry with SDMA (CSA010)), a complete blood count, and fecal analysis using zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA (CT808) for ova and parasite detection. Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days
CSA020	Superchem with SDMA, CBC The most comprehensive chemistry with SDMA (CSA010) combined with a complete blood count. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	Daily Performed each shift
CRECHECK	Superchem with SDMA, CBC Recheck Panel This is a follow up submission to a previously submitted CSA020 that consists of the most comprehensive chemistry with SDMA and a complete blood count. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	Daily Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA460	Superchem with SDMA, CBC, Feline Serology 1 Superchem with SDMA, CBC, FeLV, FIV Antibody, and FCV Titer The most comprehensive chemistry panel with SDMA (CSA010), a complete blood count, FeLV antigen and FIV antibody detection, and a feline coronavirus titer. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-4 days
CSA490	Superchem with SDMA, CBC, FeLV, FIV Superchem with SDMA, CBC, FeLV, and FIV Antibody The most comprehensive chemistry panel with SDMA (CSA010), a complete blood count, FeLV antigen and FIV antibody detection. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
CSA500	Superchem with SDMA, CBC, Lyme Superchem with SDMA, CBC, and Lyme Titer IgG The most comprehensive chemistry panel with SDMA (CSA010), a complete blood count, and Lyme IgG titer. Interferences: marked hemolysis and lipemia.	0.75 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-4 days
CSA450	Superchem with SDMA, CBC, Special Thyroid Profile Superchem with SDMA, CBC, T3, T4, Free T4 by ED, TSH, Thyroglobulin Auto Antibody (Canine Only) The most comprehensive chemistry profile (superchemistry with SDMA (CSA010)), a complete blood count, and thyroid panel (T3, T4, FT4 ED, TSH and thyroglobulin autoantibody (TGAA)). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-6 days
CSA021	Superchem with SDMA, CBC, UA Superchem with SDMA, CBC, Urinalysis The most comprehensive minimum database (superchemistry with SDMA (CSA010), complete blood count and urinalysis). Interferences: marked hemolysis and lipemia (serum and whole blood).	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	Daily Performed each shift
CSA110	Total Body Function Plus with SDMA Superchem with SDMA, CBC, T4, and Heartworm Antigen The most comprehensive chemistry panel with SDMA (CSA010), a complete blood count, a total T4, and heartworm antigen detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA788	Total Body Function Plus with SDMA, O&P, <i>Giardia</i> Superchem with SDMA, CBC, T4, Heartworm Antigen and O&P with Centrifugation and <i>Giardia</i> The most comprehensive chemistry panel with SDMA (CSA010), a complete blood count, total T4, heartworm antigen detection and fecal analysis using zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA (CT808) for ova and parasite detection. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum, 1.0 mL EDTA whole blood, and 6 grams feces Serum in red top or serum separator tube, lavender top, Antech provided fecal container	1-2 days
CSA120	Total Body Function with SDMA Superchem with SDMA, CBC, and Total T4 The most comprehensive chemistry panel with SDMA (CSA010), a complete blood count, and a total T4. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	Daily Performed each shift
CAC120	Total Body Function with SDMA, Accuplex® Superchem with SDMA, CBC, T4, and Accuplex® The most comprehensive chemistry panel with SDMA (CSA010), a complete blood count, total T4, and reflex Accuplex® for canine vector borne disease screening. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results. Accuplex® (CAC100) can be reflexed for an additional charge.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
CSA724	Total Body Function with SDMA, FeLV, FIV, Feline Heartworm Ab Superchem with SDMA, CBC, T4, FIV, FeLV, Feline Heartworm Antibody The most comprehensive chemistry panel with SDMA (CSA010), a completely blood count, a total T4, FeLV antigen detection, FIV antibody and heartworm antibody detection. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-4 days
CSA120B	Total Body Function with SDMA, T4 Immulite Superchem with SDMA, CBC, T4 (Free) by Immulite The most comprehensive chemistry panel with SDMA (CSA010), a complete blood count, and T4 (Free) run by Immulite (CT499). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CAC120B	Total Body Function with SDMA, T4 Immulite, Accuplex® Superchem with SDMA, CBC, T4 (Free) by Immulite, and Accuplex® The most comprehensive chemistry panel with SDMA (CSA010), a complete blood count, T4 (Free) run by Immulite (CT499) and Accuplex® (CAC100). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results. Accuplex® (CAC100) can be reflexed for an additional charge.	1.25 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
CSA026	Vet Screen 2 with SDMA Total Protein, Albumin, Globulin, Albumin/Globulin Ratio, AST (SGOT), ALT (SGPT), Alkaline Phosphatase, Total Bilirubin, BUN, Creatinine, BUN/Creatinine Ratio, Phosphorous, Glucose, Calcium, Sodium, Potassium, Chloride, Cholesterol, Amylase, CPK, and SDMA A comprehensive chemistry profile with SDMA (CSA025) to which amylase evaluation has been added. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CSA035	Vet Screen 2 with SDMA, CBC Vet Screen 2 with SDMA and Amylase, and CBC A comprehensive chemistry panel with SDMA (CSA026) and complete blood count. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	Daily Performed each shift
CSA025	Vet Screen with SDMA Total Protein, Albumin, Globulin, Albumin/Globulin Ratio, AST (SGOT), ALT (SGPT), Alkaline Phosphatase, GGT, Total Bilirubin, BUN, Creatinine, BUN/Creatinine Ratio, Phosphorous, Glucose, Calcium, Sodium, Potassium, Sodium/Potassium Ratio, Chloride, Cholesterol, CPK, and SDMA A comprehensive chemistry profile with SDMA. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CSA914	Vet Screen with SDMA, CBC, FeLV, FIV, FIP Vet Screen with SDMA, CBC, FeLV, FCV, and FIV A comprehensive profile with SDMA (CSA025), a complete blood count, FeLV antigen and FIV antibody detection, and feline coronavirus titer. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA908	Vet Screen with SDMA, CBC, FeLV, FIV, FIP, T4 Vet Screen with SDMA, CBC, FeLV, FCV, FIV, and T4 A comprehensive profile with SDMA (CSA025), a complete blood count, total T4, FeLV antigen and FIV antibody detection, and feline coronavirus titer. OR CSA914 to which a total T4 has been added. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-3 days
CSA910	Vet Screen with SDMA, CBC, FeLV, FIV, T4 Vet Screen with SDMA, CBC, FeLV, FIV, and T4 A comprehensive profile with SDMA (CSA025), a complete blood count, total T4, FeLV antigen and FIV antibody detection. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
CSA906	Vet Screen with SDMA, CBC, FeLV, FIV, UA Vet Screen with SDMA, CBC, FeLV, FIV, and Urinalysis A complete minimum database (comprehensive chemistry with SDMA (CSA025), complete blood count, and urinalysis), FeLV antigen and FIV antibody detection. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days
CSA912	Vet Screen with SDMA, CBC, T4, FT4ED Vet Screen with SDMA, T4, CBC, Free T4 by ED A comprehensive chemistry with SDMA (CSA025), a complete blood count, and a thyroid profile (total T4 and FT4 ED). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-6 days
CSA904	Vet Screen with SDMA, CBC, UA, O&P, Heartworm Vet Screen with SDMA, CBC, Urinalysis, O&P with Centrifugation, Heartworm Antigen A complete minimum database (comprehensive chemistry with SDMA (CSA025), complete blood count, and urinalysis), heartworm antigen detection and fecal analysis using zinc sulfate centrifugation/flotation (CT805) for ova and parasite detection. Interferences: marked hemolysis and lipemia.	0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 5 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA030	Vet Screen with SDMA, CBC A comprehensive chemistry profile with SDMA (CSA025) and a complete blood count. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	Daily Performed each shift
CRECHECK2	Vet Screen with SDMA, CBC Recheck Panel Vet Screen with SDMA, CBC This is a comprehensive chemistry with SDMA (CSA025) and complete blood count submitted within 30 days of a previous accession to assess/monitor trends in values. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Lavender top, serum in red top or serum separator	Daily Performed each shift
CSA039	Vet Screen with SDMA, CBC, Heartworm Ag Vet Screen with SDMA, CBC, Heartworm Antigen A comprehensive chemistry profile with SDMA (CSA025), a complete blood count, and heartworm antigen detection. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
CSA034	Vet Screen with SDMA, CBC, T4 Vet Screen with SDMA, CBC, Total T4 A comprehensive chemistry profile with SDMA (CSA025), complete blood count, and a total T4. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	Daily Performed each shift
CSA032	Vet Screen with SDMA, CBC, T4 (Free) Vet Screen with SDMA, CBC, T4 (Free) by Immulite A comprehensive chemistry with SDMA (CSA025), a complete blood count and T4 (Free) by Immulite. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Lavender top tube, serum separator	1-2 days
CSA038	Vet Screen with SDMA, CBC, T4 (Free), UA Vet Screen with SDMA, CBC, T4 (Free) by Immulite, Urinalysis A complete minimum database (comprehensive chemistry with SDMA (CSA025), complete blood count and urinalysis) and a T4 (Free) by Immulite. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA530	Vet Screen with SDMA, CBC, T4, Heartworm Vet Screen with SDMA, CBC, Heartworm Antigen, T4 A comprehensive chemistry panel with SDMA (CSA025), complete blood count, heartworm antigen detection, and total T4 assessment. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
CSA036	Vet Screen with SDMA, CBC, UA Vet Screen with SDMA, CBC, Urinalysis A complete minimum database including a comprehensive chemistry profile with SDMA (CSA025), a complete blood count, and urinalysis. Interferences: marked hemolysis and lipemia (serum and whole blood).	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	Daily Performed each shift
CSA922	Vet Screen with SDMA, CBC, UA, O&P, Giardia Vet Screen with SDMA, CBC, Urinalysis, O&P with Centrifugation, <i>Giardia</i>	0.5 mL serum, 1.0 mL EDTA whole blood, 6.0 mL urine, and 6 grams feces Serum in red top or serum separator tube, lavender top, urine transport tube, Antech provided fecal container	1-2 days
CSA037	Vet Screen with SDMA, CBC, UA, T4 Vet Screen with SDMA, CBC, Total T4, Urinalysis A complete minimum database (comprehensive chemistry profile with SDMA (CSA025), a complete blood count, and urinalysis) and a total T4. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.5 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator tube, lavender top, urine transport tube	Daily Performed each shift

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

HEMATOLOGY

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CADD290	Add-on PT/PTT Prothrombin Time and aPTT Prothrombin Time and Activated Partial Thromboplastin Time. Interferences: marked hemolysis or lipemia. Partially full blue top tubes may falsely increase sample's coagulation time. Note: clotting of the sample may preclude the analysis. The blue top tube needs to be filled to 2/3 or more of its capacity. Partially full blue top tubes may falsely increase coagulation times.	0.5 mL citrated plasma in blue top tube	Daily Performed each shift
CADD140	Add-on Reticulocytes Reticulocytes, Absolute Reticulocytes Includes an unadjusted reticulocyte percent and absolute reticulocyte count. Only aggregate reticulocytes are counted. Interferences: marked hemolysis.	1.0 mL EDTA whole blood in lavender top tube	Daily Performed each shift
CT315	Blood Type, Canine - DEA 1.1 Blood Type - DEA 1.1 This test determines the appropriateness of canine blood donors prior to transfusion. Dog Erythrocyte Antigen (DEA) 1.1 is the most antigenic blood group in dogs, and canine blood donors are often screened to see if they are positive or negative for this blood group. Note: keep sample refrigerated prior to transport and send with an ice pack.	1.0 mL EDTA whole blood in lavender top tube	2-4 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS16100	Blood Type, Canine - Full Panel Blood Type: DEA 1, DEA 4, DEA 5, DEA 7 This panel is used for RBC typing for identifying blood donors. It evaluates the presence of the blood group antigens; DEA 1.1, 1.2, 1.3, 4, 5, and 7. Antibody screen can be performed at no extra charge if clients submit serum with LTT.	1.0 mL EDTA whole blood in lavender top tube	2-4 days
CT331	CBC, Path Review CBC, clinical pathologist review of blood smear and results This includes WBC, RBC, HGB, HCT, MCV, MCH, MCHC, platelet count and estimate, WBC differential, RBC, and WBC morphology. A clinical pathologist will then review results and blood smears and provide an interpretative comment of the complete blood count. Interferences: marked hemolysis and lipemia.	1.0 mL EDTA whole blood in lavender top tube	1-3 days
CT337	CBC/Reticulocyte CBC, Reticulocyte Count This includes WBC, RBC, HGB, HCT, MCV, MCH, MCHC, platelet count and estimate, WBC differential, RBC, and WBC morphology (T330) and a reticulocyte count (T425). Interferences: marked hemolysis and lipemia.	1.0 mL EDTA whole blood in lavender top tube	Daily Performed each shift
CSA290	Coagulation Profile 1 CBC, PT and aPTT, Fibrinogen Quantitative, and D-Dimer Panel used to evaluate a patient exhibiting unexplained bleeding and/or thromboembolic disorders. Additionally, it can be used to evaluate risk of bleeding secondary to provocative procedures. Includes a CBC, PT/PTT, fibrinogen and D-dimer. Interferences: marked hemolysis, lipemia or clotting.	1.0 mL EDTA whole blood and 1.0 mL citrated plasma Lavender top, citrated plasma (or filled blue top)	Daily Performed each shift

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA300	Coagulation Profile 2	1.0 mL EDTA whole blood and 1.0 mL citrated plasma	Daily
	<p>PT and aPTT, Fibrinogen (Quantitativ), D-Dimer, and Platelet Count</p> <p>Panel used to evaluate a patient exhibiting unexplained bleeding and/or thromboembolic disorders. Additionally, it can be used to evaluate risk of bleeding secondary to provocative procedures. Includes PT/PTT, platelet count, fibrinogen and D- dimer.</p> <p>Interferences: marked hemolysis, lipemia or clotting.</p>	<p>Lavender top, citrated plasma (or filled blue top)</p>	<p>Performed each shift</p>
CSA305	Coagulation Profile 3	1.0 mL EDTA whole blood and 1.0 mL citrated plasma	Daily
	<p>PT and aPTT, and Platelet Count</p> <p>Panel used to evaluate a patient exhibiting unexplained bleeding. Additionally, it can be used to evaluate risk of bleeding secondary to provocative procedures. Includes PT/PTT and platelet count.</p> <p>Interferences: marked hemolysis, marked lipemia or clotting.</p>	<p>Citrated plasma (or filled blue top)</p>	<p>Performed each shift</p>
CT330	Complete Blood Count	1.0 mL EDTA whole blood in lavender top tube	Daily
	<p>Includes WBC, RBC, Hemoglobin, Hematocrit, MCV, MCH, MCHC, Platelet count and Differential: Percent and Absolute Neutrophils, Bands, Lymphocytes, Monocytes, Eosinophils, Basophils</p> <p>Includes WBC, RBC, HGB, HCT, MCV, MCH, MCHC, platelet count and estimate, WBC differential, RBC and WBC morphology.</p> <p>Interferences: marked hemolysis and lipemia.</p>		<p>Performed each shift</p>
CT540	Direct Coombs' Test (Warm)	1.0 mL of EDTA whole blood in lavender top tube	1-2 days
	<p>Used to investigate hemolytic anemia.</p> <p>Interferences: marked hemolysis and lipemia. Test should be performed within 24 hours of sample collection.</p>		

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT365	Fibrinogen This test is used to measure the concentration of functional fibrinogen in the plasma. Interferences: clotted sample precludes analysis. Citrated Whole Blood (blue top tube) or citrated plasma is the only acceptable sample. The tube should be greater than 2/3rds filled. If submitting separated citrated plasma, label it as Citrated Plasma.	0.5 mL citrated plasma collected as whole blood in blue top tube at least 2/3rds full to the fill line	Daily Performed each shift
CS86493	Lymphocyte Phenotyping A diagnostic procedure used to assess the chemical and physical characteristics of a lymphocyte population. Most commonly used to determine if there is a monoclonal subset of cells.	Submit patient history, CSU Clinical Immunology Form, and the sample that correlates to the source of disease. Choose sample based on disease source, for: A. Peripheral lymphocytosis, suspect leukemia, etc.: submit 1.0 mL whole blood or bone marrow in lavender top tube(EDTA) and CBC results that are no older than 48 hours. OR B. Multicentric lymphoma, lymphoma/ suspect neoplastic cells in lymph node(s) or organ(s): submit lymph node or other tissue/organ aspirate in non-additive red or white top tube. See Collection Method (or Special Instructions for Draw) for additional instructions. OR C. Pleural effusion and ascites: submit 0.5 mL of cavity fluid in a lavender top tube (EDTA) and 0.5 mL in a red top tube (with a few drops of patient serum added for protein samples less than 4 mg/dL as evaluated on a refractometer.) If only able to submit one effusion sample, submission of lavender top tube (EDTA) is preferred.	7-10 days Samples should be received at Antech Laboratory on Monday or Tuesday for timely processing and testing by the reference laboratory performing this test (CSU).

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CREVW	Path Review	0.5 mL EDTA whole blood in lavender top tube	1-2 days
	This code is used to Add-on a clinicopathologist review to a CBC.		
CT400	Platelet Count	1.0 mL EDTA whole blood in lavender top tube	Daily
	Interferences: marked hemolysis.		
CT410	Prothrombin Time	0.5 mL citrated plasma collected as whole blood in blue top tube at least 2/3rds full to the fill line	Daily
	Prothrombin Time measures the integrity of the extrinsic and common components of the coagulation cascade. Interferences: marked hemolysis or lipemia. Partially full blue top tube may falsely increase sample's coagulation time. Note: clotting of the sample may preclude the analysis. The blue top tube needs to be filled to 2/3 or more of its capacity.		

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT415	PT/PTT	0.5 mL citrated plasma collected as whole blood in blue top tube at least 2/3rds full to the fill line	Daily
Add-on Equivalent CADD290	<p>Prothrombin Time and Activated Partial Thromboplastin Time</p> <p>Prothrombin Time and Activated Partial Thromboplastin Time. Prothrombin time measures the integrity of the extrinsic and common components of the coagulation cascade. Partial Thromboplastin Time (PTT) measures the integrity of the intrinsic and common components of the coagulation cascade.</p> <p>Interferences: marked hemolysis or lipemia. Partially full blue top tube may falsely increase sample's coagulation time.</p> <p>Note: clotting of the sample may preclude the analysis. The blue top tube needs to be filled to 2/3 or more of its capacity. Partially full blue top tubes may falsely increase coagulation times.</p>		Performed each shift
CT395	PTT	0.5 mL citrated plasma collected as whole blood in blue top tube at least 2/3rds full to the fill line	Daily
	<p>Activated Partial Thromboplastin Time</p> <p>Partial Thromboplastin Time (PTT) measures the integrity of the intrinsic and common components of the coagulation cascade.</p> <p>Interferences: marked hemolysis or lipemia. Partially full blue top tube may falsely increase sample's coagulation time.</p> <p>Note: clotting of the sample may preclude the analysis. The blue top tube needs to be filled to 2/3 or more of its capacity.</p>		Performed each shift
CS17123	Von Willebrand Factor	1.5 mL citrated plasma collected as whole blood in blue top tube, spun and plasma placed in non-additive tube (labeled as Citrated plasma). Freeze citrated plasma.	7-10 days
	<p>To evaluate for Von Willebrands factor deficiency. Results will be expressed as a percentage. The percentage will indicate the risk of bleeding problems.</p> <p>Note: for genetic screening, do not test bitches in season, pregnant or lactating. Do not test unhealthy animals (with no evidence of primary clotting issues), those on medication for a recent illness, or those vaccinated within 14 days. Large breed pups can be tested at seven weeks of age. Small breed pups should be over 12 weeks of age prior to testing.</p>		

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CHEMISTRY

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA665	Adult Wellness Chemistry with SDMA Wellness Chemistry (ALT (SGPT), Alkaline Phosphatase, Total Protein, Albumin, Globulin, A/G Ratio, BUN, Creatinine, BUN/Creatinine Ratio, Glucose, Sodium, Potassium, Na/K Ratio, and Chloride) with Electrolytes and SDMA A comprehensive chemistry panel and SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator	Daily Performed each shift
CT223	Bile Acid Post Interferences: marked hemolysis or lipemia. Note: ursodeoxycholic acid may be detected by bile acid assay, causing falsely elevated values. Feed the patient a maintenance diet meal after a 10-12 hour fast and draw a serum sample 2 hours post-feeding. Label the tube as postprandial bile acids.	0.5 mL serum in red top or serum separator tube	1-2 days
CT220	Bile Acid Profile Fasted bile acids, 2 hour post feeding bile acids Pre and post prandial bile acids. Interferences: marked hemolysis or lipemia. Ursodeoxycholic acid may be detected by bile acid assay, causing falsely elevated values. Note: suggested protocol: 1. Obtain a fasting serum sample (1 mL) and label the tube preprandial. 2. Obtain a second serum sample (1 mL) 2 hours after feeding the animal a maintenance diet and label the tube postprandial.	0.5 mL serum in red top or serum separator tube (fasted sample labeled Pre) and 0.5 serum in red top or serum separator tube (2-hour post sample labeled Post)	1-2 days
CT225	Bile Acids Single Bile Acid Single sample labeled as either pre or post bile acids. Interferences: marked hemolysis or lipemia. Note: ursodeoxycholic acid may be detected by bile acid assay, causing falsely elevated values. 1.0 mL of serum labeled either as fasted (10-12 hours fast) or post prandial (2 hours after feeding a maintenance diet meal).	0.5 mL serum in red top or serum separator tube	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT110	Calcium Total Calcium Interferences: marked hemolysis or lipemia. Lipemia can result in falsely elevated total calcium values.	0.5 mL serum in red top tube or serum separator tube	Daily Performed each shift
CSA324	Chemistry Panel Total Protein, Albumin, Globulin, AST (SGOT), ALT (SGPT), Alkaline Phosphatase, GGTP, Total Bilirubin, BUN, Glucose, Sodium, Potassium, Chloride A liver chemistry panel that is useful when trending response of specific values to the management of a previously defined hepatobiliary disease. Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CT7008	Chemistry Renal with SDMA Profile Total Protein, Albumin, Globulin, A/G Ratio, BUN, Creatinine, BUN/Creat Ratio, Phosphorus, Calcium, Corrected Calcium, Sodium, Potassium, Na/K Ratio, Chloride, and SDMA A comprehensive profile panel that can be used to trend response to therapy when treating renal disease (glomerular or tubulointerstitial) which includes SDMA for glomerular filtration rate estimation (see CT1035). Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CT140	Electrolyte Screen Sodium, potassium, sodium/potassium ratio, chloride, bicarbonate and anion gap Evaluation of serum sodium, potassium, sodium/potassium ratio, chloride, bicarbonate, and anion gap. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CSA250	Feline Heartworm Program Miniscreen 4 Chemisty with Electrolytes, Feline Heartworm Antibody A small chemistry panel (total protein, ALT, BUN, glucose) combined with a feline heartworm antibody test. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	1-2 days
CT145	GGTP Gamma-glutamyl transferase Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS18537	Ionized Calcium Evaluation of the ionized fraction, rather than protein bound fraction, of calcium. Interferences: aerobic exposure, or severe lipemia. Note: sample must be anaerobically transferred from spun collection tube into plain red top.	2.0 mL serum anaerobically transferred from spun red top or serum separator tube into a plain, unopened red top tube. Call customer service for specimen handling protocol: 1-800-341-3440, dial 0	2-4 days
CSA321	Liver Chemistry Screen Total Protein, Albumin, Globulin, A/G Ratio, AST (SGOT), ALT (SGPT), Alkaline Phosphatase, GGT, Total Bilirubin, Direct Bilirubin, Sodium, Potassium, Chloride, and Cholesterol A chemistry panel that includes all cholestatic and hepatocellular liver enzymes, bilirubin, albumin, globulin and cholesterol. Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top tube or serum separator tube	Daily Performed each shift
CT165 Add-on Equivalent CADD90	PrecisionPSL PrecisionPSL (Pancreatic Sensitive Lipase) is used to help diagnose pancreatitis in patients with consistent clinical signs. Interferences: marked hemolysis or lipemia. Hemolysis will falsely decrease and lipemia will falsely increase PSL results.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CSA040	Pre-Op Screen with SDMA Includes Total Protein, Albumin, Globulin, A/G Ratio, ALT (SGPT), Alkaline Phosphatase, BUN, Creatinine, BUN/Creat Ratio, Sodium, Potassium, Chloride, Na/K Ratio, Glucose, and SDMA A smaller chemistry which includes total protein, albumin, globulin, ALT, ALP, BUN, creatinine, glucose and SDMA. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CT240 Add-on Equivalent CADD130	Protein Electrophoresis Total protein, Albumin, Globulin, Alpha 1, Alpha 2, Beta 1, and Gamma 1 fraction assessment with interpretation. An evaluation of the globulin fraction of the serum (alpha 1, alpha 2, beta, and gamma) to determine if the globulin fraction is monoclonal based on these components. Interferences: marked hemolysis or lipemia. Hemolysis may falsely elevate the beta globulin fraction.	0.5 mL serum in red top or serum separator tube	2-4 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA311	Renal Chemistry Screen with SDMA Total Protein, Albumin, Globulin, Albumin/globulin Ratio, BUN, Creatinine, Phosphorous, Glucose, Calcium, Sodium, Potassium, Sodium/potassium Ratio, and Chloride A chemistry panel which includes albumin, globulin, BUN, creatinine, phosphorous, glucose, calcium, sodium, potassium, chloride and SDMA. Interferences: marked hemolysis or lipemia.	0.5 mL sereum in red top or serum separator tube	Daily Performed each shift
CT1035	SDMA SDMA (symmetric dimethylarginine) is freely filtered by the kidneys, and elevation in SDMA is indicative of a reduced glomerular filtration rate.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CSA010	Superchem with SDMA Total Protein, Albumin, Globulin, Albumin/Globulin Ratio, AST (SGOT), ALT (SGPT), Alkaline Phosphatase, GGT, Total Bilirubin, BUN, Creatinine, BUN/Creatinine Ratio, Phosphorous, Glucose, Calcium, Magnesium, Sodium, Potassium, Sodium/Potassium Ratio, Chloride, Cholesterol, Triglyceride, PrecisionPSL, CPK, and SDMA Most comprehensive chemistry profile which includes SDMA for estimation of the glomerular filtration rate (CT1035) Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CS85783	Troponin I	1.0 mL frozen serum (should be fasting, non-hemolyzed sample)	
CS16016	Vitamin D	1.0 mL serum in red top or serum separator tube. Overnight fast recommended. Store and transport refrigerated or frozen (on ice packs)	4-10 days

URINE

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CADD220	Add-on Urinalysis Includes chemical analysis (Color, Appearance, Specific Gravity, pH, Protein, Glucose, Ketones, Bilirubin, Occult blood) and microscopic evaluation of urine sediment A complete urinalysis which includes urine specific gravity, urine chemistry, and full microscopic evaluation. Interferences: visible levels of hemolysis, drugs containing dyes, nitrofurantoin, or riboflavin.	6.0 mL urine in urine transport tube	Daily Performed each shift
CADD230	Add-on Urine Protein/Creat Ratio Protein/Creatinine Ratio Urine Protein/Creatinine Ratio (UPC) is used to determine the magnitude of the proteinuria.	0.5 mL urine in urine transport tube	1-2 days
CS16735	Crystallographic Stone Analysis Nidus Composition, Stone Composition and Shell Composition The report includes a comprehensive stone analysis using optical crystallography with a region-by-region stone fracture. Stone composition is provided for nidus, stone, shell, and surface. Note: do not send in formalin.	Dry stone in red top tube or plain container	5-7 days
CT764	Path Review (Urinalysis) This test can be added onto a complete urinalysis and includes sediment evaluation by a clinical pathologist.	Not applicable	1-2 days
CT760 Add-on Equivalent CADD220	Urinalysis Includes chemical analysis (Color, Appearance, Specific Gravity, pH, Protein, Glucose, Ketones, Bilirubin, Occult blood) and microscopic evaluation of urine sediment The complete urinalysis includes a physical (color, appearance, USG), chemical (pH, protein, glucose, ketones, bilirubin, occult blood), and microscopic (WBC, RBC, casts crystals, bacteria, transitional and squamous epithelial cells, fat droplets) exam of the urine. Interferences: visible levels of hemolysis, drugs containing dyes, nitrofurantoin, or riboflavin.	6.0 mL urine in urine transport tube	Daily Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT925	Urinalysis and UPC Ratio Urinalysis and Urine Protein/Creatinine Ratio A complete urinalysis (CT760) and a urine protein to creatinine ratio (CT775).	6.5 mL urine in urine transport tube	1-2 days
CT765	Urinalysis Clearance Ratio Urinary Clearance Creatinine Ratios (Phosphorous, Calcium, Sodium, Potassium, Chloride) Determination of urine fractional excretion is indicated to assess renal clearance of specific substances. Interferences: marked hemolysis or lipemia in serum.	0.5 mL urine and 0.5 mL serum Serum in red top or serum separator tube, urine transport tube	1-2 days
CM133	Urinalysis/Urine Culture with MIC Urinalysis and Urine Culture and Sensitivity A complete urinalysis (CT760), and urine culture (CM130). Culture performed using <i>FIRST</i> Tract methodology allowing for rapid growth within 5 hours and subsequent reporting of a positive or negative result in under 24 hours. A positive result will be immediately plated for bacterial identification and antimicrobial sensitivity testing. Interferences: patients should be off antibiotics for seven to ten days prior to urine culture.	6.5 mL urine in urine transport tube	1-4 days
CT830C	Urine Microalbumin Canine Measures albumin concentration in the urine.	0.5 mL urine in urine transport tube	1-2 days
CT830F	Urine Microalbumin Feline Measures albumin concentration in the urine.	0.5 mL urine in urine transport tube	1-2 days
CT775 Add-on Equivalent CADD230	Urine Protein/Creatinine Ratio Urine Protein, Urine Creatinine and Urine Protein and Creatinine Ratio Urine Protein/Creatinine Ratio (UPC) to determine the magnitude of the proteinuria.	0.5 mL urine in urine transport tube	1-2 days

ENDOCRINOLOGY

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CADD260	<p>Add-on Fructosamine</p> <p>Serum fructosamine concentrations reflect the mean blood glucose concentrations during the preceding one to two weeks. The measurement can be used to differentiate stress hyperglycemia from diabetes mellitus.</p> <p>Interferences: hemolysis and lipemia may interfere with testing.</p>	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CADD300	<p>Add-on Post T4</p> <p>Used for monitoring thyroid status in patients previously diagnosed as hypothyroid, being supplemented on synthetic thyroid hormone (4-6 hours post pill).</p> <p>Note: refrigerate and send on ice. Ideal sampling time is typically 4-6 hours post pill administration.</p>	0.5 mL serum in red top or serum separator tube	1-2 days
CADD190	<p>Add-on T4</p> <p>Total T4</p> <p>Interferences: marked hemolysis. Moderate to marked lipemia.</p>	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CADD200	<p>Add-on TSH</p> <p>Used in conjunction with other thyroid hormone values, typically T4 and FT4, to allow for the diagnosis of hypothyroidism or hyperthyroidism.</p> <p>Interferences: marked lipemia.</p>	0.5 mL serum in red top or serum separator tube	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS14410	Anti-Mullerian Hormone (C/F) A qualitative AMH test for dogs and cats is meant to distinguish between spayed and ovarian intact bitches and queens once they have reached the age of sexual maturity. This test is to be used on canine and feline species only. Interferences: marked hemolysis or lipemia. Note: serum should be separated and placed in a plain red top tube. Sample must be kept cold — freeze after collection and send on ice packs.	0.5 mL serum in non-additive tube (plastic preferred) Note: Avoid use of serum separator tubes due to possible assay interference.	7-10 days
CT445	Cortisol Cortisol Assay, Single sample Single cortisol level Interferences: marked hemolysis may affect results, and marked lipemia (if unable to be cleared by centrifugation) may falsely decrease cortisol levels.	0.5 mL serum in red top or serum separator tube	1-2 days
CACTH2	Cortisol Serial ACTH 2 Cortisol Assay, 2 Samples: Pre (Baseline) and Post ACTH Stimulation This is an ACTH stimulation test that requires two samples: one pre, one post. Results include baseline cortisol and one post-ACTH cortisol. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube (fasted sample labeled Pre) and 0.5 serum in red top or serum separator tube (2-hour post sample labeled Post)	1-2 days
CACTH3	Cortisol Serial ACTH 3 Cortisol Assay, 3 Samples: Pre (Baseline) and 2 Post ACTH Stimulation This is an ACTH stimulation test that requires three samples. One pre cortrosyn administration and then two post at established intervals Interferences: marked hemolysis may affect results, and marked lipemia (if unable to be cleared by centrifugation) may falsely decrease cortisol levels.	0.5 mL serum in red top or serum separator tube for baseline sample (labeled Pre) and 0.5 serum in red top or serum separator tube for each post sample (post samples labeled Post)	1-2 days
CACTH4	Cortisol Serial ACTH 4 Cortisol Assay, 4 Samples: Pre (Baseline) and 3 Post ACTH Stimulation This is an ACTH stimulation test that requires four samples. One pre cortrosyn administration and then three post at established intervals Interferences: marked hemolysis may affect results, and marked lipemia (if unable to be cleared by centrifugation) may falsely decrease cortisol levels.	0.5 mL serum in red top or serum separator tube for baseline sample (labeled Pre) and 0.5 serum in red top or serum separator tube for each post sample (post samples labeled Post)	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CDEX2	<p>Cortisol Serial DEX 2</p> <p>Cortisol Assay, 2 Samples: Pre (Baseline) and Post Dexamethasone Suppression</p> <p>This test measures cortisol concentrations before and after administration of exogenous dexamethasone (0 hours followed by an 8-hour sample). This is only a two cortisol assay. The low dose dexamethasone suppression test, a three cortisol sample test (mnemonic CDEX3), is preferred for diagnosing hyperadrenocorticism. The high-dose dexamethasone suppression test, a three cortisol sample test, is used to differentiate pituitary-dependent from adrenal tumor hyperadrenocorticism.</p> <p>Interferences: marked hemolysis may affect results, and marked lipemia (unable to be cleared by centrifugation) may falsely decrease cortisol levels.</p>	<p>0.5 mL serum in red top or serum separator tube for baseline sample (labeled Pre) and 0.5 serum in red top or serum separator tube for post sample (labeled Post)</p>	<p>1-2 days</p>
CDEX3	<p>Cortisol Serial DEX 3</p> <p>Cortisol Assay, 3 Samples: Pre (Baseline) and 2 Post Dexamethasone Suppression</p> <p>This test measures cortisol concentrations before and after administration of exogenous dexamethasone. This suppression test requires three samples; baseline cortisol prior to dexamethasone administration and two cortisol levels post dexamethasone administration (4 and 8 hours).</p> <p>Interferences: marked hemolysis may affect results, and marked lipemia (unable to be cleared by centrifugation) may falsely decrease cortisol levels.</p>	<p>0.5 mL serum in red top or serum separator tube for baseline sample (labeled Pre) and 0.5 serum in red top or serum separator tube for each post sample (post samples labeledPost)</p>	<p>1-2 days</p>
CDEX4	<p>Cortisol Serial DEX 4</p> <p>Cortisol Assay, 4 Samples: Pre (Baseline) and 3 Post Dexamethasone Suppression</p> <p>This test measures cortisol concentrations before and after administration of exogenous dexamethasone. This suppression test includes four samples; baseline cortisol prior to dexamethasone administration and three additional cortisol samples post dexamethasone administration (last occurring at 8 hours). Samples taken at 2 or 6 hours are interpreted similarly to those obtained at 4 hours.</p> <p>Interferences: marked hemolysis may affect results, and marked lipemia (unable to be cleared by centrifugation) may falsely decrease cortisol levels.</p>	<p>0.5 mL serum in red top or serum separator tube for baseline sample (labeled Pre) and 0.5 serum in red top or serum separator tube for each post sample (post samples labeledPost)</p>	<p>1-2 days</p>
CT770	<p>Cortisol/Creatinine Ratio</p> <p>Urine Cortisol:Creatinine Ratio</p> <p>The urine cortisol/creatinine ratio is most commonly used as a screening test for canine hyperadrenocorticism and for hypoadrenocorticism.</p> <p>Interferences: marked hemolysis and lipemia.</p>	<p>0.5 mL urine in urine transport tube</p>	<p>2-4 days</p>

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT435	Endogenous ACTH Endogenous ACTH used for small animal testing. Note: endogenous ACTH is extremely labile. Immediately separated and frozen EDTA-plasma without aprotinin may be acceptable but is NOT the preferred sample. Serum or heparinized plasma is not acceptable for ACTH testing. Call lab to obtain an aprotinin-treated lavender top tube and transfer tube.	1.0 mL aprotinin treated EDTA plasma in non-additive transport tube (labeled as AP treated plasma). Transport refrigerated (on ice packs).	2-4 days
CT460 Add-on Equivalent CADD50	Free T4 by Equilibrium Dialysis Note: this test should not be performed as an Add-on to samples older than 72 hours.	0.5 mL serum in red top or serum separator tube	2-6 days
CADD50	Free T4 Equilibrium Dialysis Add-on	0.5 mL serum in red top or serum separator tube	2-6 days
CS16345 Add-on Equivalent CADD260	Fructosamine Interferences: hemolysis and lipemia may affect results.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CT470	Insulin Insulin with Glucose Interferences: marked hemolysis. Note: Generally used to evaluate for the presence of an insulinoma. Concurrent blood glucose concentration should be below 60 mg/dl (3.3 mmol/L) when serum insulin levels are measured. Not to be used on patients receiving exogenous insulin without an appropriate withdrawal period.	0.5 mL serum in red top tube or serum separator tube	2-3 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT470F	Insulin (Feline) with Glucose Interferences: marked hemolysis. Note: Generally used to evaluate for the presence of an insulinoma. Concurrent blood glucose concentration should be below 60 mg/dl (3.3 mmol/L) when serum insulin levels are measured. Not to be used on patients receiving exogenous insulin without an appropriate withdrawal period.	0.5 mL serum in red top or serum separator tube	5-7 days
CS16520	Luteinizing Hormone	0.5 mL serum in red top or serum separator tube	2-4 days
CS16596	Parathormone Related Protein Determines the presence of the parathormone-related protein (PTHrP). Interferences: marked hemolysis or lipemia. Note: submit frozen plasma and label submission tube as plasma PTHrP.	1.0 mL frozen EDTA plasma in plain transport tube, collected as whole blood in lavender top tube (spin and separate plasma into plain transport tube). Include patient history. Frozen plasma (Call customer service for specimen handling protocol: 1-800-341-3440, dial 0)	7-14 days
CT448	Post Cortisol A single cortisol level, submitted as a post, used to complete an ACTH stimulation or dexamethasone suppression test (8 hours). Interferences: marked lipemia.	0.5 serum in red top or serum separator tube (labeledPost)	1-2 days
CT497 Add-on Equivalent CADD300	Post T4 Used for monitoring thyroid status in patients previously diagnosed as hypothyroid, being supplemented on synthetic thyroid hormone (4-6 hours post pill). Note: draw serum 4-6 hours post pill.	0.5 mL serum in red top or serum separator tube	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT475	Progesterone Progesterone level Interferences: serum separator gel may interfere with test.	0.5 mL serum in red top tube (do not use serum separator tube)	1-2 days
CS16595	PTH/Ionized Calcium Parathormone and Ionized Calcium Evaluates the total calcium, ionized calcium, and parathyroid hormone to further define the cause of abnormalities in calcium homeostasis. Note: serum samples must be anaerobically obtained (exposure to air may artifactually decrease ionized calcium). Samples should not be transported in serum separator tubes. Transfer serum into two plain red top tubes and label as serum ionized calcium and serum PTH. Samples should be frozen and sent on ice.	2.0 mL serum anaerobically transferred from spun red top or serum separator tube into a plain, unopened red top tube stored and transported frozen along with patient history	7-10 days
CT480	T3 Tri-iodothyronine (T3)	0.5 mL serum in red top or serum separator tube	2-4 days
CT495	T4 Total T4 Interferences: marked hemolysis and moderate to marked lipemia. Lipemia can falsely decrease T4 results.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CT499	T4 (Free) Free T4 T4 (Free) determined by Immulite.	0.5 mL serum in red top or serum separator tube	1-2 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS16760	Testosterone Single testosterone level	0.5 mL serum in red top or serum separator tube	3-5 days
CT505	Thyroglobulin Autoantibodies Canine Thyroglobulin Auto Antibody Level Thyroglobulin auto antibody test is used to confirm the presence of autoimmune thyroiditis	0.5 mL serum in red top or serum separator tube	3-7 days
CSA360	Thyroid Profile 1 T4 and T3 Used to evaluate T4 and T3 levels. Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	0.75 mL serum in red top or serum separator tube	1-5 days
CSA370	Thyroid Profile 2 Total T4 and Free T4 By Equilibrium Dialysis Panel includes a total T4 and free T4 by Equilibrium Dialysis (FT4 ED). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum in red top or serum separator tube	1-6 days
CSA380	Thyroid Profile 3 Total T4, Free T4 By Equilibrium Dialysis, and TSH Panel includes a total T4, free T4 by Equilibrium Dialysis (FT4 ED), and thyroid stimulating hormone (TSH). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.5 mL serum in red top or serum separator tube	1-6 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA390	Thyroid Profile 4 TSH and Free T4 by Equilibrium Dialysis Panel includes a free T4 by Equilibrium Dialysis (FT4 ED), and thyroid stimulating hormone (TSH). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	1.0 mL serum in red top or serum separator tube	1-6 days
CSA400	Thyroid Profile 5 T4, Free T4 By Equilibrium Dialysis, TSH, and Thyroglobulin Auto Antibody Panel includes a total T4, free T4 by Equilibrium Dialysis (FT4 ED), thyroid stimulating hormone (TSH), and thyroglobulin autoantibody (TGAA). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results.	2.0 mL serum in red top or serum separator tube	1-6 days
CSA410	Thyroid Profile 6 T3, T4, Free T4 By Equilibrium Dialysis, T3 Autoantibodies, T4 Autoantibodies, TSH, and Thyroglobulin Auto Antibody Panel includes total T3, T3 autoantibody, total T4, free T4 by Equilibrium Dialysis (FT4 ED), T4 autoantibody, thyroid stimulating hormone (TSH), and thyroglobulin autoantibody (TGAA). Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results. Lipemia can falsely elevate TSH results. Note: T4 and FT3 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 and FT3 results. Lipemia can falsely elevate TSH results.	2.5 mL serum in red top or serum separator tube	1-7 days
CT510	TSH Thyroid stimulating hormone (TSH) Interferences: marked lipemia. Lipemia can cause falsely elevated results.	0.5 mL serum in red top or serum separator tube	1-2 days
Add-on Equivalent CADD200			

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

INFECTIOUS

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CAC100	Accuplex® Heartworm, <i>Borrelia burgdorferi</i> , <i>E. canis</i> , <i>Anaplasma phagocytophilum</i> This is a canine vector-borne disease screening test for Heartworm (Ag), Lyme disease (includes screening for antibodies against two C6 peptides), <i>E. canis</i> , and <i>A. phagocytophilum</i> .	0.5 mL serum in red top or serum separator tube	1-2 days
CADD04	Add-on Coccidioidomycosis Screen IgM Antibody, IgG Antribody, <i>Coccidioides</i> IgG Titer Serology by agar gel immunodiffusion (AGID) used most frequently to screen for exposure to <i>Coccidioides immitis</i> . Interferences: marked lipemia.	0.5 mL serum in red top or serum separator tube	5-7 days
CADD05	Add-on Ehrlichia canis Assesses for the presence of antibodies against <i>Ehrlichia canis</i> evaluated for by IFA. Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	2-4 days
CADD07	Add-on Feline Retroviral FeLV Antigen (ELISA), FIV Antibody Detection of Feline leukemia virus (FeLV) antigen by ELISA and Feline immunodeficiency virus (FIV) antibody by IFA. Interferences: marked hemolysis or lipemia may result in false positive FeLV antigen result.	0.5 mL serum in red top or serum separator tube	1-2 days
CADD70	Add-on Heartworm Antigen Heartworm Antigen This test detects antigen from the female heartworm. A detected result (positive) indicates infection with adult female heartworms (<i>Dirofilaria immitis</i>). A not detected result (negative) indicates that at the time of submission, no antigen was detected. Adult <i>Dirofilaria immitis</i> antigens will not be detected for 5-7 months following infection with L3 larval stages.	0.5 mL serum in red top or serum separator tube	1-2 days
CADD20	Add-on Hemobartonella <i>Mycoplasma haemofelis</i> Microscopic evaluation of a fresh blood smear for hemotropic <i>Mycoplasma</i> . Interferences: marked hemolysis.	1.0 mL EDTA whole blood in lavender top tube	Daily

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS14437	<i>Anaplasma phagocytophilum</i> AB Used to detect antibodies against <i>Anaplasma phagocytophilum</i> by IFA methodology in all other species except canine and equine. Interferences: marked hemolysis and lipemia may affect results.	0.5 mL serum or CSF in red top tube	4-7 days
CS16872	<i>Anaplasma phagocytophilum</i> AB Used to detect antibodies against <i>Anaplasma phagocytophilum</i> by IFA methodology. Exposure does not necessarily indicate that clinical signs are caused by infection. Canine only. Note: use on canine samples only. Equine samples should be submitted under CS16872E.	1.0 mL serum in red top or serum separator tube	2-6 days
CT530	<i>Brucella canis</i> - Screen This test is species-specific (dogs only). This is a screening test performed at 1:50 dilution. A titer is not performed. Note: canine only. Test will not detect antibodies against <i>B. suis</i> , <i>B. melitensi</i> , and <i>B. abortus</i> . Use the <i>Brucella</i> AGID test if a dog may have been infected with <i>Brucella</i> from other species (e.g., from cattle, sheep, goats, pigs, sheep). Do not use this test for export.	0.5 mL serum in red top tube or serum separator tube	2-4 days
CS16131	<i>Brucella</i> Screen Multiplex <i>Brucella</i> Screen Multiplex. Includes confirmatory testing with <i>Brucella</i> Slide Agglutination and AGID for samples with non-negative screen results An antibody detection for <i>Babesia canis</i> antigen via the Canine <i>Brucella</i> Multiplex Assay. Samples undergo sequential testing dependent upon initial results. If the Multiplex assay is positive, further testing is performed by Canine <i>Brucella</i> Slide Agglutination and AGID testing at no additional cost. Note: canine only. This test can detect antibodies against <i>B. canis</i> , <i>B. suis</i> , <i>B. melitensi</i> and <i>B. abortus</i> .	1.0 mL serum in red top tube or serum separator tube	7-10 days
CS16135	Calici Virus - IFA A serology test evaluating for IgG antibody levels by IFA suggesting Feline calicivirus exposure. Interferences: sample submitted in tubes with separator gel.	0.5 mL serum in red top tube. Submission in serum separator tube is not recommended.	10-15 days
CT535 Add-on Equivalent CADD04	Coccidioidomycosis, Screen and Titer Serology by agar gel immunodiffusion is used most frequently to screen for exposure to <i>Coccidioides</i> . This is a titer test for antibodies directed against <i>Coccidioides immitis</i> (aka Valley Fever). Interferences: lipemia.	0.5 mL serum in red top or serum separator tube	4-6 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT550	Cryptococcus Antigen Used as a diagnostic test for Cryptococcosis. Detects the presence of Cryptococcal antigen. Interferences: marked hemolysis and lipemia. Note: serum and CSF submitted in a plain tube (white top or red top) is the only acceptable specimens. Plasma (any type) is not an acceptable specimen for this test. Body cavity fluid (e.g., CSF) cannot be submitted in a tube containing anticoagulant (e.g., LT) because the anticoagulant may invalidate the test.	0.5 mL serum or fluid in red top tube or plain transport tube (without anticoagulant)	2-4 days
CT570 Add-on Equivalent CADD05	Ehrlichia canis This test detects antibodies (IgG) directed against <i>Ehrlichia canis</i> . Interferences: marked hemolysis and lipemia. Canine only. Test is not valid for export purposes.	0.5 mL serum in red top or serum separator tube	2-4 days
CS16900	Ehrlichiosis Serology Panel (Canine) <i>Ehrlichia canis</i> , <i>Anaplasma phagocytophilum</i> , <i>Neorickettsia</i> antibody detection Evaluates for the presence of antibodies indicating exposure to <i>Ehrlichia canis</i> , <i>Anaplasma phagocytophilum</i> , and <i>Neorickettsia risticii</i> . Interferences: marked hemolysis and lipemia.	2.0 mL serum in red top or serum separator tube	2-7 days
CT595	Feline Coronavirus Titer This test determines if a cat has antibody titers against Feline Coronavirus (FCV). Titers are determined at 1:400 and 1:1600. A result of less than 1:400 is not the same as a negative titer. A titer of less than 1:400 indicates that the cat does not have high antibody concentrations against FCV. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	2-4 days
CT630	Feline Heartworm Profile Feline Heartworm Antibody and Heartworm Antigen This test evaluates for antibodies consistent with heartworm exposure and antigens consistent with the presence of female adult heartworm. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	1-2 days
CSA260 Add-on Equivalent CADD07	Feline Retroviral FeLV Antigen ELISA, FIV Antibody Evaluation for the presence of FeLV antigen (ELISA) and FIV antibodies (IFA). Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	0.5 mL serum in red top or serum separator tube	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA262	Feline Serology FeLV Antigen ELISA, FIV Antibody, <i>Toxoplasma</i> Ab - IgG/IgM Panel that assesses possible infectious etiologies including FeLV (CT580), FIV (CT610), and <i>Toxoplasma gondii</i> (CT720) for the exhibited clinical signs. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum in red top or serum separator tube	1-5 days
CSA265	Feline Serology 1 FeLV Antigen ELISA, Feline Coronavirus Titer, FIV Antibody Panel that assesses possible infectious etiologies including FeLV (CT580), FIV (CT610), and feline coronavirus (CT595) for the exhibited clinical signs. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	0.75 mL serum in red top tube or serum separator tube	1-3 days
CSA270	Feline Serology 2 FeLV Antigen ELISA, FIV Antibody, Feline Coronavirus Titer, Cryptococcal Antigen, <i>Toxoplasma</i> Ab - IgG/IgM Panel that assesses possible infectious etiologies including FeLV (CT580), FIV (CT610), feline coronavirus (CT595), <i>Cryptococcus neoformans</i> (CT550), and <i>Toxoplasma gondii</i> (CT720) for the exhibited clinical signs. Interferences: marked hemolysis or lipemia may result in false positive results on the FeLV Antigen ELISA.	1.0 mL serum in red top or serum separator tube	1-5 days
CT585	FeLV - IFA This is an immunofluorescence assay run on blood or bone marrow smears to look for evidence that Feline Leukemia Virus (FeLV) has infected the cat's bone marrow. Note: requires unstained blood or bone marrow smears.	2 freshly prepared, unstained blood smears or 2 bone marrow smears.	1-3 days
CT580 Add-on Equivalent CADD06	FeLV Antigen (ELISA) Test detects FeLV antigen by ELISA methodology. Interferences: marked hemolysis or lipemia may cause false positive results.	0.5 mL serum in red top or serum separator tube	1-2 days
CADD06	FeLV Antigen (ELISA) Add-on Test detects FeLV antigen by ELISA methodology. Interferences: marked hemolysis or lipemia. marked hemolysis or lipemia may result in false positive results.	0.5 mL serum in red top or serum separator tube	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS6234	FeLV PCR Used to detect FeLV using PCR in cases where ELISA and IFA are inconclusive. Results are reported as positive or negative. Note: samples should be refrigerated. 1-5 grams of bone marrow in a leak-proof container or 1.0 mL of whole blood in a lavender top tube.	1 to 5 grams bone marrow in leakproof container OR 1.0 mL whole blood (EDTA) in lavender top tube	7-10 days
CT610 Add-on Equivalent CADD15	FIV Antibody Evaluation for the presence of FIV antibody by IFA.	0.5 mL serum in red top or serum separator tube	1-2 days
CADD15	FIV Antibody Add-on Evaluation for the presence of FIV antibody by IFA.	0.5 mL serum in red top or serum separator tube	1-2 days
CS86454	Heartworm (Microfilaria) for Export Heartworm Antigen test and Microfilaria filtration test Heartworm antigen detection by ELISA and heartworm microfilarial filtration test.	2.0 mL EDTA whole blood in lavender top tube	7-10 days
CT613	Heartworm Ag-Heat-Treated Serum Heartworm Antigen using heat-treated serum This test is to be done on HW negative samples that are thought to be false negative because of antigen/antibody complexes. This test is NOT recommended for samples that are already HW antigen-positive or borderline positive. Interferences: marked hemolysis or lipemia. Marked hemolysis may cause false positive result. Marked lipemia may cause false negative results.	1.0 mL serum in red top or serum separator tube	2-4 days
CT625	Heartworm Antibody, Feline This test evaluates for antibodies consistent with heartworm exposure. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	2-4 days
CT615 Add-on Equivalent CADD70	Heartworm Antigen This test evaluates for the presence of heartworm antigen. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT620	Heartworm Antigen Feline This test evaluates for the presence of heartworm antigen.	0.5 mL serum in red top or serum separator tube	1-2 days
CT618	Heartworm, Ova and Parasite Heartworm Antigen, Fecal O&P with Centrifugation Heartworm antigen detection and fecal analysis via zinc sulfate centrifugation/flotation for ova and parasite detection (T805).	0.5 mL serum and 5 grams feces Serum in red top or serum separator tube, Antech provided fecal container	1-2 days
CT380 Add-on Equivalent CADD20	Hemobartonella <i>Mycoplasma haemofelis</i> Microscopic evaluation of a fresh blood smear for hemotropic <i>Mycoplasma</i> . Interferences: marked hemolysis.	1.0 mL EDTA whole blood in lavender top tube	Daily
CS86096	Influenza - Canine Acute H3N8 and H3N2 influenza antibody detection Used to detect for the presence of antibody to H3N8 and H3N2 influenza viruses.	1.0 mL serum in red top tube	10-15 days
CS16510	Leptospirosis <i>L. pomona</i> , <i>L. icterohemorrhagiae</i> , <i>L. canicola</i> , <i>L. grippityphossa</i> , <i>L. hardjo</i> , <i>L. autumnalis</i> , <i>L. bratislava</i> The results of this test include a semi-quantitative titer for serovars. Interferences: recent vaccination (within one month) may interfere with testing. Note: test is not species specific.	1.0 mL serum in red top tube or serum separator tube	5-10 days
CT670	Lyme IgG Detection of IgG consistent with exposure (natural or vaccinal) to <i>Borrelia burgdorferi</i> . Interferences: marked hemolysis or lipemia. Note: test does not distinguish between natural and vaccinal Lyme exposure.	0.5 mL serum in red top or serum separator tube	2-4 days
CT390	Microfilaria This is a concentration test involving whole blood that is more sensitive than the examination of a peripheral blood smear to determine the presence of microfilariae. Interferences: marked hemolysis.	1.0 mL EDTA whole blood in lavender top tube.	1-2 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS16560	<i>Neospora caninum</i> - IFA Results are expressed as an antibody titer consistent with exposure to <i>Neospora caninum</i> .	1.0 mL serum in red top or serum separator tube	7-10 days
CT812	O&P, <i>Giardia</i> with Heartworm Heartworm Antigen, Fecal O&P with Centrifugation, <i>Giardia</i> Heartworm antigen detection and fecal analysis via zinc sulfate centrifugation/flotation and <i>Giardia</i> Ag ELISA for ova and parasite detection (T808). Note: sample should be evaluated within 24 hours of collection. If a worm has been identified in the sample, separate the worm and place in a container labeled Worm in black. Additionally, mark Worm on the TRF being submitted with the sample.	0.5 mL serum and 6 grams feces Serum in red top or serum separator tube, Antech provided fecal container	1-2 days
CT695	Parvovirus Antigen Detects parvovirus antigen in fecal samples. Note: modified live virus vaccination for parvovirus may give false-positive results for about two weeks after vaccination. A negative test does not rule out parvovirus infection.	5 grams feces in Antech provided fecal container	2-4 days
CS1204	Rabies Diagnostic Non-Export Results are reported as an antibody titer. Note: not to be used for export. In the case of export, the individual country's guidelines for testing need to be reviewed. With regards to Rabies vaccination guidelines and the utility of titers, the clinician is recommended to refer to state/provincial guidelines for vaccination requirements.	1.5 mL serum in red top or serum separator tube	3-4 weeks
CS17108	Rabies Export - FAVN Rabies Antibody Titer, Export by FAVN Evaluation of rabies titer by fluorescent antibody virus neutralization. Interferences: marked hemolysis or lipemia. Note: complete FAVN form. For the post rabies vaccination, wait at least 3-4 weeks prior to drawing the sample.	1.5 mL serum in red top or serum separator tube (with FAVN form)	3-4 weeks
CT720	<i>Toxoplasma</i> AB - IgG/IgM Detection of IgG and IgM antibodies to <i>Toxoplasma gondii</i> . Interferences: marked hemolysis and lipemia. This code is for domestic cats only. Alternative codes are available for canines and exotic animals (including felines).	0.5 mL serum in red top or serum separator tube	2-4 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS18708	<i>Toxoplasma gondii</i> PCR Used to evaluate for the presence of <i>Toxoplasma gondii</i> in cases where toxoplasmosis is suspected. Note: preferred specimens include whole blood, bronchoalveolar lavage, transtracheal wash, or lymph nodes aspirate.	0.5 mL whole blood (EDTA) OR transtracheal wash, OR BAL, OR lymph node aspirate	9-12 days
CS85030	Toxoplasmosis IgG/IgM <i>Toxoplasma</i> IgG, <i>Toxoplasma</i> IgM Detection of IgG and IgM antibodies to <i>Toxoplasma gondii</i> . Note: when submitting CSF, it is recommended that serum be submitted as well.	1.0 mL serum (canine) in red top or serum separator tube	7-10 days
CS85819	<i>Tritrichomonas</i> PCR Used to evaluate for the presence of <i>Tritrichomonas foetus</i> in the feces of a patient with the appropriate clinical signs. Note: feces should be diarrheic and not be contaminated with cat litter. Fecal samples collected using a fecal loop are preferred.	0.3 grams feces in Antech provided container	1-2 days
CS16581	Vaccine Panel (Panleukopenia, Rhino, Calici) Feline Panleukopenia, Rhinotracheitis, and Calicivirus Vaccinal Titers This is a semi-quantitative titer that is correlated to the amount of protective humoral antibody present in relation to previous feline panleukopenia, rhinotracheitis, and calicivirus vaccination.	2.5 mL serum in red top or serum separator tube	3-15 days
CADD250	Add-on <i>Giardia</i> (ELISA) The <i>Giardia</i> antigen capture Enzyme- linked Immunosorbent Assay (ELISA) is a fecal procedure designed to detect a <i>Giardia</i> specific antigen.	4 grams feces in Antech provided fecal container	1-2 days
CT808	<i>Giardia</i> and Ova & Parasite with Centrifugation Fecal ova and parasites and <i>Giardia</i> Ag detection by ELISA. O&P samples are appropriately mixed with zinc sulfate solution, centrifuged, followed by flotation and slide evaluation. <i>Giardia</i> ELISA is considered more sensitive than a single flotation test. Note: sample should be evaluated within 24 hours of collection. If a worm has been identified in the sample, separate the worm and place in a container labeled Worm in black. Additionally, mark Worm on the TRF being submitted with the sample.	6 grams fresh fecal specimen in Antech provided fecal container	1-2 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

MICROBIOLOGY

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CM010	Acid Fast Stain - Micro Acid-fast stain applied to an air-dried smear to determine if mycobacteria present in sample provided. Note: air dried smear is the preferred sample submission.	1 air-dried smear	1-3 days
CADD210	Add-on Urine Culture MIC Includes culture and sensitivity <i>FIRST</i> Tract is a rapid urine culture analysis that allows for rapid bacterial growth (within 5 hours) and subsequent reporting of a positive or negative result in under 24 hours. A positive result will be immediately plated for bacterial identification and antimicrobial sensitivity testing. Interferences: patient should be off antibiotics for seven to ten days prior to urine culture.	0.5 mL cystocentesis, clean catch or catheterized urine in sterile red top tube or urine transport tube, Culturette, red top or white top tube with fluid, or other sterile container. Samples collected in EDTA tube are not acceptable.	1-3 days Preliminary reports available every 24 hours. Canine urine culture final result available in 48 hours. Feline urine culture final result available in 72 hours. If a fastidious organisms is observed, the listed turnaround time may be extended.
CM020	Aerobic Culture and MIC Aerobic Culture and MIC sensitivity testing This test is used when an aerobic bacterial infection is suspected in a tissue or fluid. Interferences: 1. Patient should be off antibiotics for at least 7-10 days. 2. Fluid in lavender top tube is unacceptable (EDTA inhibits growth). Note: tissue should be submitted in an RTT with a few drops of saline (to keep moist). Samples should be refrigerated prior to transportation to the laboratory. A separate sample should be submitted for each type of culture needed.	Culturette from fluid (body cavity fluids), TTW, BAL, wound, lesion, or skin. Tissue sample in saline. Culturette, red top or white top tube with fluid, or other sterile container. Samples collected in EDTA tube are not acceptable.	3-4 days Preliminary report available every 24 hours. Final culture result available in 72 hours. If a fastidious organism is observed, the listed turnaround time may be extended.

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CM040	<p>Aerobic Culture and MIC and Anaerobic Culture</p> <p>Aerobic culture with MIC sensitivity testing and anaerobic culture without sensitivities</p> <p>This test is used when a bacterial infection is suspected in a tissue or fluid, but it is uncertain whether an aerobic or an anaerobic organism is the cause.</p> <p>Interferences:</p> <ol style="list-style-type: none"> 1. Patient should be off antibiotics for at least 7-10 days. 2. Fluid in lavender top tube is unacceptable (EDTA inhibits growth). <p>Note: fluids should be submitted in a red top tube. Tissue should be submitted in an air-tight container for anaerobic culture. Samples should be refrigerated prior to submission. Sensitivities are not performed, but drugs of choice for anaerobes are provided as a guide.</p>	<p>Culturette (aerobic and anaerobic) from fluid (body cavity fluids), TTW, BAL, wound, lesion, or skin. Tissue sample in saline.</p>	<p>3-5 days</p> <p>Preliminary report available every 24 hours. Final culture result available in 72 hours.</p>
CM030	<p>Anaerobic Culture</p> <p>Anaerobic culture with identification of organism(s)</p> <p>This test is used when an anaerobic bacterial infection is suspected in a tissue or fluid.</p> <p>Interferences: sample with exposure to air (or sample that has dried) will preclude accurate testing.</p> <p>Note: fluids should be submitted in an RTT. Tissue should be submitted in an air-tight, sterile container (e.g., urine cup). No saline is necessary. Sensitivities are not performed, but drugs of choice for anaerobes are provided as a guide.</p>	<p>Anaerobic culturette of fluid (body cavity fluids), TTW or BAL, wound, lesion. Tissue sample (at least 2 cm x 2 cm) in sterile, air tight container with saline.</p>	<p>3-5 days</p> <p>Preliminary report available 3rd day. Final culture result available 4th day</p>
CM060	<p>Blood Culture - Aerobic Only</p> <p>Aerobic bacterial culture and MIC sensitivity testing</p> <p>This test is used when bacteremia caused by an aerobic bacteria is suspected. Single sample culture.</p> <p>Note: submit in BACTEC aerobic culture bottle. Clip fur and scrub venipuncture site for aseptic collection. Aim to collect 1.0 mL for cats and small dogs, 2-3 mL for larger dogs to be inoculated into 20 mL bottle. Don't unscrew caps on bottles. Remove the protective top and wipe visible parts of the rubber stopper with 70% ethanol. Allow stopper to dry or wipe with sterile gauze. Replace the drawing needle with a sterile needle before puncturing the rubber stopper, fill until vacuum stops, then gently invert the bottle to mix. Anticoagulants in the media will prevent blood from clotting.</p>	<p>Whole blood collected in BD Bactec Blood Culture Bottle(BCB)</p>	<p>5-7 days</p>

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CM061	<p>Blood Culture - Aerobic/Anaerobic</p> <p>Bacterial culture and MIC sensitivity testing</p> <p>This test is utilized when bacteremia is suspected. Submit a single sample in two separate BACTEC bottles as required for aerobic and anaerobic culture.</p> <p>Interferences: antibiotics.</p> <p>Note: submit two BACTEC culture bottles (pink label and gold label). Clip fur and scrub venipuncture site for aseptic collection. Aim to collect 1.0 mL for cats and small dogs, 2-3 mL for larger dogs to be inoculated into BacTec Peds Plus pink cap/strip bottle. Anaerobic blood culture bottles with orange cap/gold strip require 3-10 mL of blood. Don't unscrew caps on bottles. Remove the protective top and wipe visible parts of the rubber stopper with 70% ethanol. Allow stopper to dry or wipe with sterile gauze. Replace the drawing needle with a sterile needle before puncturing the rubber stopper, fill until the vacuum stops, then gently invert the bottle to mix. Anticoagulants in the media will prevent blood from clotting.</p>	<p>Whole blood collected in two separate BD Bactec Blood Culture Bottle (BCB): 1 aerobic BD Bactec BCB and 1 anaerobic BD Bactec BCB</p>	<p>5-7 days</p>
CM062	<p>Serial Blood Culture X2</p> <p>Bacterial culture and MIC Sensitivity testing</p> <p>Utilize this test when bacteremia is suspected. Aerobic/ anaerobic cultures are obtained at two separate time points. Each single time point requires two separate BACTEC bottles as indicated for aerobic and anaerobic culture.</p> <p>Interferences: antibiotics.</p> <p>Note: submit two BACTEC culture bottles (pink label and gold label) for each time point. For each collection, clip fur and scrub venipuncture site for aseptic collection. Aim to collect 1.0 mL for cats and small dogs, 2-3 mL for larger dogs to be inoculated into BacTec Peds Plus pink cap/strip bottle. Anaerobic blood culture bottles with orange cap/ gold strip require 3-10 mL of blood. Don't unscrew caps on bottles. Remove the protective top and wipe visible parts of the rubber stopper with 70% ethanol. Allow stopper to dry or wipe with sterile gauze. Replace the drawing needle with a sterile needle before puncturing the rubber stopper, fill until the vacuum stops, then gently invert the bottle to mix. Anticoagulants in the media will prevent blood from clotting.</p>	<p>2 sets of aerobic and anaerobic Bactec blood culture bottles (one set for each single time point)</p>	<p>5-7 days</p>

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CM063	<p>Serial Blood Cultures X3</p> <p>Bacterial culture and MIC Sensitivity testing</p> <p>Utilize this test when bacteremia is suspected. Aerobic/ anaerobic cultures are obtained at three separate time points. Each single time point requires two separate BACTEC bottles as indicated for aerobic and anaerobic culture.</p> <p>Interferences: antibiotics.</p> <p>Note: submit two BACTEC culture bottles (pink label and gold label) for each time point. For each collection, clip fur and scrub venipuncture site for aseptic collection. Aim to collect 1.0 mL for cats and small dogs, 2-3 mL for larger dogs to be inoculated into BacTec Peds Plus pink cap/strip bottle. Anaerobic blood culture bottles with orange cap/ gold strip require 3-10 mL of blood. Don't unscrew caps on bottles. Remove the protective top and wipe visible parts of the rubber stopper with 70% ethanol. Allow stopper to dry or wipe with sterile gauze. Replace the drawing needle with a sterile needle before puncturing the rubber stopper, fill until the vacuum stops, then gently invert the bottle to mix. Anticoagulants in the media will prevent blood from clotting.</p>	<p>3 sets of aerobic and anaerobic Bactec blood culture bottles (one set for each single time point)</p>	<p>5-7 days</p>
CM050	<p>Culture (Aerobic & Fungal)</p> <p>Aerobic culture with MIC sensitivity testing and fungal culture</p> <p>This test is used when an aerobic bacterial infection or fungal infection is suspected in a tissue or fluid.</p> <p>Note: Aerobic culture — fluid (pus, body cavity, transtracheal wash, bronchoalveolar lavage) should be submitted in a red top tube (no additive). Copan swabs from wounds, lesions, and fluids are acceptable. The tissue should be submitted in a red top tube with a few drops of saline. Fungal culture: hair, sputum, skin, and body fluids can be submitted in a sterile red top vacutainer tube or any sterile container without additives. DTM plates are acceptable for fungal ID. However, please indicate the source of the specimen plated.</p>	<p>Aerobic culture: Culturette from body fluid, TTW, BAL, wound, lesion, or skin; tissue in saline. Fungal culture: Dry hair and skin scraping collected in a sterile red top or other sterile container without additive.</p>	<p>3-21 days</p>
CDTM	<p>Culture Dermatophytes</p> <p>Fungal culture on DTM and Sabourad dextrose agar.</p> <p>Use this test to investigate ringworm (dermatophyte fungal infections) as a cause for skin lesions.</p> <p>Note: submit hair (10 hairs), skin scrapings, and nails in a container. Submit pre-inoculated DTM agar. In cases of fuzzy growth on the agar plate, use unit code CM082.</p>	<p>Hair, nails, or skin scraping in plain container</p> <p>Plain container or DTM bottle</p>	<p>14-21 days</p> <p>Preliminary report every 7 days. Final report available at 21 days (3 weeks)</p>

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CM070	Culture ID This test is for bacterial culture and identification only. No sensitivity is performed.	Culturette from fluid (body cavity fluids), TTW, BAL, wound, lesion, skin, or environmental sample. Tissue sample in saline. Culturette, red top or white top tube with fluid, or other sterile container. Samples collected in EDTA tube are not acceptable.	2-4 days If a fastidious culture organism is observed, the listed turnaround time may be extended.
CM225	Culture, <i>Campylobacter</i> Use this test when a <i>Campylobacter</i> infection is suspected to be the cause of gastroenteritis that is not responding to conservative management. Interferences: antibiotic therapy, submission of a dried swab.	2.0 grams feces in sterile container or in <i>Campylobacter</i>-Thioglycollate broth (Campy-Thio broth), or fecal swab	3-5 days
CM080	Culture, Fungal Fungal culture and identification. Does not include sensitivity to antifungal medications. This test is for fungal culture and identification. Not for suspected ringworm infections. Note: submit hair, skin, nails, and bodily fluids on swabs or in sterile containers. For cases suspected of ringworm infection, use CM240 instead or consider ringworm PCR panel. If you are submitting a DTM plate with fuzzy growth, use code CM082.	Dry hair, nails, skin scraping, body fluid, or lesion material collected in a sterile red top or other sterile container without additive, or on a culturette. DTM bottle, culturette, container	21 days Preliminary report every 7 days. Final report available at 21 days (3 weeks), or upon identification
CM110	Culture, <i>Mycoplasma</i> Use this test when <i>Mycoplasma</i> infection is suspected. Interferences: antibiotics. Note: submit uterine, cervical, vaginal, lung tissue, mucus, semen, or fluids (tracheal wash, pericardial effusion, urine) in a sterile container or on a culturette/dry swab. Do not submit milk samples under this code (use 85694). Copan swabs can also be submitted but cannot be submitted in gel. Samples can also be submitted in Aimes transport media without charcoal.	Culturette from body fluid, TTW, semen, fresh urine. Uterine, cervical, vaginal or lung tissue. Blue swab	7-10 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CM130	Culture, Urine	0.5 mL cystocentesis, clean catch or catheterized urine in sterile red top tube or urine transport tube	1-3 days
	<p><i>FIRST</i>Tract Urine Culture followed by plating in the case of a positive result for bacterial identification and sensitivity</p> <p>All urine cultures deemed positive by <i>FIRST</i>Tract, a highly accurate and rapid urine culture technique, will then undergo traditional culture for identification and susceptibility.</p> <p>Interferences: patient should be off antibiotics for seven to ten days prior to urine culture.</p>	Urine transport tube, white top with urine	Preliminary reports available every 24 hours. Canine urine culture final result available in 48 hours. Feline urine culture final result available in 72 hours. If a fastidious organism is observed, the listed turnaround time may be extended.
CM100	Culture, Acid Fast Bacilli <p>Mycobacterial culture and acid fast staining only</p> <p>Use this test when a mycobacterial infection of fluid/tissue is suspected.</p> <p>Interferences: antimicrobials.</p>	Urine, CSF, sputum, or tissue in sterile container or Copan swab (without contact with gel) <p>Call customer service for specimen handling protocol: 1-800-341-3440, dial 0</p>	8 Weeks
CM125	Fecal Culture <p>Testing for <i>Salmonella</i>, Shigella, and <i>Campylobacter</i> as well as other pathogens, reported as positive or negative. Sensitivity provided for <i>Salmonella</i> and Shigella isolates.</p> <p>Culture specifically evaluates for <i>Salmonella</i>, Shigella and <i>Campylobacter</i> spp. Fecal PCR testing is more sensitive and tests for a broader array of potential pathogens.</p>	2.0 grams feces in Antech provided fecal container	3-4 days <p>Preliminary report every 24 hours. Final report available in 72 hours</p>
CM090	Gram Stain <p>This test is used to investigate the presence and morphology of bacteria, yeast, and fungi.</p> <p>Note: place specimen on a slide or material to prepare a slide. Store at room temperature.</p>	Prepared smear, copan swab or culturette from fluid (body cavity fluids, TTW or BAL), feces, urine transport tube, or wound or lesion material. <p>Culturette, slide holder</p>	1-2 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

PATHOLOGY

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CBMBC	Bone Marrow Core Biopsy The histopathological evaluation of a bone marrow core biopsy. A CBC submitted at the time of the bone marrow sampling is utilized in the interpretation of the bone marrow evaluation to allow for a more complete assessment of the hematopoietic system. A detailed history and time course of events will further aid the pathologist in the interpretation of the findings.	Bone marrow tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin.	3-5 business days
CBONE	Bone Marrow Cytology This test includes a microscopic evaluation of bone marrow. The report includes descriptions of cellularity, assessment of cell lineage, and pathologist comments. Note: for maximal diagnostic utility, a CBC should be submitted at the time of bone marrow sampling. Additionally, submission of a core biopsy along with an aspirate is necessary for a complete bone marrow evaluation. If submitting bone marrow aspirate slides and a core biopsy sample (formalin), ensure that the formalin sample is in a separate zip lock bag.	4 to 5 air dried bone marrow slides or bone marrow sample in lavender top tube Lavender top, slide holder	1-3 days
CT325	Buffy Coat for Mast Cells Buffy Coat Cytology for Mast Cells Used to assess for the presence of mast cells, evaluate for hemoparasites/infectious agents, or to do a differential on a complete blood count with low white cells.	1.0 mL EDTA whole blood in lavender top tube	1-3 days
CS86793	Canine Melanoma Diagnostic Panel	Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin.	2 Weeks

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS86792	Canine Melanoma Prognostic Panel	Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin.	2 Weeks
CCYTO	Cytology Includes preparation of submitted sample and microscopic interpretation by clinical pathologist Testing includes a microscopic evaluation of cells. The report includes cytologic interpretation, diagnosis, and comments regarding the etiology and biological behavior where applicable. Note: single source.	2 or more air-dried, unstained slides with clinical history Lavender top with fluid, red top with fluid, air-dried unstained slides	1-3 days
CDERM	Dermatopathology Consult Includes preparation of sample and microscopic interpretation by an anatomical pathologist This specialty biopsy is supported by a team of board-certified anatomic pathologists with a specialization in dermatopathology.	Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin.	5-7 days
CFLUA	Fluid Analysis Includes preparation of submitted sample, Cell Count (WBC & RBC), Specific Gravity, Protein measurement, and microscopic interpretation by clinical pathologist This is the appropriate test for the evaluation of fluids from a cavity or tissue space. The fluid will be characterized based on protein, specific gravity, white and red cells, and a pathologist will evaluate the fluid and provide a microscopic description and, when possible, diagnosis. Note: a comprehensive history should be included with the fluid analysis submission.	1.0 mL body cavity fluid in lavender or red top tube with 2 unstained smears prepared from fluid Red top with fluid, lavender top, air dried unstained slides	1-3 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CFBX	<p>Histopathology</p> <p>Preparation of submitted sample and microscopic interpretation by anatomic pathologist</p> <p>A boarded pathologist will evaluate the tissue submitted for Histopathology. A full written biopsy report will be provided, including source, history, description of submitted tissue, microscopic description, and diagnosis where possible. The pathologist's comments will include, where applicable, margin evaluation, grading, interpretation of special stains, and recommendations for further testing if needed.</p> <p>Note: firmly and evenly tighten the lid of the formalin jar and check for any leaks prior to placing the sample in a ziplock bag with the test requisition form. A comprehensive history should be included with each biopsy submission.</p>	<p>Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid</p> <p>To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin.</p>	<p>3-5 business days</p>
CS86569	<p>Histoplasma capsulatum Ag, EIA</p>	<p>2.0 mL urine in urine transport tube</p>	<p>7-9 days</p>
IHC1	<p>Immunohistochemistry 1 Stain</p> <p>This test includes a single immunohistochemistry stain. It's an antibody-based method to detect a specific protein.</p> <p>Note: typically an Add-on test after histopathological evaluation of submitted tissue based on pathologist recommendations.</p>	<p>Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid</p> <p>To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin.</p>	<p>10-14 days</p>
IHC2	<p>Immunohistochemistry 2 Stains</p> <p>This test includes two immunohistochemistry stains. It's an antibody-based method to detect specific proteins.</p> <p>Note: typically an Add-on test after histopathological evaluation of submitted tissue based on pathologist recommendations.</p>	<p>Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid</p> <p>To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin.</p>	<p>10-14 days</p>

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
IHC3	<p>Immunohistochemistry 3 Stains</p> <p>This test includes three immunohistochemistry stains. It's an antibody-based method to detect specific proteins.</p> <p>Note: typically an Add-on test after histopathological evaluation of submitted tissue based on pathologist recommendations.</p>	<p>Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid</p> <p>To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin.</p>	10-14 days
IHC4	<p>Immunohistochemistry 4 Stains</p> <p>This test includes four immunohistochemistry stains. It's an antibody-based method to detect specific proteins.</p> <p>Note: typically an Add-on test after histopathological evaluation of submitted tissue based on pathologist recommendations.</p>	<p>Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid</p> <p>To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin.</p>	10-14 days
IHC5	<p>Immunohistochemistry 5 Stains</p> <p>This test includes five immunohistochemistry stains. It's an antibody-based method to detect specific proteins.</p> <p>Note: typically an Add-on test after histopathological evaluation of submitted tissue based on pathologist recommendations.</p>	<p>Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid</p> <p>To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin.</p>	10-14 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CFBXORL	<p>Oral Path Biopsy</p> <p>Preparation of each submitted oral tissue sample and microscopic interpretation by anatomic pathologist</p> <p>A boarded pathologist with a special interest in oral cavity pathology will evaluate the tissue. A full written biopsy report will be provided, including source, history, description of submitted tissue, microscopic description, and diagnosis where possible. Additionally, a dentist will provide a treatment plan based on the histopathology or recommendations for further diagnostic evaluation.</p> <p>Note: firmly and evenly tighten the lid of the formalin jar and check for any leaks prior to placing the sample in a ziplock/ sealable bag with the test requisition form. A comprehensive history should be included with each biopsy submission.</p>	<p>Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid</p> <p>To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin.</p>	<p>3-5 business days</p>
CFBXNEO	<p>Pet Cancer Specialty Biopsy</p> <p>This specialty biopsy service is comprised of a team of six board-certified anatomic pathologists with specialization in cancer. The biopsy is reviewed by three of these pathologists including one report written with a consensus opinion.</p> <p>Note: firmly and evenly tighten the lid of the formalin jar and check for any leaks prior to placing the sample in a ziplock/ sealable bag with the test requisition form. A comprehensive history should be included with each biopsy submission.</p>	<p>Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid</p> <p>To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin.</p>	<p>3-5 business days</p>
CFBXLSP	<p>Superliver (Liver Biopsy with Liver Staining Panel)</p> <p>Liver Biopsy with Liver Staining Panel</p> <p>The tissue will be evaluated by a boarded pathologist with a special interest in hepatic disease. A full written biopsy report will be provided, including source, history, description of submitted tissue, microscopic description, and diagnosis where possible. Trichrome, rhodanine, and reticulum staining and interpretation are included. The report will include the staging of fibrosis, evaluation of the degree of parenchymal collapse, and a qualitative copper grade. Internal pathologist to pathologist opinions are sought when needed. There is an option for further discussion with an internal medicine consultant with a special interest in hepatic disease.</p> <p>Note: firmly and evenly tighten the lid of the formalin jar and check for any leaks prior to placing the sample in a ziplock/ sealable bag with the test requisition form. A comprehensive history should be included with each biopsy submission, including the breed, age, and time when liver values were first noted to be elevated. Additional diagnostics performed as well as responses to any specific treatments should be included in the history.</p>	<p>Liver tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid</p> <p>To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin.</p>	<p>3-5 business days</p>

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CFBXTBR	Tumor Board Review Biopsy Preparation of each submitted sample and microscopic interpretation by anatomic pathologist A second opinion review of an Antech or other lab biopsy case performed by three pet cancer specialty pathologists (Tumor Board), with one report containing consensus opinion of all three specialists. Note: firmly and evenly tighten the lid of the formalin jar and check for any leaks prior to placing the sample in a ziplock/ sealable bag with the test requisition form. A comprehensive history should be included with each biopsy submission.	Slides prepared from a previous histopathology sample submission. Formalin container	3-5 business days
CS11068 - SearchLight DNA® Histology Add-on CS10058 - SearchLight DNA® Cytology Add-on	SearchLight DNA® SearchLight DNA® is an Add-on test to any biopsy or cytology where cancer is suspected or diagnosed. This test provides precision diagnostic, prognostic, and therapeutic guidance for multiple types of cancers. This canine cancer genomic panel identifies mutations in 120 relevant cancer genes and provides insights into cancer's origin, behavior, and the optimal approach for treatment. Report includes: <ul style="list-style-type: none">• Information about mutations in 120 cancer genes including diagnostic, prognostic, and therapeutic biomarker associations with supporting evidence from peer-reviewed literature, clinical consensus, and inference from human FDA guidance.• Pharmacogenomic marker (MDR1) mutation status.• A list of targeted therapeutic drug(s) if the patient's tumor mutations are a match for an existing drug.• Information on clinical trials by tumor type curated from public databases and individual academic centers by Vidium Animal Health. To order SearchLight DNA®, please contact customer service. Note: SearchLight DNA® can only be added to a biopsy or cytology performed at Antech where neoplasia has been diagnosed or suspected.	Fine-Needle Aspirate (FNA) Samples: FNA Slides - ≥2 unstained or stained slides + 1 Diffquick FNA Liquid - ≥0.5 mL in sterile vial/tube with no additives Biopsy Samples: Formalin-Fixed, Paraffin-Embedded (FFPE) Scrolls – 10 x 10- micron sections with 1 adjacent hematoxylin and eosin (H&E) slide Unstained FFPE Slides – 10 x 10-micron sections with 1 adjacent H&E slide Bone Tumor – Non-decalcified sample	12-17 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

MOLECULAR

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT1025	<p>CADET® BRAF</p> <p>CADET <i>BRAF</i>. Includes CADET <i>BRAF</i>-PLUS when reflexed based on CADET <i>BRAF</i> result</p> <p>A non-invasive, highly sensitive assay that evaluates free-catch urine samples from dogs for the presence of cells harboring the b-raf mutation or specific copy number variations associated with Transitional Cell Carcinoma/ Urothelial Carcinoma/ Prostatic carcinoma. The assay identifies 95% of TCC/UC cases. It can be used to diagnose TCC/UC and monitor dogs undergoing treatment for remission and relapse.</p> <p>Interferences: cystocentesis drawn urine.</p> <p>Radiation and Chemotherapy are not interferences for the Cadet <i>BRAF</i>/<i>BRAF</i> Plus test assay; however, they may affect the tumor itself which may influence results.</p> <p>Note: urine sample should be free catch. Do not freeze.</p>	<p>40 mL urine collected in CADET BRAF urine container</p> <p>CADET <i>BRAF</i> container (urine must be put in <i>BRAF</i> container within 15 minutes of collection and can be collected over multiple days)</p>	3-6 days
CT980	<p>Canine Ehrlichia/Anaplasma PCR</p> <p><i>Anaplasma phagocytophilum</i>, <i>Anaplasma platys</i>, <i>Ehrlichia canis</i>, <i>Ehrlichia chaffeensis</i>, <i>Ehrlichia ewingii</i></p> <p>A comprehensive PCR panel used when the clinical signs are suggestive of Anaplasmosis or Erhlichiosis.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results.</p>	<p>1.0 mL whole blood, synovial fluid or CSF</p> <p>Lavender top, white top with fluid, CSF</p>	3-5 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT950	<p>Canine GI PCR Panel</p> <p><i>Campylobacter jejuni/coli</i>, canine enteric coronavirus, canine parvovirus, <i>Clostridium difficile</i> toxins A/B, <i>Clostridium perfringens</i> enterotoxin, <i>Cryptosporidium</i> spp., <i>Giardia</i> spp., <i>Salmonella</i> spp.</p> <p>Highly sensitive and specific panel used to evaluate for viral (canine enteric coronavirus, canine parvovirus, and canine rotavirus), bacterial (<i>Campylobacter</i> spp., <i>Clostridium</i> spp. and <i>Salmonella</i> spp.) and protozoal (<i>Giardia</i> spp. and <i>Cryptosporidium</i> spp.) causes for the gastrointestinal signs exhibited by the patient.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results for protozoal or bacterial organisms.</p> <p>Note: follow-up cultures are immediately performed on submissions found to be PCR-positive for <i>Salmonella</i> (when a Copan fecal swab is received by the laboratory). Additionally, samples found to be PCR-positive for DNA of <i>C. perfringens</i> enterotoxin or <i>C. difficile</i> toxins are immediately tested for these toxins by ELISA.</p>	<p>0.5 grams feces and 1 copan fecal swab</p> <p>Culturette, Antech provided fecal container</p>	3-5 days
CT953	<p>Canine GI PCR with SARS CoV-2 PCR Panel</p> <p>Canine GI PCR Panel, SARS-CoV-2 PCR Add-on Panel</p> <p>Canine GI PCR panel to which evaluation for SARS-CoV-2 PCR has been added to determine whether the gastrointestinal signs being exhibited could be related to the typical etiologic agents or potentially SARS-CoV-2.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results for protozoal or bacterial organisms.</p> <p>Note: follow-up cultures are immediately performed on submissions found to be PCR-positive for <i>Salmonella</i> (when a Copan fecal swab is received by the laboratory). Additionally, samples found to be PCR-positive for DNA of <i>C. perfringens</i> enterotoxin or <i>C. difficile</i> toxins are immediately tested for these toxins by ELISA.</p>	<p>0.5 grams feces and 1 copan fecal swab, conjunctival or nasal swab</p> <p>Culturette, Antech provided fecal container</p>	3-11 days
CS8710	<p>Canine Parvovirus PCR</p> <p>Detection of canine parvovirus in whole blood, feces or tissues.</p> <p>Note: vaccination for canine parvovirus within the previous 2-3 weeks may result in false positive results.</p>	<p>2.0 - 3.0 mL EDTA whole blood in lavender top tube OR 1-2 grams feces OR small tissue in red top</p>	7-10 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT995	<p>Canine Respiratory PCR Panel</p> <p>Canine influenza virus (H3N8), Pan Influenza (detects all Influenza Type A, including H3N2; when positive, confirmation by H3N2 specific qPCR), H5N1 influenza virus, canine adenovirus type 2, canine distemper virus, canine herpesvirus, canine parainfluenza virus, canine resp. coronavirus, <i>Bordetella bronchiseptica</i>, <i>Mycoplasma cynos</i>, <i>Streptococcus equi</i> subsp. <i>zooepidemicus</i></p> <p>A highly sensitive and specific panel covering the more common etiological agents causing acute respiratory signs in dogs. Includes Distemper, Parainfluenza, Respiratory coronavirus, Canine influenza (H3N8), Influenza (H5N1), PanInfluenza, Herpes, Adenovirus Type 2, <i>Bordetella bronchiseptica</i>, <i>Mycoplasma cynos</i>, and <i>Streptococcus zooepidemicus</i>.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results.</p>	<p>2 Sterile Swabs (nasal, deep pharyngeal or conjunctival) for PCR (submitted dry without transport media) and 1 Copan Swabin Gel for Bordatella culture. All swabs sampled from samesite.</p> <p>Culturette and 2 swabs in white tube</p>	3-5 days
CT998	<p>Canine Respiratory PCR with SARS CoV-2 PCR Panel</p> <p>Canine Respiratory PCR Panel(Canine influenza virus (H3N8), Pan Influenza (detects all Influenza Type A, including H3N2; when positive, confirmation by H3N2 specific qPCR), H5N1 influenza virus, canine adenovirus type 2, canine distemper virus, canine herpesvirus, canine parainfluenza virus, canine resp. coronavirus, <i>Bordetella bronchiseptica</i>, <i>Mycoplasma cynos</i>, <i>Streptococcus equi</i> subsp. <i>zooepidemicus</i>), SARS-CoV-2 PCR Add-on Panel</p> <p>Canine Respiratory PCR panel to which evaluation for SARS-CoV-2 PCR has been added to determine whether the respiratory signs being exhibited could be related to the typical etiologic agents (see T995) or potentially SARS-CoV-2.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results.</p>	<p>2 Sterile Swabs (nasal, deep pharyngeal or conjunctival) for PCR (submitted dry without transport media) and 1 Copan Swab in Gel for Bordatella culture. All swabs sampled from same site.</p> <p>Culturette and 2 swabs in white tube</p>	3-11 days
CT960	<p>Canine Tick Borne PCR Panel</p> <p><i>Anaplasma phagocytophilum</i>, <i>Anaplasma platys</i>, <i>Babesia canis</i>, <i>Babesia</i> sp. (Coco), <i>Babesia conradae</i>, <i>Babesia gibsoni</i>, <i>Bartonella henselae</i>, <i>Bartonella vinsonii</i>, <i>Ehrlichia canis</i>, <i>Ehrlichia chaffeensis</i>, <i>Ehrlichia ewingii</i>, <i>Mycoplasma haemocanis</i>, Cand. <i>M. haematoparuvum</i>, <i>Neorickettsia risticii</i></p> <p>A highly sensitive and specific panel used to evaluate for the evidence of infection in cases suspected of having a vector-borne disease. Includes <i>Anaplasma phagocytophilum</i>, <i>Anaplasma platys</i>, <i>Babesia canis</i>, <i>Babesia</i> spp. (non <i>canis</i>), <i>Bartonella henselae</i>, <i>Bartonella vinsonii</i>, <i>Ehrlichia canis</i>, <i>Ehrlichia</i> spp. (non-<i>canis</i>), <i>Mycoplasma heamocanis/</i> heamatoparvum, <i>Neorickettsia risticii</i>, and <i>Rickettsia rickettsii</i>.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results.</p>	<p>1.0 mL of EDTA whole blood in lavender top tube</p>	3-5 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT961	Canine Tick Borne PCR with Lyme Canine Tick Borne PCR Panel, Lyme Titer IgG Canine Tick Borne PCR panel (T960) to which Lyme IgG has been added. Interferences: previous or current antibiotic usage may interfere with PCR test results.	1.0 mL of EDTA whole blood and 0.2 mL of serum in red top tube Lavender top, red top or green top	3-5 days
CT982	Canine/Feline Ringworm PCR Panel <i>Microsporum</i> spp., <i>M. canis</i> , <i>M. gypseum</i> (<i>Arthroderma gypseum</i> , <i>A. fulvum</i> , <i>A. incurvatum</i>), <i>Trichophyton</i> spp., <i>T. mentagrophytes</i> (<i>Arthroderma benhamiae</i> , <i>A. vanbreusegheimii</i>) A highly sensitive and specific test used when clinical signs are suggestive of ringworm infection.	Minimum of 10 plucked hair with roots, skin scraping, or tooth brush sample in sterile, dry container free of liquids or preservative Hair, toothbrush, DTM	5-7 days
CT986	Canine/Feline Ringworm PCR with Dermatophyte Culture PCR for <i>Microsporum</i> spp., <i>M. canis</i> , <i>M. gypseum</i> (<i>Arthroderma gypseum</i> , <i>A. fulvum</i> , <i>A. incurvatum</i>), <i>Trichophyton</i> spp., <i>T. mentagrophytes</i> (<i>Arthroderma benhamiae</i> , <i>A. vanbreusegheimii</i>), and dermatophyte culture Evaluation for ringworm infection via both PCR and standard culture.	Minimum of 12 plucked hair with roots, skin scraping, or tooth brush sample in sterile, dry container free of liquids or preservative Hair, toothbrush, DTM	5-21 days
CT965	Feline Flea and Tick Borne PCR <i>Anaplasma phagocytophilum</i> , <i>Bartonella clarridgeiae</i> , <i>Bartonella henselae</i> , <i>Bartonella quintana</i> , <i>Ehrlichia</i> spp., <i>Mycoplasma haemofelis</i> , <i>Candidatus M. haemominutum</i> , <i>Candidatus M. turicensis</i> , <i>Rickettsia felis</i> , <i>Rickettsia rickettsii</i> A vector borne disease panel that evaluates for the presence of <i>Anaplasma phagocytophilum</i> , <i>Bartonella henselae</i> , <i>Bartonella clarridgeai</i> , <i>Bartonella quintana</i> , <i>Ehrlichia</i> spp, <i>Mycoplasma hemofelis</i> , <i>Mycoplasma haemominutum</i> , <i>Mycoplasma turicensis</i> , <i>Rickettsia rickettsii</i> and <i>Rickettsia felis</i> by PCR Interferences: previous or current antibiotic usage may interfere with PCR test results.	1.0 mL of whole blood in lavender top tube	3-5 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT955	<p>Feline GI PCR Panel</p> <p><i>Campylobacter coli</i>, <i>Campylobacter jejuni</i>, <i>Clostridium difficile</i> toxins A/B, <i>Clostridium perfringens</i> enterotoxin, <i>Cryptosporidium</i> spp. and <i>C. felis</i>, feline parvovirus, <i>Giardia</i> spp. <i>Salmonella</i> spp. <i>Tritrichomonas foetus</i></p> <p>Highly sensitive and specific panel used to evaluate for viral (feline panleukopenia), bacterial (<i>Campylobacter</i> spp, <i>Clostridium</i> spp toxins, and <i>Salmonella</i> spp.) and protozoal (<i>Giardia</i> spp., <i>Tritrichomonas foetus</i> and <i>Cryptosporidium</i> spp.) causes for the gastrointestinal signs exhibited by the patient by PCR.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results for protozoal or bacterial organisms.</p> <p>Note: follow-up cultures are immediately performed on submissions found to be PCR-positive for <i>Salmonella</i> (when a Copan fecal swab is received by the laboratory). Additionally, samples found to be PCR-positive for DNA of <i>C. perfringens</i> enterotoxin or <i>C. difficile</i> toxins are immediately tested for these toxins by ELISA.</p>	<p>0.5 grams feces and 1 copan fecal swab</p> <p>Culturette, Antech provided fecal container</p>	3-5 days
CT958	<p>Feline GI PCR with SARS CoV-2 PCR Panel</p> <p>Feline GI PCR Panel, SARS-CoV-2 PCR Add-on Panel</p> <p>Feline GI PCR panel (T955) to which evaluation for SARS-CoV-2 PCR has been added to determine whether the gastrointestinal signs being exhibited could be related to the typical etiologic agents or potentially SARS-CoV-2.</p> <p>Note: follow-up cultures are immediately performed on submissions found to be PCR-positive for <i>Salmonella</i> (when a Copan fecal swab is received by the laboratory). Additionally, samples found to be PCR-positive for DNA of <i>C. perfringens</i> enterotoxin or <i>C. difficile</i> toxins are immediately tested for these toxins by ELISA.</p>	<p>0.5 grams feces and 1 copan fecal swab</p> <p>Culturettes, Antech provided fecal container</p>	3-11 days
CT985	<p>Feline Hemoplasma PCR Panel</p> <p><i>Mycoplasma haemofelis</i>, <i>Candidatus M. haemominutum</i>, <i>Candidatus M. turicensis</i></p> <p>Used in the evaluation of a regenerative anemia where hemotropic <i>Mycoplasma</i> is a differential. PCR testing for the presence of <i>Mycoplasma haemofelis</i>, <i>Mycoplasma haemominutum</i>, and <i>Mycoplasma turicensis</i>.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results.</p>	<p>1.0 mL EDTA whole blood in lavender top tube</p>	3-5 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT990	<p>Feline Respiratory PCR Panel</p> <p>Feline calicivirus (FCV), Feline herpesvirus-1 (FHV-1), <i>Bordetella bronchiseptica</i>, <i>Chlamydomphila felis</i>, <i>Mycoplasma felis</i>, H1N1 influenza virus.</p> <p>A high sensitive and specific panel covering the more common etiological agents associated with upper respiratory signs in cats. Includes feline calicivirus, feline herpesvirus 1, <i>Chlamydia felis</i>, <i>Mycoplasma felis</i>, <i>Bordetella bronchiseptica</i>, Influenza H1N1.</p> <p>Interferences: previous or current antibiotic/antiviral usage may interfere with PCR test results.</p> <p>Note: complimentary culture/susceptibility testing on all submissions PCR-positive for <i>Bordetella bronchiseptica</i> (when the laboratory receives a Copan swab).</p>	<p>2 Sterile Swabs (nasal, deep pharyngeal or conjunctival) for PCR (submitted dry without transport media) and 1 Copan Swabin Gel for <i>Bordetella</i> culture. All swabs sampled from same site.</p> <p>Culturette and 2 swabs in white tube</p>	3-5 days
CT993	<p>Feline Respiratory PCR with SARS CoV-2 PCR Panel</p> <p>Feline Respiratory PCR Panel, SARS-CoV-2 PCR Add-on Panel</p> <p>Feline Respiratory PCR panel (T990) to which evaluation for SARS-CoV-2 PCR has been added to determine whether the respiratory signs being exhibited could be related to the typical etiologic agents or potentially SARS-CoV-2.</p> <p>Interferences: previous or current antibiotic/ antiviral usage may interfere with PCR test results.</p>	<p>2 Sterile Swabs (nasal, deep pharyngeal or conjunctival) for PCR (submitted dry without transport media) and 1 Copan Swabin Gel for <i>Bordetella</i> culture. All swabs sampled from same site.</p> <p>Culturette and 2 swabs in white tube</p>	3-11 days
CT600	<p>FIP mRNA PCR</p> <p>This PCR test detects mRNA of the M gene of all known feline coronavirus strains in any sample. However, for diagnosis of FIP, only the detection of mRNA outside of the intestinal tract is indicative since active replication of the virus in circulating mononuclear cells is typical for FIP.</p> <p>Note: most common sample submission is peritoneal effusion that on fluid analysis is consistent with FIP. Other sample options include CSF, tissue, and tissue aspirates (lymph node, etc.).</p>	<p>1.0 mL fresh fluid (body cavity, CSF, joint, milk, sputum, or urine) in EDTA (lavender top tube), EDTA whole blood, tissue samples in plain, non-additive tube (+/- saline), or unstained slides from aspirated tissue (i.e., lymph node, liver or kidney)</p>	5-10 days
CADD350	<p>KeyScreen® GI Parasite PCR Add-on</p> <p>Detects the following parasites by PCR: <i>Ancylostoma</i> spp. (if positive evaluated for benzimidazole resistance by PCR), <i>Uncinaria stenocephala</i>, <i>Toxocara/Toxascaris</i> spp. (if <i>Toxocara</i> spp. positive further speciated into <i>T. canis</i>, <i>T. cati</i>, <i>T. leonina</i>), <i>Baylisascaris procyonis</i>, <i>Trichuris vulpis</i>, <i>Dipylidium caninum</i>, <i>Echinococcus granulosus</i>, <i>Echinococcus multilocularis</i>, <i>Taenia</i> spp., <i>Giardia duodenalis</i>, <i>Giardia</i> zoonotic (if positive evaluated for potentially zoonotic strains A and B), <i>Cystoisospora</i> spp., <i>Eimeria</i> spp., <i>Cryptosporidium canis</i>, <i>Cryptosporidium felis</i>, <i>Toxoplasma gondii</i> (feline specific), <i>Neospora caninum</i>, <i>Tritrichomonas blagburni</i> (feline specific).</p>	<p>0.3 - 0.5 grams of fresh feces in Antech provided fecal container</p>	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT991	KeyScreen® GI Parasite PCR Panel	0.3 grams of fresh feces in Antech provided fecal container (minimum 0.15 grams)	1-2 days
Add-on Equivalent CADD350	Detects the following parasites by PCR: <i>Ancylostoma</i> spp. (if positive evaluated for benzimidazole resistance by PCR), <i>Uncinaria stenocephala</i> , <i>Toxocara/Toxascaris</i> spp. (if <i>Toxocara</i> spp. positive further speciated into <i>T. canis</i> , <i>T. cati</i> , <i>T. leonina</i>), <i>Baylisascaris procyonis</i> , <i>Trichuris vulpis</i> , <i>Dipylidium caninum</i> , <i>Echinococcus granulosus</i> , <i>Echinococcus multilocularis</i> , <i>Taenia</i> spp., <i>Giardia duodenalis</i> , <i>Giardia</i> zoonotic (if positive evaluated for potentially zoonotic strains A and B), <i>Cystoisospora</i> spp., <i>Eimeria</i> spp., <i>Cryptosporidium canis</i> , <i>Cryptosporidium felis</i> , <i>Toxoplasma gondii</i> (feline specific), <i>Neospora caninum</i> , <i>Tritrichomonas blagburni</i> (feline specific).		
	Detects 20 individual parasite species.		
	Hookworms:		
	<ul style="list-style-type: none"><i>Ancylostoma caninum</i><i>Uncinaria stenocephala</i><i>Ancylostoma</i> benzimidazole resistance		
	Roundworms:		
	<ul style="list-style-type: none"><i>Toxocara</i> spp.<i>Toxocara canis</i><i>Toxocara cati</i><i>Toxascaris leonina</i><i>Baylisascaris procyonis</i>		
	Whipworms:		
	<ul style="list-style-type: none"><i>Trichuris vulpis</i>		
	Tapeworms:		
	<ul style="list-style-type: none"><i>Dipylidium caninum</i><i>Echinococcus granulosus</i><i>Echinococcus multilocularis</i><i>Taenia</i> spp.		
	<i>Giardia</i>		
	<ul style="list-style-type: none"><i>Giardia duodenalis</i><i>Giardia</i> zoonotic (strain A or B)		
	Coccidia		
	<ul style="list-style-type: none"><i>Cystoisospora</i> spp.<i>Eimeria</i> spp.		
	Additional Protozoa:		
	<ul style="list-style-type: none"><i>Cryptosporidium canis</i><i>Cryptosporidium felis</i><i>Toxoplasma gondii</i> (feline specific)<i>Neospora caninum</i><i>Tritrichomonas blagburni</i> (feline specific)		
	If detected, <i>Toxocara</i> spp. will be identified at the species level as <i>T. cati</i> , <i>T. canis</i> , and <i>T. leonina</i> . <i>Ancylostoma caninum</i> will be evaluated for benzimidazole resistance, and <i>Giardia</i> evaluated for the potentially zoonotic strains A or B.		

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT974	<p>Leptospira PCR Blood</p> <p>Canine <i>Leptospira</i> PCR - Blood: <i>Leptospira interrogans</i> serovars <i>icterohaemorrhagiae</i>, <i>canicola</i>, <i>pomona</i>, <i>australis</i>, <i>bratislava</i>, <i>autumnalis</i>, <i>ballum</i> and <i>pyrogenes</i>, <i>L. kirschneri</i> serovar <i>grippotyphosa</i>, <i>L. interrogans/borgpetersenii</i> serovar <i>sejroe</i></p> <p>A highly sensitive and specific test used to evaluate patients suspected of having acute leptospirosis. <i>Leptospira</i> bacteremia is brief and is followed by bacteriuria. In cases where a patient has been ill anywhere from five days to two weeks, submission of blood and urine PCR will increase diagnostic sensitivity.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results.</p>	0.5 mL EDTA whole blood in lavender top tube	3-5 days
CT978	<p>Leptospira PCR Blood/Urine</p> <p>Canine <i>Leptospira</i> PCR Blood/Urine</p> <p>A highly sensitive and specific test used to evaluate patients suspected of having leptospirosis. <i>Leptospira</i> bacteremia is followed by bacteriuria. In cases where a patient has been ill anywhere from five days to two weeks, submission of blood and urine PCR will increase diagnostic sensitivity.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results.</p>	<p>0.5 mL EDTA whole blood and 2.0 mL urine</p> <p>Lavender top or blue top or green top, urine</p>	3-5 days
CT976	<p>Leptospira PCR Urine</p> <p>Canine <i>Leptospira</i> PCR - Urine: <i>Leptospira interrogans</i> serovars <i>icterohaemorrhagiae</i>, <i>canicola</i>, <i>pomona</i>, <i>australis</i>, <i>bratislava</i>, <i>autumnalis</i>, <i>ballum</i> and <i>pyrogenes</i>, <i>L. kirschneri</i> serovar <i>grippotyphosa</i>, <i>L. interrogans/borgpetersenii</i> serovar <i>sejroe</i></p> <p>A highly sensitive and specific test used to evaluate patients suspected of having leptospirosis. <i>Leptospira</i> bacteremia is followed by bacteriuria. Urine usually becomes positive day 7-10 post-exposure in clinically ill patients.</p> <p>Interferences: previous or current antibiotic usage may interfere with PCR test results.</p>	2.0 mL urine in urine transport tube,	3-5 days
CS14493	<p>OncoK9®</p> <p>OncoK9®- Cancel Signal</p> <p>OncoK9® is a non-invasive multi-cancer early detection (MCED) liquid biopsy test for the detection and characterization of cancer-associated genomic alterations in DNA isolated from canine blood.</p> <p>Interferences: refrigeration, freezing; pregnancy; clots; hemolysis; trauma/surgery.</p>	Specific collection kit for this test is required and can be ordered on Antech Online. Two cell-free DNA collection tubes, each with a minimum fill of 7 mL of whole blood, total of 14-17 mL whole blood.	12-17 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS85503	<p>PARR</p> <p>PARR (PCR for Antigen Receptor Rearrangement) is an assay that amplifies DNA. This test helps differentiate between monoclonal cells (most consistent with neoplasia) or polyclonal cells (most consistent with a reactive process).</p> <p>Note: PARR cannot be run on formalin-fixed, paraffin-embedded samples/slides or slides that have coverslips attached to them. Non- diagnostic samples are obtained when there is too little DNA or in the presence of an inhibitor of the PCR reaction.</p>	<ol style="list-style-type: none">1. Blood and bone marrow: 0.5 mL EDTA whole blood or bone marrow sample in lavendertop tube. 4-5 very cellular bone marrow slides can also be submitted. Peripheral blood smears are not acceptable. If available, please provide a complete blood count withpathology review or cytology/histology report of the bonemarrow.2. Lymph node and other organ aspirates: 4-5 cellular slides (smears or cytopspins) or fluid in EDTA(lavender top tube). No additional fluid (i.e.,saline/serum) needs to be added to aspirated fluid. Gluedcover-slipped slides cannot be used. Stained cytology slidesare acceptable.If available, please provide a copy of relevant cytology and/or fluid analysis with submission.3. Body cavity fluid: 4-5 cellular slides (smears or cytopspins) or fluid in EDTA(lavender top tube).Stained cytology slides are acceptable. If available, please provide a copy of relevant cytologyand/or fluid analysis with submission.4. CSF:Multiple cytopspin preparations are preferred (stained orunstained) or fluid in EDTA (lavender top tube). For the sample to be diagnostic, it is estimated that 50,000 lymphoid cells are needed. If the CSF has a lymphocyte count of 100/uL, at least 0.5 mL of fluid is needed, or the cells from 0.5 mL spun onto slides. If available, please provide a copy of relevant cytologyand/or fluid analysis with submission.5. Formalin-Fixed Paraffin Embedded (FFPE) Tissue:For this sample type, use code S86965 (Feline IntestinalLymphoma Panel performed at MSU) or S85562 (Lymphoma PanelPARR, performed at MSU). If specifically requesting thatPARR testing be performed on FFPE tissue at CSU, please use test code S85503 and specify FFPE tissue for CSU PARR.	7-9 business days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT984	Rabbit/Small Rodent Ringworm PCR Panel <i>Microsporum</i> spp., <i>M. equinum</i> / <i>canis</i> , <i>M. gypseum</i> (<i>Arthroderma gypseum</i> , <i>A. fulvum</i> , <i>A. incurvatum</i>), <i>Trichophyton</i> spp., <i>T. benhamiae</i> / <i>mentagrophytes</i> (<i>Arthroderma benhamiae</i> , <i>A. vanbreuseghemii</i>) A highly sensitive and specific test used when clinical signs are suggestive of ringworm infection.	Minimum of 10 plucked hair with roots, skin scraping, or tooth brush sample in sterile, dry container free of liquids or preservative	5-7 days
CT988	Rabbit/Small Rodent Ringworm PCR with Dermatophyte Culture <i>Microsporum</i> spp., <i>M. equinum</i> / <i>canis</i> , <i>M. gypseum</i> (<i>Arthroderma gypseum</i> , <i>A. fulvum</i> , <i>A. incurvatum</i>), <i>Trichophyton</i> spp., <i>T. benhamiae</i> / <i>mentagrophytes</i> (<i>Arthroderma benhamiae</i> , <i>A. vanbreuseghemii</i>), Culture, Dermatophytes Evaluation for ringworm infection via both PCR and standard culture.	Minimum of 12 plucked hair with roots, skin scraping, or tooth brush sample. For dermatophyte culture can also submit preinoculated DTM bottle.	5-21 days
CT996	SARS-CoV-2 PCR Add-on Panel This test may be added to either a GI PCR or Respiratory PCR panel, where history may suggest that SARS-CoV-2 may be playing an etiological role. Note: not a standalone test. Must be an Add-on to Canine/Feline Respiratory or GI PCR.	Conjunctival and deep pharyngeal swabs, nasal swab when combined with 2nd swab, nasopharyngeal or oropharyngeal aspirates or washes, bronchoalveolar lavage, tracheal aspirates, or fecal sample. Culturette, fecal sample	3-11 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS11068 - SearchLight DNA® Histology Add-on CS10058 - SearchLight DNA® Cytology Add-on	<p>SearchLight DNA®</p> <p>SearchLight DNA® is an Add-on test to any biopsy or cytology where cancer is suspected or diagnosed. This test provides precision diagnostic, prognostic, and therapeutic guidance for multiple types of cancers.</p> <p>This canine cancer genomic panel identifies mutations in 120 relevant cancer genes and provides insights into cancer’s origin, behavior, and the optimal approach for treatment.</p> <p>Report includes:</p> <ul style="list-style-type: none">• Information about mutations in 120 cancer genes including diagnostic, prognostic, and therapeutic biomarker associations with supporting evidence from peer-reviewed literature, clinical consensus, and inference from human FDA guidance.• Pharmacogenomic marker (MDR1) mutation status.• A list of targeted therapeutic drug(s) if the patient’s tumor mutations are a match for an existing drug.• Information on clinical trials by tumor type curated from public databases and individual academic centers by Vidium Animal Health. <p>To order SearchLight DNA®, please contact customer service.</p> <p>Note: SearchLight DNA® can only be added to a biopsy or cytology performed at Antech where neoplasia has been diagnosed or suspected.</p>	<p>Fine-Needle Aspirate (FNA) Samples:</p> <p>FNA Slides - ≥2 unstained or stained slides + 1 Diffquick FNA Liquid - ≥0.5 mL in sterile vial/tube with no additives</p> <p>Biopsy Samples:</p> <p>Formalin-Fixed, Paraffin-Embedded (FFPE) Scrolls – 10 x 10- micron sections with 1 adjacent hematoxylin and eosin (H&E) slide Unstained FFPE Slides – 10 x 10-micron sections with 1 adjacent H&E slide</p> <p>Bone Tumor – Non-decalcified sample</p>	12-17 days
CS14497 Add-on Equivalent CS14515	<p>Canine Wisdom Panel™ Premium</p> <p>Canine Wisdom Panel™ Premium DNA panel including:</p> <ul style="list-style-type: none">• 268 genetic health-associated variants (inc. MDR1 mutation, vWD)• 50+ physical trait variants• Breed background detection (350+ breeds)• Genetic diversity (heterozygosity) scoring• Genetic relatives matching with other dogs in the database <p>Support with interpretation available from Wisdom Panel</p>	1 Canine Wisdom swab kit (2 swabs per kit)	14-21 days
CS14498 Add-on Equivalent CS14516	<p>Feline Wisdom Panel™ Complete</p> <p>Feline Wisdom Panel™ Complete DNA panel including:</p> <ul style="list-style-type: none">• 45 genetic health condition-related variants (inc. MDR1 mutation, Polycystic Kidney Disease)• 25 physical traits variants• Identifies genetic blood type• Provides breed (mix) identification (70+ breeds)• Genetic diversity (heterozygosity) scoring, and genetic relatives matching with other cats in the database <p>Support with interpretation available from Wisdom veterinarians.</p>	1 Feline Wisdom swab kit (2 swabs per kit)	14-21 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS14515	Canine Wisdom Panel™ Premium, Add-on Canine Wisdom Panel™ Premium DNA panel including: <ul style="list-style-type: none">• 268 genetic health-associated variants (inc. MDR1 mutation, vWD)• 50+ physical trait variants• Breed background detection (350+ breeds)• Genetic diversity (heterozygosity) scoring• Genetic relatives matching with other dogs in the database Support with interpretation available from Wisdom Panel	1 Canine Wisdom swab kit (2 swabs per kit)	14-21 days
CS14516	Feline Wisdom Panel™ Complete, Add-on Feline Wisdom Panel™ Complete DNA panel including: <ul style="list-style-type: none">• 45 genetic health condition-related variants (inc. MDR1 mutation, Polycystic Kidney Disease)• 25 physical traits variants• Identifies genetic blood type• Provides breed (mix) identification (70+ breeds)• Genetic diversity (heterozygosity) scoring, and genetic relatives matching with other cats in the database Support with interpretation available from Wisdom veterinarians.	1 Feline Wisdom swab kit (2 swabs per kit)	14-21 days
CKS14497	Canine Wisdom Panel™ with Keyscreen® Canine Wisdom Panel™ Premium and 2x KeyScreen GI Parasite PCR Panel For Canine only. Wisdom Panel and first Keyscreen must be submitted at the same time. Follow up submissions must be for the same dog and include matching patient information. The follow up Keyscreen can be submitted within 365 days of the original. Call customer service for submission details: 1-800-341-3440, dial 0.	1 Canine Wisdom swab kit (2 swabs per kit), 0.3 grams of fresh feces in Antech provided fecal container	1-21 days 1-3 days for KeyScreen, 14-21 days for Wisdom Panel

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

GENERAL

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS16070	<i>Babesia canis</i> Detects the presence of antibody to <i>Babesia canis</i> . Performed by IFA.	0.5 mL serum in red top or serum separator tube	7-10 days
CS16075	<i>Babesia gibsoni</i> Detects the presence of antibody to <i>Babesia gibsoni</i> . Performed by IFA.	0.5 mL serum in red top or serum separator tube	7-10 days
CS16502	<i>Babesia gibsoni</i> IFA, Export This test evaluates the presence of antibodies to <i>Babesia gibsoni</i> , as required for export to specific countries. Ensure export test requirements prior to test submission.	1.0 mL serum in red top or serum separator tube	7-10 days
CS85889	<i>Bartonella henselae</i>, ELISA Evaluation of antibody indicating previous exposure to <i>Bartonella henselae</i> or <i>Bartonella clarridgeiae</i> (ELISA methodology). Note: recently infected cats may not have seroconverted. Titer magnitude does not prove current infection or relevant clinical disease. Active infection is suggested by a rising titer (a four-fold increase over two weeks).	0.5 mL serum in red top or serum separator tube	5-7 days
CS16005	Acetylcholine Receptor Antibody Assessment for the presence of antibodies against the acetylcholine receptor. A positive result supports a diagnosis of generalized acquired Myasthenia Gravis. Note: immunosuppressive therapy with corticosteroids instituted prior to submission of results may lower antibody levels and give a false negative result.	1.0 mL serum in red top or serum separator tube	9-11 days
CADD270	Add-on Distemper/Parvo Vacc Titer Distemper Vaccinal Titer, Parvovirus Vaccinal Titer This test detects IgG antibodies against canine parvovirus and canine distemper virus using ELISA. Results are reported as positive (protective) or negative (insufficient neutralizing antibody detected). Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	3-5 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CADD280	Add-on Panleukopenia Panleukopenia Vaccine Titer The test detects IgG antibodies against Panleukopenia via Immunofluorescent antibody assessment (IFA).	0.5 mL serum in red top or serum separator tube	3-5 days
CADD130	Add-on Protein Electrophoresis Total protein, Albumin, Globulin, Alpha 1, Alpha 2, Beta 1, and Gamma 1 fraction assessment with interpretation. An evaluation of the globulin fraction of the serum (alpha 1, alpha 2, beta, and gamma) to determine if the globulin fraction is monoclonal based on these components. Interferences: hemolysis may interfere with results, particularly the Beta globulin fraction.	0.5 mL serum in red top or serum separator tube	2-4 days
CT010	Albumin Interferences: marked hemolysis or lipemia. marked lipemia can cause false elevations. Oxalate anticoagulants interfere with albumin assessment.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CT020	Alkaline Phosphatase Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CT215	Alkaline Phosphatase Isoenzymes Alkaline phosphatase, steroid induced alkaline phosphatase, percent steroid induced alkaline phosphatase This evaluates the percentage of steroid-induced alkaline phosphatase in the serum relative to the total serum alkaline phosphatase. Elevation of steroid-induced alkaline phosphatase is not specific for hyperadrenocorticism and is of limited diagnostic utility in the diagnosis of this disease.	0.5 mL serum in red top or serum separator tube	2-4 days
CT030	ALT (SGPT) Alanine aminotransferase Interferences: marked hemolysis or lipemia. Serum or plasma should be separated within 1 hour of draw.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CT040	Amylase Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CT050	Amylase and Lipase (PSL) Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT515	ANA (Antinuclear Antibodies) Titer Then Antinuclear Antibody test is a titer test used as part of the Systemic Lupus Erythematosus (SLE) evaluation in patients with appropriate history, clinical signs, and other laboratory findings. Positive ANA titers are most commonly seen in pets with Systemic Lupus Erythematosus. However, positive results can be seen sometimes in pets with other auto-immune disorders (e.g., Rheumatoid arthritis), inflammatory/infectious disease, or neoplasia. Many medications are also known to cause positive ANA results. Canine and feline samples only. Interferences: marked hemolysis and lipemia. Note: canine and feline only.	0.5 mL serum in red top or serum separator tube	2-4 days
CS16040	Anticoagulant Screen Use test when ingestion of warfarin, bromadiolone, coumachlor, brodifacoum, diphacinone, chlorophacinone, or difethialone is suspected. Note: send out test requiring serum, stomach contents, bait, or frozen liver tissue.	5.0 mL serum, 20 grams stomach contents or bait, or 10 grams of frozen unfixed liver	7-14 days
CT060	AST (SGOT) Aspartate aminotransferase Interferences: marked hemolysis and lipemia. Serum or plasma should be separated within 1 hour of draw and centrifuged.	0.5 mL serum in red top or spun serum separator tube	Daily Performed each shift
CSA170	Autoimmune Profile 1 CBC, Antinuclear Antibodies (ANA), and Direct Coombs' Test Warm An autoimmune panel evaluating a complete blood count, Direct Coombs test warm, and Antinuclear antibodies. Interferences: marked hemolysis and lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	2-5 days
CT840	B-12 Cobalamin level. Generally used as part of a larger panel when evaluating for malabsorption/maldigestion. Note: ideally freeze and ship on ice. Otherwise submit fresh on ice.	0.5 mL serum in red top or serum separator tube	2-4 days
CS16195	B12, Folate Cobalamin and folate level. Use this test to evaluate intestinal diseases in dogs and cats. Additional measurement of TLI is recommended to ensure the most accurate interpretation of results. Note: pets should fast overnight prior to sampling. Hemolysis can affect folate results. Previous vitamin supplementation will cause elevated concentrations. Send/ship test on ice packs.	0.5 mL serum in red top or serum separator tube	2-5 days
CT115	Bicarbonate Interferences: marked hemolysis or lipemia. Serum stable for 1 hour in an ice bath prior. Due to the instability of CO ₂ in specimens, low levels may not be accurate. Most accurately done within minutes after drawing and on blood gas analyzer.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT070	Bilirubin, Direct Interferences: marked hemolysis or lipemia. Lipemia can falsely increase results.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CT090	Bilirubin, Total Interferences: marked hemolysis or lipemia. Lipemia can falsely increase results.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CT520	Bladder Tumor Analytes This test may be used as an adjunct in the diagnosis of transitional cell carcinomas (TCC). The test is a reasonably sensitive (~90%) but not a specific test for TCC. Consider CADET BRAF (T1025). Interferences: false-positive results may occur in the presence of proteinuria, aciduria, urinary tract infection, and hematuria. An alternative test to consider is the CADET BRAF if looking for a non-invasive diagnostic test with a far higher sensitivity (95%) and specificity (99%).	2.0 mL urine in urine transport tube,	2-4 days
CT525	Blastomyces Antibody To investigate the possibility of blastomycosis as the cause of clinical signs. Interferences: marked hemolysis or lipemia. Note: test is reported as positive or negative. No titer is given.	0.5 mL serum in red top or serum separator tube	2-4 days
CS86293	Blastomyces Quantitative Ag Assay Blastomyces dermatitidis Ag, EIA An antigen test used to diagnose blastomycosis.	2.0 mL urine in urine transport tube,	7-9 days
CS16003	Brucella Titer For Export (KSU) Evaluation for antibody suggestive of <i>Brucella canis</i> exposure. <i>Brucella canis</i> antibody evaluation is performed by serum agglutination. The test is typically used when required for export. Canine only. Note: Clearly write microchip number on the test requisition form (TRF) and on submitted sample tube.	1.0 mL serum in red top tube or serum separator tube	7-10 days
CT105	BUN/Creatinine BUN, Creatinine with SDMA and BUN/Creat Ratio BUN, creatinine, and SDMA for glomerular filtration rate estimation (see T1035). Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top tube or serum separator tube	Daily Performed each shift
CS86344	C Reactive Protein CRP Note: valid in dogs only.	1.0 mL serum in red top tube or serum separator tube	5-15 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA160	Canine Maldigestion Profile Cobalamin B12/Folate, TLI - Canine Cobalamin, folate, TLI Interferences: see individual components. Note: pet should fast overnight. Hemolysis of sample may affect folate results. Previous vitamin supplementation will cause elevated concentrations.	1.0 mL serum in red top or serum separator tube	2-5 days
CT1010	Cardio BNP-Canine Antech Cardio BNP - Canine This test provides a quantitative assessment of canine heart health by measuring the concentration of the carboxy terminal of the BNP peptide (cBNP), which is released by cardiac myocytes in response to stretch and stress. Note: draw 2.0 mL of whole blood into an LT tube, centrifuge for plasma, draw 0.5 mL plasma from the tube, and inject into closed BNP tube. Keep sample cold. It requires a specific BNP tube for submission.	0.5 mL EDTA plasma injected into BNP tube	2-4 days
CT120	Chloride Interferences: marked lipemia. Serum should be separated within 1 hour after draw.	0.5 mL serum in red top or serum separator	Daily Performed each shift
CT125	Cholesterol Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator	Daily Performed each shift
CT127	Cholesterol, Triglycerides Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator	Daily Performed each shift
CT235	Cholinesterase-Serum Test used to diagnose organophosphate toxicity.	1.0 mL serum in red top or serum separator	7-10 days
CSA275	Cobalimine, Folate, TLI (Feline) Cobalamin B12/Folate, TLI - Feline Used to further evaluate gastrointestinal signs where malassimilation and maldigestion are suspected. Note: pet should fast overnight. Hemolysis of sample may affect folate results. Previous vitamin supplementation will cause elevated concentrations.	1.5 mL serum in red top or serum separator tube	1-10 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS16210	Copper Level This is used to detect exposure to excessive copper. Note: submit sample on ice. In the case of a primary copper, hepatopathy only a quantitative copper on liver tissue combined with histopathology can confirm the diagnosis.	1.0 mL serum or plasma, 5 grams fresh liver, 50 mg fresh liver biopsy (approx. 3 Tru-Cut samples), 3.0 mL urine, or 500 g Feed. For birds: 0.3 mL serum or plasma is adequate. For ruminants and camelids: 5 grams fresh kidney is required.	10-15 days
CS16215	Copper Storage Disease Tissue Copper Quantification To be used when a quantitative copper level is required on hepatic tissue. Note: a minimum of 5 grams of fresh liver in saline is the ideal sample.	5 grams fresh liver tissue placed in saline	10-15 days
S16225	Corona Virus IgG and IgM Evaluation for IgM and IgG antibody to canine coronavirus.	0.5 mL serum in red top or serum separator tube	3-5 days
CT130	CPK Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CS7592	CPK Isoenzymes CPK Isoenzymes: CK- MM, CK- MB, CK- BB Evaluation of CPK isoenzymes. Note: clinical utility of this test has not been validated in dogs and cats.	0.5 mL serum in red top tube or plain transport tube	2-7 days
CT135	Creatinine Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT340	<p>Crossmatch</p> <p>Donor ID, Major and Minor Cross Match</p> <p>Screen for pre-existing antibodies directed against red cell antigens to determine serological compatibility prior to transfusion. If serum, as well as whole blood, is submitted. Both minor and major cross-match will be performed. This test is not intended for use in determining pre-breeding compatibility testing, neonatal isoerythrolysis, or the prevention of neonatal isoerythrolysis.</p> <p>Interferences: marked hemolysis and lipemia. Clotting or freezing of sample may preclude analysis. Use 1.0 mL whole blood and 1.0 mL of serum from the patient and prospective donor. Ensure all tubes are appropriately labeled as patient and donor.</p>	1.0 mL EDTA whole blood in lavender top tube and 1.0 mL serum in red top tube for recipient and each donor	2-4 days
CT345	<p>Crossmatch Additional Donor 2</p> <p>Donor ID, Major and Minor Cross Match</p> <p>Screens for pre-existing antibodies directed against red cell antigens to determine serological compatibility prior to transfusion. If serum, as well as whole blood, is submitted both minor and major cross-match will be performed. This test is not intended for use in determining pre-breeding compatibility testing, neonatal isoerythrolysis, or the prevention of neonatal isoerythrolysis.</p> <p>Interferences: marked hemolysis and lipemia, clotting or freezing of sample may preclude analysis.</p> <p>Note: this test must be used in addition to Cross Match - CT340.</p>	1.0 mL EDTA whole blood in lavender top tube and 1.0 mL serum in red top tube for recipient and each donor	2-4 days
CSA350	<p>Crypto, Giardia, Clostridium Enterotoxin</p> <p><i>Clostridium</i> Perfringens Enterotoxin, <i>Giardia</i> (FA), <i>Cryptosporidium</i> (FA), <i>Giardia</i> ELISA</p> <p>A panel that can be considered when suspecting <i>Clostridium perfringens</i>, <i>Cryptosporidium</i> spp., or <i>Giardia</i> spp. as the cause for gastroenteritis. Submission includes <i>Cryptosporidium</i> spp. <i>Giardia</i> spp. evaluated by FA, <i>Clostridium</i> perfringens enterotoxin assessment and a <i>Giardia</i> ELISA Antigen detection.</p>	10 grams feces in Antech provided fecal container	2-7 days
CT350	<p>D-Dimer</p> <p>This test is used to measure the concentration of D dimers.</p> <p>Interferences: marked hemolysis or lipemia.</p> <p>Note: moderately to markedly lipemic samples can cause D dimer results to be spuriously decreased. Citrated Whole Blood (Blue Top Tube) or citrated plasma is the only acceptable sample. The tube should be >2/3rds filled. If submitting separated citrated plasma, label it as Citrated Plasma.</p>	0.5 mL citrated plasma or citrated whole blood in blue toptube at least 2/3rds full to the fill line	2-4 days Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS16250	Distemper Smear, FA This test is helpful in making a diagnosis of distemper, in the absence of recent vaccination, when used in conjunction with clinical signs.	EDTA whole blood in lavender top tube, or unstained peripheral blood, buffy coat or urine sediment smears Call customer service for specimen handling protocol: 1-800-341-3440, dial 0	3-7 days
CT560	Distemper Vaccinal Titer This test detects IgG antibodies against canine distemper virus. The result is reported as positive or negative. A titer is not reported. This test is not to be used in patients suspected to be ill from distemper virus infection. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	2-4 days
CT675	Distemper/Parvo End Point Distemper IgG End PT Titer, Parvo IgG End Point This test provides an endpoint IgG titer for Canine Distemper and Canine Parvovirus exposure. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	2-4 days
CT565 Add-on Equivalent CADD270	Distemper/Parvo Vaccine Titer This test detects IgG antibodies against canine parvovirus and canine distemper virus using ELISA. Results are reported as positive (protective) or negative (insufficient neutralizing antibody detected). This test is not to be used in patient suspected of having distemper or parvovirus infection. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	2-4 days
CSA280	Feline Autoimmune Profile CBC, Antinuclear Antibodies (ANA), Direct Coombs' Test Warm An autoimmune panel evaluating a complete blood count, Direct Coombs test warm, and Antinuclear antibodies. Interferences: marked hemolysis and lipemia.	1.0 mL serum and 1.5 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-4 days
CS16865	FIV Antibody, Western Blot This is a western blot test to detect antibodies directed against the Feline Immunodeficiency Virus. It may be used as a confirmatory test for the presence of antibodies against FIV. Note: does not differentiate cats infected with FIV from cats vaccinated against or kittens with passive transfer of antibodies against FIV from the queen.	0.5 mL serum in red top or serum separator tube	2-5 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSA340	Fungal Serology with Cocci Fungal Serology (<i>Histoplasma</i> , <i>Blastomyces</i> , <i>Aspergillus</i>), <i>Coccidioidomycosis</i> Screen and Titer Detects the presence of antibodies to <i>Histoplasma</i> , <i>Blastomyces</i> , <i>Aspergillus</i> , and <i>Coccidioides</i> spp. Interferences: marked lipemia. A titer is only reported for <i>Coccidioides</i> . All other results are reported as positive or negative for antibody detection. Both the <i>Blastomyces</i> and <i>Histoplasma</i> antigen tests are preferred over antibody testing as a diagnostic test.	1.0 mL serum in red top or serum separator tube	3-8 days
CT790	Giardia ELISA/Crypto FA <i>Giardia</i> ELISA, <i>Giardia</i> IFA, and <i>Cryptosporidium</i> IFA Evaluation for <i>Giardia</i> Ag (ELISA) and <i>Cryptosporidium</i> (FA) presence in stool. Note: <i>Giardia</i> FA is less sensitive than the <i>Giardia</i> ELISA and O&P tests. <i>Cryptosporidium</i> FA is likely more sensitive than an O&P test. The GI PCR panel would be preferred to either FA test with regards to sensitivity and specificity.	5 grams fresh fecal specimen in Antech provided fecal container	2-5 days
CT150	Glucose Interferences: hemolysis or lipemia. Hemolysis may cause false decreases in glucose. Lipemia may cause false increases in glucose.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CSA140	Heartworm Program 1 Miniscreen 4 Chem with Electrolytes, Heartworm Antigen Miniscreen chemistry and heartworm antigen detection. Interferences: marked hemolysis and lipemia.	1.0 mL serum Serum in red top or serum separator	1-2 days
CSA130	Heartworm Program 2 Miniscreen 4 Chem, CBC, Heartworm Antigen Miniscreen chemistry, a complete blood count, and heartworm antigen detection. Interferences: marked hemolysis and lipemia.	1.0 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator tube, lavender top	1-2 days
CS16400	Herpes - Conjunctival Smear Conjunctival smears are used to evaluate for the presence of herpes virus antigen in epithelial cells via immunofluorescence. Note: ensure slides have adequate cellularity.	2 conjunctival smears	9-12 days
CS86022	Herpes Antibody, IFA Detection of Herpes antibody by IFA.	0.3 mL serum in red top or serum separator tube	5-7 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS17029	Heska ALLERCEPT Environmental & Food Panel	3.0 mL serum in red top or serum separator tube	4-5 business days
CS17026	Heska ALLERCEPT Environmental Panel	2.0 mL serum in red top or serum separator tube	4-5 business days
CS17028	Heska ALLERCEPT Food Panel	1.0 mL serum in red top or serum separator tube	4-5 business days
CS17027	Heska Equine ALLERCEPT Panel	3.0 mL serum in red top or serum separator tube	4-5 business days
CS16405	Histoplasma Titer Detection of antibody indicating exposure to <i>Histoplasma</i> spp. Titer provided.	1.0 mL serum in red top or serum separator tube	9-12 days
CT660	Immunoglobulins A, G, M This test quantifies IgG, IgM and IgA levels by radial immunodiffusion in dogs. Interferences: marked hemolysis or lipemia. Note: reference intervals for adult dogs and cats have been established. However, if submitting serum for a young animal, age-matched controls should be submitted concurrently. References intervals for other species are provided by the manufacturer and are age-dependent. It is imperative to provide the entire signalment, including age, when submitting samples.	0.5 mL serum in red top tube	7-9 days
CS86102	Influenza - Canine Convalescent H3N8 and H3N2 influenza antibody detection Used to detect for the presence of antibody to H3N8 and H3N2 influenza viruses.	1.0 mL serum in red top tube	10-15 days
CT155	Iron Serum Interferences: marked hemolysis. Hemolysis can result in falsely increased iron levels.	0.5 mL serum in red top tube or serum separator tube	2-4 days
CT160	Lactic Dehydrogenase Lactate dehydrogenase (LDH) Interferences: hemolysis. Hemolysis or delayed separation can falsely elevate results.	0.5 mL serum in red top tube or serum separator tube	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS5017	<i>Leptospira canicola</i>, Export Serum assessed for antibodies to <i>Leptospira canicola</i> using microscopic agglutination test (MAT). Note: clearly write microchip number on the requisition form, sample tube, and KSU submission form.	0.3 mL serum in red top tube or serum separator tube	10-15 days
CT170	Magnesium Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CS16535	Masticatory Muscle Myositis Masticatory Muscle Myositis - 2M Antibody This test is used to diagnose masticatory myositis by detecting autoantibodies specifically directed against masticatory muscle proteins. Note: this assay may be negative if immunosuppressive dosages of cortisocosteroids have been given for longer than 7-10 days prior to test submission.	1.0 mL serum in red top or serum separator tube	12-15 days
CT810	Occult Blood, Feces This is a test for the detection of occult blood in fecal samples. The test is not highly sensitive and is dependent on the amount of heme present and how evenly the heme is spread throughout the sample. Note: diets of red meat, chicken, or fresh/uncooked vegetables may cause false-positive results and should not be eaten for 48 hours prior to testing.	2 grams feces in Antech provided fecal container	1-2 days
CS16575	Osmolality - Serum This may be of utility in the evaluation of the causes of PU/PD.	1.0 mL serum in red top or serum separator tube	7-10 days
CS85364	Pancreatic Lipase Immunoreactivity Measures concentrations of pancreatic lipase in the serum. Interferences: hemolysis. Note: a 12 hour fast is recommended prior to sample submission.	0.5 mL serum in red top or serum separator tube	7-9 business days
CS86468	Pancreatitis Profile - Canine Cobalamin B12/Folate, Canine cTLI, Pancreatic Lipase Immunoreactivity Includes cobalamin (B12), folate, trypsin like immunoreactivity (TLI), and pancreatic lipase immunoreactivity (cPLI). Note: pet should fast overnight. Hemolysis of sample may affect folate results. Previous vitamin supplementation will cause elevated concentrations.	1.5 mL serum in red top or serum separator tube	2-10 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS86288	Pancreatitis Profile - Feline Cobalamin B12/Folate, Feline fTLI, Pancreatic Lipase Immunoreactivity Includes cobalamin (B12), folate, trypsin like immunoreactivity (TLI), and pancreatic lipase immunoreactivity (fPLI).	1.5 mL serum in red top or serum separator tube	2-10 days
CS16580 Add-on Equivalent CADD280	Panleukopenia Titer IgG, IgM This test evaluates for the presence of IgM and IgG antibodies to the feline panleukopenia virus. It will be positive either due to infection or vaccination against panleukopenia. Kittens may test positive due to transmammary transfer of antibodies. This test will also detect antibodies directed against canine parvovirus in cats infected with CPV-2a and CPV-2b. Interferences: marked hemolysis or lipemia (if unable to clear with centrifugation).	0.5 mL serum in red top or serum separator tube	3-7 days
CS16053 Add-on Equivalent CADD280	Panleukopenia Vaccinal Titer The test detects IgG antibodies against Panleukopenia via Immunofluorescent antibody assessment (IFA). A titer of 1:5 or greater, in the absence of clinical signs of disease, indicates an immunological response to vaccination. A vaccine titer of less than 1:5 indicates a low level of circulating antibodies. Interferences: marked hemolysis or lipemia (if unable to clear with centrifugation).	0.5 mL serum in red top or serum separator tube	3-5 days
CT825	Parasite Identification Note: document where on the animal the parasite was noted.	Parasite sent in 70% ethanol. Include host species and body location where parasite was collected. If submitting slides with worms, place slides in a leak proof container with 70% ethanol. Other representative sample types include feather, intestine, liver, fresh tissue, urine, wash/lavage, feces or vomitus)	10-12 days
CT700	Parvovirus Antibody/Antigen Parvovirus Antigen, Canine Parvovirus Antibody Titer (IgM and IgG) Detects parvovirus antigen in fecal samples and evaluates a concurrent serum sample for IgM and IgG titers consistent with Canine Parvovirus exposure. Note: modified live virus vaccination for parvovirus may give false-positive results for about two weeks after vaccination. A negative test does not rule out parvovirus infection.	0.5 mL serum and 5 grams feces Serum in red top or serum separator tube, Antech provided fecal container	1-4 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT705	<p>Parvovirus Vaccinal Titer</p> <p>Used to evaluate current IgG levels to Parvovirus in light of previous vaccination.</p> <p>Interferences: marked hemolysis or lipemia.</p> <p>Note: not to be used in dogs with current or recent clinical signs of parvovirus infection as a positive may indicate response to infection rather than protection.</p>	0.5 mL serum in red top or serum separator tube	2-4 days
CT180	<p>Phosphorus</p> <p>Interferences: marked hemolysis or lipemia. Hemolysis or delayed separation will falsely increase results.</p>	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CT185	<p>Potassium</p> <p>Interferences: hemolysis and marked lipemia. Hemolysis will falsely elevate results.</p>	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CT9810	<p>Relaxin</p> <p>Relaxin is a hormone produced by the placenta, and detection of relaxin in serum or plasma is a sensitive and specific test for pregnancy in dogs and cats.</p> <p>Interferences: marked hemolysis.</p>	0.2 mL serum in red top or serum separator tube	2-3 days
CS16702	<p>Rhinotrachitis Feline</p> <p>Viral Neutralization Antibody Titer</p> <p>This is a semi-quantitative titer that is correlated to the amount of protective humoral antibody present in relation to previous feline rhinotracheitis vaccination (feline herpesvirus).</p>	1.0 mL serum in red top or serum separator tube	10-12 days
CS16730	<p>Selenium Level</p> <p>Interferences: marked hemolysis or lipemia.</p>	1.0 mL serum in red top or serum separator tube (spun), 1.0 mL EDTA whole blood in lavender top tube, or 10.0 grams offeed (representative sample) Call customer service for specimen handling protocol: 1-800-341-3440, dial 0	7-10 days
CT195	<p>Sodium</p> <p>Interferences: marked hemolysis or lipemia. Hemolysis falsely elevates results.</p>	0.5 mL serum in red top or serum separator tube	Daily Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT200	Sodium and Potassium Sodium and potassium evaluation. Interferences: hemolysis and marked lipemia. Hemolysis will falsely elevate results.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CT485	T3 Autoantibodies Detection of T3 autoantibody is used to assess the presence of autoimmune thyroid disease. This test is less sensitive to autoimmune thyroid disease than the thyroglobulin autoantibody (TGAA) test.	0.5 mL serum in red top or serum separator tube	2-4 days
CT500	T4 Autoantibodies Detection of T4 autoantibody is used to assess the presence of autoimmune thyroid disease. This test is less sensitive to autoimmune thyroid disease than the thyroglobulin autoantibody (TGAA) test.	0.5 mL serum in red top or serum separator tube	2-4 days
CS16755	Taurine Measurement of taurine concentration in whole blood via HPLC. Note: do not fast prior to sample submission. Whole blood submitted in a green-top refrigerated is preferred. EDTA may also be used but taurine concentrations in EDTA anticoagulated blood will be slightly lower than in a GTT.	1.0 mL heparinized whole blood in green top tube	10-15 days
CSA330	Tick Serology 1 <i>Ehrlichia canis</i> , Lyme Titer IgG, and Rocky Mountain Spotted Fever IFA Detects the presence of antibodies to <i>Ehrlichia canis</i> , <i>Borrelia burgdorferi</i> (Lyme), and Rocky Mountain Spotted fever. Interferences: marked hemolysis or lipemia.	1.0 mL serum in red top or serum separator tube	2-5 days
CS16800	TLI, Feline Trypsin-Like Immunoreactivity Use to confirm feline exocrine pancreatic insufficiency (EPI). Note: recommend fasting for 8-12 hours prior to drawing sample.	1.0 mL serum in red top or serum separator tube	7-10 days
CT190	Total Protein Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CT205	Triglycerides Interferences: marked hemolysis.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CT230	Trypsin-Like Immunoreac Canine cTLI Canine Trypsin- Like Immunoreactivity (cTLI is) used to diagnose exocrine pancreatic insufficiency. Note: a 12-hour fast is required. Canine only.	0.5 mL serum in red top or serum separator tube	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT227	Ur Bile Acid (UBA):Ur Creat (UCR) Urine Bile Acid and Urine Creatinine Ratio Evaluation of the amount of urinary bile acid present with respect to creatinine. Note: urine containing ascorbic acid may result in false positive results.	1.0 mL urine in urine transport tube	2-5 days
CT100	Urea Nitrogen Blood urea nitrogen. Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CS16870	Zinc Evaluation of blood zinc level. Note: plasma or serum needs to be separated from cells and transferred to a plastic vial. Do not use tubes with rubber stoppers.	0.3 mL of serum in royal blue top tube or plastic transfertube Call customer service for specimen handling protocol: 1-800-341-3440, dial 0	7-10 days
CT730 Add-on Equivalent CADD320	Bromide This test is used to monitor bromide therapy. Therapeutic serum bromide concentrations should be measured at three weeks after initiating therapy, but steady-state concentrations may fluctuate among dogs due to differences in drug clearance and bioavailability. Interferences: gel may interfere with test. Do not use a serum separator tube for sample draw or submission.	0.5 mL serum in red top tube. Submission in serum separator tube is not recommended.	1-2 days
CS18702	Cyclosporine Used to assess the cyclosporine drug level being achieved, at the current drug dosing, in a given patient. Note: peak drug level occurs 2 hours post medication. Trough drug level will occur just prior to the next 12-hour dose.	2.0 mL EDTA whole blood in lavender top tube	7-10 days
CT735	Digoxin Digoxin level Interferences: hemolysis, gel may interfere with test. Do not use a serum separator tube for sample draw or submission.	0.5 mL serum in red top tube collected 8 hours post pill	1-2 days
CS86541	Keppra Level Interferences: hemolysis and gel in serum separator tube may interfere with test.	2.0 mL serum in red top or non-additive tube. Submission in serum separator tube is not recommended	5-7 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT750	Phenobarbital Phenobarbital level Interferences: hemolysis and gel in serum separator tube may interfere with test. Do not use a serum separator tube for sample draw or submission.	0.5 mL serum in red top tube. Submission in serum separator tube is not recommended.	1-2 days
CSA830	Phenobarbital Panel Plus Liver Chemistry with Electrolytes, CBC, Bile Acids, Phenobarbital Includes a chemistry (CSA324), complete blood count, a single bile acids, and a phenobarbital level. Interferences: marked hemolysis and lipemia. All therapeutic drug monitoring should be submitted as serum in a red top tube. A serum separator tube should not be used in sample collection as the drug concentration being measured may be erroneously decreased by as much as 30% due to binding with gel. Once steady-state drug levels have been reached, phenobarbital levels can be evaluated at any time of the day in relation to medication administration. Ursodeoxycholic acid may be detected by bile acid assay, causing falsely elevated values.	1.5 mL serum and 1.0 mL EDTA whole blood Serum in red top tube, lavender top	1-2 days
CT755	Pre & Post Phenobarbital Pre and Post Pill Phenobarbital Levels Submission of a peak and trough phenobarbital level. Interferences: hemolysis and gel in serum separator tube may interfere with test.	0.5 mL serum in red top tube (labeled Pre) and 0.5 serum in red top tube (labeled Post). Submission in serum separator tube is not recommended.	1-2 days
CS86480	Zonisamide Zonegran Zonisamide concentration. Interferences: hemolysis, gel may interfere with test.	0.5 mL serum in red top tube. Submission in serum separator tube is not recommended.	3-5 days

AVIAN • EXOTICS

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS17116	Adrenal Androgen Panel For Ferrets Estradiol, 17-OH Progesterone, Androstenedione	0.5 mL heparinized plasma spun in green microvial with gel separator	10-15 days
CS16011	Aspergillus AB Avian (AFMP1P ELISA) Aspergillus Antibody, Avian	0.5 mL serum in microtainer or regular serum separator tube(spun)	7-10 days
CS85358	Aspergillus Ag, Galactomannan Assay Aspergillus Antigen, Galactomannan Assay (Avian) Antigen detection test (EIA immunoenzymatic sandwich microplate assay) for galactomannan.	Canine or Feline: 0.8 mL serum in red top or serum separator tube Avian: 0.1 mL serum in red top or serum separator tube	7-10 days
CAE021	Avian Comprehensive Profile with Bile Acid Comprehensive Avian Profile with Bile Acids Avian/Exotic	0.75 mL of heparinized plasma, 0.5 mL of heparinized whole blood and 2 freshly prepared blood smears 2 green top micro vials (with gel), green top micro vial (no gel) or 2 hematocrit tubes, slide holders	Daily Performed each shift
CAE025	Avian Comprehensive Profile with Bile Acid & EPH Comprehensive Avian Profile with Bile Acids (Avian/Exotic), and Avian Protein Electrophoresis	0.75 mL of heparinized plasma, 0.5 mL of heparinized whole blood and 2 freshly prepared blood smears 2 green top micro vials (gel), green top micro vial (no gel) or 2 hematocrit tubes, slide holders	1-4 days
CAE300	Avian Protein Electrophoresis Includes Total Protein, Pre-Albumin, Albumin, Alpha 1, Alpha 2, Beta and Gamma	0.5 mL of heparinized plasma in green micro vial (with gel)	2-5 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CAE051	Avian Standard Profile with Bile Acid Standard Avian Profile with Electrolytes, Avian/Exotic CBC, Bile Acids	0.75 mL of heparinized plasma, 0.5 mL of heparinized whole blood and 2 freshly prepared blood smears 2 green micro vials (gel), 2 hematocrit tubes, 2 blood film slides	1-2 days
CAE055	Avian Standard Profile with Bile Acid & EPH Standard Avian Profile with Electrolytes, Bile Acids Avian/Exotic, Avian Protein Electrophoresis	0.75 mL of heparinized plasma, 0.5 mL of heparinized whole blood and 2 freshly prepared blood smears 2 green micro vials (gel), 2 hematocrit tubes, 2 blood film slides	1-5 days
CS16012	Avian Zinc Assay	0.5 mL serum in yellow micro vial with gel	7-10 days
CAE270	Avian/Exotic CBC Includes WBC Estimate, Thrombocyte Estimate, Hematocrit, Blood Parasites, and Differential: Percent and Absolute Heterophils, Bands, Lymphocytes, Monocytes, Eosinophils, Basophils and Azurophilic Monocytes Interferences: marked hemolysis.	0.5 mL of heparinized whole blood in green top tube with 2 freshly prepared blood smears	Daily Performed each shift
CAE260	Bile Acids (Avian/Exotic) Single Bile Acids, Avian/Exotic Interferences: marked hemolysis or lipemia. Ursodeoxycholic acid may be detected by bile acid assay, causing falsely elevated values.	0.5 mL serum in red top or serum separator tube	1-2 days
CS16671	Chlamydia Antibody Titer IFA	0.5 mL serum in red top or serum separator tube	7-10 days
CS16788	Chlamydia PCR Blood	0.5 mL heparinized whole blood in green top tube	7-10 days
CS16672	Chlamydia PCR Swab	Avian: Swab in red top tube (combined choanal & cloacal) Feline & other species: conjunctival swab in red top	7-10 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS85206	Chlamydia Profile <i>Chlamydophila</i> PCR Blood, <i>Chlamydophila</i> Antibody Titer IFA, <i>Chlamydophila</i> Titer EBA, <i>Chlamydophila</i> PCR Swab, <i>Chlamydophila</i> Titer EBA, <i>Chlamydophila</i> PCR Swab	0.4 mL heparinized whole blood, 0.4 mL serum, swab in red top tube (combined choanal & cloacal) Green micro vial (gel), Green microvial (no gel), Swab in white tube (combined choanal & cloacal)	7-14 days
CS16670	Chlamydia Titer EBA Psittacosis	0.5 mL serum in serum separator tube	7-14 days
CAE010	Comprehensive Avian Chemistries Total Protein, Albumin, Globulin, SGOT (AST), Phosphorus, Glucose, Calcium, Sodium, Potassium, Chloride, Cholesterol, CPK, Uric Acid	0.5 mL heparinized plasma in green microtainer with gel	Daily Performed each shift
CAE030	Comprehensive Avian Post Purchase Avian Comprehensive Chemistries, <i>Chlamydophila</i> Titer EBA, <i>Giardia</i> ELISA, Gram Stain, Psittacine Beak and Feather Disease, Polyoma PCR Swab, Avian Protein Electrophoresis, Avian/Exotic CBC	1.0 mL heparinized plasma, 1.0 mL heparinized whole blood, 0.5 grams feces, and 2 fecal swabs Antech provided fecal container, 3 green micro vials (gel), green micro vial (no gel), 2 hematocrit tubes, 2 slides, swabs in white tube (2 fecals)	1-10 days
CAE020	Comprehensive Avian Profile Avian/Exotic CBC, Avian Comprehensive Chemistries	0.5 mL of heparinized plasma, 0.5 mL of heparinized whole blood and 2 freshly prepared blood smears Green top micro vial (with gel), green top micro vial (no gel) or 2 hematocrit tubes, slide holder	Daily Performed each shift
CAE200	Comprehensive Mammalian Profile CBC Small Mammalian, Mammalian Comprehensive Chemistries	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator or green micro vial (gel), lavender or green vial (no gel)	Daily Performed each shift

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CAE160	Comprehensive Reptilian Profile Avian/Exotic CBC, Reptilian Comprehensive Chemistries	0.5 mL of heparinized plasma, 0.5 mL of heparinized whole blood, and 2 freshly prepared blood smears Green top micro vial (with gel), green top micro vial (no gel) or 2 hematocrit tubes, slide holder	Daily Performed each shift
CRECHECKREP	Comprehensive Reptilian Reckeck Panel Avian/Exotic CBC, Reptilian Comprehensive Chemistries	0.5 mL of heparinized plasma, 0.5 mL of heparinized whole blood and 2 freshly prepared blood smears Green top micro vial (with gel), green top micro vial (no gel) or 2 hematocrit tubes, slide holder	Daily Performed each shift
CAE070	Diarrhea Profile (Avian) Comprehensive Avian Profile, Gram Stain, Culture, Aerobic Culture, <i>Giardia</i> ELISA Interferences: marked hemolysis and lipemia.	0.5 mL heparinized plasma in green microtainer with gel, 0.5 mL whole blood in green microtainer, 2 fecal slides, 0.5 grams feces, and culturette Culturette, Antech provided fecal container, green micro vial (gel), 2 hematocrit tubes, 2 slides each (blood and fecal)	1-2 days
CS16501	Distemper PCR	2.0 mL urine in urine transport tube, whole blood (EDTA) in lavender top tube, CSF OR other body cavity fluid in red top or lavender top tube	5-9 days
CS16107	Distemper Titer-Virus Neutraliz.	0.5 mL serum in red top or serum separator tube	14-21 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS14484	Eastern Encephalitis Virus PCR	Avian and Equine: Brain in leak-proof container Call customer service for specimen handling protocol: 1-800-341-3440, dial 0	5-7 days
CS16877	Encephalitozoon cuniculi IgG Ab Interferences: marked hemolysis.	0.10 mL spun serum or spun heparinized plasma in plain, non-additive transport tube	7-10 days
CAE080	Feather Picker Profile Comprehensive Avian Profile, Culture (Aerobic), <i>Giardia</i> ELISA, Gram Stain, Psittacine Beak and Feather Disease, Avian Protein Electrophoresis	1.0 mL heparinized plasma in green micro vial with separator gel, 1.0 mL heparinized whole blood in green top tube, culturette, 2 fecal slides and 0.5 grams feces Culturette, Antech provided fecal container, 2 green top micro vial, green top micro vial (no gel), 2 hematocrit tubes, 2 slides, (blood and fecal)	1-7 days
CAE230	Geriatric/Weak Ferret Profile Mammalian Comprehensive Chemistries, Insulin Ferret, Urinalysis, CBC Small Mammalian	1.0 mL serum, 1.0 mL EDTA whole blood, and 6.0 mL urine Serum in red top or serum separator or 2 green micro vial (gel), lavender or green vial (no gel), urine transport tube	1-5 days
CT820 Add-on Equivalent CADD250	Giardia (ELISA) The <i>Giardia</i> antigen capture Enzyme- linked Immunosorbent Assay (ELISA) is a fecal procedure designed to detect a <i>Giardia</i> specific antigen.	5 grams fresh fecal specimen in Antech provided fecal container	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CAE090	Hepatic Profile Comprehensive Avian Profile, Bile Acids Avian/Exotic, <i>Chlamydomphila</i> Titer EBA, Avian Protein Electrophoresis	1.0 mL heparinized plasma (in green top micro via with gel separator) and 0.5 mL heparinized whole blood 3 green micro vial (gel), 2 hematocrit tubes, slides	1-14 days Performed each shift
CAE190	Mammalian Comprehensive Chems Total Protein, Albumin, Globulin, SGOT (AST), SGPT (ALT), Alk Phos, T. Bilirubin, Bun, Creatinine, Phosphorus, Glucose, Calcium, Sodium, Potassium, Chloride, Cholesterol, CPK	0.5 mL serum in red top or serum separator OR green top micro vial (gel)	Daily Performed each shift
CAE210	Mammalian Standard Chemistries Total Protein, SGPT (ALT), Alk Phos, T. Bilirubin, Bun, Creatinine, Sodium, Potassium, Chloride, Phosphorus, Glucose, Calcium	0.5 mL serum in red top or serum separator OR green top micro vial (gel)	Daily Performed each shift
CAE060	Mini Avian Post Purchase Comprehensive Avian Profile, Gram Stain, Avian Protein Electrophoresis	0.75 mL heparinized plasma, 0.5 mL heparinized whole blood and 2 freshly prepared blood smears, and fecal slides Green top micro vial (with gel), green top micro vial (no gel) or 2 hematocrit tubes, slide holder, culturette	1-5 days
CS16789	Mycoplasma PCR	Tissue (fresh) OR fluid (ocular, conjunctival, BAL) OR nasal swab	10-15 days
CT805	Ova & Parasite O&P with Centrifugation Fecal ova and parasites. Samples are appropriately mixed with zinc sulfate solution, centrifuged, followed by flotation and slide evaluation. Note: sample should be evaluated within 24 hours of collection. If a worm has been identified in the sample, separate the worm and place it in a container labeled Worm in black, additionally mark Worm on the TRF being submitted with the sample.	5 grams feces in Antech provided fecal container	1-2 days
CS16085	Psittacine Beak and Feather Disease	0.5 mL EDTA or heparinized whole blood in lavender or green top tube	7-10 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CAE110	PU/PD Profile Avian Comprehensive Chemistries, Avian/Exotic CBC, Avian Protein Electrophoresis, Urinalysis	0.5 mL heparinized plasma, 0.5 mL heparinized whole blood, and 1 mL urine 2 green top micro vials (gel), green top micro vial (no gel) or 2 hematocrit tubes, slides, urine transport tube	1-5 days
CAE240	Rabbit Neurologic Profile CBC Small Mammalian, Mammalian Comprehensive Chemistries, <i>Encephalitozoon cuniculi</i> IgG Ab, Pasteurella Antibody Titer	1.0 mL of serum or heparinized plasma and 1.0 mL EDTA whole blood 2 green top micro vials (gel), lavender top or green top micro vial (no gel)	1-10 days
CAE150	Reptilian Comprehensive Chems Total Protein, Albumin, Globulin, SGOT (AST), Bun, Phosphorus, Glucose, Calcium, Sodium, Potassium, Chloride, CPK, and Uric Acid Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator OR green top micro vial (gel)	Daily Performed each shift
CAE170	Reptilian Standard Chemistries Total Protein, SGOT (AST), Sodium, Potassium, Chloride, Phosphorus, Glucose, Calcium, CPK, and Uric Acid Interferences: marked hemolysis and lipemia.	0.5 mL serum in red top or serum separator OR green top micro vial (gel)	Daily Performed each shift
CT425	Reticulocytes Reticulocytes, Absolute Reticulocytes Includes an unadjusted reticulocyte percent and absolute reticulocyte count. Only aggregate reticulocytes are counted. Interferences: marked hemolysis.	1.0 mL EDTA whole blood in lavender top tube	Daily Performed each shift
Add-on Equivalent CADD140			
CAE050	Standard Avian Profile Avian/Exotic CBC, Avian Standard Chemistries with Electrolytes Interferences: marked hemolysis or lipemia.	0.5 mL heparinized plasma and 0.5 mL heparinized whole blood Green top micro vial (gel), green top micro vial (no gel) or 2 hematocrit tubes, slides	Daily

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CAE220	Standard Mammalian Profile CBC Small Mammalian, Mammalian Standard Chemistries with Electrolytes Interferences: marked hemolysis or lipemia.	0.5 mL serum and 1.0 mL EDTA whole blood Serum in red top or serum separator or green micro vial (gel), lavender or green vial (no gel)	Daily
CAE180	Standard Reptilian Profile Avian/Exotic CBC, Reptilian Standard Chemistries with Electrolytes Interferences: marked hemolysis or lipemia.	0.5 mL heparinized plasma and 0.5 mL heparinized whole blood Green micro vial (gel), green top micro vial (no gel) or 2 hematocrit tubes, slides	Daily
CS16792	Toxoplasma antibody - Exotics	1.0 mL serum in red top or serum separator tube	7-14 days
CS85448	West Nile Titer (PRNT & IgM ELISA)	1.0 mL serum in red top or serum separator tube (mammals) OR 0.1 mL heparinized plasma (avians)	7-14 days
CS85449	West Nile Virus PCR	Equine: 1.0 mL EDTA whole blood in lavender top tube, CSF, brain or spinal cord. Avian: Kidney in sterile container, EDTA or heparinized whole blood, or tissue (brain or heart). Other mammals: EDTA whole blood, CSF, brain or spinal cord.	7-10 days
CS14485	West Nile Virus Plaque Reduction SN	1.0 mL serum in red top or serum separator tube	7-17 days

EQUINE • LARGE ANIMAL ESSENTIAL PANELS

	EQUINE CHEMISTRY CL010	EQUINE PERFORMANCE HORSE CHEMISTRY CL601	EQUINE HEPATIC SCREEN CL225	EQUINE RENAL SCREEN CL240	RUMINANT CHEMISTRY CL1001
A/G Ratio	●	●			●
Albumin	●	●	●	●	●
Alkaline Phosphatase	●	●	●		●
AST (SGOT)	●	●	●		●
BUN	●	●	●	●	●
BUN/Creatinine Ratio	●	●			●
Calcium	●			●	●
Chloride	●	●	●	●	●
Cholesterol	●				●
CPK	●	●			●
Creatinine	●	●		●	●
Direct Bilirubin	●		●		●
GGT	●	●	●		●
Globulin	●	●	●		●
Glucose	●			●	●
LDH	●		●		●
Magnesium					●
NA/K Ratio	●				●
Phosphorus	●			●	●
Potassium	●	●	●	●	●
Sodium	●	●	●	●	●
Total Bilirubin	●		●		●
Total Protein	●	●	●	●	●
Triglycerides	●				●

EQUINE

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
EQUINE HEALTH PROGRAM (EHP)			
CL511E	Equine Health Program GVL (ELISA) Equine Chemistry, CBC, Fibrinogen, EIA (ELISA) by GVL, FEC (MST)	2.0 mL serum, 1.0 mL EDTA whole blood, 0.5 mL of citrated plasma, and 10 grams feces	1-4 days
CL510NE	Equine Health Program (no EIA) Equine Chemistry, CBC, Fibrinogen, FEC (MST)	1.0 mL serum, 1.0 mL EDTA whole blood, 0.5 mL of citrated plasma, and 10 grams feces	1-2 days
CL611E	Equine Performance Horse Health Program GVL (ELISA) Equine Performance Horse Chemistry, CBC, Fibrinogen, EIA (ELISA) by GVL, FEC (MST)	2.0 mL of serum, 1.0 mL EDTA whole blood, 0.5 mL of citrated plasma, and 10 grams feces	1-4 days
CL611NE	Equine Performance Horse Health Program (no EIA) Equine Performance Horse Chemistry, CBC, Fibrinogen, FEC (MST)	1.0 mL serum, 1.0 mL EDTA whole blood, 0.5 mL citrated plasma, and 10 grams feces	1-2 days
CL616E	Equine Senior Health Program GVL (ELISA) Equine Chemistry, CBC, Fibrinogen, Endogenous ACTH, Insulin, EIA (ELISA) by GVL, FEC (MST) Note: see CL525 and CT470E for additional drawing instructions.	3.0 mL of serum, 1.0 mL EDTA whole blood, 0.5 mL of citrated plasma, 10 grams feces, and 1.0 mL EDTA plasma	1-4 days
CL615NE	Equine Senior Health Program (no EIA) Equine Chemistry, CBC, Fibrinogen, Endogenous ACTH, Insulin, FEC (MST) Note: see CL525 and CT470E for additional drawing instructions.	3.0 mL of serum, 1.0 mL EDTA whole blood, 0.5 mL of citrated plasma, 10 grams feces, and 1.0 mL EDTA plasma	1-3 days
CL616T	Equine Senior Health Program GVL with TRH Equine Chemistry, CBC, Fibrinogen, TRH Stimulation Test for PPID, Insulin, EIA (ELISA) by GVL, FEC (MST) Note: see CL535 and CT470E for additional drawing instructions.	3.0 mL of serum, 1.0 mL EDTA whole blood, 0.5 mL of citrated plasma, 10 grams feces, 1.0 mL Pre EDTA plasma and 1.0 mL Post EDTA plasma	1-4 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CL121	Equine Infectious Anemia (EIA) ELISA by GVL Note: <ul style="list-style-type: none">• Equine Infectious Anemia (EIA) testing submission must include GlobalVetLINK (GVL) form• GVL form must include the following to be processed (District information, Owner's Name and Address, and Location of animal at time of the blood draw)• Valid veterinary license for the submitting veterinarian will be verified at time of submission• Sample tube labelling must match official form by either horse's name or tube number EXACTLY	1.0 mL serum in red top or serum separator tube	1-4 days
CBC AND CHEMISTRY			
CT332	CBC and Fibrinogen Interferences: marked hemolysis or lipemia.	1.0 mL EDTA whole blood and 0.5 mL citrated plasma (preferred for Fibrinogen)	Daily Performed each shift
CL365	Fibrinogen	0.5 mL citrated plasma collected as whole blood in blue top tube at least 2/3rds full to the fill line	Daily Performed each shift
CL070	Equine Inflammatory Profile CBC, Fibrinogen, Plasma Protein Interferences: marked hemolysis or lipemia.	1.0 mL EDTA whole blood and 0.5 mL citrated plasma	1-2 days
CT415	PT and aPTT Prothrombin Time, Activated Partial Thromboplastin Time Prothrombin Time and Activated Partial Thromboplastin Time. Prothrombin time measures the integrity of the extrinsic and common components of the coagulation cascade. Partial Thromboplastin Time (PTT) measures the integrity of the intrinsic and common components of the coagulation cascade. Interferences: marked hemolysis or lipemia. Partially full blue top tube may falsely increase sample's coagulation time. Note: clotting of the sample may preclude the analysis. The blue top tube needs to be filled to 2/3 or more of its capacity. Partially full blue top tubes may falsely increase coagulation times.	0.5 mL citrated plasma or citrated whole blood in blue top tube at least 2/3rds full to the fill line	Daily Performed each shift
Add-on Equivalent CADD290			

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CL010	Equine Chemistry Total Protein, Albumin, Globulin, A/G Ratio, AST (SGOT), Alk Phos, GGT, T. Bilirubin, D. Bilirubin, BUN, Creatinine, BUN/Creatinine Ratio, Phosphorus, Glucose, Calcium, Sodium, Potassium, NA/K Ratio, Chloride, Cholesterol, Triglyceride, CPK, LDH Panel is recommended for all large animal patients except for Bovine (use code CL1001). Sample Handling: serum separator samples should be spun prior to submission. Red top tubes should be spun and serum transferred to a non-additive tube and marked as SERUM. If the serum is not separated from the red blood cells, chemistry values can be affected.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CL601	Equine Performance Horse Chemistry Total Protein, Albumin, Globulin, A/G Ratio, AST (SGOT), Alk Phos, GGT, BUN, Creatinine, BUN/Creatinine Ratio, Sodium, Potassium, Chloride, CPK Sample Handling: serum separator samples should be spun prior to submission. Red top tubes should be spun and serum transferred to a non-additive tube and marked as SERUM. If the serum is not separated from the red blood cells, chemistry values can be affected.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CL050	Equine Chemistry, CBC	0.5 mL serum in red top or serum separator tube and 1.0 mL EDTA whole blood	Daily Performed each shift
CL040	Equine Chemistry, CBC, Fib Panel is recommended for all large animal patients except for Bovine (use code CL080).	0.5 mL serum in red top or serum separator tube, 1.0 mL EDTA whole blood, and 0.5 mL citrated plasma	Daily Performed each shift
CL040R	Equine Chemistry, CBC, Fib (Recheck) Note: resubmission must be within 30 days of original accession and previous accession number provided on the new request form.	0.5 mL serum in red top or serum separator tube, 1.0 mL EDTA whole blood, and 0.5 mL citrated plasma	Daily Performed each shift
CL640	Equine Performance Horse Panel Equine Performance Horse Chemistry, CBC, Fibrinogen	1.0 mL of serum in red top or serum separator tube, 1.0 mL EDTA whole blood, and 0.5 mL of citrated plasma	Daily Performed each shift
CL565	Equine Chemistry, CBC, Fib, T4, fT4 (ED) Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results.	2.0 mL serum in red top or serum separator tube, 1.0 mL EDTA whole blood, and 0.5 mL citrated plasma	1-4 days
CL035	Equine Chemistry, CBC, Fib, fT4 (ED) *Plasma is NOT an acceptable sample	1.5 mL serum in red top or serum separator tube, 1.0 mL EDTA whole blood, and 0.5 mL citrated plasma	1-4 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CL030	Equine Chemistry, CBC, Fib, T4 Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results.	0.5 mL serum in red top or serum separator tube, 1.0 mL EDTA whole blood, and 0.5 mL citrated plasma	Daily Performed each shift
CL290	Neonatal Foal Panel Equine Chemistry, CBC, Fibrinogen, Equine IgG Total	1.0 mL serum in red top or serum separator tube, 1.0 mL EDTA whole blood, and 0.5 mL citrated plasma	1-3 days
CL090	Equine IgG Total IgG is generally used to assess immunoglobulin concentration and passive transfer of immunity in foals or immunoglobulin status in adult horses. This test is species specific for horses. Note: IgG testing is species specific (see Large Animal Section for all non-equine species) Methodology: Immunoturbidimetric method	0.5 mL serum in red top or serum separator tube, 1.0 mL EDTA whole blood, or plasma from a EDTA whole blood sample *If submitting milk transfer to a non-additive tube and clearly marked as MILK	1-3 days
CL225	Equine Hepatic Screen Total Protein, Albumin, Globulin, AST (SGOT), Alk Phos, GGT, T. Bilirubin, D. Bilirubin, I. Bilirubin, BUN, Sodium, Potassium, Chloride, LDH Sample Handling: serum separator samples should be spun prior to submission. Red top tubes should be spun and serum transferred to a non-additive tube and marked as SERUM. If the serum is not separated from the red blood cells, chemistry values can be affected.	0.5 mL serum in red top or serum separator tube	Daily Perfomed each shift
CL230	Equine Hepatic Panel Equine Hepatic Screen with Electrolytes, Bile Acids	1.0 mL serum in red top or serum separator tube	1-2 days
CT225	Bile Acids Interferences: marked hemolysis or lipemia.	0.5 mL serum in red top or serum separator tube	1-2 days
CL1025	Equine Hepatic Profile Equine Hepatic Screen, Bile Acids, SDH Note: SDH sample must be kept cold or frozen and submitted within 24 hrs.	1.5 mL serum in red top or serum separator tube AND 0.5 mL separated serum in non-additive tube marked SERUM (Cold or Frozen) for SDH testing	1-5 days
CT250	Sorbital Dehydrogenase (SDH) Note: sample must be submitted within 24 hrs.	0.5 mL separated serum in non-additive tube marked SERUM (Cold or Frozen)	2-5 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CL240	Equine Renal Screen Total Protein, Albumin, BUN, Creatinine, Phosphorus, Glucose, Calcium, Sodium, Potassium, Chloride Sample Handling: serum separator samples should be spun prior to submission. Red top tubes should be spun and serum transferred to a non-additive tube and marked as SERUM. If the serum is not separated from the red blood cells, chemistry values can be affected.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CL275	Equine Muscle Enzyme Screen AST (SGOT), CPK Sample Handling: serum separator samples should be spun prior to submission. Red top tubes should be spun and serum transferred to a non-additive tube and marked as SERUM. If the serum is not separated from the red blood cells, chemistry values can be affected.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CL280	Equine Muscle Enzyme Screen, CBC, FIB CBC, AST (SGOT), CPK, Fibrinogen	1.0 mL of serum in red top or serum separator tube, 1.0 mL EDTA whole blood, and 0.5 mL of citrated plasma	Daily Performed each shift
CL190	Equine Rhabdomyolysis Screen AST (SGOT), BUN, Creatinine, Phosphorus, Calcium, Sodium, Potassium, Chloride, CPK, LDH Interferences: marked hemolysis and lipemia. Sample Handling: serum separator samples should be spun prior to submission. Red top tubes should be spun and serum transferred to a non-additive tube and marked as SERUM. If the serum is not separated from the red blood cells, chemistry values can be affected.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CT240 Add-on Equivalent CADD130	Protein Electrophoresis (Serum) Total Protein, Albumin, Globulin, Alpha 1, Alpha 2, Beta 1, Gamma 1, fractional assessment with interpretation. An evaluation of the globulin fraction of the serum (alpha 1, alpha 2, beta, and gamma) to determine if the globulin fraction is monoclonal based on these components. Interferences: marked hemolysis or lipemia. Hemolysis may falsely elevate the beta globulin fraction.	0.5 mL serum in red top or serum separator tube	2-4 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS1680	Equine Drug Screen Drug Screen (performed at TVMDL) includes the following: Corticosteroids (LC/MS): Betamethasone, dexamethasone, methylprednisolone, prednisolone, prednisone and triamcinolone acetonide. NSAIDs: Acetaminophen, acetylsalicylic acid, carprofen, celcoxib, deracoxib, diclofenac, diflunisal, eltenac, ethacrynic acid, etodolac, fenbufen, fenoprofen, firocoxib, flufenamic acid, flunixin, flurbiprofen, ibuprofen, indomethacin, indoprofen, ketoprofen, ketorolac, meclofenamic acid, mefanamic acid, meloxicam, nabumetone, naproxen, oxyphenbutazone, phenylbutazone, piroxicam, salicylic acid, tenoxicam, tolfenamic acid and tolmetin. Tranquilizers: Acepromazine, fluphenazine, fluoxetine, guanabenz, reserpine, romifidine, xylazine, and detomidine Muscle relaxant: Methocarbamol. Note: urine is the preferred specimen but serum is acceptable.	10 mL urine in urine transport tube OR 5.0 mL serum in red top or serum separator (Mark accordingly)	7-10 days
URINE			
CT760	Urinalysis-Complete The complete urinalysis includes a physical (color, appearance, USG), chemical (pH, protein, glucose, ketones, bilirubin, occult blood), and microscopic (WBC, RBC, casts crystals, bacteria, transitional and squamous epithelial cells, fat droplets) exam of the urine. Interferences: visible levels of hemolysis, drugs containing dyes, nitrofurantoin, or riboflavin.	6.0 mL urine in urine transport tube	Daily Performed each shift
CL340	Fractional Excretion of Electrolytes Serum: Calcium, Chloride, Creatinine, Phosphorus, Sodium, Potassium Urine: Calcium, Chloride, Creatinine, Phosphorus, Sodium, Potassium Note: serum and urine must be collected within 24 hours of each other and submitted together.	1.0 mL serum in red top or serum separator tube AND 1.0 mL urine in urine transport tube (mark each sample as URINE or SERUM)	1-2 days
FECAL			
CT826	FEC:Modified Stoll's Technique (MST) <i>Strongyle</i> sp. and <i>Parascaris</i> sp. ova (EPG) performed by Modified Stoll's Technique (MST) and qualitative ID of all other parasites Fecal Egg Count by Modified Stoll's Technique test provides quantitative eggs counts for <i>Strongyle</i> & <i>Parascaris</i> ova reported as eggs per gram (EPG). The lowest detection limit is 2 EPG. Other eggs or parasites seen will be reported. Note: fresh fecal samples should be kept refrigerated until submission. Specimens should be tested within 72 hours of collection.	10 grams feces in Antech provided fecal container	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CT828	FEC: McMaster's Method (MM) - Equine <i>Strongyle</i> sp. and <i>Parascaris</i> sp. ova (EPG) performed by McMaster's Method (MM) and qualitative ID of all other parasites Fecal Egg Count by McMaster's Method test provides an estimated egg count for <i>Strongyle</i> & <i>Parascaris</i> ova reported as Eggs per gram (EPG). The lowest detection limit is 100 EPG. Other eggs or parasites seen will be reported. Note: fresh fecal samples should be kept refrigerated until submission. Specimens should be tested within 72 hours of collection.	10 grams feces in Antech provided fecal container	1-2 days
CL86181	<i>Clostridium difficile</i> toxins A/B Note: see CL950 Equine PCR Gastrointestinal (Diarrhea) Panel for a broader array of potential pathogens.	5 grams feces in Antech provided fecal container Send on ice to keep cold	2-5 days
CT16007	<i>Clostridium perfringens</i> enterotoxin Note: see CL950 Equine PCR Gastrointestinal (Diarrhea) Panel for a broader array of potential pathogens.	5 grams feces in Antech provided fecal container Send on ice to keep cold	2-4 days
CM160	Culture, Feces Culture specifically evaluates for <i>Salmonella</i> , <i>Shigella</i> and <i>Campylobacter</i> spp. Note: see CL950 Equine PCR Gastrointestinal (Diarrhea) Panel for a broader array of potential pathogens.	5 grams feces in Antech provided fecal container.	3-4 days Preliminary report every 24 hours. Final report available in 72 hours.
CM121	Culture, <i>Salmonella</i> Negative <i>Salmonella</i> spp. culture result does not rule out intermittent or lower levels of fecal shedding of <i>Salmonella</i> spp. Five-day serial submissions are recommended.	5 grams feces in Antech provided fecal container	3-4 days Preliminary report every 24 hours. Final report available in 72 hours.
CL496	Acute Diarrhea Panel Fecal Culture, <i>Clostridium perfringens</i> enterotoxin, <i>Clostridium difficile</i> toxins A/B	10 grams feces in Antech provided fecal container Send on ice to keep cold	2-4 days
CL492	Foal Diarrhea Panel FEC (MST), Fecal Culture, <i>Clostridium perfringens</i> enterotoxin, <i>Clostridium difficile</i> toxins A/B	10 grams feces in Antech provided fecal container Send on ice to keep cold	2-4 days
CL420	Chronic Diarrhea Panel Equine Chemistry, CBC, Fibrinogen, FEC (MST), Fecal Culture, <i>Clostridium perfringens</i> enterotoxin, <i>Clostridium difficile</i> toxins A/B	1.0 mL serum in red top or serum separator tube, 1.0 mL EDTA whole blood, 0.5 mL citrated plasma, and 10 grams feces Send on ice to keep cold	1-4 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
ENDOCRINOLOGY			
CT495	Total T4	0.5 mL serum in red top or serum separator tube	Daily
Add-on Equivalent CADD190	Interferences: marked hemolysis and moderate to marked lipemia. Lipemia can falsely decrease T4 results.		Performed each shift
CT460	Free T4 by Equilibrium Dialysis	0.5 mL serum in red top or serum separator tube	2-6 days
Add-on Equivalent CADD50	Note: this test should not be performed as an Add-on to samples older than 5 days.		
CSA370	Thyroid Profile 2	1.0 mL serum in red top or serum separator tube	1-6 days
	Total T4, Free T4 By Equilibrium Dialysis		
	Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results.		
CL590	Equine TRH Stimulation (Thyroid Function)	0.5 mL serum labeled Pre, 0.5 mL serum labeled Post T3, and 0.5 mL serum labeled Post T4. All samples should be in red top or serum separator tube.	4-5 days
	Baseline T3, Baseline T4, 2-hr Post T3, 4-hr Post T4		
	Drawing Instructions: <ul style="list-style-type: none">• Draw a baseline serum sample for T3 and T4• Inject 1 mg of TRH IV• Draw a 2 hour post for T3• Draw a 4 hour post for T4• MARK ALL SAMPLES ACCORDINGLY		
CL500	Pituitary Pars Intermedia Dysfunction (PPID) Panel	1.0 mL serum in red top or serum separator tube, 1.0 mL EDTA whole blood, 0.5 mL citrated plasma, and 1.0 mL EDTA plasma	1-3 days
	Equine Chemistry, CBC, Fibrinogen, Endogenous ACTH (Equine), Insulin (Equine), T4		
	Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results.		
	Note: see CL525 and CT470E for additional drawing instructions.		
CL500TRH	Pituitary Pars Intermedia Dysfunction (PPID) with TRH Stim	1.0 mL serum in red top or serum separator tube, 1.0 mL EDTA whole blood, 0.5 mL citrated plasma, 1.0 mL Pre EDTA plasma and 1.0 mL Post EDTA plasma	1-3 days
	Equine Chemistry, CBC, Fibrinogen, TRH Stimulation for PPID, Insulin (Equine), T4		
	Interferences: T4 evaluation can be affected in a hemolyzed or moderately lipemic sample. Lipemia can falsely decrease T4 results.		
	Note: see CL535 and CT470E for additional drawing instructions.		
CL560	Pituitary Pars Intermedia Dysfunction (PPID) with fT4ED	2.0 mL serum in red top or serum separator tube, 1.0 mL EDTA whole blood, 0.5 mL citrated plasma, and 1.0 mL EDTA plasma	1-6 days
	Equine Chemistry, CBC, Fibrinogen, Endogenous ACTH (Equine), Insulin (Equine), T4, Free T4 by Equilibrium Dialysis		
	Note: see CL525 and CT470E for additional drawing instructions.		

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CL525	Endogenous ACTH (Baseline) Equine Plasma ACTH levels show a seasonal variation for horses in the northern hemisphere. Seasonally elevated ACTH levels are usually observed between Mid-July to Mid-November. Horses with PPID may have exaggerated levels during this time. Interpretation range values will vary by season. Sample Handling: <ul style="list-style-type: none"> • 12 hour fast with one flake of hay left in the stall • Draw sample using EDTA whole blood tube • Spin sample within 8 hours of draw (lab will not run the test on unspun samples) • Transfer the plasma to new non-additive tube and mark as PLASMA • Send sample on ice pack 	1.0 mL EDTA plasma in non-additive transport tube (labeled as EDTA plasma)	1-2 days
CL540	Endogenous ACTH, Insulin, Glucose Note: see CL525 and CT470E for additional drawing instructions.	1.0 mL serum in red top or spun serum separator tube and 1.0 mL EDTA plasma (mark accordingly)	1-2 days
CL575	Endogenous ACTH, Insulin, Glucose, Leptin (Equine) Note: see CL525 and CT470E for additional drawing instructions.	2.0 mL serum in red top or spun serum separator tube and 1.0 mL EDTA plasma (mark accordingly)	1-10 days
CS14402	Leptin (Equine) Note: transport refrigerated (on ice packs)	2.0 mL serum in red top or non-additive tube	7-10 days
CL535	TRH Stimulation for PPID Endogeneous ACTH (Pre and 10 minute Post) Tests for Cushings Disease (PPID). If Endogenous ACTH results (CL525) are normal but disease is strongly suspected, this panel is recommended. Interferences: serum or whole blood submissions will preclude analysis and will NOT be accepted. Note: see CL525 for additional drawing instructions.	1.0 mL EDTA plasma (labeled as pre plasma) and 1.0 mL EDTA plasma (labeled as post plasma) Send on ice to keep cold	1-2 days
CT470E	Equine Insulin and Glucose In horses, a high baseline insulin concentration is indicative of Insulin Dysregulation (ID). However, many horses will have a normal baseline insulin with an abnormal response to sugar ingestion. If baseline results are normal but disease is strongly suspected, an oral sugar test (OST) is recommended see code CL545 or CL550. Interferences: marked hemolysis. Sample Handling: <ul style="list-style-type: none"> • 12 hour fast with one flake of hay left in the stall • Draw sample in a red top or serum separator tube • Spin sample within 8 hours of draw • Transfer serum into new non-additive tube and mark as SERUM 	0.5 mL serum in red top or spun serum separator tube	1-2 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CL580	TRH Stimulation for PPID with Insulin and Glucose Note: see CL525 and CT470E for additional drawing instructions.	1.0 mL serum in red top or spun serum separator tube, 1.0 mL Pre EDTA plasma and 1.0 mL Post EDTA plasma (mark accordingly)	1-2 days
CL545	Oral Sugar Test (2 samples) Insulin (Equine)/Glucose (60 and 90 minute Post) Oral sugar test assesses the insulin response to ingested sugars for a diagnosis of Insulin Dysregulation (ID). Note: see CT470E for additional drawing instructions.	1.0 mL serum in red top or spun serum separator tube (labeled 60 min Post) AND 1.0 mL serum in red top or spun serum separator tube (labeled 90 min Post)	1-2 days
CL550	Oral Sugar Test (3 samples) Insulin (Equine)/Glucose (Pre, 60 and 90 minute Post) Oral sugar test assesses the insulin response to ingested sugars for a diagnosis of Insulin Dysregulation (ID). Note: see CT470E for additional drawing instructions.	1.0 mL serum in red top or spun serum separator tube (labeled Baseline), 1.0 mL serum in red top or spun serum separator tube (labeled Serum 60 min Post), AND 1.0 mL serum in red top or spun serum separator tube (labeled Serum 90 min Post)	1-2 days
NEUROLOGY			
CS14388	EPM SAG 2,4/3 ELISA The ELISA format contains surface antigen proteins (SAG 2,4/3) of the <i>S. neurona</i> parasite. A negative serum result is a good rule out for EPM caused by <i>S. neurona</i> . Note: if submitting BOTH serum and CSF, should be submitted separately as two different requests in two different bags.	1.0 mL serum in red top or serum separator OR 1.0 mL CSF in red top tube	4-7 days
CS14390	EPM SAG 2,4/3 ELISA Serum/CSF Ratio The specific antibody index is performed on paired serum and CSF samples in conjunction with the <i>S. neurona</i> SAG 2,4/3 ELISA. This test can help discern an intrathecal IgG response from contamination by utilizing a ratio of serum and CSF albumin and IgG levels. Note: both serum and CSF MUST be submitted together in the same bag.	1.0 mL serum in red top or serum separator AND 1.0 mL CSF in red top tube	4-7 days
CS14392	EPM SAG 2,4/3 ELISA, N. hughesi ELISA The <i>Neospora hughesi</i> ELISA utilizes a specific protein, SAG 1, found on the <i>N. hughesi</i> parasite. Serum and/or CSF samples can be tested for an IgG response to the recombinant NhSAG 1 protein. The ELISA generates an endpoint titer. Note: if submitting BOTH serum and CSF, should be submitted separately as two different requests in two different bags.	1.0 mL serum in red top or serum separator OR 1.0 mL CSF in red top tube	4-7 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS16275	Equine Encephalitis Viral Panel EEE, WEE, VEE (PRNT) and EEE (IgM Capture ELISA)	1.0 mL serum in red top or serum separator tube	10-15 days
CS17500	Equine Encephalitis Viral Panel Plus Equine Encephalitis Viral Panel, West Nile titer (PRNT/IgM), EHV-I Ab	3.0 mL serum in red top or serum separator tube	10-15 days
CS14477	Equine Comprehensive Neurological Panel Equine Encephalitis Viral Panel, EPM SAG 2, 4/3 ELISA, West Nile Titer (PRNT/IgM), EHV-I PCR	3.0 mL serum in red top or serum separator tube. For EHV-1 PCR component: 5.0 - 7.0 mL EDTA whole blood in lavender top tube, 5.0 mL nasal wash in red top tube, and/or nasal swab in non-additive tube with no media marked NASAL	5-15 days
CS7591	Equine Rhinopneumonitis virus Ab (EHV-1)	1.0 mL serum in red top or serum separator tube	10-15 days
CS85448	West Nile Titer (PRNT/IgM ELISA)	1.0 mL serum in red top or serum separator tube	7-15 days
PCR			
CL950	Equine PCR Gastrointestinal (Diarrhea) Panel <i>Clostridium difficile</i> A & B, <i>Clostridium perfringens</i> , <i>C. perfringens</i> toxin A, <i>C. perfringens</i> toxin B, <i>C. perfringens</i> toxin E, <i>C. perfringens</i> toxin NetF, Equine Coronavirus, <i>N. risticii</i> (Potomac Horse Fever), Equine Rotavirus, and <i>Salmonella</i> spp. Interferences: wooden swabs can interfere with testing and plastic swabs are preferred. Note: A single negative <i>Salmonella</i> spp. PCR result does NOT rule out intermittent or lower levels of fecal shedding of <i>Salmonella</i> spp. Three (3) serial submissions are recommended (use code CS14416 for sample 2 and/or sample 3).	5 grams feces in Antech provided fecal container	4-7 days
CL955	Equine PCR Reproductive Panel EVA, EHV-1, <i>Strep equi</i> subsp. <i>zooepidemicus</i> , <i>Leptospira</i> spp., <i>Klebsiella pneumoniae</i> , <i>E. coli</i> Interferences: wooden swabs can interfere with testing and plastic swabs are preferred.	Uterine lavage fluid, semen, endometrial or cervical swab in red top tube (without media), OR placenta	4-7 days
CL960	Equine PCR Respiratory Panel <i>Streptococcus equi</i> , <i>Streptococcus zooepidemicus</i> , <i>Rhodococcus equi</i> , Equine Herpes Virus 1 (EHV1), Equine Herpes Virus 4 (EHV4), Equine Influenza Virus (EIV), Equine Rhinitis Virus A and B (ERAV & ERBV) Interferences: wooden swabs can interfere with testing and plastic swabs are preferred.	Transtracheal wash, bronchoalveolar lavage, guttural pouch wash, nasal swab (in non-additive tube without media) or respiratory tract tissue	4-7 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CL965	Equine PCR Fever of Unknown Origin (FUO) Panel Swab: <i>Streptococcus equi</i> , Equine Herpes Virus 1, Equine Herpes Virus 4, Equine Influenza Virus, Equine Rhinitis Virus A & B EDTA whole blood: Equine Herpes Virus-1, <i>Anaplasma phagocytophilum</i> , <i>Neorickettsia risticii</i> (PHF) Fecal: <i>Neorickettsia risticii</i> (PHF), Equine Coronavirus (ECoV) Interferences: wooden swabs can interfere with testing and plastic swabs are preferred.	1.0 mL EDTA whole blood, nasal swab in non-additive tube (without media), AND 5 grams feces in Antech provided fecal container	4-7 days
CL970	Equine PCR Fever of Unknown Origin (FUO) Mini EDTA whole blood: <i>Anaplasma phagocytophilum</i> , <i>Neorickettsia risticii</i> (PHF) Fecal: <i>Neorickettsia risticii</i> (PHF), Equine Coronavirus (ECoV)	1.0 mL EDTA whole blood AND 5 grams feces in Antech provided fecal container	4-7 days
CS14421	<i>Anaplasma phagocytophilum</i> PCR	1.0 mL EDTA whole blood	4-7 days
CS14414	Equine Coronavirus (ECoV) PCR	5 grams feces in Antech provided fecal container	4-7 days
CS14394	Equine Herpes Virus 1 (EHV1) PCR Viremia can be detected in the blood as soon as four days post-infection. Viral shedding in nasal secretions can be detected as soon as six days post-infection and for a longer time than in blood. Sample Handling: submission of BOTH nasal wash/nasal swab AND EDTA whole blood is recommended.	5.0 mL EDTA whole blood, 5.0 mL nasal wash in red top tube, and/or nasal swab in non-additive tube without media	4-7 days
CL974	Equine <i>Leptospira</i> PCR Blood	0.5 mL EDTA whole blood	3-5 days
CL976	Equine <i>Leptospira</i> PCR Urine	2.0 mL urine in urine transport tube	3-5 days
CL978	Equine <i>Leptospira</i> PCR Blood/Urine	0.5 mL EDTA whole blood and 2.0 mL of urine in urine transport tube	3-5 days
CL973	Equine <i>Leptospira</i> PCR Aqueous or Vitreous Humor	1.5 mL aqueous or vitreous humor fluid in lavender top tube	3-5 days
CS14479	<i>Neorickettsia risticii</i> (PHF) PCR PCR testing will be done on both whole blood and fecal samples for no additional charge when both samples are submitted.	1.0 mL EDTA whole blood in lavender top tube and/or 5 grams feces in Antech provided fecal container	4-7 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS14396	<i>Rhodococcus equi</i> PCR Identifies strains of <i>R. equi</i> carrying virulence plasmid gene vapA. Transtracheal wash sample is the recommended specimen. Interferences: wooden swabs can interfere with testing and plastic swabs are preferred. A fecal sample can be tested but negative results will be unreliable.	Transtracheal wash, bronchoalveolar lavage, and/or nasal swab in non-additive tube without media.	4-7 days
CT983	Equine Ringworm PCR Panel <i>Microsporum</i> spp., <i>M. equinum</i> / <i>canis</i> , <i>M. gypseum</i> (<i>Arthroderma gypseum</i> , <i>A. fulvum</i> , <i>A. incurvatum</i>), <i>Trichophyton</i> spp., <i>T. equinum</i> / <i>mentagrophytes</i> (<i>Arthroderma benhamiae</i> , <i>A. vanbreuseghemii</i>)	Minimum of 10 plucked hair with roots, skin scraping, or tooth brush sample in sterile, dry container free of liquids or preservatives	5-7 days
CT987	Equine Ringworm PCR with Dermatophyte Culture <i>Microsporum</i> spp., <i>M. equinum</i> / <i>canis</i> , <i>M. gypseum</i> (<i>Arthroderma gypseum</i> , <i>A. fulvum</i> , <i>A. incurvatum</i>), <i>Trichophyton</i> spp., <i>T. equinum</i> / <i>mentagrophytes</i> (<i>Arthroderma benhamiae</i> , <i>A. vanbreuseghemii</i>), and Dermatophyte Culture	Minimum of 15 plucked hair with roots, skin scraping, or tooth brush sample in sterile, dry container free of liquids or preservatives	5-21 days
CS14416	<i>Salmonella</i> spp. PCR Detects DNA from pathogenic species of <i>Salmonella</i> . <i>Salmonella</i> PCR testing can also be used as a bio-surveillance tool to identify asymptomatic/shedding horses and monitor in-patients in a hospital setting. Note: a single negative <i>Salmonella</i> spp. PCR result does not rule out intermittent or lower levels of fecal shedding of <i>Salmonella</i> spp. Three (3) serial submissions are recommended.	5 grams feces in Antech provided container	4-7 days
CS86308	<i>Streptococcus equi</i> PCR Detects <i>S. equi</i> bacterial DNA and is used primarily to identify asymptomatic carriers. Nasal washes, nasal swabs, or guttural pouch washes are appropriate samples for the <i>S. equi</i> PCR. Back-up culture is also available. Interferences: wooden swabs can interfere with testing and plastic swabs are preferred.	5.0 mL nasal wash, OR 5.0 mL guttural pouch wash, OR 1 pharyngeal swab in non-additive tube without media	4-7 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
MARE REPRODUCTION			
CS16295	Estradiol Recommended for testing if greater than 120 days post-breeding in conjunction with progesterone (CL140) at the same time using code CL470.	1.0 mL serum in red top or serum separator tube	7-10 days
CS16300	Estrone Sulfate Recommended for testing if greater than 100 days post breeding.	1.0 mL serum in red top or serum separator tube	7-10 days
CS16635	Pregnant Mare Serum Gonadotropin (PMSG) Recommended for testing at day 45-120 post breeding in conjunction with progesterone (CL140) using test code CL460.	2.0 mL serum in red top or serum separator tube	7-9 days
CL140	Progesterone	0.5 mL serum in red top tube. Submission of serum separator tube is not recommended.	1-2 days
CL460	Progesterone/PMSG	2.5 mL serum in red top tube. Submission of serum separator tube is not recommended.	1-9 days
CL470	Progesterone/Estradiol Panel (Equine)	1.5 mL serum in red top tube. Submission of serum separator tube is not recommended.	1-10 days
CL480	Progesterone/PMSG/Estradiol	3.0 mL serum in red top tube. Submission of serum separator tube is not recommended.	1-10 days
GRANULOSA CELL TUMOR TESTING			
CS85857	Equine Granulosa Cell Tumor Progesterone, Testosterone, Inhibin	3.0 mL serum in red top tube. Submission of serum separator tube is not recommended.	2-21 days
CS4131	Inhibin Inhibin is a heat-sensitive protein. Interferences: marked hemolysis or lipemia.	1.0 mL serum in red top tube or serum separator tube	2-3 Weeks
CS14320	Anti-Mullerian Hormone (Equine)	1.0 mL serum in red top or serum separator tube	5-10 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
STALLION/GELDING REPRODUCTION			
CS16760	Testosterone Single testosterone level.	0.5 mL serum in red top or serum separator tube	3-5 days
CSEMEN	Semen Analysis Sperm Count (including Total Volume and pH), with Cytologic Evaluation	2.0 mL semen (fresh or extended) in lavender top tube	1-2 days
CRYPTORCHID TESTING			
CS16300	Estrone Sulfate Valid for testing if greater than 3 years of age	1.0 mL serum in red top or serum separator tube	7-10 days
CS14320	Anti-Mullerian Hormone (Equine) Valid for testing if less than 3 years of age	1.0 mL serum in red top or serum separator tube	5-10 days
CS85530	Testosterone Panel Pre, 20 min post, 1 hr post, 2 hr post HCG Used to evaluate horses for cryptorchids and granulosa cell tumors. Equine HCG Stimulation Testing Protocol: 1. Take baseline serum sample for testosterone. 2. Administer 6000 IU (total dose) of HCG IV. 3. Take post-HCG testosterone at 20 minutes, 1 hour and 2 hour post administration. 4. If testicular tissue is present, the post values should be greater than 0.4 nmol/L and/or double the baseline level.	0.5 mL serum for baseline and each post sample (mark each sample with time)	3-5 days
MICROBIOLOGY			
CM020	Culture, Aerobic This test is used when an aerobic bacterial infection is suspected in a tissue or fluid. Interferences: 1. Patient should be off antibiotics for at least 7-10 days. 2. Fluid in lavender top tube is unacceptable (EDTA inhibits growth). Note: tissue should be submitted in an RTT with a few drops of saline (to keep moist). Samples should be refrigerated prior to transportation to the laboratory. A separate sample should be submitted for each type of culture needed.	Fluid (body cavity fluids), TTW or BAL, sample from wound or lesion in red top tube, or tissue sample in saline. Submission in lavender top tube (with EDTA) is not acceptable.	3-4 days Preliminary report available every 24 hours. Final culture result available in 72 hours. If a fastidious organism is observed, the listed turnaround time may be extended.

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CM030	<p>Culture, Anaerobic</p> <p>This test is used when an anaerobic bacterial infection is suspected in a tissue or fluid.</p> <p>Interferences: sample with exposure to air (or sample that has dried) will preclude accurate testing.</p> <p>Note: fluids should be submitted in an RTT. Tissue should be submitted in an air-tight, sterile container (e.g., urine cup). No saline is necessary. Sensitivities are not performed, but drugs of choice for anaerobes are provided as a guide.</p>	<p>Fluid (body cavity fluids), TTW or BAL, or sample from wound or lesion in red top tube. Tissue sample (at least 2 cm x 2 cm) in sterile, air tight container. Submission in lavender top tube (with EDTA) is not acceptable.</p>	<p>3-5 days</p> <p>Preliminary report available 3rd day. Final culture result available 4th day</p>
CM040	<p>Culture, Aerobic and Anaerobic</p> <p>This test is used when a bacterial infection is suspected in a tissue or fluid, but it is uncertain whether an aerobic or an anaerobic organism is the cause.</p> <p>Interferences:</p> <ol style="list-style-type: none">1. Patient should be off antibiotics for at least 7-10 days.2. Fluid in lavender top tube is unacceptable (EDTA inhibits growth). <p>Note: fluids should be submitted in an RTT. In the case of an aerobic culture, a copan swab can be submitted. In the case of an anaerobic culture, an anaerobic culturette swab can be submitted. Tissue should be submitted in an air-tight container for anaerobic culture. Tissue should be submitted in a container with a few drops of saline in the case of an aerobic culture. Samples should be refrigerated prior to submission. Sensitivities are not performed, but drugs of choice for anaerobes are provided as a guide.</p>	<p>Fluid (body cavity fluids), TTW or BAL, sample from wound or lesion in red top tube, or tissue sample (2 x 2 cm) in saline. Submission in lavender top tube (with EDTA) is not acceptable. Sample submitted for anaerobic culture should not be dry or exposed to air</p>	<p>3-5 days</p>
CM060	<p>Culture, Blood Aerobic Only</p> <p>This test is used when bacteremia caused by an aerobic bacteria is suspected. Submit a single sample in a BACTEC bottle.</p> <p>Interferences: antibiotics.</p> <p>Note: submit in BACTEC aerobic culture bottle. Clip fur and scrub venipuncture site for aseptic collection. Aim to collect 1.0 mL for cats and small dogs, 2-3 mL for larger dogs to be inoculated into 20 mL bottle. Don't unscrew caps on bottles. Remove the protective top and wipe visible parts of the rubber stopper with 70% ethanol. Allow stopper to dry or wipe with sterile gauze. Replace the drawing needle with a sterile needle before puncturing the rubber stopper, fill until vacuum stops, then gently invert the bottle to mix. Anticoagulants in the media will prevent blood from clotting.</p>	<p>Whole blood collected in BD Bactec Blood Culture Bottle (BCB)</p>	<p>3-5 days</p>
CM220	<p>Culture, Enhanced Fluid Aerobic</p> <p>Interferences: patient should be off antibiotics seven to ten days prior to culture submission.</p> <p>Note: if submitting multiple joints, must be submitted separately in different bags and forms.</p>	<p>Fluid in BD Bactec Blood Culture Bottle (BCB)</p>	<p>3-4 days</p>

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CM080	Culture, Fungal This test is for fungal culture and identification. Not for suspected ringworm infections. Note: submit hair, skin, nails, and bodily fluids on swabs or in sterile containers. For cases suspected of ringworm infection, use CM240 instead or consider ringworm PCR panel. If you are submitting a DTM plate with fuzzy growth, use code CM082.	Dry hair, nails, skin scraping, body fluid, or lesion material collected in a sterile red top or other sterile container without additive, or on a culturette. DTM bottle, culturette, container	3 weeks Preliminary report every 7 days. Final report available at 3 weeks (or upon identification)
CM130 Add-on Equivalent CADD210	Culture, Urine with <i>FIRST</i>Tract <i>FIRST</i> Tract Urine Culture followed by plating in the case of a positive result for bacterial identification and sensitivity All urine cultures deemed positive by <i>FIRST</i> Tract, a highly accurate and rapid urine culture technique, will then undergo traditional culture for identification and susceptibility. Interferences: patient should be off antibiotics for seven to ten days prior to urine culture.	0.5 mL cystocentesis, clean catch or catheterized urine in sterile red top tube or urine transport tube	3-4 days Preliminary reports available every 24 hours. If a fastidious organism is observed, the listed turnaround time may be extended.
CYTOLOGY			
CCYTO	Cytology Includes preparation of submitted sample and microscopic interpretation by clinical pathologist Testing includes a microscopic evaluation of cells. The report includes cytologic interpretation, diagnosis, and comments regarding the etiology and biological behavior where applicable. Note: single source.	Fluid in red top or lavender top tube, or air-dried, unstained slides with clinical history	1-3 days
CFLUA	Fluid Analysis Includes preparation of submitted sample, Cell Count (WBC & RBC), Specific Gravity, Protein measurement, and microscopic interpretation by clinical pathologist This is the appropriate test for the evaluation of fluids from a cavity or tissue space. The fluid will be characterized based on protein, specific gravity, white and red cells, and a pathologist will evaluate the fluid and provide a microscopic description and, when possible, diagnosis. Note: a comprehensive history should be included with the fluid analysis submission.	1.0 mL body fluid in lavender top or red top with 2 unstained smears prepared from fluid	1-3 days
CCSF	Fluid Analysis, CSF Includes preparation of submitted sample, CSF Cell Count (WBC & RBC), Specific Gravity, Protein & Glucose measurement, and microscopic interpretation by clinical pathologist	1.0 mL cerebrospinal fluid in lavender or red top tube	1-3 days

*SEE ANTECHDIAGNOSTICS.COM/TERMS-SERVICE-CA FOR MORE DETAILS ON TURNAROUND TIMES

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CSYFLUA	Fluid Analysis, Synovial Synovial Fluid Analysis (Color, Clarity, Specific Gravity, WBC, RBC), Total Protein, and Cytologic Evaluation Note: single source.	1.0 mL synovial fluid in lavender top tube	1-3 days
CBAL	Brochoalveolar Lavage Cytology Includes preparation of submitted sample and microscopic interpretation by clinical pathologist Note: if culture is also being submitted, use Equine Test Form 6 to send in a separate sample in non-additive tube marked CM020.	1.0 mL airway wash fluid in red top or lavender top tube, or 2 unstained air-dried smears prepared from fluid	1-3 days
CTTW	Transtracheal Wash Cytology Includes preparation of submitted sample and microscopic interpretation by clinical pathologist Note: if culture is also being submitted, use Equine Test Form 6 to send in a separate sample in non-additive tube marked CM020.	1.0 mL airway wash fluid in red top or lavender top tube, or 2 unstained air-dried smears prepared from fluid	1-3 days
CFBX	Histopathology Preparation of submitted sample and microscopic interpretation by anatomic pathologist A boarded pathologist will evaluate the tissue submitted for Histopathology. A full written biopsy report will be provided, including source, history, description of submitted tissue, microscopic description, and diagnosis where possible. The pathologist's comments will include, where applicable, margin evaluation, grading, interpretation of special stains, and recommendations for further testing if needed. Note: firmly and evenly tighten the lid of the formalin jar and check for any leaks prior to placing the sample in a ziplock bag with the test requisition form. A comprehensive history should be included with each biopsy submission.	Tissue in 10% neutral buffered formalin in Antech approved container with screw-on lid To prevent severe biopsy damage when temperatures are below freezing, we recommend adding 1:10 ratio of isopropyl alcohol (70% or greater) to 10% neutral buffered formalin	3-5 business days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
INDIVIDUAL TESTS			
CS14466	Equine 11 Metal & Mineral Panel Serum: Cobalt, Copper, Iron, Magnesium, Molybdenum, Selenium, Zinc EDTA Whole Blood: Arsenic, Cadmium, Lead, Thallium Serum selenium provides current-day levels, while whole blood selenium indicates selenium status over a longer period of time.	One or more of the following: 1. 1.0 mL serum in non-additive tube AND EDTA whole blood 2. 5 grams fresh liver sample in sealed, leak-proof container 3. 50 mg fresh liver biopsy (3 Tru-Cut samples, no fluid added) in sealed, leak-proof container 4. 500 grams of feed (representative sample) in sealed container	5-10 days
CS16285	Equine Viral Arteritis (EVA)	1.0 mL serum in red top or serum separator tube	7-10 days
CS16510	Leptospirosis <i>L. pomona</i> , <i>L. icterohemorrhagiae</i> , <i>L. canicola</i> , <i>L. grippotyphossa</i> , <i>L. hardjo</i> , <i>L. autumnalis</i> , <i>L. bratislava</i> The results of this test include a semi-quantitative titer for serovars. Interferences: recent vaccination (within one month) may interfere with testing. Note: test is not species specific.	1.0 mL serum in red top tube or serum separator tube	5-10 days
CT672	Lyme Equine Multiplex	0.5 mL serum in red top or serum separator tube	5-10 days
CS16848	Piroplasmosis, Equine <i>Theileria (Babesia) equi</i> , <i>Babesia caballi</i> .	2.0 mL serum in red top or serum separator tube	10-12 days
CS16740	Streptococcus Equi ELISA The SeM ELISA test assesses the IgG titer against the <i>S. equi</i> specific M protein. It is helpful for making vaccination decisions on horses with existing titers, identifying horses at risk of complications due to elevated titers or for horses with aberrant abscesses (bastard strangles). The SeM ELISA should not be used to determine the infection status of horses with clinical signs of strangles.	5.0 mL serum in red top or serum separator tube	7-10 days
CS16850	Vitamin E Note: cover sample to prevent exposure to light	1.5 mL serum in red top or serum separator tube (spun)	7-10 days
CS17505	Vitamin E/Selenium	1.5 mL serum in red top or serum separator tube (spun)	7-10 days

LARGE ANIMAL

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CBC AND CHEMISTRY			
CT332	CBC and Fibrinogen Interferences: marked hemolysis or lipemia.	1.0 mL EDTA whole blood and 0.5 mL citrated plasma	Daily Performed each shift
CT365	Fibrinogen This test is used to measure the concentration of functional fibrinogen in the plasma. Interferences: clotted sample precludes analysis. Citrated Whole Blood (blue top tube) or citrated plasma is the only acceptable sample. The tube should be greater than 2/3rds filled. If submitting separated citrated plasma, label it as Citrated Plasma.	0.5 mL citrated plasma collected as whole blood in blue top tube at least 2/3rds full to the fill line	Daily Performed each shift
CL1001	Ruminant Chemistry Total Protein, Albumin, Globulin, A/G Ratio, AST (SGOT), Alk Phos, GGT, T. Bilirubin, D. Bilirubin, BUN, Creatinine, BUN/Creat Ratio, Phosphorus, Glucose, Calcium, Magnesium, Sodium, Potassium, NA/K Ratio, Chloride, Cholesterol, Triglyceride, CPK, LDH Sample Handling: serum separator samples should be spun prior to submission. Red top tubes should be spun and serum transferred to a non-additive tube and marked as SERUM. If the serum is not separated from the red blood cells, chemistry values can be affected.	0.5 mL serum in red top or serum separator tube	Daily Performed each shift
CL080	Ruminant Chemistry/CBC/Fib	0.5 mL serum, 1.0 mL EDTA whole blood, and 0.5 mL citrated plasma	Daily Performed each shift
INDIVIDUAL TESTS			
CS16035	Anaplasma CF, Bovine Interferences: severe hemolysis or lipemia will preclude testing.	2.0 mL serum in red top or serum separator tube	5-10 days
CS16124	Babesia Bovis	1.0 mL serum in red top or serum separator tube	7-10 days
CT785	Baermann Baermann test is to detect lungworm larvae in fecal samples. Sample Handling: specimen should be tested within 24 hours of collection.	10 grams feces in Antech provided fecal container	2-4 days
CS14456	Bovine Herpes Virus	1.0 mL serum in red top or serum separator tube	5-10 days

CODE	TEST NAME DESCRIPTION COMPONENTS	SPECIMEN	TAT*
CS16844	Bovine Viral Diarrhea (BVD) ELISA	1.0 mL of serum from precolostral newborn calves or calves older than three months are suitable for testing. Ear notches from animals of all ages may also be tested	7-10 days
CS86551	Bovine Viral Diarrhea (BVD) PCR	2.0 mL EDTA whole blood in lavender top tube, fresh tissue, semen, or milk	7-10 days
CS16145	Caprine Arthritis Encephalitis (CAE)	1.0 mL serum in red top or serum separator tube	7-10 days
CS16425	IgG Bovine Test is species specific for bovine. Methodology is the RID Method.	1.0 mL serum in red top or serum separator tube	7-10 days
CS16430	IgG Llama Test is species specific for llama. Methodology is the RID Method.	1.0 mL serum in red top or serum separator tube	3-7 days
CS16302	Johne's Disease Antibody Whole herd screening.	1.0 mL serum in red top or serum separator tube	5-10 days
CT805	Ova and Parasite Fecal ova and parasites. Samples are appropriately mixed with zinc sulfate solution, centrifuged, followed by flotation and slide evaluation. Note: sample should be evaluated within 24 hours of collection. If a worm has been identified in the sample, separate the worm and place it in a container labeled Worm in black, additionally mark Worm on the TRF being submitted with the sample.	5 grams feces in Antech provided fecal container	1-2 days

TEST INDEX

TEST NAME	CODE	PAGE
Accuplex®	CAC100	20, 79
Accuplex®, KeyScreen®	CKAC100	20
Acetylcholine Receptor Antibody	CS16005	111
Acid Fast Stain - Micro	CM010	87
Acute Diarrhea Panel	CL496	142
Add-on Coccidioidomycosis Screen	CADD04	79
Add-on Distemper/Parvo Vacc Titer	CADD270	111
Add-on <i>Ehrlichia canis</i>	CADD05	79
Add-on Feline Retroviral	CADD07	79
Add-on Fructosamine	CADD260	71
Add-on <i>Giardia</i> (ELISA)	CADD250	86
Add-on Heartworm Antigen	CADD70	79
Add-on <i>Hemobartonella</i>	CADD20	79
Add-on Panleukopenia	CADD280	112
Add-on Post T4	CADD300	71
Add-on Protein Electrophoresis	CADD130	112
Add-on PT/PTT	CADD290	59
Add-on Reticulocytes	CADD140	59
Add-on T4	CADD190	71
Add-on TSH	CADD200	71
Add-on Urinalysis	CADD220	69
Add-on Urine Culture MIC	CADD210	87
Add-on Urine Protein/Creat Ratio	CADD230	69
Adrenal Androgen Panel For Ferrets	CS17116	127
Adult Chem with Lytes, SDMA, CBC, Accuplex®, KeyScreen®	CKA535	20
Adult Chem with Lytes, SDMA, CBC, Accuplex®, T4, KeyScreen®	CKA583	20
Adult Chem with Lytes, SDMA, CBC, Feline Heartworm, FIV, FeLV, KeyScreen®	CKS519	20
Adult Chem with Lytes, SDMA, CBC, FIV, FeLV, KeyScreen®	CKS575	21
Adult Chem with Lytes, SDMA, CBC, Heartworm, KeyScreen®	CKS534	21
Adult Chem with Lytes, SDMA, CBC, Heartworm, T4, KeyScreen®	CKS586	21

TEST NAME	CODE	PAGE
Adult Chem with Lytes, SDMA, CBC, KeyScreen®	CKS535	21
Adult Chem with Lytes, SDMA, CBC, T4, Feline Heartworm, FIV, FeLV, KeyScreen®	CKS591	21
Adult Chem with Lytes, SDMA, CBC, T4, FIV, FeLV, KeyScreen®	CKS597	22
Adult Chem with Lytes, SDMA, CBC, T4, KeyScreen®	CKS583	22
Adult Chem with Lytes, SDMA, CBC, T4, UA, Feline Heartworm, FIV, FeLV, KeyScreen®	CKS589	22
Adult Chem with Lytes, SDMA, CBC, T4, UA, FIV, FeLV, KeyScreen®	CKS619	22
Adult Chem with Lytes, SDMA, CBC, UA, Accuplex®, T4, KeyScreen®	CKA673	23
Adult Chem with Lytes, SDMA, CBC, UA, Feline Heartworm, FIV, FeLV, KeyScreen®	CKS623	23
Adult Chem with Lytes, SDMA, CBC, UA, FIV, FeLV, KeyScreen®	CKS634	23
Adult Chem with Lytes, SDMA, CBC, UA, Heartworm, KeyScreen®	CKS592	23
Adult Chem with Lytes, SDMA, CBC, UA, Heartworm, T4, KeyScreen®	CKS587	24
Adult Chem with Lytes, SDMA, CBC, UA, KeyScreen®	CKS053	24
Adult Chem with SDMA, CBC, Accuplex®, KeyScreen®	CKA600	24
Adult Chem with SDMA, CBC, Accuplex®, T4, KeyScreen®	CKA670	24
Adult Chem with SDMA, CBC, Feline Heartworm, FIV, FeLV, KeyScreen®	CKS622	25
Adult Chem with SDMA, CBC, FIV, FeLV, KeyScreen®	CKS675	25
Adult Chem with SDMA, CBC, Heartworm, KeyScreen®	CKS605	25
Adult Chem with SDMA, CBC, Heartworm, T4, KeyScreen®	CKS685	25
Adult Chem with SDMA, CBC, KeyScreen®	CKS600	25
Adult Chem with SDMA, CBC, T4, Feline Heartworm, FIV, FeLV, KeyScreen®	CKS590	26
Adult chem with SDMA, CBC, T4, FIV, FeLV, KeyScreen®	CKS697	26
Adult Chem with SDMA, CBC, T4, KeyScreen®	CKS670	26
Adult Chem with SDMA, CBC, T4, UA, Feline Heartworm, FIV, FeLV, KeyScreen®	CKS588	26
Adult Chem with SDMA, CBC, T4, UA, FIV, FeLV, KeyScreen®	CKS621	27
Adult Chem with SDMA, CBC, T4, UA, KeyScreen®	CKS672	27
Adult Chem with SDMA, CBC, UA, Accuplex®, KeyScreen®	CKA607	27
Adult Chem with SDMA, CBC, UA, Accuplex®, T4, KeyScreen®	CKA672	27
Adult Chem with SDMA, CBC, UA, Feline Heartworm, FIV, FeLV, KeyScreen®	CKS624	28
Adult Chem with SDMA, CBC, UA, FIV, FeLV, KeyScreen®	CKS631	28
Adult Chem with SDMA, CBC, UA, Heartworm, KeyScreen®	CKS625	28
Adult Chem with SDMA, CBC, UA, Heartworm, T4, KeyScreen®	CKS687	28
Adult Chem with SDMA, CBC, UA, KeyScreen®	CKS607	29
Adult Chem with SDMA, Lytes, CBC, T4, UA, KeyScreen®	CKS673	29
Adult Chem with SDMA, Lytes, CBC, UA, Accuplex®, KeyScreen®	CKA053	29
Adult Wellness Chemistry with SDMA	CSA665	65
Adult Wellness with SDMA (No Heartworm)	CSA600	35

TEST NAME	CODE	PAGE
Adult Wellness with SDMA, Fecal O&P	CSA615	35
Adult Wellness with SDMA, Fecal Combo	CSA655	35
Adult Wellness with SDMA, Fecal Combo with Accuplex®	CAC655	35
Adult Wellness with SDMA, FeLV, FIV, O&P	CSA682	36
Adult Wellness with SDMA, FeLV, FIV, O&P, <i>Giardia</i>	CSA684	36
Adult Wellness with SDMA, Heartworm, UA, O&P, <i>Giardia</i>	CSA651	36
Adult Wellness with SDMA, O&P, UA	CSA647	36
Adult Wellness with SDMA, T4	CSA670	36
Adult Wellness with SDMA, T4, UA	CSA672	37
Adult Wellness with SDMA, UA, Accuplex®	CAC607	37
Adult Wellness with SDMA, UA, O&P, <i>Giardia</i>	CSA649	37
Adult Wellness with SDMA, UA, O&P, <i>Giardia</i> , Accuplex®	CAC651	37
Adult Wellness with SDMA, Urinalysis	CSA607	37
Aerobic Culture and MIC	CM020	87, 150
Aerobic Culture and MIC and Anaerobic Culture	CM040	88, 151
Albumin	CT010	112
Alkaline Phosphatase	CT020	112
Alkaline Phosphatase Isoenzymes	CT215	112
ALT (SGPT)	CT030	112
Amylase	CT040	112
Amylase and Lipase (PSL)	CT050	112
ANA (Antinuclear Antibodies) Titer	CT515	113
Anaerobic Culture	CM030	88, 151
<i>Anaplasma</i> CF, Bovine	CS16035	155
<i>Anaplasma phagocytophilum</i> AB	CS14437	80
<i>Anaplasma phagocytophilum</i> PCR	CS14421	147
<i>Anaplasma phagocytophylium</i> AB	CS16872	80
Anticoagulant Screen	CS16040	113
Anti-Mullerian Hormone (Equine)	CS14320	149
Anti-Mullerian Hormone(C/F)	CS14410	72
<i>Aspergillus</i> AB Avian (AFMP1P ELISA)	CS16011	127
<i>Aspergillus</i> Ag, Galactomannan Assay	CS85358	127
AST (SGOT)	CT060	113
Autoimmune Profile 1	CSA170	113
Avian Comprehensive Profile with Bile Acid	CAE021	127
Avian Comprehensive Profile with Bile Acid & EPH	CAE025	127

TEST NAME	CODE	PAGE
Avian Protein Electrophoresis	CAE300	127
Avian Standard Profile with Bile Acid	CAE051	128
Avian Standard Profile with Bile Acid & EPH	CAE055	128
Avian Zinc Assay	CS16012	128
Avian/Exotic CBC	CAE270	128
B-12	CT840	113
B12, Folate	CS16195	113
<i>Babesia Bovis</i>	CS16124	155
<i>Babesia canis</i>	CS16070	111
<i>Babesia gibsoni</i>	CS16075	111
<i>Babesia gibsoni</i> IFA, Export	CS16502	111
Baermann	CT785	155
<i>Bartonella henselae</i> , ELISA	CS85889	111
Basic Wellness Screen with SDMA	CSA710	38
Bicarbonate	CT115	113
Bile Acid Post	CT223	65
Bile Acid Profile	CT220	65
Bile Acids	CT225	65, 139
Bile Acids (Avian/Exotic)	CAE260	128
Bilirubin, Direct	CT070	114
Bilirubin, Total	CT090	114
Bladder Tumor Analytes	CT520	114
Blastomyces Antibody	CT525	114
Blastomyces Quantitative Ag Assay	CS86293	114
Blood Culture - Aerobic Only	CM060	88, 151
Blood Culture - Aerobic/Anaerobic	CM061	89
Blood Type, Canine - DEA 1.1	CT315	59
Blood Type, Canine - Full Panel	CS16100	60
Bone Marrow Core Biopsy	CBMBC	93
Bone Marrow Cytology	CBONE	93
Bovine Herpes Virus	CS14456	155
Bovine Viral Diarrhea (BVD) ELISA	CS16844	156
Bovine Viral Diarrhea (BVD) PCR	CS86551	156
Brochoalveolar Lavage Cytology	CBAL	153
Bromide	CT730	125
<i>Brucella canis</i> - Screen	CT530	80

TEST NAME	CODE	PAGE
<i>Brucella</i> Screen Multiplex	CS16131	80
<i>Brucella</i> Titer For Export (KSU)	CS16003	114
Buffy Coat for Mast Cells	CT325	93
BUN/Creatinine	CT105	114
C Reactive Protein CRP	CS86344	114
CADET® <i>BRAF</i>	CT1025	99
Calcium	CT110	66
Calici Virus - IFA	CS16135	80
Canine Adult Wellness with SDMA, Heartworm, <i>Ehrlichia</i> & Lyme IgG	CSA620	38
Canine Adult Wellness with SDMA, Heartworm, Fecal Combo	CSA650	38
Canine Adult Wellness with SDMA, Heartworm, Ova & Parasite	CSA635	38
Canine Adult Wellness with SDMA, T4	CSA685	39
Canine Adult Wellness with SDMA, UA, O&P	CSA680	39
Canine Adult Wellness with SDMA, UA, O&P, Vaccine Titer	CSA645	39
Canine Adult Wellness with SDMA, Urinalysis	CSA625	38
Canine <i>Ehrlichia</i> / <i>Anaplasma</i> PCR	CT980	99
Canine GI PCR Panel	CT950	100
Canine GI PCR with SARS CoV-2 PCR Panel	CT953	100
Canine Maldigestion Profile	CSA160	115
Canine Melanoma Diagnostic Panel	CS86793	93
Canine Melanoma Prognostic Panel	CS86792	94
Canine Parvovirus PCR	CS8710	100
Canine Respiratory PCR Panel	CT995	101
Canine Respiratory PCR with SARS CoV-2 PCR Panel	CT998	101
Canine Senior Profile with SDMA, Fecal Combo	CSA765	39
Canine Tick Borne PCR Panel	CT960	101
Canine Tick Borne PCR with Lyme	CT961	102
Canine Wellness Profile with SDMA	CSA605	40
Canine Wisdom Panel™ Premium	CS14497	109
Canine Wisdom Panel™ Premium, Add-on	CS14515	110
Canine Wisdom Panel™ with Keyscreen®	CKS14497	110
Canine/Feline Ringworm PCR Panel	CT982	102
Canine/Feline Ringworm PCR with Dermatophyte Culture	CT986	102
Caprine Arthritis Encephalitis (CAE)	CS16145	156
Cardio BNP-Canine	CT1010	115
Cat Scan Plus with SDMA	CSA220	40

TEST NAME	CODE	PAGE
Cat Scan with SDMA	CSA230	40
CBC and Fibrinogen	CT332	137, 155
CBC, Path Review	CT331	60
CBC/Reticulocyte	CT337	60
Chemistry Panel	CSA324	66
Chemistry Renal with SDMA Profile	CT7008	66
<i>Chlamydia</i> Antibody Titer IFA	CS16671	128
<i>Chlamydia</i> PCR Blood	CS16788	128
<i>Chlamydia</i> PCR Swab	CS16672	128
<i>Chlamydia</i> Profile	CS85206	129
<i>Chlamydia</i> Titer EBA	CS16670	129
Chloride	CT120	115
Cholesterol	CT125	115
Cholesterol, Triglycerides	CT127	115
Cholinesterase-Serum	CT235	115
Chronic Diarrhea Panel	CL420	142
<i>Clostridium difficile</i> toxins A/B	CL86181	142
<i>Clostridium perfringens</i> enterotoxin	CT16007	142
Coagulation Profile 1	CSA290	60
Coagulation Profile 2	CSA300	61
Coagulation Profile 3	CSA305	61
Cobalimine, Folate, TLI (Feline)	CSA275	115
Coccidioidomycosis, Screen and Titer	CT535	80
Complete Blood Count	CT330	61
Comprehensive Avian Chemistries	CAE010	129
Comprehensive Avian Post Purchase	CAE030	129
Comprehensive Avian Profile	CAE020	129
Comprehensive Mammalian Profile	CAE200	129
Comprehensive Reptilian Profile	CAE160	130
Comprehensive Reptilian Reckech Panel	CRECHECKREP	130
Copper Level	CS16210	116
Copper Storage Disease	CS16215	116
Corona Virus IgG and IgM	S16225	116
Cortisol	CT445	72
Cortisol Serial ACTH 2	CACTH2	72
Cortisol Serial ACTH 3	CACTH3	72

TEST NAME	CODE	PAGE
Cortisol Serial ACTH 4	CACTH4	72
Cortisol Serial DEX 2	CDEX2	73
Cortisol Serial DEX 3	CDEX3	73
Cortisol Serial DEX 4	CDEX4	73
Cortisol/Creatinine Ratio	CT770	73
CPK	CT130	116
CPK Isoenzymes	CS7592	116
Creatinine	CT135	116
Crossmatch	CT340	117
Crossmatch Additional Donor 2	CT345	117
Crypto, <i>Giardia</i> , <i>Clostridium</i> Enterotoxin	CSA350	117
<i>Cryptococcus</i> Antigen	CT550	81
Crystallographic Stone Analysis	CS16735	69
Culture (Aerobic & Fungal)	CM050	90
Culture Dermatophytes	CDTM	90
Culture ID	CM070	91
Culture, Acid Fast Bacilli	CM100	92
Culture, <i>Campylobacter</i>	CM225	91
Culture, Enhanced Fluid Aerobic	CM220	151
Culture, Feces	CM160	142
Culture, Fungal	CM080	91, 152
Culture, <i>Mycoplasma</i>	CM110	91
Culture, <i>Salmonella</i>	CM121	142
Culture, Urine	CM130	92, 152
Cyclosporine	CS18702	125
Cytology	CCYTO	94, 152
D-Dimer	CT350	117
Dermatopathology Consult	CDERM	94
Diabetes Monitoring Panel with SDMA	CSA800	40
Diarrhea Profile (Avian)	CAE070	130
Digoxin	CT735	125
Direct Coombs' Test (Warm)	CT540	61
Distemper PCR	CS16501	130
Distemper Smear, FA	CS16250	118
Distemper Titer-Virus Neutraliz.	CS16107	130
Distemper Vaccinal Titer	CT560	118

TEST NAME	CODE	PAGE
Distemper/Parvo End Point	CT675	118
Distemper/Parvo Vaccine Titer	CT565	118
Eastern Encephalitis Virus PCR	CS14484	131
<i>Ehrlichia canis</i>	CT570	81
Ehrlichiosis Serology Panel (Canine)	CS16900	81
Electrolyte Screen	CT140	66
<i>Encephalitozoon cuniculi</i> IgG Ab	CS16877	131
Endogenous ACTH	CT435	74
Endogenous ACTH (Baseline) Equine	CL525	144
Endogenous ACTH, Insulin, Glucose	CL540	144
Endogenous ACTH, Insulin, Glucose, Leptin (Equine)	CL575	144
EPM SAG 2,4/3 ELISA	CS14388	145
EPM SAG 2,4/3 ELISA Serum/CSF Ratio	CS14390	145
EPM SAG 2,4/3 ELISA, <i>N. hughesi</i> ELISA	CS14392	145
Equine 11 Metal & Mineral Panel	CS14466	154
Equine Chemistry	CL010	138
Equine Chemistry, CBC	CL050	138
Equine Chemistry, CBC, Fib	CL040	138
Equine Chemistry, CBC, Fib (Recheck)	CL040R	138
Equine Chemistry, CBC, Fib, fT4 (ED)	CL035	138
Equine Chemistry, CBC, Fib, T4	CL030	139
Equine Chemistry, CBC, Fib, T4, fT4 (ED)	CL565	138
Equine Comprehensive Neurological Panel	CS14477	146
Equine Coronavirus (ECoV) PCR	CS14414	147
Equine Drug Screen	CS1680	141
Equine Encephalitis Viral Panel	CS16275	146
Equine Encephalitis Viral Panel Plus	CS17500	146
Equine Granulosa Cell Tumor	CS85857	149
Equine Health Program GVL (ELISA)	CL511E	136
Equine Health Program (no EIA)	CL510NE	136
Equine Hepatic Panel	CL230	139
Equine Hepatic Profile	CL1025	139
Equine Hepatic Screen	CL225	139
Equine Herpes Virus 1 (EHV1) PCR	CS14394	147
Equine IgG Total	CL090	139
Equine Infectious Anemia (EIA) ELISA by GVL	CL121	137

TEST NAME	CODE	PAGE
Equine Inflammatory Profile	CL070	137
Equine Insulin and Glucose	CT470E	144
Equine <i>Leptospira</i> PCR Aqueous or Vitreous Humor	CL973	147
Equine <i>Leptospira</i> PCR Blood	CL974	147
Equine <i>Leptospira</i> PCR Blood/Urine	CL978	147
Equine <i>Leptospira</i> PCR Urine	CL976	147
Equine Muscle Enzyme Screen	CL275	140
Equine Muscle Enzyme Screen, CBC, FIB	CL280	140
Equine PCR Fever of Unknown Origin (FUO) Mini	CL970	147
Equine PCR Fever of Unknown Origin (FUO) Panel	CL965	147
Equine PCR Gastrointestinal (Diarrhea) Panel	CL950	146
Equine PCR Reproductive Panel	CL955	146
Equine PCR Respiratory Panel	CL960	146
Equine Performance Horse Chemistry	CL601	138
Equine Performance Horse Health Program (no EIA)	CL611NE	136
Equine Performance Horse Health Program GVL (ELISA)	CL611E	136
Equine Performance Horse Panel	CL640	138
Equine Renal Screen	CL240	140
Equine Rhabdomyolysis Screen	CL190	140
Equine Rhinopneumonitis virus Ab (EHV-1)	CS7591	146
Equine Ringworm PCR Panel	CT983	148
Equine Ringworm PCR with Dermatophyte Culture	CT987	148
Equine Senior Health Program (no EIA)	CL615NE	136
Equine Senior Health Program GVL (ELISA)	CL616E	136
Equine Senior Health Program GVL with TRH	CL616T	136
Equine TRH Stimulation (Thyroid Function)	CL590	143
Equine Viral Arteritis (EVA)	CS16285	154
Estradiol	CS16295	149
Estrone Sulfate	CS16300	149, 150
Feather Picker Profile	CAE080	131
FEC: McMaster's Method (MM) - Equine	CT828	142
FEC:Modified Stoll's Technique (MST)	CT826	141
Fecal Culture	CM125	92
Feline Adult Wellness Profile with SDMA	CSA610	40
Feline Adult Wellness with SDMA, FeLV, FIV	CSA622	41
Feline Adult Wellness with SDMA, Heartworm, Fecal Combo	CSA660	41

TEST NAME	CODE	PAGE
Feline Adult Wellness with SDMA, UA, O&P, Panleukopenia	CSA640	41
Feline Adult Wellness with SDMA, Urinalysis	CSA630	41
Feline Autoimmune Profile	CSA280	118
Feline Comprehensive Plus with SDMA, FT4ED	CSA205	41
Feline Comprehensive Wellness Screen with SDMA	CSA715	42
Feline Comprehensive with SDMA	CSA210	42
Feline Coronavirus Titer	CT595	81
Feline Flea and Tick Borne PCR	CT965	102
Feline GI PCR Panel	CT955	103
Feline GI PCR with SARS CoV-2 PCR Panel	CT958	103
Feline Heartworm Profile	CT630	81
Feline Heartworm Program	CSA250	66
Feline Heartworm Program Plus	CSA240	42
Feline Hemoplasma PCR Panel	CT985	103
Feline Respiratory PCR Panel	CT990	104
Feline Respiratory PCR with SARS CoV-2 PCR Panel	CT993	104
Feline Retroviral	CSA260	81
Feline Serology	CSA262	82
Feline Serology 1	CSA265	82
Feline Serology 2	CSA270	82
Feline Total Health Check with SDMA	CSA190	42
Feline Total Health Check with SDMA, UA	CSA786	42
Feline Total Health Plus with SDMA	CSA180	43
Feline Wellness Profile with SDMA	CSA700	43
Feline Wellness with SDMA, FeLV, FIV, O&P	CSA683	43
Feline Wisdom Panel™ Complete	CS14498	109
Feline Wisdom Panel™ Complete, Add-on	CS14516	110
FeLV - IFA	CT585	82
FeLV Antigen (ELISA)	CT580	82
FeLV Antigen (ELISA) Add-on	CADD06	82
FeLV PCR	CS6234	83
FeLV, FIV, KeyScreen®	CKS260	34
Fibrinogen	CT365	62, 155
Fibrinogen - Equine	CL365	137
FIP mRNA PCR	CT600	104
FIV Antibody	CT610	83

TEST NAME	CODE	PAGE
FIV Antibody Add-on	CADD15	83
FIV Antibody, Western Blot	CS16865	118
Fluid Analysis	CFLUA	94, 152
Fluid Analysis, CSF	CCSF	152
Fluid Analysis, Synovial	CSYFLUA	153
Foal Diarrhea Panel	CL492	142
Fractional Excretion of Electrolytes	CL340	141
Free T4 by Equilibrium Dialysis	CT460	74, 143
Free T4 Equilibrium Dialysis Add-on	CADD50	74
Fructosamine	CS16345	74
Fungal Serology with Cocci	CSA340	119
General Senior Profile with SDMA	CSA705	43
Geriatric/Weak Ferret Profile	CAE230	131
GGTP	CT145	66
<i>Giardia</i> (ELISA)	CT820	131
<i>Giardia</i> and Ova & Parasite with Centrifugation	CT808	86
<i>Giardia</i> ELISA/Crypto FA	CT790	119
Glucose	CT150	119
Gram Stain	CM090	92
Heartworm (Microfilaria) for Export	CS86454	83
Heartworm Ag-Heat-Treated Serum	CT613	83
Heartworm Antibody, Feline	CT625	83
Heartworm Antigen	CT615	83
Heartworm Antigen Feline	CT620	84
Heartworm Program 1	CSA140	119
Heartworm Program 2	CSA130	119
Heartworm, KeyScreen®	CKT615	34
Heartworm, Ova and Parasite	CT618	84
<i>Hemobartonella</i>	CT380	84
Hepatic Profile	CAE090	132
Herpes - Conjunctival Smear	CS16400	119
Herpes Antibody, IFA	CS86022	119
Heska ALLERCEPT Environmental & Food Panel	CS17029	120
Heska ALLERCEPT Environmental Panel	CS17026	120
Heska ALLERCEPT Food Panel	CS17028	120
Heska Equine ALLERCEPT Panel	CS17027	120

TEST NAME	CODE	PAGE
Histopathology	CFBX	95, 153
<i>Histoplasma capsulatum</i> Ag, EIA	CS86569	95
<i>Histoplasma</i> Titer	CS16405	120
Hyperthyroid Monitoring Profile with SDMA	CSA440	43
Hyperthyroid Panel with SDMA	CSA805	44
IgG Bovine	CS16425	156
IgG Llama	CS16430	156
Immunoglobulins A, G, M	CT660	120
Immunohistochemistry 1 Stain	IHC1	95
Immunohistochemistry 2 Stains	IHC2	95
Immunohistochemistry 3 Stains	IHC3	96
Immunohistochemistry 4 Stains	IHC4	96
Immunohistochemistry 5 Stains	IHC5	96
Influenza - Canine Acute	CS86096	84
Influenza - Canine Convalescent	CS86102	120
Inhibin	CS4131	149
Insulin	CT470	74
Insulin (Feline) with Glucose	CT470F	75
Ionized Calcium	CS18537	67
Iron Serum	CT155	120
Johne's Disease Antibody	CS16302	156
Keppra Level	CS86541	125
KeyScreen® GI Parasite PCR Add-on	CADD350	104
KeyScreen® GI Parasite PCR Panel	CT991	105
Lactic Dehydrogenase	CT160	120
Leptin (Equine)	CS14402	144
<i>Leptospira canicola</i> , Export	CS5017	121
<i>Leptospira</i> PCR Blood	CT974	106
<i>Leptospira</i> PCR Blood/Urine	CT978	106
<i>Leptospira</i> PCR Urine	CT976	106
Leptospirosis	CS16510	84, 154
Liver Chemistry Screen	CSA321	67
Liver Profile	CSA320	44
Luteinizing Hormone	CS16520	75
Lyme Equine Multiplex	CT672	154
Lyme IgG	CT670	84

TEST NAME	CODE	PAGE
Lymphocyte Phenotyping	CS86493	62
Magnesium	CT170	121
Mammalian Comprehensive Chems	CAE190	132
Mammalian Standard Chemistries	CAE210	132
Masticatory Muscle Myositis	CS16535	121
Microfilaria	CT390	84
Mini Avian Post Purchase	CAE060	132
Miniscreen 4 Chem, CBC	CSA070	44
<i>Mycoplasma</i> PCR	CS16789	132
Neonatal Foal Panel	CL290	139
<i>Neorickettsia risticii</i> (PHF) PCR	CS14479	147
<i>Neospora caninum</i> - IFA	CS16560	85
NSAID 1 with SDMA	CSA810	44
NSAID 2 with SDMA	CSA815	44
NSAID 3 with SDMA	CSA820	45
O&P, <i>Giardia</i> with Heartworm	CT812	85
Occult Blood, Feces	CT810	121
OncoK9®	CS14493	106
Oral Path Biopsy	CFBXORL	97
Oral Sugar Test (2 samples)	CL545	145
Oral Sugar Test (3 samples)	CL550	145
Osmolality - Serum	CS16575	121
Ova & Parasite	CT805	132, 156
Pancreatic Lipase Immunoreactivity	CS85364	121
Pancreatitis Profile - Canine	CS86468	121
Pancreatitis Profile - Feline	CS86288	122
Panleukopenia Titer IgG, IgM	CS16580	122
Panleukopenia Vaccinal Titer	CS16053	122
Parasite Identification	CT825	122
Parathormone Related Protein	CS16596	75
PARR	CS85503	107
Parvovirus Antibody/Antigen	CT700	122
Parvovirus Antigen	CT695	85
Parvovirus Vaccinal Titer	CT705	123
Path Review	CREVW	63
Path Review (Urinalysis)	CT764	69

TEST NAME	CODE	PAGE
Pet Cancer Specialty Biopsy	CFBXNEO	97
Phenobarbital	CT750	126
Phenobarbital Panel Plus	CSA830	126
Phosphorus	CT180	123
Piroplasmosis, Equine	CS16848	154
Pituitary Pars Intermedia Dysfunction (PPID) Panel	CL500	143
Pituitary Pars Intermedia Dysfunction (PPID) with fT4ED	CL560	143
Pituitary Pars Intermedia Dysfunction (PPID) with TRH Stim	CL500TRH	143
Platelet Count	CT400	63
Post Cortisol	CT448	75
Post T4	CT497	75
Potassium	CT185	123
Pre & Post Phenobarbital	CT755	126
PrecisionPSL	CT165	67
Pregnant Mare Serum Gonadotropin (PMSG)	CS16635	149
Pre-Op Chem with Electrolytes, SDMA, CBC	CSA516	45
Pre-Op Chem with Electrolytes, SDMA, CBC, Urinalysis	CSA053	45
Pre-Op Panel Plus with SDMA	CSA050	45
Pre-Op Panel Plus with SDMA, Heartworm	CSA512	45
Pre-Op Panel with SDMA, T4, UA	CSA920	45
Pre-Op Screen with SDMA	CSA040	67
Pre-Op Screen with SDMA, CBC	CSA055	46
Pre-Op Screen with SDMA, CBC, Heartworm	CSA056	46
Pre-Op Screen with SDMA, CBC, PT, PTT, FeLV, FIV	CSA510	46
Pre-Op Screen with SDMA, CBC, T4, FeLV, FIV	CSA514	46
Pre-Op Screen with SDMA, CBC, Urinalysis	CSA052	46
Pre-Op with SDMA, CBC with Accuplex®	CAC055	46
Pre-Op with SDMA, CBC, FeLV, FIV	CSA051	47
Pre-Op with SDMA, CBC, T4	CSA508	47
Progesterone	CT475	76
Progesterone - Equine	CL140	149
Progesterone/Estradiol Panel (Equine)	CL470	149
Progesterone/PMSG	CL460	149
Progesterone/PMSG/Estradiol	CL480	149
Protein Electrophoresis	CT240	67, 140
Prothrombin Time	CT410	63

TEST NAME	CODE	PAGE
Psittacine Beak and Feather Disease	CS16085	132
PT/PTT	CT415	64, 137
PTH/Ionized Calcium	CS16595	76
PTT	CT395	64
PU/PD Profile	CAE110	133
Rabbit Neurologic Profile	CAE240	133
Rabbit/Small Rodent Ringworm PCR Panel	CT984	108
Rabbit/Small Rodent Ringworm PCR with Dermatophyte Culture	CT988	108
Rabies Diagnostic Non-Export	CS1204	85
Rabies Export - FAVN	CS17108	85
Relaxin	CT9810	123
Renal Chemistry Screen with SDMA	CSA311	68
Renal Profile with SDMA	CSA310	47
Reptilian Comprehensive Chems	CAE150	133
Reptilian Standard Chemistries	CAE170	133
Reticulocytes	CT425	133
Rhinotrachitis Feline	CS16702	123
<i>Rhodococcus equi</i> PCR	CS14396	148
Ruminant Chemistry	CL1001	155
Ruminant Chemistry/CBC/Fib	CL080	155
<i>Salmonella</i> spp. PCR	CS14416	148
SARS-CoV-2 PCR Add-on Panel	CT996	108
SDMA	CT1035	68
SearchLight DNA® Cytology Add-on	CS10058	98, 109
SearchLight DNA® Histology Add-on	CS11068	98, 109
Selenium Level	CS16730	123
Semen Analysis	CSEMEN	150
Senior Comprehensive Plus with SDMA	CSA080	47
Senior Comprehensive Profile with SDMA	CSA090	48
Senior Comprehensive with SDMA, T3	CSA105	48
Senior Feline Plus with SDMA	CSA740	48
Senior Profile 1 with SDMA, Accuplex®	CAC705	49
Senior Profile 1 with SDMA, Accuplex®, KeyScreen®	CKA705	29
Senior Profile 1 with SDMA, Fecal Combo	CSA770	49
Senior Profile 1 with SDMA, Fecal Combo, Accuplex®	CAC770	49
Senior Profile 1 with SDMA, Heartworm, KeyScreen®	CKS710	30

TEST NAME	CODE	PAGE
Senior Profile 1 with SDMA, Heartworm, Ova & Parasite	CSA745	49
Senior Profile 1 with SDMA, KeyScreen®	CKS705	30
Senior Profile 1 with SDMA, Ova & Parasite	CSA750	50
Senior Profile 1 with SDMA, UMIC	CSA566	48
Senior Profile 1 with SDMA, Vaccine Titers	CSA730	50
Senior Profile 2 with SDMA	CSA720	50
Senior Profile 2 with SDMA, Fecal Combo	CSA775	51
Senior Profile 2 with SDMA, Feline Heartworm Ab, O&P	CSA790	51
Senior Profile 2 with SDMA, FeLV, FIV	CSA780	50
Senior Profile 2 with SDMA, Ova & Parasite	CSA755	51
Senior Profile 2 with SDMA, UMIC	CSA746	51
Senior Profile 2 with SDMA, Vaccine Titers	CSA760	52
Senior Profile with SDMA, Feline Heartworm Ab, Fecal Combo	CSA768	52
Serial Blood Culture X2	CM062	89
Serial Blood Cultures X3	CM063	90
Sodium	CT195	123
Sodium and Potassium	CT200	124
Sorbtal Dehydrogenase (SDH)	CT250	139
Standard Avian Profile	CAE050	133
Standard Mammalian Profile	CAE220	134
Standard Reptilian Profile	CAE180	134
<i>Streptococcus</i> Equi ELISA	CS16740	154
<i>Streptococcus equi</i> PCR	CS86308	148
Superchem with SDMA	CSA010	68
Superchem with SDMA, CBC	CSA020	52
Superchem with SDMA, CBC Recheck Panel	CRECHECK	52
Superchem with SDMA, CBC, Accuplex®, KeyScreen®	CKA020	30
Superchem with SDMA, CBC, Accuplex®, T4, KeyScreen®	CKA114	30
Superchem with SDMA, CBC, Feline Heartworm, FIV, FeLV, KeyScreen®	CKS737	31
Superchem with SDMA, CBC, Feline Serology 1	CSA460	53
Superchem with SDMA, CBC, FeLV, FIV	CSA490	53
Superchem with SDMA, CBC, FIV, FeLV, KeyScreen®	CKS490	31
Superchem with SDMA, CBC, Heartworm, KeyScreen®	CKS028	31
Superchem with SDMA, CBC, Heartworm, T4 Free by Immulite	CSA111	47
Superchem with SDMA, CBC, Heartworm, T4, KeyScreen®	CKS110	31
Superchem with SDMA, CBC, KeyScreen®	CKS020	31

TEST NAME	CODE	PAGE
Superchem with SDMA, CBC, Lyme	CSA500	53
Superchem with SDMA, CBC, O&P, <i>Giardia</i>	CSA955	52
Superchem with SDMA, CBC, Special Thyroid Profile	CSA450	53
Superchem with SDMA, CBC, T4, Feline Heartworm, FIV, FeLV, KeyScreen®	CKS724	32
Superchem with SDMA, CBC, T4, FIV, FeLV, KeyScreen®	CKS220	32
Superchem with SDMA, CBC, T4, KeyScreen®	CKS120	32
Superchem with SDMA, CBC, T4, UA, Feline Heartworm, FIV, FeLV, KeyScreen®	CKS715	32
Superchem with SDMA, CBC, T4, UA, FIV, FeLV, KeyScreen®	CKS700	32
Superchem with SDMA, CBC, UA	CSA021	53
Superchem with SDMA, CBC, UA, Accuplex®, KeyScreen®	CKA021	33
Superchem with SDMA, CBC, UA, Feline Heartworm, FIV, FeLV, KeyScreen®	CKS494	33
Superchem with SDMA, CBC, UA, FIV, FeLV, KeyScreen®	CKS702	33
Superchem with SDMA, CBC, UA, Heartworm, KeyScreen®	CKS112	34
Superchem with SDMA, CBC, UA, KeyScreen®	CKS021	34
Superliver (Liver Biopsy with Liver Staining Panel)	CFBXLSP	97
T3	CT480	76
T3 Autoantibodies	CT485	124
T4	CT495	76, 143
T4 (Free)	CT499	76
T4 Autoantibodies	CT500	124
Taurine	CS16755	124
Testosterone	CS16760	77, 150
Testosterone Panel	CS85530	150
Thyroglobulin Autoantibodies	CT505	77
Thyroid Profile 1	CSA360	77
Thyroid Profile 2	CSA370	77, 143
Thyroid Profile 3	CSA380	77
Thyroid Profile 4	CSA390	78
Thyroid Profile 5	CSA400	78
Thyroid Profile 6	CSA410	78
Tick Serology 1	CSA330	124
TLI, Feline	CS16800	124
Total Body Function Plus with SDMA	CSA110	53
Total Body Function Plus with SDMA, O&P, <i>Giardia</i>	CSA788	54
Total Body Function with SDMA	CSA120	54
Total Body Function with SDMA, Accuplex®	CAC120	54

TEST NAME	CODE	PAGE
Total Body Function with SDMA, FeLV, FIV, Feline Heartworm Ab	CSA724	54
Total Body Function with SDMA, T4 Immulite	CSA120B	54
Total Body Function with SDMA, T4 Immulite, Accuplex®	CAC120B	55
Total Protein	CT190	124
<i>Toxoplasma</i> AB - IgG/IgM	CT720	85
<i>Toxoplasma</i> antibody - Exotics	CS16792	134
<i>Toxoplasma gondii</i> PCR	CS18708	86
Toxoplasmosis IgG/IgM	CS85030	86
Transtracheal Wash Cytology	CTTW	153
TRH Stimulation for PPID	CL535	144
TRH Stimulation for PPID with Insulin and Glucose	CL580	145
Triglycerides	CT205	124
<i>Tritrichomonas</i> PCR	CS85819	86
Troponin I	CS85783	68
Trypsin-Like Immunoreac Canine cTLI	CT230	124
TSH	CT510	78
Tumor Board Review Biopsy	CFBXTBR	98
Ur Bile Acid (UBA):Ur Creat (UCR)	CT227	125
Urea Nitrogen	CT100	125
Urinalysis	CT760	69, 141
Urinalysis and UPC Ratio	CT925	70
Urinalysis Clearance Ratio	CT765	70
Urinalysis/Urine Culture with MIC	CM133	70
Urine Microalbumin Canine	CT830C	70
Urine Microalbumin Feline	CT830F	70
Urine Protein/Creatinine Ratio	CT775	70
Vaccine Panel (Panleukopenia, Rhino, Calici)	CS16581	86
Vet Screen 2 with SDMA	CSA026	55
Vet Screen 2 with SDMA, CBC	CSA035	55
Vet Screen with SDMA	CSA025	55
Vet Screen with SDMA, CBC	CSA030	57
Vet Screen with SDMA, CBC Recheck Panel	CRECHECK2	57
Vet Screen with SDMA, CBC, FeLV, FIV, FIP	CSA914	55
Vet Screen with SDMA, CBC, FeLV, FIV, FIP, T4	CSA908	56
Vet Screen with SDMA, CBC, FeLV, FIV, T4	CSA910	56
Vet Screen with SDMA, CBC, FeLV, FIV, UA	CSA906	56

TEST NAME	CODE	PAGE
Vet Screen with SDMA, CBC, Heartworm Ag	CSA039	57
Vet Screen with SDMA, CBC, T4	CSA034	57
Vet Screen with SDMA, CBC, T4 (Free)	CSA032	57
Vet Screen with SDMA, CBC, T4 (Free), UA	CSA038	57
Vet Screen with SDMA, CBC, T4, FT4ED	CSA912	56
Vet Screen with SDMA, CBC, T4, Heartworm	CSA530	58
Vet Screen with SDMA, CBC, UA	CSA036	58
Vet Screen with SDMA, CBC, UA, O&P, <i>Giardia</i>	CSA922	58
Vet Screen with SDMA, CBC, UA, O&P, Heartworm	CSA904	56
Vet Screen with SDMA, CBC, UA, T4	CSA037	58
Vitamin D	CS16016	68
Vitamin E	CS16850	154
Vitamin E/Selenium	CS17505	154
Von Willebrand Factor	CS17123	64
West Nile Titer (PRNT & IgM ELISA)	CS85448	134, 146
West Nile Virus PCR	CS85449	134
West Nile Virus Plaque Reduction SN	CS14485	134
Zinc	CS16870	125
Zonisamide Zonegran	CS86480	126

**BETTER DIAGNOSTICS.
BETTER MEDICINE.**



ANTECH[®]

THE

RIGHT DIAGNOSTICS

AT THE RIGHT TIME

We are proud to deliver the most comprehensive portfolio of veterinary products and services to you and your patients, now all offered by Antech.



ANTECH

1-800-341-3440



HESKA

1-866-382-6937



sound

1-800-268-5354



AIS

1-877-727-6800

EXPERIENCE A COMPLETE DIAGNOSTIC SOLUTION

Antech partners with veterinary care teams by providing a flexible, complete diagnostic portfolio that empowers you to practice medicine your way — because better diagnostics, means better medicine.

LEARN HOW TODAY



[ANTECHDIAGNOSTICS.COM/LEARN](https://AntechDiagnostics.com/Learn)

© 2023 Antech Diagnostics, Inc. All rights reserved. Antech, the Antech logo, and all other trademarks used herein are the registered trademarks of Antech Diagnostics, Inc. or its affiliates. 113023